**ZuZi**

SERIE SERIES SÉRIE 122

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Thank you for choosing this equipment. We sincerely wish that you enjoy your professional microscope Zuzi series 122. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user. The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

**TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

Please bear in mind the following:

- ◆ This manual is inseparable from the professional microscope Zuzi series 122, so it should be available for all the users of this equipment.
- ◆ Establish all actions required according current regulations in subject of health and safety at work (electrical safety devices, flammable, toxic and / or pathological products, etc).
- ◆ You should carefully handle the microscope avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the microscope to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the microscope.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Zuzi Technical Assistance Department by email to [asistencia@auxilab.com](mailto:asistencia@auxilab.com).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the microscope warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the microscope warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

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Reference: 90122430      Composed of: Filtro analizador y filtro polarizador

**Phase contrast equipments**

Reference	Description	Phase contrast objectives	For models
90122022	Simple model	Achromatics: 4x, 10x, 40x, 100x	122/107 and 122/108
90122122	urret model	Plan achromatics: 4x, 10x, 40x, 100x	122/147 and 122/148

**Dark field Condenser**

Reference	Description	Numerical aperture
90122024	Dry	0.90
90122026	Immersion	1.25

**Eyepiece video cameras**

Code	Resolution	Output
59140052	640x480 pixels, 0.3 Mp	Video composite
59140060	640x480 pixels, 0.3 Mp	USB 2.0
59150060	1280x1024 pixels, 1.3 Mp	USB 2.0
59150063	2048x1536 pixels, 3 Mp	USB 2.0
59150065	2592x1944, 5 Mp	USB 2.0
59150068	3264x2884, 8 Mp	USB 2.0

**Putting the eyepiece video camera**

In Zuzi microscope model 122/107 (binocular), take away one of the eyepieces and put the video camera.

In Zuzi microscope model 122/108 and 122/148 (trinoculars) insert the video camera into the vertical tube and pull the sliding rod to open the path of the light beam to the vertical tube.

A wide range of eyepieces of different magnifications and features (wide field, plan, micrometric, etc.) are available. Consult your wholesaler or the web page [www.auxilab.com](http://www.auxilab.com)

**Objectives**

There is also a wide range of objectives of different magnifications: achromatic, semi-plan, plan, etc. Consult your wholesaler or the web page [www.auxilab.com](http://www.auxilab.com)

**INSTRUCTIONS ON ENVIRONMENTAL PROTECTION**

Do not dispose of this product in the usual household garbage at the end of its life cycle; hand it over at a collection point for the recycling of electrical and electronic appliances. The symbol on the product, the instructions for use or the packing will inform about the methods for disposal. The materials are recyclable as mentioned in its marking. By recycling, material recycling or other forms or re-utilization of old appliances, you are making an important contribution to protect our environment. Please inquire at the community administration for the authorized disposal location.

**ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

**INDEX OF LANGUAGES**

Spanish ..... 2-19  
 English ..... 20-36  
 Francais ..... 37-54

**INDEX OF CONTENTS**

1. USES OF THE INSTRUMENT ..... 21  
 2. DESCRIPTION ..... 22  
 3. TECHNICAL SPECIFICATIONS ..... 23  
 4. INSTALLATION / SETTING UP ..... 26  
 5. MAINTENANCE AND CLEANING ..... 31  
 6. CHOOSING OBJECTIVES AND EYEPIECES ..... 33  
 7. CAUSES OF A DEFECTIVE IMAGE ..... 33  
 8. TROUBLESHOOTING ..... 34  
 9. USEFUL RECOMMENDATIONS ..... 35  
 10. ACCESSORIES ..... 36  
 ANNEX I: CE CERTIFICATE ..... 55

**1. USES OF THE INSTRUMENT**


Zuzi professional microscopes series 122 are very hard equipments designed to provide the highest satisfaction in a wide range of demands. They are equipped with high quality optics and a wide range of accessories to complete the microscope.

Zuzi professional microscopes series 122 are suitable for educational purposes as well as for being used in hospitals, veterinary, etc.

## 2. DESCRIPTION

- |   |  |
|---|--|
| 1.1 Head  | 1.9 Stage clip   |
| 1.2 Eyepieces   | 1.10 Stage movement knobs  |
| 1.3 Dioptre compensation ring                         | 1.11 Fine adjustment knob  |
| 1.4 Head holding screw                                | 1.12 Coarse adjustment knob                                      |
| 1.5 Sliding rod to direct light (122/108 and 122/148) | 1.13 Tube for eyepiece/photographic camera (122/108 and 122/148) |
| 1.6 Nosepiece   | 1.14 On/off switch   |
| 1.7 Objectives  | 1.15 Light intensity control                                     |
| 1.8 Stage   | 1.16 Köhler illumination system                                  |

- ◆ Check the prepared slide is not the other way round (being the slide over the cover glass) and place it properly in case it is appropriate to do so.
- ◆ Check also if the cover glass has the suitable thickness (0.17mm).

## 9. PRACTICAL RECOMMENDATIONS

- ◆ As the head is rotary you should rather observe by the stage's frontal part, as it facilitates access to the mechanical control knobs.

### Use of the immersion objective:

- ◆ To use all the numerical aperture of the immersion objective (100x), add a drop of immersion oil in the space between the objective and the slide. Proceed as follows:
  - Focus the sample with a lower magnifying objective (40x) and put it in the centre of the field of view. This is essential because it is very difficult to find structures with the immersion objective.
  - Turn slightly the nosepiece so as no objective stays over the slide.
  - Add a drop of immersion oil over the region of the sample to be observed.
  - Turn completely the nosepiece and put the immersion objective over the slide (both the sample and the objective lens must be in contact with the immersion oil). Adjust focusing with the fine adjustment knob.
- ◆ Avoid the appearance of bubbles in the oil film, since it spoils the resulting image.
- ◆ After each observation, clean the objective lens with a special lens tissue to remove the rests of immersion oil. Never let the oil dry on the lens.

**NOTE:** never submerge the objective in any liquid nor use solvents such as xylol or toluene for cleaning, since it may provoke the detachment of the lens.

### Use of the coarse and fine adjustment knobs

- ◆ Depending on the observant the tension of the coarse adjustment knob can be modified, since it is provided with a friction control ring that is placed next to the right coarse adjustment knob of the microscope (*figure 8*). Avoid loosening it excessively since this could produce the fall of the stage together with problems on focusing accuracy.

Figura. 8



Figura. 9

- ◆ There is also a device for the control of the upper focus limit located next to the left coarse adjustment knob (*figure 9*). This device avoids that the plate raises beyond a point, that can be determined once the sample has been focussed, and this fact can avoid possible accidental collisions during the observation. This mechanism also allows the possibility of getting a quick focusing, accelerating the observation process when working with different samples.

any diffuser filter.

- ◆ Eyepieces are well fit and objectives well screwed.
- ◆ Eyepieces and objectives are perfectly cleaned and free of particles. To check cleaning of the whole optical system proceed as follows:

- Make the eyepieces turn checking the little specks are moving; if so, clean them. Make the head turn. You should never disassemble the head, but you can clean it with delicacy by blowing the accessible surfaces of prisms with a plastic bulb.

- Turn the objective and, in case the parasite images also turn, clean it with the help of a dry brush to remove the dust. Watch the front part of the microscope with a magnifying glass or an inverted eyepiece.

- ◆ There is enough immersion oil and it does not contain bubbles or impurities.
- ◆ The thickness of microscope slide, microscope cover glass and assembly medium is not too high so as to avoid the focusing at medium or high magnification (there are standard dimensions for both microscope slides and cover glasses). Both microscope slide and cover glass should be clean. Check that the slide is properly placed and that there are not two superimposed cover glasses.

## 8. TROUBLESHOOTING

Before sending the microscope back to Zuzi Technical Assistance Department, check the following:

### *If the lamp does not switch on:*

- ◆ Check it is properly connected to the right current intake.
- ◆ Check both fuse and lamp are in good conditions. If not, replace the fuse or the lamp.

### *If the field of vision is cut:*

- ◆ Check the nosepiece is properly fit. To do that, turn it slightly back and forth until it fits.
- ◆ Check the condenser is well centred. This operation is performed through the centring screws located at the condenser support.

### *If there is dust or dirt in the field of view:*

- ◆ Check there is dust on the precondenser lens, upper lens of the condenser or the eyepiece, as well as the cleanliness of the sample. Once you find out where the dirt is, clean it as explained before.

### *If focus is lost while observation:*

- ◆ Check the friction control ring located at the right coarse focus knob, since it is probably too loose and makes the stage to fall.

### *If an area in the field of view is out of focus:*

- ◆ Check the objective is properly located on the luminous beam's path. If not, turn it until it is in an appropriate position.
- ◆ Check the sample is properly located on the stage, making sure it is held by the stage clips.

*If the frontal lens of the objective touches the prepared slide when focusing or when passing from a low-magnification objective to a superior one:*

## 3. TECHNICAL SPECIFICATIONS

	122/107	122/108	122/147	122/0148
<b>Optics</b>				
Universal infinity corrected			•	•
Standard	•	•		
<b>Head</b>				
Binocular	•		•	
Trinocular		•		•
<b>Eyepieces</b>				
WF/10X (18 mm)	•	•		
WF/10X (20 mm)			•	•
WF/16X (11 mm)	•	•		
<b>Eyepiece diameter</b>				
23 mm	•	•		
30 mm			•	•
<b>Objectives</b>				
Achromatic (antifungal)	•	•		
Plan achromatic (antifungal)			•	•
<b>Nosepiece</b>				
Quadruple				
4X A.N.: 0.10	•	•	•	•
10X A.N.: 0.25	•	•	•	•
40X A.N.: 0.65	•	•	•	•
100X (R) (I) A.N.: 1.25	•	•	•	•
<b>Stage</b>				
	Double layer mechanical stage, 180x150 mm			
	Simple movement 75(X) x 50(Y) mm			
	Clamp to hold samples			
<b>Illumination system</b>				
	Type Köhler. Condenser Abbe (N.A. 1.25) with aperture diaphragm and filter holder.			
	Precondenser With field diaphragm			
	Lamp Halogen, low tension, 20 W			

### **Base:**

Metallic, very stable, with coaxial knobs at both sides for fine and coarse adjustment. The right knob is provided with a friction control ring and the left knob presents an upper limit blocking lever.

### Head

The Zuzi series 122 is provided with different heads depending on the model:

- ◆ Binocular, 30° inclined and rotary 360°. Interpupillary distance adjustment through a free mechanism with graduations between 55 and 75 mm; the dioptre compensation ring is at the left eyepiece holder.
- ◆ Trinocular, 30° inclined and rotary 360°. Interpupillary distance adjustment through a free mechanism with graduations between 55 and 75 mm; the dioptre compensation ring is at the left eyepiece holder. With vertical tube to fit a photographic camera or an eyepiece video camera.

### Stage

- ◆ Movement control through vertical coaxial knobs
- ◆ Dimensions: 180x150 mm
- ◆ Movement range: 75 (X) x 50 (Y) mm
- ◆ Graduated scales at 0.1 mm

### Illumination system

High luminosity Köhler system:

- ◆ Double lens Abbe condenser (N.A.: 1.25) with aperture diaphragm and movable filter holder (32 mm green and blue filters included). Height adjustable by rack and pinion mechanism.
- ◆ Precondenser with field diaphragm
- ◆ Low tension halogen lamp 20 W with intensity control. Built-in power input (220-230 V, 50-60 Hz  $\pm$  10%).

### Optical specifications

- ◆ Achromatic objectives: They correct optical aberrations. This optical system is able to transmit white light without separating constituent colours.
- ◆ Plan Achromatic objectives: They correct optical aberrations and also the spherical aberrations.

Zuzi microscopes series 122 are provided with a quadruple nosepiece with three “dry” achromatic objectives (4x, 10x and 40x) and an immersion objective (100x). 40x and 100x objectives are retractile (R).

The features of objectives are inscribed on the barrel of each objective in the following way (*Figure 2*)

- 40: Objective magnification
- 0.65: Numerical aperture
- 160: Tube length
- 0.17: Thickness of cover glass

### Cleaning

**NOTE: Zuzi recommends using the microscope cleaner kit (code 89000001).**

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the microscope and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.
- ◆ Lenses must not be disassembled by the user. Were there any dust or dirt to be cleaned, clean it with a natural horse hair brush or a smooth piece of cloth, fluff-free, dampened with a bit of xylol or toluene.
- ◆ To remove dust from lenses, blow them with a plastic bulb or clean them with a natural horsehair brush.
- ◆ Use non-corrosive lubricants in metallic parts, being careful of not touching optical parts.

**ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

## 6. CHOOSING OBJECTIVES AND EYEPIECES

The Zuzi microscopes series 122 are supplied with WF16x (only in 122/107 and 122/108 models) and WF10x eyepieces and with 4x, 10x, 40x and 100x achromatic objectives, but can be used with other eyepieces and objectives. To know our offer of eyepieces and objectives look up our catalogue or visit our web page ([www.auxilab.com](http://www.auxilab.com)).

*Note: The microscopes 122/147 and 122/148 only can be used with eyepieces WF/10x-20 mm*

The observed image loses surface and sharpness as magnification increases. The mentioned increase should be done by changing objectives and putting each time a more powerful one and not by changing eyepieces to a higher magnification, as eyepieces only magnify the image obtained with the objective. Thus, the more magnification the eyepiece has, the higher will be the loss of sharpness, clarity and surface of the resulting image.

For routine observations use the eyepiece with lesser magnification with more powerful objectives. The eyepieces with the higher magnification should be kept back for particular occasions, bearing in mind that they decrease definition and do not increase resolution.

## 7. CAUSES OF A DEFECTIVE IMAGE

In case of a defective image check the following:


- ◆ Illumination is well done and luminous intensity is neither excessive nor too weak. Never adjust it with the condenser's diaphragm. Both the condenser and the lamp should be well centred. Check that between the field diaphragm and the aperture diaphragm there is not

**IMPORTANT:** before changing the lamp bulb or the fuse, make sure that the microscope is disconnected from the net. Check you are using the proper lamps, as other types can provoke malfunction.

**Changing the lamp**

- ◆ In case the lamp blows up, replace it for another one, making sure it is a 20 W Zuzi original one. Please take care and do not touch the lamp with bare hands. Do not use lamps with a higher power, as they could provoke over heating or any malfunction.
- ◆ To change the lamp, remove the precondenser (Figure 8)
- ◆ Replace the lamp, holding it firmly. Centre the lamp holder again, if necessary.


Figura 8

**Changing the fuse**

- ◆ If it is necessary to change the fuse loosen the protective lid that is at the back of the microscope over the socket; change the fuse and place it in the same position, putting the protective lid again. (Figure 9)

Figura 9: Cambio de fusible

Magnification	4x	10x	20x*	40x (R)	60x (R)*	100x (R) (I)
Working distance (mm)	17.50	7.32	1.26	0.63	0.41	0.19
Numerical aperture	0.10	0.25	0.40	0.65	0.85	1.25

(R) Retractable; (I) Immersion

\* Optional

Eyepieces: made of two lenses separated by a diaphragm; they project the image from the objectives to the eyes.

**Note:** The microscopes model 122/107 and 122/108 have eyepiece of 23 mm diameter, and the 122/147 and 122/148 models have 30 mm diameter eyepieces.

Type / magnification	WF/10x	WF/16x
Focal distance	24.99 mm	15.58 mm
Field diameter	18 ó 20 mm	11 mm

**Total magnification**

It is the result of multiplying the eyepiece magnification and the objective magnification.

	Objectives	4x	10x	20x*	40x (R)	60x (R)*	100x (R) (I)
<b>Eyepieces</b>							
WF/10x	Total magnification	40	100	200	400	600	1000
	Field of view (mm)	4.50	1.80	0.90	0.45	0.30	0.18
WF/16x	Total magnification	64	160	320	640	960	1600
	Field of view (mm)	2.75	1.10	0.55	0.27	0.18	0.11

(R) Retractable; (I) Immersion

\* Optional

**Note:** The microscopes 122/147 and 122/148 only can be used with eyepieces WF/10x-20 mm

**Numerical aperture**

It determines the properties of the objectives. A higher numerical aperture makes the image brighter and better resolved and permits to visualize more details.

**Working distance**

It is the distance (mm) between the sample specimen and the frontal lens of the objective when the microscope is focused.

**Focal distance**

It is the distance, in mm, from the main image plane to its image focus.

**Resolution**

It is the reciprocal value of the separating power, that is to say, the smallest distance bet-

ween two points on a sample that can still be distinguished as two separate entities.

#### Field of view

Dimensions of the real field that is being observed, expressed in mm.

## 4. INSTALLATION / SETTING UP

### Preliminary inspection

Unwrap the microscope, take off the involving plastic and take off the polyspan protection in which it comes fitted. Take off all the protective items and, without connecting the microscope to the net, make sure that it does not present any damage because of the shipment. In case the microscope presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- Base
- Head (depending on the model supplied)
- Achromatic objectives: 4x, 10x, 40x (R) and 100x (R) (I)
- Eyepieces: WF/16x (only in 122/107 and 122/108) and WF/10x
- Green and blue filters
- Halogen lamp, 20 W
- Fuse
- Standard Schuko wire
- Plastic cover
- User's manual
- Warranty certificate

We will only accept the equipment return within 15 days after delivery provided it comes in its original wrapping.

### Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

**PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT**

◆ Please put the microscope on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the microscope near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.

se it directly to the sun. Avoid vibrations, dust and dry environments.

- ◆ During operation dangerous materials such as flammable or pathological substances must be out of the working area.
- ◆ In case of using the immersion objective, be careful with the immersion oil; if it comes into contact with the eyes or skin, wash them very well with water.
- ◆ When you are not using the microscope for a long period of time please make sure it is unplugged from the net in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

## 5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

*Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.*

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The microscope is supplied with a Schuko standard wire and it should be connected to a current wire provided with an earth wire, it should be handy to be disconnected in case of emergency.
- ◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the microscope (burns, hurts...) or damages in nearby equipments.
- ◆ Please use the plastic cover provided whenever the microscope is not used so as to avoid dust laying on the optics.
- ◆ Please keep the original packaging to transport the equipment, either when it is not being used for a long time or in case you send it for an overhaul.
- ◆ In the event of breakdown please contact your distributor to overhaul through Zuzi Technical Assistance Department.

diaphragm to check that it is perfectly centred in the field of view.

#### Use of the field diaphragm

The field diaphragm is used to adjust the diameter of the light beam that reaches the objective; this way it excludes the external light increasing the image contrast and the field depth and decreasing the resolving power.

- ◆ For a correct use of the field diaphragm, proceed as follows with each objective:
- ◆ Focus the sample
- ◆ Close the field diaphragm by turning clockwise the ring located at the precondenser, so as the image of the diaphragm is projected at the same plane of the sample and appears in the field of view.
- ◆ Move vertically the condenser until focusing the image of the diaphragm.
- ◆ Open the diaphragm until it completely circumscribes the field of view.
- ◆ The filament of the lamp must be perfectly centred with regard to the condenser.

#### Use of the aperture diaphragm

The aperture diaphragm determines the numerical aperture of the illumination system. Its main effect is adjusting the resolution and contrast of the image; if the aperture diaphragm is closed the focusing depth is increased. Since the contrast of samples is usually low, it is recommended to adjust the aperture diaphragm of the condenser to the 70-80% of the numerical aperture of the objective used.

In order to carry out the adjustment, start with maximum aperture. As the contrast increases because the diaphragm is closer, diffracting phenomena appear which limit the resolution. The correct adjustment will be obtained when the maximum contrast is get without changing resolution.

For the correct use of the aperture diaphragm, proceed as follows:


Figura 7: Diafragma iris

intensity when this is not adequate.

#### Security

- ◆ The microscope must be used by previously qualified staff that knows how to operate it thanks to the user manual.
- ◆ Put the microscope in a horizontal plane stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the microscope near any warm supply (burners, blowlamps, etc), nor expo-

◆ **VERY IMPORTANT:** Never hold the microscope by the stage nor the tube, as all the equipment's weight would fall on the fine adjustment screw and mechanical and precision parts would be slowly eroded. Hold it by the slot placed at the back and by the base.

- ◆ Avoid flammable or toxic substances in the working area.
- ◆ The operator has to assume a comfortable, upright position with the back straight. It is convenient to work in a dark table in order to avoid an awkward fatigue and parasite light, as it may dazzle the operator or decrease the image sharpness.
- ◆ The microscope is supplied with a Schuko standard wire. Please insert the wire that feeds the AC electric current in the base of current 220 V, 50 Hz  $\pm 10\%$  provided with earth wire and the other end to the microscope connector. The tension should be 220 V, 50 Hz  $\pm 10\%$ .

*Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz  $\pm 10\%$ .*

◆ Please use the plastic cover provided whenever the microscope is not used so as to avoid dust laying on the optics.

◆ If you are not using the microscope for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

#### Setting up

- ◆ Lower the stage with the help of the coarse adjustment knob. Screw the objectives on the nosepiece so as when it is turned clockwise the objectives follow an ascendant way (4x-10x-40x-100x) (Figure 3).


Figura.3

- ◆ Assemble the head by adjusting it with the built-in screw. Insert the eyepiece/s (Figure 4). The head is normally placed facing the front of the microscope, but if it is necessary it can be placed in any other direction.


Figura. 4

◆ Models provided with binocular observation system must have the eyepiece tubes properly adjusted. This will depend on:


Figura. 5

- Interpupillary distance (*Figure 5*): Adjust the interpupillary distance by separating or joining the eyepieces until there is a total fusion of the two images.

- Dioptre compensation: a more or less marked difference in the focusing is appreciated by closing alternatively one eye, then the other; the dioptre compensation corrects this effect. The microscope is provided with a dioptre compensation ring at the left eyepiece holder. First of all, adjust the dioptre ring to zero and looking through the right eyepiece, focus the sample specimen using the coarse and fine adjustment knobs. Once the sample is focused, look through the left eyepiece and adjust the dioptre compensation ring until obtaining a clear, sharp image.

**NOTE:** Once the values of interpupillary distance and dioptre compensation are found it will be very useful to learn them by heart, mostly if the microscope is shared with other users. This is highly advisable to avoid repeating the adjustment each time the microscope is used.

◆ Lower the stage according to the appropriate working distance for the objective with lesser magnification before placing the sample on it.

◆ Put the slide on the stage holding it with the clip and, by looking directly at it (not by the eyepiece), place the sample centred on the stage opening.

◆ Looking through the eyepieces using the objective with lesser magnification, raise the stage by using the coarse adjustment knob until the image of the sample appears. Then, move the fine adjustment knob until the image is focused.

**IMPORTANT:** operation should always start by using the objective with lesser magnification, as it facilitates focusing and makes impossible to ruin the slides or get the objectives dirty. Furthermore, it is essential for adjusting the lighting previous to using the higher magnification, and it helps to see the whole topographic structure of the sample so as later concentrate on the more appealing parts to observe them at higher magnifications.

◆ By turning the nosepiece, pass to an objective with higher magnification. As they are parfocal objectives it will only be necessary to adjust focusing with the fine adjustment knob.

◆ Adjust the condenser's height to obtain a homogeneous lighting of the field of view. The higher magnifying objective is used, the smaller will be the field of view, and thus, the condenser's precision will be higher with high magnification values. At 40x (objective 4x and eyepiece 10x), descend the condenser to obtain a wider field of view and avoid having excessive lighting. The condenser will have to be raised progressively using higher objectives to concentrate the light beam and get a better illumination and contrast. With the immersion objective the condenser will have to be raised to its maximum height.

### Use of the Köhler illumination system

◆ Zuzi microscopes series 122 are provided with a Köhler illumination system that consists of a precondenser placed at the level of the light source and equipped with a field diaphragm, and a condenser equipped with an aperture diaphragm (*Figure 6*)

6.1 Support of the condenser

6.2 Condenser height adjustment knob

6.3 Condenser


6.4 Aperture diaphragm

6.5 Condenser holding screw

6.6 Condenser centring screws

6.7 Filter holder

6.8 Precondenser with field diaphragm

### Assembling the condenser

◆ Raise the stage to its uppermost position by using the coarse adjustment knob.

◆ Turn the condenser height adjustment knob to lower the condenser support at its lowest position.

◆ Loosen the condenser holding screw.

◆ Place the condenser so as the scale of the aperture diaphragm is located at the front. Slide gently the condenser through the dovetail joint until it stops.

◆ Tighten the condenser holding screw and raise the support to the uppermost position by turning the height adjustment knob.

### Centring the condenser

◆ Raise the condenser support to its uppermost position by using the condenser height adjustment knob.

◆ Focus the sample with the 10x objective.

◆ Close the field diaphragm by turning clockwise the ring located at the precondenser, so as the image of the diaphragm appears in the field of view.

◆ Adjust the height of the condenser until focusing the image of the diaphragm.

◆ Centre the diaphragm by using the condenser centring screws. Open and close the field