AUXILAB S.L.

Material de laboratorio
Laboratory supplies

DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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Parte 1: Requisitos generales

Safety requirements for electrical equipments for measurement, control and laboratory use

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Partie 1: Prescriptions générales.

Fdo: Alfonso Ainciburu Sanz
DIRECTOR | GERENTE


A handwritten signature in black ink, appearing to be 'A. Ainciburu Sanz', written over a white background.

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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Embedding rings

Code	Description
40203001	Embedding rings, white, B/1000
40203002	Embedding rings, yellow, B/1000
40203003	Embedding rings, green, B/1000
40203004	Embedding rings, pink, B /1000
40203005	Embedding rings, blue, B/1000

Universal clamp for cassettes

Code	Description
90500001	Universal clamp for cassettes

**INSTRUCTIONS ON ENVIRONMENT PROTECTION**

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

- Lock the manual handwheel
- Remove the knife from the knife holder
- Remove the sample waste with a dry brush
- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the microtome and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

6. TROUBLESHOOTING

TROUBLE	CAUSE	SOLUTION
Alternately thick and thin sections are produced In extreme cases, no section is obtained.	Knife insufficiently clamped Blunt cutting edge Knife clearance angle too narrow	Re-clamp knife Adjust knife holder laterally or use new blade Systematically try wider clearance angle until optimum angle width has been found.
Sections are extremely compressed, wrinkled or jammed together	Blunt cutting edge Sample too warm Sectioning speed too high	Use a different zone of cutting blade or use new blade Cool sample prior to sectioning Work at lower sectioning speed
Noise during sectioning Knife "sing." when sectioning hard samples. Sections show scratches	Sectioning speed too high Knife clearance angle too wide Insufficient clamping of sample and/or knife holder	Work at lower sectioning speed Systematically reduce clearance angle until optimum angle width has been found Check all screws and clamping mechanisms at sample holder system and knife holder. Retighten if necessary

microtome 508. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita rotary microtome 508, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the equipment avoiding sudden movements, knocks, free fall of heavy / sharp objects on it.
- ◆ Never dismantle the different pieces of the equipment to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the microtome.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Accessories (including their loss) are not covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish.....	2-11
English.....	12-20
French.....	21-30

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT.....	13
2. DESCRIPTION.....	13
3. TECHNICAL SPECIFICATIONS.....	15
4. INSTALLATION / SETTING UP.....	15
5. MAINTENANCE AND CLEANING.....	18
6. TROUBLESHOOTING.....	19
7. ACCESSORIES.....	20
ANNEX I: CE CERTIFICATE	31

1. USES OF THE INSTRUMENT

Nahita rotary microtome model 508 is designed to slice paraffin embedded plant or animal tissues in histology routine applications in hospital and research laboratories.

2. DESCRIPTION

Microtome is a very precise instrument used to slice paraffin embedded tissues. The obtained sections have a thickness of micrometer and this fact allows its observation under the microscope.

There are several types of microtomes (hand, balance, rotary, sliding microtomes, cryostats, ultramicrotomes, etc.) but the most used is the rotary microtome. Some of its advantages are high precision and the possibility of producing very thin serial sections thanks to the desmultiplication process that changes a rotation movement into an advance movement.

- the blade dispenser
- ◆ Loosen the sample holding screw and remove the sample
- ◆ Remove all section waste
- ◆ Clean the equipment (see Cleaning)

Seguridad

- ◆ The microtome must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ Never try to catch a falling knife.
- ◆ Never place a knife anywhere with the cutting edge facing upwards.
- ◆ Always place the sample before placing the knife.
- ◆ Wear glasses in case of working with samples that can produce splinters.
- ◆ Always lock the handwheel prior to equipment cleaning (see Cleaning), prior to any manipulation of knife or sample and during all work breaks.
- ◆ When working with cleaning detergents, bear in mind the manufacturer safety instructions as well as laboratory regulations in force.
- ◆ No fluid must enter the interior of the equipment during cleaning.
- ◆ Never try to repair the microtome by yourself, since you will lose the warranty and may provoke damages to the general operating system.
- ◆ Made under the European regulations for security on machines.

5. MAINTENANCE AND CLEANING

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

To get the best results and a higher duration of the microtome it is essential to follow the processes of use.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ Never try to repair the microtome by yourself, since you will lose the warranty and may provoke damages to the general operating system.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Cleaning

Before each cleaning:

- ◆ Remove the sample and place the sample holder at the uppermost position

Placing and orienting the knife

- 3.1 Knife holding screws
- 3.2 Cutting angle fixing screws
- 3.3 Cutting angle adjusting lever
- 3.4 Cutting angle scale

The microtome is supplied with a versatile knife in a wedge very useful when tissue hardness varies notably and with a blade holder for disposable blades.

- ◆ Very carefully insert the knife in the knife holder (1.1)
- ◆ Hold the knife firmly with the knife holding screws (3.1)
- ◆ Adjust the cutting angle of the knife; for this purpose, loose the cutting angle fixing screw (3.2) and then adjust the blade to the desired angle (see scale (3.4)) by moving the adjusting lever (3.3). Then hold the position with the fixing screw (3.2).
- ◆ Adjust the position of the knife-holder; for this purpose, turn rightwards the manual advance locking lever (1.6) and then adjust manually the position of the knife-holder approaching it to the sample. Then, turn again the lever (1.6) to hold this position. By turning the coarse feed knob (1.5) (graduated from 0 to 500 µm/turn) adjust the final position.

Trimming the sample

Sample trimming consists of removing the upper layers of paraffin until the whole surface of the tissue is obtained in each slice. For this:

- ◆ Select a thick slice thickness (max. 25 µm) by turning the coarse feed knob (1.5).
- ◆ Turn the handwheel (1.3) until the sample has been trimmed down to the desired level.

Tissue sectioning

Never cut the sample with the same region of the knife edge used for trimming:

- ◆ Select the desired slice thickness (in microns) by turning the corresponding knob (1.7); take into account that sections of paraffin embedded tissues are very difficult to obtain when thickness is lower than 5 µm.
- ◆ Pull the locking device (1.4) to release the handwheel.
- ◆ For sectioning, turn the manual handwheel (1.3) uniformly and clockwise at a constant speed.
- ◆ Pick up the obtained sections and place them on microscope slides.

NOTE: always rotate the handwheel evenly. Rotating speed must be adjusted to the sample hardness; harder samples should be cut at slower speed.

Finishing operation

- ◆ Move sample clamp to uppermost position and lock the handwheel
- ◆ Turn the lever (1.6) rightwards and manually move backwards the knife holder
- ◆ Release the knife and carefully insert it into the receptacle for used blades in the base of

The rotary microtome presents the following different parts:

Sample holder: it is the device in which the paraffin embedded tissue is placed. Thanks to the screws the paraffin block surface can be oriented in a position parallel to the blade.


Knife holder: it is the device in which the knife is placed; the approach of the knife towards the sample can be adjusted through a coarse feed knob.

Knife: in a wedge type knife or universal disposable blades

Manual handwheel: it transmits to the advance mechanism the impulse necessary for its correct operation; each complete turn of the handwheel produces a cut sequence. The handwheel presents a safety device that locks it for its transport or when it is not being used. So, before using the microtome it is essential to pull the blocking mechanism and give a half-turn to the handwheel to place it on cutting position.

Advance mechanism: the sample holder moves up and down while the knife holder moves towards the sample. A complex mechanical system is necessary to perform these movements that finally allow obtaining serial tissue sections from the paraffin block.

- | | |
|-------------------------------|-------------------------------------|
| 1.1. Knife holder | 1.5. Coarse feed knob |
| 1.2. Sample holder | 1.6. Manual advance locking lever |
| 1.3. Manual handwheel | 1.7. Slice thickness adjusting knob |
| 1.4. Handwheel locking device | 1.8. Slice collecting tray |

3. TECHNICAL SPECIFICATIONS

<i>Slice thickness</i>	1-25 µm
<i>Min. adjustment</i>	1 µm
<i>Accuracy</i>	±10%
<i>Max. specimen size</i>	40x30 mm
<i>Dimensions (LxWxH)</i>	42x30x32 cm
<i>Weight</i>	28 Kg (approx.)

4. INSTALLATION / SETTING UP

Preliminary inspection

- ◆ Unwrap the microtome and take off all the protections in which it comes fitted.
- ◆ Make sure that the equipment does not present any damage because of the shipment. In case it presents any damage, tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping with all the accessories and documents supplied.

- ◆ Check all the accessories that must be supplied together with the equipment:
 - Knife in a wedge and blade holder for universal disposable blades
 - Lubricating oil

Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.


TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please put the microtome on top of a horizontal, plane and stable table making the enough free space around it.

Setting up

Placing and orienting the sample

- 2.1 Sample holding screw
- 2.2 Sample orienting screws
- 2.3 Sample position locking lever

NOTE: Always, place the sample before placing the knife.

The manual handwheel must be always locked before changing the sample as well as during all work breaks.

- ◆ Put the sample holder (1.2) at the uppermost position; for this, turn the manual handwheel (1.3) and lock it with the safety device located at the lower side (1.4).
- ◆ Place the paraffin block in the sample holder (1.2) and hold it firmly with the sample holding screw (2.1).
- ◆ Turn counterclockwise the locking lever (2.3) and use the sample orienting screws (2.2) to orient the paraffin block surface so as it is perfectly parallel to the edge of the knife.