ESPECTROFOTÓMETROS SPECTROPHOTOMETERS SPECTROPHOTOMÈTRE

Modelos / Models/ Modèle 4201/50

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Zuzi spectrophotometer model 4201/50. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Zuzi spectrophotometer model 4201/50, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the spectrophotometer avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the spectrophotometer to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the spectrophotometer.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Zuzi Technical Assistance Department by email to asistencia@auxilab.es.
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

ENGLISH

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish	2-12
English	13-23
Français	24-36

INDEX OF CONTAINS

1. USES OF THE INSTRUMENT	14
2. DESCRIPTION	14
3. TECHNICAL SPECIFICATIONS	15
4. INSTALLATION / SETTING UP	15
5. MAINTENANCE AND CLEANING	21
6. TROUBLESHOOTING	23
ANNEX I: CE CERTIFICATE	36

1. USES OF THE INSTRUMENT

Zuzi spectrophotometer model 4201/50 is a versatile instrument that performs precise quantitative and qualitative analysis in the ultraviolet spectral range. It is very useful in conventional laboratories and in a wide range of fields such as petrochemical and mining industry, agricultural engineering, environmental preservation, biochemistry, hospitals and research laboratories.

2. DESCRIPTION

Figure 1

- 1.1 LCD display and control panel
- 1.2 Sample compartment
- 1.3 Cell changer manual knob
- 1.4 Power ON/OFF
- 1.5 Fuse holder
- 1.6 Socket
- 1.7 USB and parallel ports
- 1.8 Manual power switch 110V/220V



Figure 2

- 2.1 GOTO λ : selection of wavelength
- 2.2 ZERO: blank
- 2.3 SET: selection of the different functions
- 2.4 PRINT: Printing of results
- 2.5 (-) left (F1): confirmation of the left command on display
- 2.6 (-) right (F2): confirmation of the right command on display
- 2.7 \blacktriangledown : parameter decrease
- 2.8 \blacktriangle : parameter increase

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<i>Code</i>	54201050
<i>Wavelength range</i>	200-1000 nm
<i>Spectral bandwidth</i>	4 nm
<i>Optical system</i>	Simple beam, grating 1200 lines/mm
<i>Wavelength accuracy</i>	± 2nm
<i>Wavelength repeatability</i>	1 nm
<i>Wavelength resolution</i>	± 1 nm
<i>Photometrical range</i>	-0.097/1.999 A, 0/125% T
<i>Photometrical accuracy</i>	± 0.5 %T
<i>Photometrical repeatability</i>	± 0.3 %T
<i>Stray light</i>	0.3% T
<i>Stability</i>	± 0.004 A/h to 500 nm
<i>Detector</i>	Silicon photodiode
<i>Sample compartment</i>	4 cells of 10 mm
<i>Lamps</i>	Tungsten and deuterium
<i>Output</i>	USB port, parallel port (printer)
<i>Power</i>	AC 220 V / 50 Hz or AC 110V/50Hz
<i>Dimensions</i>	480x360x160 mm

4. INSTALLTION / SETTING UP

Unwrap the spectrophotometer, take off the involving plastic and take off the protection in which it comes fitted. Take off all the protective items and, without connecting the spectrophotometer to the net, make sure that it does not present any damage because of the shipment. In case the spectrophotometer presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied. Please check that all the accessories are enclosed with the equipment:

- 4 glass cells
- 2 quartz cells
- Cable USB
- Software MWave Basic
- Schuko wire
- User's manual
- Warranty

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.

Please put the spectrophotometer on top of a horizontal, plane and stable table making a free space of at least 30 cm per side. Do not put the spectrophotometer near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc. Avoid inflammable or toxic substances in the working area.

The spectrophotometer is supplied with a Schuko standard wire. Please insert the wire that feeds the AC electric current in the base of current 220 V, 50 Hz ± 10% provided with earth wire and to the other end to the spectrophotometer connector.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz ± 10%.

Setting up

- ◆ Press On/Off (1.4).
- ◆ The equipment initializes and makes an autotest. This process can spend a few minutes; during this time *Zuzi self test...* will be displayed. After this time, *Warm up 20 minutes any key to skip* will be displayed; wait for 20 minutes before operating with the equipment.


NOTE: Do not open the sample compartment during autotest.

BASIC FUNCTIONS

1. Wavelength adjustment (Fig. 3)

- ◆ Press GOTOλ (2.1) to enter the wavelength adjustment display. Select the desired value of wavelength by pressing ▼ ▲ (2.7) (2.8).
- ◆ Select OK with F1 (2.5) to save the wavelength value.
- ◆ The equipment automatically returns to the main menu.
- ◆ If there is no need of changing the wavelength value, select Return with F2 (2.6) to go to the main menu.


Figure 3

2. Blank (Fig. 4)

- ◆ Put the cell with the reference sample into the cell holder compartment.
- ◆ Press ZERO (2.2); the equipment will make the blank automatically.


Figure 4

3. Utilities (Fig. 5)

- ◆ Press SET (2.3) to enter the function selection display.
- ◆ Use ▼ ▲ (2.7) (2.8) to move along the different functions and select OK with F1 (2.5) to enter the selected function.
- ◆ If no function needs to be selected, select Cancel with F2 (2.6) to return to the main menu.

Figure 5


3.1 Test Mode (Fig. 6)

To select the measuring mode, absorbance or transmittance, proceed as follows:

- ◆ At the function selection display, move the cursor with ▼ ▲ (2.7) (2.8) up to the Test Mode function.
- ◆ Confirm OK with F1 (2.5) to enter the test mode selection display.
- ◆ Select the desired measuring mode with ▼ ▲ (2.7) (2.8) : Abs (absorbance), T% (transmittance) or Energy. Confirm OK with F1 (2.5).
- ◆ The symbol ✓ will be displayed beside the selected measuring mode.
- ◆ Confirm Cancel with F2 (2.6) to return to the function selection display.

Figure 6

3.2 Deuterium lamp on/off (D2 Lamp On/Off) (Fig. 7).


Note: Zuzi recommends having the deuterium lamp switched off when working out of the ultraviolet range.

- ◆ At the function selection display, move the cursor with ▼ ▲ (2.7) (2.8) up to the D2 Lamp function.
- ◆ Confirm OK with F1 (2.5) to enter the deuterium lamp state selection display.
- ◆ With ▼ ▲ (2.7) (2.8) select On or Off and confirm OK with F1 (2.5). Beside the selected state the symbol ✓ will be displayed.
- ◆ Confirm Cancel with F2 (2.6) to return to the function selection display.

Note: When the deuterium lamp is switched on, it will take a few minutes to warm up.

Figure 7


3.3 Tungsten lamp on/off (W Lamp On/Off) (Fig. 8)

Note: Zuzi recommends having the tungsten lamp switched off when working out of the visible range

- ◆ At the function selection display, move the cursor with ▼ ▲ (2.7) (2.8) up to the W Lamp function.
- ◆ Confirm OK with F1 (2.5) to enter the tungsten lamp state selection display.
- ◆ With ▼ ▲ (2.7) (2.8) select On or Off and confirm OK with F1 (2.5). Beside the selected state the symbol ✓ will be displayed.
- ◆ Confirm Cancel with F2 (2.6) to return to the function selection display.

Figure 8


3.4 Dark Current (0%T) (Fig. 9)

Zuzi recommends to carry out this function whenever changing the environmental conditions (changes in temperature or equipment location, peaks in electrical current, etc.). For this purpose:

- ◆ At the function selection display, move the cursor with ▼ ▲ (2.7) (2.8) up to the Dark current function.
- ◆ Confirm OK with F1 (2.5); Get dark current will be displayed.
- ◆ After a few seconds, the equipment will automatically return to the function selection display and will be ready to operate.


Figure 9

3.5 Wavelength calibration (Reset WL) (Fig. 10)

- ◆ At the function selection display, move the cursor with ▼ ▲ (2.7) (2.8) up to the Reset WL function.
- ◆ Confirm OK with F1 (2.5); Calibrating WL... will be displayed.
- ◆ After a few seconds, the equipment will automatically return to the function selection display and will be ready to operate.

Figure 10

4. Sample analysis


Before proceeding to analysis clean the cell with distilled water. Wipe the cell with a cellulose paper to remove drops of water and fingerprints.

4.1 Basic mode

To determine the absorbance or transmittance of different samples at a certain wavelength:

- ◆ Put the cell with the blank reference into the cell holder and place it into the light path.
- ◆ Select Basic with F1 (2.5). The equipment will make the blank automatically and then will enter the measuring mode display (Fig. 11).

Figure 11


- ◆ On the upper part, the selected wavelength (WL) and 100% T or 0.0000 Abs depending on the working measuring mode will be displayed. Immediately below a chart headline will appear:

No.	WL	T%	Abs
Number assigned in order to each result	Wavelength	Transmittance value	Absorbance value

The results of the successive analysis will be shown on this chart.

- ◆ Put the cell with solution to be analyzed into the light path. Select Test with F1 (2.5); the result of the measurement will be displayed (figure 12)

Figure 12

- ◆ Repeat this process with the different samples. The spectrophotometer can save up to 200 data.
- ◆ To return to the main menu select Cancel with F2 (2.6).

Note: At the basic mode, the wavelength can be adjusted and the blank can be made as many times as needed by pressing GOTO λ (2.1) and ZERO (2.2) respectively.


4.2 Quantitative mode

Zuzi spectrophotometer model 4201/50 allows performing quantitative analysis by the single point standard curve method. The standard curve is obtained through a sample of known concentration. The absorbance of this sample is measured and the equipment calculates the equation of the corresponding curve, which will be used for the analysis of samples of unknown concentration.

Note: the curve will not be displayed but it will be used to measure the concentration of the samples.

- ◆ Select Quantitative with F2 (2.6). First the system will ask to insert the blank (Fig. 13).

Figure 13


- ◆ Insert the blank, the standard sample and the sample of unknown concentration into the sample compartment. Put the blank into the light path.

- ◆ Press GOTO λ (2.1) and adjust the working wavelength value.
- ◆ Select OK with F1 (2.5) to make the blank.

Then the system will ask to insert the standard sample (Fig 14).

- ◆ Put the standard sample into the light path and input the value of its concentration. For this, adjust the value digit

Figure 14

by digit with $\blacktriangle/\blacktriangledown$ (2.7) (2.8); to go from one digit to the next one select OK with F1 (2.5). When the value of the last digit is adjusted select OK with F1 (2.5) (Fig. 15).

- ◆ Put again the blank into the light path and press ZERO (2.2) to make the blank.
- ◆ Put the sample of unknown concentration into the light path and the concentration value will be displayed.

4.3 Deleting and printing results

At the result chart of any measuring mode, press PRINT (2.4) and the following options will be displayed (Fig. 16):

1. Print and clear: for printing and then, deleting data
2. Only clear, not print: for deleting data but not printing
3. Cancel

- ◆ Select the desired option with $\blacktriangle/\blacktriangledown$ (2.7) (2.8) and then confirm OK with F1 (2.5).

If you are not using the spectrophotometer for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Figure 15



Figure 16

Security

- ◆ The spectrophotometer must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ You should put the spectrophotometer on top of a horizontal, plane, stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the spectrophotometer near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During operation dangerous materials such as flammable or pathological substances must be out of the safety area.
- ◆ When you are not using the spectrophotometer for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the spectrophotometer by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the spectrophotometer (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The spectrophotometer is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the spectrophotometer by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the spectrophotometer (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Zuzi Technical Assistance Department.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the spectrophotometer and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

6. TROUBLESHOOTING

TROUBLE	CAUSE	SOLUTION
The equipment does not turn on	The wire is not connected	Connect the wire
	The power source is broken	Change the power source
	Voltage error	Check the voltage of the electric installation
	Fuse blown or defect in an electric component	Contact the Zuzi Technical Assistance Department
100% T (0.000A) can not be selected	The sample is not well placed in the cell holder	Put correctly the sample into the cell holder
	The lamp is worn out or defective	Change the lamp
	The lamp is turned off	Look up the user's manual, lamp management section
	Defect in an electric component	Contact the Zuzi Technical Assistance Department
The equipment strays and there is noise	Not enough period of warm up	Leave 20 minutes to warm up
	Use of glass cell in the UV range	Use a quartz cell
	The lamp is worn out or defective	Change the lamp
	Instable power source	Stabilize the voltage
	Significant change of temperature The lamp is not adjusted	Adjust correctly the lamp Contact the Zuzi Technical Assistance Department
	Defect or dirtiness in detector or defect in an electric component	Contact the Zuzi Technical Assistance Department
Incorrect readings	Not enough volume of sample	Fill the cell with sample
	Equipment out of calibration	Contact the Zuzi Technical Assistance Department
	Bubbles or particles in sample	Check sample preparation, analytical procedure
	Lost of vapours of sample	Adapt the sample to the equipment Use the adequate ventilation
	Wrong blanking	Check the wavelength and the analytical procedure Make Dark current and calibrate the equipment

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