DECLARACIÓN DE CONFORMIDAD CE  
CE DECLARATION OF CONFORMITY  
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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## 8. TROUBLE SHOOTING

PROBLEM	CAUSE	SOLUTION
Digits do not appear on the screen	<ol style="list-style-type: none"> <li>1. Check that the electrical power gets into the balance. The adapter could be damaged</li> <li>2. The fuse does not work.</li> </ol>	<ol style="list-style-type: none"> <li>1. Adjust the power plugs. Change the current adapter.</li> <li>2. Replace the fuse.</li> </ol>
The value on the screen is changing all the time	<ol style="list-style-type: none"> <li>1. Your balance is in a place with vibrations.</li> <li>2. The balance's door is not completely closed.</li> <li>3. An odd object can be between the pan and its placing, causing the destabilization.</li> <li>4. The underhook where the thread is hitched in hydrostatic measuring, is open.</li> </ol>	<ol style="list-style-type: none"> <li>1. Move the balance to a place far from vibrations.</li> <li>2. Close every balance's door.</li> <li>3. Check the inner part of the balance is clean and there is no odd object in it.</li> </ol>
The value on the screen is wrong	<ol style="list-style-type: none"> <li>1. The balance is not calibrated</li> <li>2. The balance has not been tared before weighing</li> <li>3. The balance is not leveled</li> </ol>	<ol style="list-style-type: none"> <li>1. Calibrate</li> <li>2. Tare the balance before each weighing.</li> <li>3. Level the balance until the bubble in the back is centered.</li> </ol>

**NOTE:** According to the in force legislation regarding "Non-automatic weighing instruments" in which balances are included, by means of writ dating from 22nd October, 1994 (BOE 3/1/95), these balances must not be used for:

- Commercial Transactions.
- Calculating taxes, tariffs, rates, indemnities and other similar canons.
- Judicial surveys.
- Pharmaceutical medicine preparations, as well as analysis made in medical or pharmaceutical laboratories.
- Determining the price or total amount in sale price or in prepackaged preparations.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your analytical balance Series 5034. We highly recommend looking after this equipment according to what is stated in this manual.

NAHITA develops its products being remarkable in the ergonomics and security of its user. The incorrect use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages. Please read the point of Maintenance, where we expose the security notes.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

**TO GET THE BEST RESULTS AND HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

Please bear in mind the following:

- ◆ This manual should be available for all the users of these equipments
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ As an accuracy instrument, we should manipulate carefully avoiding the sudden movements, knocks, free fall of heavy objects, manipulation of sharp objects, and the pouring of liquids in the interior of the balance.
- ◆ Never dismantle the different pieces of the balance to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty.
- ◆ To prevent fire or shock hazard, do not expose the balance to rain or humid environments.
- ◆ If you have any doubts do not hesitate in contacting your distributor or our Nahita Technical Assistant Department by email [asistencia@auxilab.es](mailto:asistencia@auxilab.es).
- ◆ This equipment is protected under the warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the balance warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the balance warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

## 7. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is necessary to follow the processes of use.

*Note: All the processes of use mentioned before will not have any value unless you keep a continued and careful maintenance.*

- ◆ Please follow the processes of use of this manual
- ◆ This manual should be available for all users of these equipments
- ◆ Always use original components and supplies. Other dispositives can be similar but they can damage the equipment.
- ◆ Preserve the balance from heavy shocks and look for a dry place, away from dust, direct sun beams or air flows. Handle the balance gently, since it is an accurate instrument.
- ◆ Do not press the display or any key with a pointed object such as a ballpoint pen or a screwdriver, as this may scratch or damage the key. Use your only your fingers.
- ◆ Do not let objects fall on the pan in order to avoid damages to the weighing system.
- ◆ Never place into the balance an object heavier than the maximum weight of the balance's range, the sensor could be damaged.
- ◆ Never disconnect it by pulling the cable for security reasons.
- ◆ Place the unit in a location with adequate ventilation to prevent heat build-up the unit.
- ◆ Do not submerge the balance in liquids.
- ◆ The balance has a Schuko plug. This plug must be connected to a switch connected to ground. This switch must be close in order to be able to disconnect it in case of emergency.
- ◆ In case of damage or doubt it can be solved getting in touch with the NAHITA Technical Service

### Cleaning

- ◆ For the cleaning of the inox steam parts use always a special cloth that does not give off fluff. Do not use any cloth that can produce damages to the equipment.
- ◆ For the cleaning of the different pieces of the equipment, we highly recommend to use neutral products that do not have a PH superior than 6 (soaped water without abrasives)

**ATTENTION! AUXILAB DOES NOT ADMIT REPARATIONS IF THE EQUIPMENT IS NOT CLEAN AND DISINFECTED.**

### INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances.

It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking.

By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

**Command mode**

Data will be transmitted from the balance to an external device when they receive an order. Balance can be controlled from the external device by the following orders in the same way as by the frontal panel of the balance.

**Orders (Commands)**

<O><CR><LF> Orden ON/OFF (4F 0D 0A) the same as ON/OFF function from the frontal panel.


<P><CR><LF> Printing order (50 0D 0A) the same as PRINT function from the frontal panel.

<T><CR><LF> Tare order (54 0D 0A) the same as TARE function from the frontal panel.


<M><CR><LF> Mode order (4D 0D 0A) the same as MODE function from the frontal panel.

<C><CR><LF> Calibration order (43 0D 0A) the same as CAL function from the frontal panel.

**Connection with external services: 'Connection with PC'**

**Connection with external services: 'Connection with printer'**

**ATTENTION! AUXILAB DOES NOT ADMIT REPARATIONS IF THE EQUIPMENT IS NOT CLEAN AND DISINFECTED.**

**INDEX OF LANGUAGES**

Spanish ..... 2-16  
 English ..... 17-32  
 French ..... 33-47

**INDEX OF CONTENTS**

1. APLICATIONS..... 19  
 2. DESCRIPTION ..... 20  
 3. TECHNICAL SPECIFICATIONS ..... 21  
 4. SETTING UP..... 21  
 5. FUNCTIONS OF THE BALANCE ..... 23  
 6. SOFTWARE ..... 27  
 7. MANTEINANCE AND CLEANING ..... 31  
 8. TROUBLESHOOTING..... 32  
 ANNEX I: CE CERTIFICATE ..... 48


**1. USES OF THE EQUIPMENT**

The Nahita analytical balance is an essential instrument in laboratory work when weighing tasks need to be carried out. It is a reliable and long-lasting instrument, which has been developed with the latest technology. Its simple, sturdy, ergonomic design makes this balance suitable for any industrial use, and its smooth surface makes the cleaning easier.


**2. DESCRIPTION**

- 1.1 Glass chamber
- 1.2 Front panel
- 1.3 Pan
- 1.4 Screwed legs
- 1.5 Bubble level
- 1.6 Adapter connection
- 1.7 RS232 output

**Front Panel**

- 2.1 Display
- 2.2 ON/OFF
- 2.3 TARE: Tare
- 2.4 CAL: Calibration
- 2.5 PRINT: Data printing
- 2.6 MODE: Unit and measurement mode selection
- 2.7 STAB: Stability indicator
- 2.8 %: percentage mode
- 2.9 COUNT: piece counting mode

**RS-232C OUTPUT. SPECIFICATIONS.**

Transmission way: double, asynchronous and bi-directional  
 Data format: digital transmission speed: 1200, 2400, 4800 and 9600.  
 Informative bit: 8  
 Bit of parity: none  
 Bit of stop: 1  
 Code: ASCII  
 Data output: Auto-print, PRINT key, data emission

Data output format:

```
ST, +001.000xxg <CR><LF><LF>
US, +001.000xxg <CR><LF><LF>
```

**Note:**

**ST = Stable (53H, 54H) (Estable)**  
**US = Unstable (55H, 53H) (Inestable)**  
**X = Space (20H) (Espacio)**  
**<CR> = Carriage return (0DH)**  
**<LF> = Line Feed <0AH>**

**DATA OUTPUT.**

**Auto-print mode**

There are three different ways in Auto Print mode:

◆ **AUTO PRINT A**

**Note: the pan must be clean and screen must show Zero after each weighing in AUTO PRINT A mode.**

Data transmission starts when balance begins to stabilize (stability indicator appears (2.7)).

◆ **AUTO PRINT B**

Data transmission happens whenever balance stabilizes.

◆ **AUTO PRINT C**

Data are transmitted continuously when plugging the balance without giving any additional order to neither the balance nor the computer or printer.

**Print mode**

Data will only be transmitted if PRINT key (2.5) is pressed.

Table 2:

Function	Parameters (CX-Y)	Name	Meaning
C1	0~3	Display actualization	
C2	0~3	RS232 data output	
C3	C3-0* C3-1	Zero and tare	ON OFF
C4	C4-0* C4-1 C4-2 C4-3	Digital transmission speed	2400 1200 4800 9600
C5	C5-0* C5-1 C5-2 C5-3	Data exit	Zero stabilization Balance stabilization Print Continuous exit
C6	C6-0* C6-1	Back illumination	ON OFF
C7	C7-0* C7-1 C7-2 C7-3	Stability band's amplitude allowed	1 2 3 4

\*Factory settings

#### DATA RECORDING

All the above features (*table 2*) can be modifying directly from the screen.

- ◆ Press TARE (2.3) and ON/OFF (2.2) keys simultaneously with the balance turned off.
- ◆ Release the key and "88888888" will be displayed.
- ◆ Press PRINT key (2.5).
- ◆ "C1-0" (CX-Y) will be displayed.
- ◆ Press MODE key (2.6) to change the X digit between values from 1 to 7 (*see table 2*).
- ◆ By pressing PRINT key (2.5) change the Y value. This digit can take different values depending on X digit value (*see table 2*).
- ◆ Press ON/OFF key (2.2) to save the new parameters in the memory.
- ◆ The balance is ready to work with the new parameters, that will be recorded when you turn on/off the balance.

### 3. TECHNICAL SPECIFICATIONS

Among its various characteristics we would like to emphasize the following:

- Fast setting up and stabilization (indicator)
- Easy-to-use, light and simple design
- Subtractive tare
- Connection by adapter (included)
- Automatic calibration by external weight
- Weighing in grams, miligrams, ounces and carats
- Weighing modes, percentage, piece counting and below weighing
- Bubble level and screwed legs
- RS-232 output for connection to computer or printer

Code	55034001	55034002	55034003
Capacity	60 g	110 g	200 g
Readability		0.0001 g	
Repeatability		± 0.0001 g	
Linearity		± 0.0002 g	
Operating temperature		5°C ~ 40°C	
Pan size		90 mm	
Chamber size		170 x 190 x 190 mm	
Externals dimensions		330 x 200 x 260 mm	
Weight		6 Kg	
Power		Adapter AC 100-240V / 50-60Hz	
Calibration mass	50 g	100 g	200 g

### 4. INSTALLATION

#### Preliminary inspection

- ◆ Unwrap the balance, take off the involving plastic and take off the poliespan protection in which it comes fitted.
- ◆ Please make sure that it does not present any damage because of the shipment. If it does please communicate to your transport agent or suministrator immediately, so that they can make the claims in the correct time limit.

**We will only accept any equipment return within 15 days after delivery and provided it comes in its original package and with all accessories and documents supplied**

- ◆ Please check that all the accessories are enclosed with the equipment:
  - Pan
  - Calibration weight (depending on capacity)

- Adapter
- Cover
- Cloth
- User manual
- Warranty certificate

### Installation


Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

**PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT**

- ◆ Please put the balance on top of a horizontal, plane, stable table making a free space at least at 30 cm per side of the balance.
- ◆ Do not put the balance near any warm supply (burners, blowlamp...), do not expose directly to the sun, etc.
- ◆ The balance is suminstrated with a Schuko standard wire. Please insert the wire that feeds the A/C in the base of current provided with earth wire.

*The manufacturer and the distributor will not assume any responsibility for the damages produced to the equipment in its installing or damages to persons suffered by the wrong use of the electric connection. The tension should be 100-240V ±10% 50-60Hz.*

- ◆ Position the pieces on the balance as shown in the following sequence: **Picture 3**  
dust guard, draft ring, pan support and pan. (**Picture 3**)
- ◆ Check that the level (**1.5**) is in the right position (centred), if not, centre the bubble with the knobs situated at the frontal legs of the balance (**1.4**).
- ◆ If your balance is not used for a long period, please turn the display OFF then unplug the AC adaptor. This way you will save electricity and your equipment will have a longer life, also you will avoid accidents and an unnecessary use of energy.

### Setting up

**NOTE:** Before using the balance, it must be calibrated according to "Calibration" point of this manual. Recalibration must be performed whenever observing important variations or mistakes in weightings.

- ◆ Make sure that there is nothing on the pan and switch the balance on (**2.2**). The display will light up (**2.1**) and the balance will make and autotest that finishes when "0.0000" is visualized on display.
- ◆ Select the desired function: weighing, percentage or pieces count.
- ◆ By pressing MODE key (**2.6**) select the measuring unit of weight (grams (g), milligrams (mg), ounces (oz) or carats (ct)), percentage or piece count mode.

6. Press TARE key (**2.3**) and after a few seconds, the corresponding number of pieces will be displayed.
7. Remove the pieces from the pan.
8. "0" will be displayed.
9. Place on the pan the amount of pieces to be counted and the number of pieces will be displayed.

**Note:** To set a new reference, press CAL key (**2.4**) and proceed as previously described (from point 4).

### PERCENTAGE MODE

In percentage mode the balance perform measurements based on a reference value.

1. By pressing MODE key (**2.6**), select percentage mode (%).
2. "0.00" will be displayed.
3. Press CAL key (**2.4**).
4. By default, "--100--" will be displayed; this value indicates the weight percentage corresponding to the reference sample.

**NOTE:** To modify this value, press MODE (**2.6**) key (to increase) or PRINT (**2.5**) key (to decrease). You can choose between 50, 100, 150, 200, 250, 300, 350 and 400 %.

5. Place on the pan the sample whose weight is going to be considered as the 100 % (or the selected value).
6. Press TARE key (**2.3**), and after a few seconds, the corresponding percentage will be displayed.
7. Removed the sample from the pan.
8. "0.00 %" will be displayed.
9. Place the sample to be compared on the pan. As a result of the comparison between the reference sample and the problem sample, the deviation in percentage will be displayed.

**Note:** To set a new reference, press CAL (**2.4**) and proceed as previously described (from point 4).

## 6. SOFTWARE

### PARAMETERS.

These balances are supplied with Software to adapt their features, as best as possible, to users' necessities:

Fig. 4





Fig. 5

The balance is provided with a built-in hook to perform below weighings; this technique is widely used for weighing magnetic material and its use can be extended to industry.

For below weighing, proceed as follows:

- ◆ Remove the 2 screws that hold the rhombus-shaped lid placed at the bottom of the balance (see fig.4)
- ◆ Remove the lid to see the hook to hold weighing device (not included with the balance)
- ◆ Put the thread of the weighing device through the hook
- ◆ Place the balance on a table with a hole to let the thread hanging through freely
- ◆ Hang the weighing device at the other end of the thread (see fig.5)
- ◆ Proceed in the same way as for normal weighing

#### PIECE COUNTING MODE (MODE "COUNT")

In "COUNT" mode the balance calculates the number of pieces placed on the pan. Proceed as follows:

1. By pressing "MODE" key (2.6), select "COUNT" mode (2.9).
2. "0" will be displayed.
3. Press CAL key (2.4).
4. By default, "10" will be displayed; this number indicates the number of pieces that should be placed on the pan and that will be taken as a reference for a later count.

**NOTE: To modify this value, press MODE (2.6) key (to increase) or PRINT (2.5) key (to decrease). You can choose between 5, 10, 15, 20, 25, 30, 35 and 40 pieces.**

5. Place 10 pieces of those to be counted, or the selected number of pieces, on the pan.

- ◆ The balance is ready to weigh.
- ◆ Switch the balance off by pressing On/Off key (2.2) once operation has finished.
- ◆ The following table shows the conversion of these units into grams.

TABLE 1: Units conversion

Abbreviation	Name	Conversion
Kg	Kilogram	1000 g
g	Gram	1 g
mg	Milligram	0.001 g
oz	Ounce	28.3495231 g
ct	Carat	0.1999694 g

**NOTE: To avoid possible oscillations in weighing, it is important to avoid air flow and vibrations on the table in which it is situated the balance while weighing.**

#### Security

- ◆ The balance has to be used by qualified personal so it gets its correct use.
- ◆ You should put the balance in a horizontal plane stable table, having a free space of at least 30 cm per side.
- ◆ Do not put the balance near any warm supply (burners, blowlamp...), do not expose directly to the sun, avoid vibrations and dust.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ If you are not going to use the equipment during a long period of time please make sure that it is unplugged from the net
- ◆ Never dismantle the different pieces of the balance to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product's warranty
- ◆ Made under European regulations.

## 5. FUNCTIONS OF BALANCE

To increase your balance's accuracy while weighing, please check that the place where you are going to take your weights has an average temperature of  $\pm 20$  °C and a relative humidity of 50%. Before using your balance, let it 10 or 15 minutes for its stabilizing.

If your balance is used almost every day, just turn the display off by pressing ON/OFF (2.2) key with power connected. In long periods of no usage, keep the balance unplugged.

#### CALIBRATION

The balance should be calibrated the first time it is used. It is necessary to recalibrate your balance when it is moved to another place or after a period of use (about 30 days).

Proceed as following:

- ◆ Turn the balance on by pressing On/Off key (2.2); the pan should be empty.
  - ◆ The balance will perform an autotest and finally stabilizes on "0.0000" g. If "0.0000" g does not appear on display, press TARE key (2.3) until it appears and the stabilization indicator (2.7) lights up.
  - ◆ Press CAL key (2.4) to enter calibration mode.
  - ◆ "CAL--0" will be displayed.
  - ◆ Press TARE key (2.3).
  - ◆ "CAL-F" will be displayed.
  - ◆ Place the calibration weight on the pan.
  - ◆ Press TARE key (2.3) again.
  - ◆ "CA--End" and the weight of the standard calibration weight used will be displayed.
  - ◆ Removed the calibration weight from the pan.
- 0.000 g will be displayed, indicating that balance calibration has been well performed.

It is necessary to check the results of calibration:

Place the calibration weight on the pan when "0.0000" g is displayed, and check that the displayed value coincides with the weight value with a maximum deviation of  $\pm 0.0001$  g; otherwise, balance is not correctly calibrated, and calibration should be performed again.

*Note: If "CAL no" is displayed during calibration, it means that calibration weights are not the correct ones because their weight is up  $\pm 10\%$  on the maximum capacity of the balance.*

Change the weights for ones below these limits and do the calibration again.

If "CAL no" message remains in spite of using the appropriate weights, please ask Nahita authorised technical service.

#### **BALANCE ADJUSTMENT.**

If you need more accuracy and precision in measures and you also know the uncertainty of calibration weights, you can adjust manually the weighing values of the balance. Taking as an example that the calibration weight value is 50.0005, calibration should be performed as follows:

- ◆ Turn the balance on by pressing ON/OFF key (2.2) with the pan empty.
- ◆ Press CAL key (2.4) when "0.0000" g is displayed.
- ◆ "CAL 0" will be displayed.
- ◆ Press CAL key (2.4) again.
- ◆ The nominal value of the calibration weight will be displayed, "50.0000 g" on the example.
- ◆ Adjust this value until reaching the exact value of the calibration weight to be used, "50.0005" g in the example. For this purpose, press MODE (2.6) or PRINT (2.5) key either to increase or decrease the value initially displayed.

- ◆ Once the calibration weight value has been adjusted, press TARE" (2.3) key.
- ◆ "CAL 0" will be displayed.
- ◆ Press TARE key (2.3) again.
- ◆ "CAL F" will be displayed.
- ◆ Place the calibration weight on the pan.
- ◆ Press TARE key (2.3).
- ◆ "CAL End" will be displayed.
- ◆ The adjusted value "50.0005" g will be displayed.
- ◆ Remove the calibration weight from the pan and "0.0000" g will be displayed; the balance is correctly calibrated and ready to use.

#### **SIMPLE WEIGHING**

Take in mind the balance's temperature if you have moved it. In this case, wait for at least an hour to obtain the plain temperature's stabilization. Calibrate the balance if you consider it is necessary based on the CALIBRATION point exposed before.

- ◆ Turn on the display by pressing ON/OFF key (2.2) with the pan empty
- ◆ When the value "0.0000" g is shown in the screen with stabilized indicator (2.7) the balance is ready to take readings.
- ◆ You can tare the balance if you are going to use a recipient by putting it on the balance, then pressing TARE key (2.3) and waiting until the screen shows "0.0000" g. After the value is stabilized, you can weight.

#### **Note**

**1. The value of subtractive tare is kept in the balance's memory until TARE key (2.3) is pressed once again.**

**2. If you want to weight more than one substance in the same container and you require having their separated weights, press TARE key (2.3) again each time you reach your wanted weight and "0.0000" g will be shown on the screen, so you will obtain the independent/separate weights.**

- ◆ If you need to weight a substance that is in the interior of a recipient, put the full recipient on the balance, then when the weight gets stabilized, press TARE (2.3) and you will read value "0.0000" g on the screen. With the help of a spatula or a similar tool, move the substance you want to weight to another recipient. When you finish moving the substance to the other recipient you will obtain its weight with negative sign.

#### **BELOW WEIGHING**

These balances are provided with a system to measure density and magnetic material's weighing.