MICROSCOPIOS INVESTIGACION RESEARCH MICROSCOPES MICROSCOPES D'INVESTIGACION

Serie/ Series/Série 137

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Tel(948) 310 513 - Fax(948) 310 500 - 312 071

e-mail: asistencia@auxilab.com <http://www.auxilab.es>

Thanks for choosing a professional achromatic optical microscope ZUZI, its high mechanical and optical quality will allow you to have this equipment available for many years. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

In order to use this equipment, it is advisable to follow these points:

- ◆ This manual should be available for all the users of these equipments.
- ◆ Any doubt can be solved getting into touch with the ZUZI Technical Service.
- ◆ As an accuracy instrument, it should be used carefully avoiding the sudden movements and knocks.
- ◆ Never dismantle the different pieces of the microscope to repair yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty.
- ◆ Do not use the equipment with fluids that can give off steams or produce explosive or inflammable mixtures.

ATTENTION!!
IF AN EQUIPMENT IS NOT PROPERLY CLEANED AND DISINFECTED IT WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

CONTENTS

1- USES OF THE INSTRUMENT	17
2- DESCRIPTION	17
3- TECHNICAL SPECIFICATIONS	17
4- INSTALLATION	21
5- INSTRUCTIONS	22
6-USEFUL ADVICES	24
7- MAINTENANCE.....	25
8- TROUBLESHOOTING	26
9- ACCESORIES	27
ANNEX I: WARRANTY INFORMATION	29


1. USES OF THE INSTRUMENT

ZUZI Professional Microscopes Series 137 are very firm equipments designed to provide the highest satisfaction in a wide range of demands. They are equipped with a high quality optics and a wide range of accessories to complete the microscope.

ZUZI research microscopes offers you a wide top of the range accessories specially designed for the most exigent profesional and to make tasks more confortable and efficient. Their advanced and up to date design, as well as their modular configuration, makes of this microscope, an essential tool in long usual tasks, research and teaching.

2. DESCRIPTION

1. Eyepiece
2. Head
3. Objective
4. Mechanical stage
5. Condenser
6. Stage Displacement
7. Collector, illumination system
8. On/Off
9. Light regulator
10. Macrometric Adjusting control
11. Micrometric Adjusting control

Nota: The high range of components of the Microscopes Series 137 allows to form the microscope according yours necessities, we have three heads different, and a lot of variety of objectives, eyepiece, condensers...

3. TECHNICAL SPECIFICATIONS

	Head			Nosepiece		Objectives	Eyepieces
	Monocular	Binocular	Triocular	Quadruple	Quintuple		
137/6	X			X		Achromatic	1WF10X
137/7		X		X		Achromatic	2WF10X
137/8			X	X		Achromatic	2WF10X
137/16	X			X		Semiplan achromatic	1WF10X-1P16X
137/17		X		X		Semiplan achromatic	2WF10X-2P16X
137/18			X	X		Semiplan achromatic	2WF10X-2P16X
137/46	X				X	Plan achromatic	1WF10X-1P16X
137/47		X			X	Plan achromatic	2WF10X-2P16X
137/48			X		X	Plan achromatic	2WF10X-2P16X

MICROSCOPE STAND:

- ◆ Metal, very firm, with coaxial controls in both sides for macro and micrometric focusing. The right control is provided with friction control ring, the left one with a device of quick focusing.

HEAD:

According with your model:

- ◆ Monocular, inclined 45° and rotary 360°. It can be fixed in any position.
- ◆ Binocular, inclined 45° and rotary 360°. Symmetrical interpupillary graduation with scale between 55 and 70mm; and two dioptic equalizing rings.
- ◆ Triocular, it is provided with a binocular head, inclined 45° and rotary 360°. Symmetrical interpupillary graduation with scale between 55 and 70mm, two dioptic equalizing rings and a vertical tube to be adapted to a photographic camera or a video.

FOCUSING PATH:

- ◆ Macro: 25 mm with final stop.
- ◆ Micro: 25 mm, 200 microns in each graduated turn in 2 microns.

NOSEPIECE:

- ◆ Inverted.
- ◆ According with model, quadruple or quintuple

MECHANICAL STAGE:

- ◆ Double lecho
- ◆ Coaxial control
- ◆ Vertical controls
- ◆ Dimensions: 140 x140 mm
- ◆ Traversing area: 70(X) x 40(Y) mm
- ◆ Graduated scales: 0,1 mm

ABBE CONDENSER: double lens.

- ◆ Centrabable focusing, vertically variable by pinion and rack (25 mm)
- ◆ Numerical aperture: 1.25.
- ◆ The Iris Diaphragm is located under the condenser and the stage, and it is charge of regulating the light entrance to the condenser.
- ◆ Movable filter holder. It includes a blue and a green filter (32 mm diameter).

ILLUMINATION:

- ◆ Halogen lamp 6V 30W, power input (220-230V 50 Hz.) incorporated, interruptor and intensity control potentiometer. The electrical components are installed at the rear of the microscope. The lamp replacement is easily done through the lid at the bottom. A little knob helps to access to the compartment where it is the alojen lamp.
- ◆ Precondenser that avoids masive dispersion light

OBJECTIVES:

They consist of a system of lens. We can talk about two kind of objectives:

- ◆ Dry objectives: There is only air between the objective and the preparation.
- ◆ Immersion objectives: It is necessary to place a liquid between the lens and the preparation that allows a higher luminosity

Dry objectives: 4X, 10X, 20X, 40X, 60X

Immersion objective: 100x.


The 40x, 60x and 100x objectives are retractil, they are provided with a device that avoids its break in case that they get in contact with the sample, avoiding the objective gets damaged.

- Achromatic: 4x, 10x, 40x, 100x

- Semiplan achromatic: 4x, 10x, 40x, 100x

- Plan achromatic: 4x, 10x, 20x, 40x, 100x

The characteristics of each objective are codified according to the following numbers:

40: Objective magnification.

0.65: Numerical aperture.

160: Tube length.

0.17: Thickness of the microscope slide.

Oil: indicator of immersion objective.

OBJECTIVES						
MAGNIFICATIONS	4X	10X	20X	40X	60X	100X
WORKING DISTANCE	17.50 mm	7.32 mm	1.26 mm	0.63 mm	0.41 mm	0.19 mm
NUMERICAL APERTURE	0.10	0.25	0.40	0.65	0.85	1.25

We can supply you with Achromatic, Semi plan and Plan Achromatic objectives. The difference between them is that both Semi plan and Plan corrects the spherical aberration caused by the effect of a simple lens' prism as the origin, with a punctual luminose focus, of a caustic cape of beams, not a punctual image. It is corrected by sticking two or three lens of different refraction or dispersión, in a higher grade plan than semiplan objectives.

You can choose one or another objective based on your requirements.

KELLNER EYEPIECES:

They are composed of two lens that are separated by a diaphragm. Its function is to take the image from the objective to the eye.

TYPE / MAGNIFICATION	W.F. / 10X	P / 16X
FOCUS DISTANCE	24.99 mm	15.58 mm
NUMBER OF FIELD	18	11

TOTAL MAGNIFICATIONS:

It is the result of multiplying the eyepiece magnification and the objective magnification.

	Eyepieces	Objectives					
		4x	10x	20x	40x (R)	60x (R)	100x (R)(I)
Total magnification	WF 10x	40	100	200	400	600	1000
Field of view		4.50 mm	1.80 mm	0.90 mm	0.45 mm	0.30 mm	0.18 mm
Total magnification	P 16x	64	160	320	640	960	1600
Field of view		2.75 mm	1.10 mm	0.55 mm	0.27 mm	0.18 mm	0.11 mm

W.F: Wide Field

P: Plan achromatic

R: Retractil

I: Immersion objective

N.A.: It determines the characteristics of the objective. The larger digital aperture makes a brighter and better resolved image.

Working distance: It is the distance (mm) between the preparation and the frontal lens of the objective when the microscope is focussed.

Focal distance: It is the distance (mm) between the main plan image of the system and its focus image.

Resolution: It is the reciprocal value of the separator power, which expresses the minimum distance in which two small particles under the lens can be separated. It will be calculated in lines/mm.

Number of field: It expresses the diameter (mm) of the diaphragm of field that is composed of the eyepiece.


Field of view: It is the size (mm) of the real field we are observing.

4. INSTALLATION

- ◆ After taking out the microscope of the packing and remove the protections, open the lock and take out the inner packing of the wooden box, check that the sign ^UP points up before opening it.

- ◆ Put the microscope stand on a steady table carefully (never hold it by the mechanical stage), keeping it apart from humidity, heat, dust and vibrations. It is advisable to put it on a felt, rubber rug or filter paper in order to make sure of its stability and avoid that it slides.

Picture 1


- ◆ The observant should be in a comfortable position, well sit down and with the back in a straight way. It is advisable to work on a dark table in order to avoid the light that dazzles and reduces the good definition of the images. Pay attention to the light of the enviroment, install the microscope far from the windows.

- ◆ Take down the mechanical stage by means of the macro focusing controls. After taking out the protective caps, screw the objectives into the nosepiece from low magnifying power to higher magnifying power (4x-10x-40x-100x) in clockwise direction (picture 1).

IMPORTANT: Never hold the microscope by the mechanical stage or the tube, since in this way all the weight is supported by the micrometric screw whose notches are worn out slowly.

- ◆ Mount the head adjusting it with the screw that is incorporated. Insert the eyepieces (10x,16x) (Pic. 2). The binocular tube is normally located in the direction of the microscope front, but it can be turned to any other direction, if necessary.

- ◆ These microscopes have a Schuko plug. This plug must be connected to a plugging device connected to ground. The grounded device must be at hand to be able to disconnect the plug quickly in case of emergency.

Picture 2


Neither the manufacturer nor the supplier will assume any responsibility for the damages caused to the microscope or for any other damage suffered by any person due to the neglect of the right process of electric connection. The alternating voltage must be 220-230 v, 50 Hz.

- ◆ Use always the microscope dust cover when it is out of use in order to avoid that the dust settles on the optical parts.

- ◆ Keep the accesories box.

5. INSTRUCTIONS


Picture 3

In those models provided with binocular observation systems, the two observation tubes must be regulated correctly. This regulation will depend on:

- ◆ The interpupillary distance (pic. 3): regulate the interpupillary distance separating or joining the eyepieces. It is achieved when there is a total fusion of the two images.

- ◆ Dioptric compensation of the tubes: You will appreciate a difference in the focusing more or less marked, closing alternatively one eye after the other. To correct this defect you should adjust the dioptric ring to zero. Looking through the right eyepiece you have to focus the sample using the macro and micrometric controls, once you have focused the sample, look through the left eyepiece and then adjust the dioptric adjusting ring until you can obtain a clear bright image. If the observer needs glasses, he must keep his glasses on, taking into account previously that the observation distance is far enough from the lens to avoid that it rubs on the glasses.

When the values are found out it will be very easy to learn them by heart, above all if the microscope is shared with other users. This is advisable in order not to repeat the location of the suitable values each time the microscope is used.

Picture 4

Before placing the preparation on the slide, the slide must get down to a higher distance than the working distance of the lower magnifying objective. Then the preparation is placed on the slide, fitting it in the square and fastening it with the moveable clip. The sample must be centered with respect to the slide aperture.

Looking with the eyepiece and the lower magnifying objective, raise the slide with the macrometric control until the image appears, at this moment operate the micrometric control until obtaining a correct focusing (Pic.4).

IMPORTANT: Starting always the observation with the lower magnifying objective makes the focusing easy, makes impossible that the preparations break or the objectives become dirty. Besides this is essential to regulate the illumination previous to a higher magnifying observation.

Operating the nosepiece it will pass to another higher magnifying objective, and you will only have to retouch the focusing with the micrometric control.


Regulate the height of the condenser to achieve a steady illumination of the field of view (the more magnification the less field), in this way the regulation of the condenser will be more accurate depending on the magnification: At 40x (objective 4x and eyepiece 10x) the condenser will get down to cover a larger field and to avoid an intense illumination. The condenser will have to be raised progressively with the higher objectives to concentrate the light beam and get a better illumination and contrast. With the immersion objective the condenser will have to be raised to its highest height.

The iris diaphragm allows to reduce the aperture to a similar value as the objective, and in this way it can avoid a marginal illumination that decreases the contrast. Therefore it will be closer when a low magnifying objective with small numerical aperture is used. In fact, the total aperture will be only used with the immersion objective.

In order to carry out the adjustment you must start with the maximum aperture, as the contrast increases because the diaphragm is closer, diffracting phenomena appear which limit the resolution. The correct adjustment will be obtained when we have the maximum contrast without any change in the resolution.

An useful way of making the adjustment is the following one:

Take away the eyepiece and observe the objective's exit pupil through the empty tube, adjusting the diaphragm aperture, generally it is convenient to keep the diaphragm aperture about approx. 70-80% of the objective numerical aperture. (Pic.5)

Picture 5

The closing of the diaphragm, reducing the numerical aperture of the optical system, increases the focal depth or penetration power. The diffraction limit can not be surpassed under pretext of increasing the focal depth.

The diaphragm must not be used to reduce the intensity of the illumination.

6. USEFUL ADVICES

The head of all models of microscopes serie 137 is rotary, then is is advisable the obser-
vation through the frontal part of the plate, since it makes easy the approach to the mechan-
ical controls.


When you use the immersion objective you have to take into account the following points

◆ To use the numerical aperture of the immersion objective, both the objective and the pre-
paration will be immersed in oil, then you have to follow this process:

- focus the preparation with a low magnifying objective (4x).
- put a drop of immersion oil in the cover glass.

rotate the nosepiece and place the immersion objective, adjust the focus with the micro-
metric control.


◆ You should avoid the appearance of bubbles in the oil film, since it produces a deteriora-
tion in the resulting image.

◆ It is advisable to dry the frontal lens after each observa-
tion. **Never use a yellow duster!!** It is better to dry it with a
special paper for the cleaning of glasses or with a piece of
linen cloth. The objective must not be immersed in any kind
of liquid, since it can bring about the release of the lens.

Picture 6

Depending on the observant the tension of the macrome-
tric focal controls can be modified, since it is provided with a
tension adjustment ring that is placed next to the right macrometric control of the microscop-
e (Pic.6). However, you must avoid to loosen it excessively since this could produce the fall
of the plate together with the problems of focusing accuracy.

There is also a device for the control of the focal limit. All
models of this serie are provided with a ring that avoids that
the plate raises beyond a point, that can be determined once
the sample has been focussed, and this fact can avoid possible
accidental collisions during the observation (pic. 7). This
mechanism also allows the possibility of getting a quick focus-
ing, acelerating the observation process when you are work-
ing with different samples.

Picture 7

Causes of a defective image

In case of a defective image you must check that:

- ◆ The illumination is well made, the light intensity is not excessive enough or too weak.
- ◆ Never regulate it with the diaphragm of the condenser. The condenser and the lamp must be well centered.
- ◆ The eyepieces must fit well and the objectives must be well tightened.
- ◆ The cleaning of the whole optical system, following this process: make the eyepieces rotate and if the specks move, clean them. Make the whole of the head rotate. *Never dismatle the head.* However, you can clean carefully those accesible surfaces of the prisms blowing with a rubber pear. Rotate the objective, if the parasitic images turn around clean it with a small dry brush trying to remove the dust, observe the frontal surface with a magnifying glass or an inverted eyepiece.
- ◆ The immersion oil must be enough and it must not contain bubbles or impurities.
- ◆ The microscope slides, microscope cover glasses and assembly medium must not be too thick so that they do not prevent the focusing at medium or high magnifications. There are standard dimensions for both microscope slides and cover glasses.

7. MAINTENANCE

- ◆ All the operating methods we have exposed previously will not have any value if you do not carry out a continuous maintenance.
- ◆ The maintenance of these microscopes is the same as any other optical instrument. They must be always clean and covered with the dust cover to avoid that the dust settles on them.
- ◆ The user must never dismantle the lens, if there would be any dirtiness on the outside surfaces of the lens, clean them with a soft cloth without fluff.
- ◆ To eliminate the dust settled on the lens, blow with a rubber pear or clean it with a small soft brush.
- ◆ For the cleaning of the mechanical pieces use non-corrosive lubricants, so that the paint does not become damaged. Be careful not to touch the optical pieces. Use always a neutral detergent.
- ◆ If the lamp burns out, change it with an URA halogen lamp 6V 30W, being careful not to touch the bulb with your hands.

IMPORTANT: Before changing the lamp bulb or the fuse, make sure that the microscope is disconnected. Make also sure of using lamps 6V 30W, since the use of any other kind of lamps will produce a defective functioning of the instrument.

- ◆ For the change of lamps these microscopes are provided with a device at the bottom of its base; loosen the screw and draw it out, loosen the screws of the lamp holder and change it, holding it firmly again. Center the lamp holder again, if necessary.
- ◆ If it is necessary to change the fuse you must follow this process: Loosen the protective lid that is on the base of the microscope (*fuse*), change the fuse and place it in the same position, putting the protective lid again.


8. TROUBLESHOOTING

<i>PROBLEM</i>	<i>SOLUTION</i>
The lamp doesn't switch on	-Check that there is a good power cord connection... -Check that the fuse and the lamp are in good state.
The viewfield is uneven.	- Check that when you revolve the nosepiece the objective is centered in the optical path correctly. -Check the correct centering of condenser. This will be done with three small screws located round the condenser harness.
Dirt or dust in the view-field.	- Check dirt or dust on the lens (pre-condenser eyepiece, field lens.) Once the problem is located,clean the sample.
Image shifts during focus.	-Check the friction ring tension, it may be too weak and creates an involuntary fall down of the slide.
Part of a viewfield is out of focus.	- Check the objective is correctly located in the optical path. If not, revolve to clic-stop position..Check the sample is correctly located-in the slide, make sure it's firmly held by the clip.
If the frontal lens of the objective is touching the slide during focus or when changing from a low power objective to a high power objective.	- Check if the slide is upsidedown (the microscope slide is over the coverglass). Proceed to the correct position if necessary. The coverglass has a specific thickness (0.17)
Image not focusable when changing slides.	-Check that focus control device is not on. If it's on, it won't let move macrometrical adjustment as needed.

9. ACCESSORIES

The research series can be completed with several Zuzi accessories that increase the field of applications.

ENGLISH

<i>Code</i>	<i>Description</i>
90137901	Quadruple nosepiece 20W.
90137902	Quadruple nosepiece 30 W.
90137903	Quintuple nosepiece 20 W.
90137904	Quintuple nosepiece 30 W.
90137911	Monocular head
90137912	Compensation free Binocular head
90137913	Compensation free Triocular head
90100590	Double layer mechanical
90100125	Wide Field eyepiece WF10x
90100135	Wide Field eyepiece WF 10x with pointer
90100175	Micrometer eyepiece WF 10x
90100126	Eyepiece P 16x
90137167	High point eyepiece WF 10x 18.5
90137115	Wide Field eyepiece WF10x 21 mm.
90100193	Photography eyepiece 4x.
90100201	Achromatic objectives 4x.
90100202	Achromatic objectives 10x.
90100203	Achromatic objectives 20x.
90100205	Achromatic objectives 40x(R).
90100206	Achromatic objectives 60x(R).
90100209	Achromatic objectives 100x(R)(I)
90100251	Semi- plan achromatic objectives 4x.
90100252	Semi- plan achromatic objectives 10x.
90100255	Semi- plan achromatic objectives 40x (R).
90100259	Semi- plan achromatic objectives 100x (R) (I).
90100241	Plan achromatic objectives 4x.
90100242	Plan achromatic objectives 10x.
90100243	Plan achromatic objectives 20x (R).
90100245	Plan achromatic objectives 40x(R).
90100249	Plan achromatic objectives 100x(R)(I)
90137210	Objective 100x with dark field diaphragm
90100026	Condensador de campo claro con diafragma iris
90137024	Dark field condenser Dry (A:N:=0.9)
90137026	Dark field condenser Oil (A.N.= 1.25-1.36)
90137421-2	Green or blue filters for condenser
90137001	Collector
90137011	Collector, Iris diaphragm and Koehler illumination
90137417	Photography attachment
90137419	Video attachment without lens
90137420	Video attachment 0.5x
90137430	Basic polarization set
90100022	Single annular phase contrast kit
90137022	Turret annular phase contrast kit
90137030	Drawing attachment
90137035	Epi Fluorescent attachment
90100981	Wooden box

WARRANTY

AUXILAB, S.L. warrants this equipment be free from defects in material and workmanship for a period of one year from date of purchase from AUXILAB S.L. or an authorised AUXILAB S.L. dealer. This warranty is limited to such purchaser and is not transferable.

During the year warranty period, AUXILAB S.L. will repair or replace, at its option, any defective equipments or pieces at no additional charge, provided that the product is returned, shipping prepaid, to AUXILAB S.L., an authorised dealer or an authorised service location. The purchaser is responsible for insuring any returned product and assumes the risk of loss during shipping. All replaced parts become the property of AUXILAB.

This limited warranty does not extend to any products which have been damaged as a result of accident, misuse, abuse or as a result of service or modification by anyone other than AUXILAB, an authorised dealer or an authorised service location.

Except as expressly set forth above, no other warranties are expressed or implied, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. AUXILAB expressly disclaims all expressed or implied warranties not stated herein.

In order to obtain a free reparation or a free piece-replacement by our Technical Service over the above-mentioned period of time, it is absolutely necessary to submit the Warranty Certificate. This warranty doesn't apply for the maintenance of the equipment or for damaged caused by a natural, everyday use.

Some states do not allow the exclusion or limitation of incidental or consequential damages for consumer products, and some states do not allow limitations on how long an implied warranty lasts, so the above limitations or exclusions may not apply to you. This limited warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

Dealer:			
Equipment:		Model:	
Serial Number:			
Date:			
Brief failure description:			

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