ASPI-PODO 1 E 2

ASPI-PODO 1 AND 2

ASPI-PODO 1 ET 2

ASPI-PODO 1 UND 2

ASPI-PODO 1 Y 2

- MANUALE ISTRUZIONI
- INSTRUCTIONS HANDBOOK
- MANUEL D'INSTRUCTIONS
- GEBRAUCHSANWEISUNG
- MANUAL DE INSTRUCCIONES

CATTANI S.p.A.

ASPI-PODO 1 E 2

INDICE

Pagina

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ASPI-PODO 1 AND 2

INDEX

	<i>Page</i>
— GENERAL RUNNING DATA	8
— LIST OF COMPONENTS	9
— GENERAL FEATURES	10
— STARTING THE ASPIRATION UNITS	10
— SIGNALS	10
— RECOMMENDATION	11
— FUNCTIONING	11
— ORDINARY MAINTENANCE OPERATIONS	11
— IMPORTANT NOTICE	11
— TRANSPORT AND STORAGE	12

General running data

Podiatric aspirator

Model:	Aspi-Podo 1 and 2
Rated voltage:	220 V~ ±5%
Rated frequency:	50 Hz
Rated current:	3.1 A
Insulation class:	I
Type of appliance:	B
Use:	continuous service
Protection against liquids:	common
Level of protection against direct or indirect contact type:	B 

	Alternating current	IEC 417-5032
	Earthing	IEC 417-5019
	Type B appliance	IEC 878-02-02
	Off	IEC 417-5008
	On	IEC 417-5007

This appliance cannot work in contact with a flammable anaesthetic mixture with air, oxygen or nitrous oxide.

Motor thermal protector.

Output power: 0.4 kW

Maximum flow: 1250 l/min

Maximum operating head for continuous service: 1300 mm H₂O

Aspi-Podo 1 sound pressure level: 64.5 dB (A)*

Aspi-Podo 2 sound pressure level: 67.5 dB (A)*

Other available voltages:

240 V~ 50 Hz 2.95 A

220 V~ 60 Hz 3.5 A

120 V~ 60 Hz 6.0 A

110 V~ 60 Hz 7.0 A

100 V~ 60 Hz 7.0 A

**Sound pressure level tested (with canalized air) according to ISO 3746-1979 (E) regulation.
Parameters: r or d=1 - background noise ≤38 dB (A) - instrument: Brüel & Kjær Type 2232.*

Manufactured by CATTANI S.P.A. - Parma - Italy

LIST OF COMPONENTS (drw. pages 30-31-32-33-34-35)

- 1 - Aspiration arm**
- 2 - Flexible tube (only for Aspi-Podo 2)**
- 3 - Micro-motor seat (only for Aspi-Podo 2)**
- 4 - Aspiration inlet for micro-motor**
- 5 - Clogged bag bip**
- 6 - Clogged bag indicator**
- 7 - Dust filter**
- 8 - Basket filter container**
- 9 - Bag container**
- 10 - Slide-bracket**
- 11 - Support with hinge and slide-bracket
(only for Aspi-Podo 1)**
- 12 - Uni-Jet 75 aspiration unit**
- 13 - Safety switch**
- 14 - Electrical circuit**
- 15 - Ventilation grid**
- 16 - On/Off switch**
- 17 - Silencer for exhausted air**
- 18 - Active carbon**
- 19* - Active carbon container (only for Aspi-Podo 2)**
- 20* - Active carbon container (only for Aspi-Podo 1)**

General features (fig. 5 page 32)

Aspi Podo 1 and 2 cabinets are fitted with a reliable and silent asynchronous aspiration unit. Aspiration can be performed both with the micromotor (det. A fig. 5) and the arm (det. B fig. 5) at the same time.

Once the job is finished, before the patient leaves the chair, the hose enables the podiatrist to aspirate any dust remaining on the floor (det. C fig. 5).

These appliances are mobile and can easily be drawn close to the chair or moved away from it.

You can use plain paper filter bags or long-life BIA-certified material (det. D fig. 5) filter bags. The appliance is not suitable for aspiration of liquids or damp substances.

Starting the aspiration unit (fig. 5 page 32)

This appliance has been designed according to C.E.I. reg. 62-5 1991 (2nd edition), and must be connected to the feeding line according to C.E.I. reg. 64-8 for class I equipment. The mains must be protected against direct and indirect contact, against overloading and overcurrent according to C.E.I. reg. 64-8.

The winding of the motor has a thermal device which cuts off the power supply when electrical absorption and temperature are above supported level.

If you open the door when the appliance is in use, the safety switch (det. 13 fig. 5) cuts off the power supply and aspiration stops; the appliance starts working again as you close the door.

Signals

Electrical shock risk: also 230 V ~ can be lethal

Biological danger, risk of infections from contagious diseases (fungi, warts, etc...)

High temperature

Compulsory direction of flow or of rotation

Signs cannot always fully express danger warnings, therefore it is necessary that the user reads the warnings and keeps them in due consideration.

Failure to respect a danger sign or warning may harm operator or patients.

Do not remove protections; do not tamper with machines or their operation.

Despite all our efforts, it is still possible that danger warnings are not exhaustive: we apologise with the users and kindly request them to care for all danger sources that might have pass unnoticed and to inform us accordingly.

On the web site: **www.cattani.it** you can find our **up-dated** instructions manuals.

We recommend reading them, particularly for the **safety** updating.

Recommendations

The retailer or the installation engineer will take care to train the surgery staff with trials on a brand new, non-contaminated appliance.

The installation of the aspirator is reserved to engineers authorized by the manufacturer. Aspirated debris is always contaminated and infected: for this reason, we stress that the greatest care must be used to prevent contamination of operators or the environment. Contamination may also result from an appliance in bad working order so we recommend to contact only dental engineers whose teaching and training is duly certified by the manufacturer. Any modification of the appliance must be agreed upon with the manufacturer.

Functioning (fig. 1-2-3-4 pag. 31-32-33-34)

To start aspiration press main switch (connection to the power supply) and unit switch (det. 16 fig. 1) (I=ON / 0=OFF).

Aspirated particles go through the arm (det. 1 fig. 1) and hose of the micro-motor (det. 3 fig. 1) and are collected in the bag filter (det. 9 fig. 2), while the air passes through the active carbon (det. 18 fig. 3-4) and is conveyed then to the aspiration unit (det. 12 fig. 2) where it is exhausted outside (det. 17 fig. 1).

Ordinary maintenance operations

In order to ensure correct functioning of the aspiration unit, it is advisable to check the plug and the power supply cable periodically: any damage can be dangerous.

The bag filter must be replaced every day in order to prevent aspirated particles from putrefying and giving off bad smells inside the appliance.

In case the filter should not need to be emptied, in order to avoid adulteration of aspirated particles you can spray some Eco-Jet (surface disinfectant) directly into the bag; the disinfecting action of Eco-Jet will freshen the exhausted air, making it more pleasant.

The active carbon will be replaced every three months, proportionally to the work carried out and as soon as a decrease of aspiration power is noticed.

Important notice

- The appliance is guaranteed for one year from date of sale, provided that guarantee card addressed to manufacturer is returned to the manufacturer reporting date of sale, retailer stamp and customer's name.
- Guarantee and manufacturer liability cease in case appliances and/or plants are treated with unsuitable products, not suggested by the manufacturer and improperly used and in case they are found tampered by any kind of action performed by unable and thus unauthorised people.
- The manufacturer, together with distributors, agents, authorised retailers and service engineers, is willing to supply spare parts, literature, instructions and any other useful information.
- Technical requirements, improvement of the product, legislation and logistic problems, availability of products or semi-manufactured products can lead to modifications without notice.
- On the web site: www.cattani.it you can find our **up-dated** instructions manuals. We recommend reading them, particularly for the **safety** updating.

Transport and storage

- *Packed equipment can be transported and stored at a temperature ranging from - 10 °C to + 60 °C.*
- *Packages must be kept away from water and splashing and cannot tolerate humidity >70%.*
- *Packages with the same weight can be stored in piles of three only.*

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