**IGF-I (Insulin-Like Growth Factor I)**  
**<sup>125</sup>I RIA Kit**

For the quantitative determination of IGF-I  
(Somatomedin C) in human serum

**Instruction Manual**

Manuel d'Instructions  
Testanleitung  
Manual de instruções  
Bruksanvisning

**Catalog No./REF./KAT.-NR./ N° de Catálogo/  
Katalognummer: 53065**

**DiaSorin**

**Stillwater, Minnesota 55082-0285, U.S.A.**

**TABLE OF CONTENTS**

English .....	Page 1
Français .....	13
Deutsch.....	26
Português.....	39
Svenska .....	52

## INSULIN-LIKE GROWTH FACTOR I RIA KIT

### 1. INTENDED USE

FOR IN VITRO DIAGNOSTIC USE.

This kit contains instructions and materials for the quantitative determination of IGF-I in human serum by radioimmunoassay (RIA). The procedure is to be used as an adjunct procedure to hGH determination in the diagnosis of hypopituitarism and acromegaly.

### 2. SUMMARY AND EXPLANATION

Somatomedin C is a 70 amino acid peptide which is part of a family of growth factors which includes insulin-like growth factors I and II, somatomedins A and B, and MSA. Somatomedin C was recently found to be identical in structure to IGF-I<sup>4</sup> whose sequence was published in 1978 by Rinderknecht and Humbel.<sup>6</sup> There are similarities between proinsulin and somatomedin in that both contain A and B peptides with a C-peptide linked and both have internal disulfides. IGF-I is highly basic in nature due to its amino acid composition and circulates bound to high molecular weight proteins of approximately 140,000 daltons.

IGF-I mediates the growth promoting actions of growth hormone. It is released from a variety of tissues in the body in response to growth hormone and is under the control of a feedback mechanism forming a tightly controlled loop between the hypothalamus and the pituitary gland.<sup>8</sup>

Other than growth hormone, two of the most important factors influencing IGF-I concentration in the serum are age and sex. IGF-I levels are lowest at birth and rise throughout infancy and childhood. The highest levels are attained during adolescence and decline slowly during adulthood. Females exhibit slightly higher IGF-I levels than males during adolescence and adult years.<sup>1,10</sup>

IGF-I measurement has been advocated as a screening and management tool in growth hormone deficient children.<sup>10</sup> Its use in diagnosis along with growth hormone measurements<sup>9</sup> or as a tool to assess a child's response to administered growth hormone<sup>4</sup> has led it to a prominent place in the endocrine laboratory, particularly when dealing with growth disorders.<sup>2,7</sup>

A second major use of IGF-I measurement is in the diagnosis and treatment of acromegaly.<sup>5</sup> IGF-I levels may be helpful to assess the results of bromocryptine treatment of acromegaly.<sup>9</sup>

### 3. PRINCIPLE OF THE ASSAY

The DiaSorin IGF-I RIA is a double antibody disequilibrium assay which includes an ODS-silica extraction procedure for serum samples. After the extraction procedure, the RIA is performed employing addition of sample and rabbit anti-IGF-I, followed by a 2 hour incubation at 2-8°C. Iodine-125 IGF-I is then added followed by a second incubation for 20 hours at 2-8°C. Pre-precipitated carrier, second antibody and polyethylene glycol are added in a single step. The assay can be centrifuged and decanted after the 2 hour second antibody incubation at 2-8°C. IGF-I from many other mammalian sources can also be measured with this procedure.

### 4. REAGENTS PROVIDED IN THE KIT

IgF-1 Calibrator 0	1 vial/20 mL
IgF-1 Calibrators 1-5	5 vials/0.5 mL
IgF-1 Antiserum	1 vial/14 mL
<sup>125</sup> I IgF-1 Tracer	1 vial/14 mL
IgF-1 Controls	2 vials/ 1.0 mL
IgF-1 Precipitating Complex	1 vial/ 35 mL
ODS-Silica	1 pkg./25 columns
Number of tests	65

STORAGE: Upon receipt, the kit should be stored at -15°C or lower. After reconstitution store each reagent at -15°C or lower until the expiration date on the label. Reagents should not be used past the expiration date.

When reconstituting the contents of the vials, mix gently to avoid foaming. Reagents from different batches must not be mixed.

**4.1 IGF-I calibrator 0:** lyophilized reagent

BSA-borate buffer contains 0.25% sodium azide as a preservative. Reconstitute the vial with 20 mL of purified water and allow it to stand for 15 -20 minutes until the contents are completely dissolved; mix thoroughly before using.

**4.2 IGF-I Calibrators (1-5):** lyophilized reagent

Five human intact IGF-I calibrators at nominal concentrations of 3.0-75 nmol/L (0.12-3.0 units/mL) are prediluted in BSA-borate buffer containing 0.1% sodium azide as a preservative. Exact values are assigned with each lot. Reconstitute each vial with 0.5 mL of purified water and allow the vials to stand for 15 - 20 minutes until the contents are completely dissolved; mix thoroughly before using. The DiaSorin IGF-I calibrators have been calibrated against the World Health Organization (W.H.O.) standard NBSB 87/518. Any comparisons with other products or procedures should be done with this reference preparation. To approximate the units used by many commercial laboratories, a value of 25 nmol/L equal to 1 unit/mL has been assigned. Nanomoles/liter of IGF-I is converted to ng/mL by the multiplier 7.7. The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as recommended.

**4.3 IGF-I Antiserum:** lyophilized reagent

Rabbit anti-IGF-I serum is diluted in BSA-borate buffer with heparin containing 0.002% thimerosal. Reconstitute the vial with 14 mL of purified water and allow it to stand for 15 - 20 minutes until the contents are completely dissolved; mix thoroughly before using.

**4.4 <sup>125</sup>I IGF-I:** lyophilized reagent

Human synthetic IGF-I (53-70) is labeled with iodine-125 and is diluted in BSA-borate-EDTA buffer which contains 0.008% thimerosal as a preservative. Reconstitute the vial with 14 mL of purified water and allow it to stand for 15 - 20 minutes until the contents are completely dissolved; mix thoroughly before using.

**4.5 IGF-I Controls (Level 1 and 2):** lyophilized reagent

Human serum is spiked, if necessary, with the appropriate amount of intact human IGF-I to obtain a concentration within a specified range. 0.02% sodium azide is added as a preservative. Reconstitute the vials with 1.0 mL of purified water and allow it to stand for 15 - 20 minutes until the contents are completely dissolved; mix thoroughly and treat the quality control serum as an unknown sample. Ranges are printed on the control vials.

**4.6 IGF-I Precipitating Complex:** lyophilized reagent

Normal rabbit serum pre-precipitated with goat anti-rabbit serum and polyethylene glycol (PEG), is diluted in BSA-borate buffer with 0.03% thimerosal added. Reconstitute the vial with 35 mL of purified water; mix THOROUGHLY until the suspension appears homogeneous and then allow the vial to stand for a minimum of 30 minutes at room temperature with occasional mixing.

**4.7 ODS-Silica:** ready to use reagent

Plastic columns contain octadecasilyl-silica (Sep-Paks). Columns must be "primed" before applying sample (see Extraction Procedure section). The ODS-silica columns may be stored at room temperature. Columns are stable for 3 years after package is opened and indefinitely if unopened.

**5. WARNINGS AND PRECAUTIONS**

**FOR *IN VITRO* DIAGNOSTIC USE.**

Not for internal or external use in humans or animals. The procedure is to be used as an adjunct procedure to hGH determination in the diagnosis of hypopituitarism and acromegaly.

**REAGENTS CONTAINING HUMAN SOURCE MATERIAL**

**Treat as potentially infectious.**

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV, and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing

human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4<sup>th</sup> ed., May 1999 or current edition.

#### **REAGENTS CONTAINING SODIUM AZIDE**

**CAUTION:** Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, U.S. 1976.

#### **European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)**

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

#### **REAGENTS CONTAINING THIMEROSAL**

Some reagents in this kit contain thimerosal which contains a mercury compound. Disposal of elemental mercury, inorganic mercury, mercury oxides and mercury compounds should be done in strict compliance with all local, state, and federal regulations.

**WARNING:** This product contains a chemical known to the State of California to cause birth defects or other reproductive harm.

#### **REAGENTS CONTAINING IODINE-125**

This kit contains radioactive material which does not exceed 2  $\mu$ Ci (74 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

**WARNING:** This product contains a chemical known to the State of California to cause cancer.

**ATTENTION:** Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

#### **6. INDICATIONS OF POSSIBLE DETERIORATION OF KIT REAGENTS**

- 6.1 The presence of abnormal particulate matter in any of the reagents.
- 6.2 A shift in the slope or position of the calibrator curve from what is normally obtained.
- 6.3 A decrease in maximum binding.
- 6.4 A high nonspecific binding.

## 7. SPECIMEN COLLECTION AND PREPARATION

Two hundred and fifty microliters of serum are required for the assay.

Collect blood by venipuncture in 5 or 10 mL evacuated glass tube. Allow the blood to clot at room temperature. Centrifuge for 15 minutes using  $760 \times g^*$  to obtain hemolysis free serum. No additives or preservatives are required to maintain the integrity of the sample. All plastics, glassware or other material coming into contact with the specimen should be entirely free of any contamination.

EDTA or heparinized plasma may cause slightly lower results.

The serum should be promptly separated from the cells and stored at  $-15^{\circ}\text{C}$  or lower. For long term storage, samples should be frozen at  $-15^{\circ}\text{C}$  or lower. Specimens can be frozen and thawed several times.

## 8. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- 8.1 Disposable borosilicate glass tubes, 12 x 75 mm.
- 8.2 Temperature controlled centrifuge to accommodate 12 x 75 mm tubes.
- 8.3 Gamma scintillation counter capable of counting iodine-125.
- 8.4 Vortex.
- 8.5 Pipetting devices:
  - a. Micropipettors calibrated to deliver 50  $\mu\text{L}$ , 250  $\mu\text{L}$  and 500  $\mu\text{L}$ .
  - b. Repeating dispensers, calibrated to deliver 200  $\mu\text{L}$  and 500  $\mu\text{L}$ .

### Materials and Reagents Required for ODS-silica

#### Column Extraction

- 8.6 12 mL disposable syringes.
- 8.7 Isopropyl alcohol.
- 8.8 Methanol, HPLC grade or better.
- 8.9 4% acetic acid.
- 8.10 0.5 N hydrochloric acid.
- 8.11 Compressed air supply.
- 8.12  $37^{\circ}\text{C}$  water bath.
- 8.13 16 x 100 mm glass tubes.

## 9. EXTRACTION OF SERUM WITH ODS-SILICA COLUMNS

- 9.1 Pipette 250  $\mu\text{L}$  of controls and unknown samples into a 12 x 75 mm glass test tube and add 1 mL of 0.5 N HCl. Vortex gently until well mixed.
- 9.2 Attach each ODS-silica column to a syringe hub.
- 9.3 Add 5 mL of reagent grade isopropyl alcohol to each syringe.
- 9.4 Press through the syringe using the syringe plunger.
- 9.5 Remove the ODS-silica column. (Remove the ODS-silica column from the end of the syringe before removing the plunger to avoid contaminating the syringe contents with the fluid in the ODS-silica column.)
- 9.6 Remove the syringe plunger from the syringes.
- 9.7 Replace the ODS-silica column on the syringes.
- 9.8 Add 5 mL of HPLC grade methanol.
- 9.9 Repeat steps 4 through 7.
- 9.10 Add 5 mL of 4% acetic acid.
- 9.11 Repeat steps 4 through 7.
- 9.12 Repeat steps 10 and 11.
- 9.13 Apply the sample from step 1. Push slowly through the cartridge at a flow rate of 1 mL per 2-3 minutes.
- 9.14 Repeat steps 5 through 7.
- 9.15 Wash the ODS-silica column 2 times with 10 mL of 4% acetic acid repeating steps 5 through 7 in between washes.

$$*g = (1118 \times 10^{-6}) (\text{radius in cm}) (\text{rpm})^2$$

- 9.16 Collect the eluate in 16 x 100 mm glass tubes by eluting the peptide with 4 mL of HPLC grade methanol. Push slowly so that the methanol remains in contact with the ODS-silica column for a minimum of 3 minutes.
- 9.17 Evaporate the methanol to dryness in a 37°C water bath using compressed air or nitrogen.
- 9.18 Reconstitute the dried samples with 500 µL of calibrator 0 and incubate at 37°C for 10 minutes to ensure complete reconstitution (vortex intermittently during incubation). A sample that is expected to be elevated should be reconstituted with 1.0 mL of calibrator 0 at this time. The final result should be corrected by the appropriate dilution factor (see Results section for calculation).
- 9.19 Assay 50 µL of the serum extracts in duplicate in the radioimmunoassay.

## 10. ASSAY PROCEDURE

- 10.1 Reconstitute the lyophilized reagents and allow any frozen reagents to thaw completely. Do not allow reagents to reach temperatures above 20-25°C. Mix all reagents gently before using.
- 10.2 Set up labeled 12 x 75 mm tubes in duplicate according to the Scheme of the Assay, on the last page. All volumes are in microliters.
- 10.3 Add reagents to the tubes as follows:
- Total count tubes**  
Set aside until step 5
  - Nonspecific binding (NSB)**  
50 µL of calibrator 0
  - IGF-I calibrator 0**  
50 µL of calibrator 0  
200 µL of IGF-I Antiserum
  - IGF-I calibrators (1-5)**  
50 µL of IGF-I calibrator  
200 µL of IGF-I Antiserum
  - Controls and unknown samples**  
50 µL of serum extracts  
200 µL of IGF-I Antiserum
- 10.4 Incubate for 2 hours ( $\pm$ 15 minutes) at 2-8°C.
- 10.5 Add 200 µL of 125I IGF-I to all tubes.
- 10.6 Vortex the tubes gently without foaming and incubate for 20 hours ( $\pm$ 1 hour) at 2-8°C.
- 10.7 Vigorously mix the precipitating complex; add 500 µL to all tubes except the total count tubes.
- 10.8 Vortex the tubes gently without foaming and incubate for 2 hours ( $\pm$ 15 minutes) at 2-8°C.
- 10.9 Centrifuge the tubes for 20 minutes using a minimum of 760 x g\* at 20-25°C.
- 10.10 Immediately decant the supernatant from all the tubes except the total count tubes by inverting them for a minimum of 2 minutes. Blot the tubes with absorbent paper to remove any drops of supernatant that may be remaining on the rims before turning the tubes upright.
- 10.11 Using a gamma scintillation counter, count the precipitate of each tube and the total count tubes for a sufficient time to achieve statistical accuracy (see Limitations of the Procedure section).

## 11. PROCEDURAL COMMENTS

- 11.1 Add each aliquot of reagent to the lower third of the assay tube to ensure complete mixture of reagents.
- 11.2 Some manufacturer's disposable tubes yield elevated nonspecific bindings.
- 11.3 If you choose to aspirate the supernatant from the precipitate, be careful not to disturb the precipitate.
- 11.4 To completely monitor the consistent performance of an RIA there are additional factors which may be checked. DiaSorin suggests a check of the following

parameters to assure consistent kit performance.

**a. Total Counts**

**b. Maximum Binding**

Average counts per minute (CPM) of calibrator 0 Tube / Average CPM of Total Count Tubes.

**c. Nonspecific Binding**

Average CPM of NSB Tube / Average CPM of Total Count Tubes.

**d. Slope of Calibrator Curve**

For example, monitor the 80, 50 and 20% suppression points of the calibrator line.

**12. QUALITY CONTROL**

Each laboratory should include at least two control samples in every assay to ensure the validity of each assay's results. A mean and standard deviation should then be determined for each control using a minimum of ten (10) assays. An acceptable range of values may then be obtained for these controls using  $\pm 2$  standard deviations of the values previously determined. The DiaSorin Quality Control Laboratory has determined a range for the controls included in this kit.

**13. CALCULATION OF RESULTS**

There are many methods in existence for calculating results of RIAs. Each is based on obtaining a calibration curve by plotting the extent of binding against stated concentrations of the calibration calibrators. This graph may be either linear or logarithmic scale.

Each of these methods gives essentially the same values for controls and samples, although certain assays may "fit" better into one particular method versus another. The calculation method for the DiaSorin Quality Control Laboratory is % B/B<sub>0</sub> versus log concentration.

**13.1** Calculate the average CPM for each calibrator, control and unknown sample.

**13.2** Subtract the average CPM of the NSB tubes from all counts.

**13.3** Divide the corrected CPM of each calibrator, control or unknown sample by the corrected CPM of the calibrator 0.

$$B/B_0 (\%) = \frac{\text{CPM of Calibrator or Unknown Sample} - \text{CPM of NSB}}{\text{CPM of calibrator 0} - \text{CPM of NSB}} \times 100$$

**13.4** Using 2 cycle semi-log or log-logit graph paper, plot percent B/B<sub>0</sub> for the IGF-I calibrators (vertical axis) versus the concentration (horizontal axis).

**13.5** Draw a best-fit line through the points.

**13.6** Interpolate the levels of IGF-I in the unknown samples from the plot.

**13.7** Correct for appropriate reconstitution volume. For example: If 250  $\mu$ L of serum was extracted and reconstituted with 500  $\mu$ L of calibrator 0, multiply the value interpolated from the curve by 2. If 250  $\mu$ L of serum was extracted and reconstituted with 1.0 mL of calibrator 0, multiply the value interpolated from the curve by 4.

**13.8** Calculate maximum binding by dividing CPM of calibrator 0 by the average total counts obtained in the total count tubes.


\*g = (1118 x 10<sup>-6</sup>) (radius in cm) (rpm)<sup>2</sup>

**TABLE II**  
DiaSorin IGF-I Sample Data

Tube	Duplicate CPM	Average CPM	Corrected CPM	% Bound (B/T)	% (B/B <sub>0</sub> )	Graph Conc. (nmol/L)	Final Conc. (nmol/L)
Total	10,830 10,991	10,991					
NSB	351 338	345		3.2			
Calibrator 0	4,060 4,055	4,058	3,713	37.2			
Calibrators (nmol/L)							
1 (2.9)	3,332 3,264	3,298	2,953		79.5		
2 (7.3)	2,852 2,683	2,673	2,328		62.7		
3 (16.4)	2,032 2,059	2,046	1,701		45.8		
4 (29.9)	1,685 1,634	1,660	1,315		35.4		
5 (82.5)	991 964	978	633		17.0		
Patient Unknowns							
1	2,564 2,590	2,577	2,232		60.1	8.3	16.6
2	1,701 1,723	1,712	1,367		36.8	27.6	55.2

Typical sample data and a calibrator curve are shown in TABLE II and FIGURE 1; this curve is for reference only and should not be used for the calculation of any value.

**FIGURE 1**

## REDUCTION DATA

The DiaSorin QC lab uses a smoothed spline curve fit.

### 14. LIMITATIONS OF THE PROCEDURE

- 14.1** Any unknown sample whose value reads greater than the highest kit calibrator must be appropriately diluted from the ODS-silica column extract with calibrator 0 and assayed again. Correct result using the appropriate dilution factor.
- 14.2** Counting times should be sufficient to prevent the introduction of error due to statistical error in counting (for example, accumulation of 2,000 CPM will yield 5% error; 10,000 CPM will yield 1% error).
- 14.3** IGF-I is elevated in pregnancy and the values were found to rise with gestation age. Other conditions are known to reduce levels of IGF-I, including poor nutritional status, fasting or hypothyroidism. Hepatitis can also influence the levels of IGF-I. In the diagnosis of hypopituitarism, provocative tests for human growth hormone should also be used when IGF-I results are low. A normal IGF-I result is very unlikely in the diagnosis of hypopituitarism except in the case of hyperprolactinemia or craniopharyngioma. No drugs at this time are found to interfere with the IGF-I assay.<sup>9</sup>

### 15. EXPECTED VALUES

#### Normal Range

Normal ranges should be verified by each laboratory. A large number of normal samples (geographic origin: Minnesota, Israel, Maryland and Australia), were studied in the development of our IGF-I assay and we have found that IGF-I values vary with age and in some age groups with sex; therefore, age and sex match controls should be used where appropriate. TABLE III and FIGURES 2 and 3 show the results of normal values assembled from several laboratories and collated by DiaSorin. Because there is pronounced skewing of the distribution during puberty, it is necessary to use a log transformation statistical model to develop 95% confidence limits in setting children's ranges. This is not necessary for treatment of adult values.

**TABLE III**  
Normal Range Study According to Age and Sex  
**Nanomoles/L**

(years) Age Group	Male			Female		
	Mean	(n)	Range*	Mean	(n)	Range*
1 mo.-1.99	7.55	(16)	4.32-13.18	7.35	(12)	3.42-15.79
2.00-5.99	12.27	(17)	5.70-26.40	11.10	(12)	4.43-27.85
6.00-8.00	15.56	(15)	6.39-36.13	11.89	(13)	6.01-23.51
8.01-10.00	15.34	(13)	5.92-39.75	17.29	(11)	9.16-32.64
10.01-13.99	21.90	(24)	8.65-55.46	32.15	(14)	9.30-111.10
14.00-18.00	42.84	(15)	17.33-105.83	42.68	(15)	12.94-140.80
Adult —	27.55	(36)	9.09-46.01	30.03	(41)	11.63-48.43

\*95% confidence limits log transformation model except for adult.

#### Units/mL

(years) Age Group	Male			Female		
	Mean	(n)	Range*	Mean	(n)	Range*
1 mo.-1.99	0.30	(16)	0.17-0.53	0.29	(12)	0.14-0.63
2.00-5.99	0.49	(17)	0.23-1.06	0.44	(12)	0.17-1.11
6.00-8.00	0.62	(15)	0.26-1.44	0.48	(13)	0.24-0.94
8.01-10.00	0.61	(13)	0.24-1.59	0.69	(11)	0.37-1.31
10.01-13.99	0.88	(24)	0.35-2.22	1.28	(14)	0.37-4.44
14.00-18.00	1.71	(15)	0.69-4.23	1.71	(15)	0.52-5.63
Adult —	1.10	(36)	0.36-1.84	1.20	(41)	0.47-1.94

\* 95% confidence limits log transformation model except for adult.


FIGURE 2


FIGURE 3

**16. SPECIFIC PERFORMANCE CHARACTERISTICS**

**16.1 Precision**

Within Assay Variation (values = nmol/L)

	Mean Value	S.D.	% C.V.	n
Low	7.1	0.6	8.4	15
Medium	13.4	1.3	10.1	15
High	24.8	2.2	9.1	15

Between Assay Variation (values = nmol/L)

	Mean Value	S.D.	% C.V.	n
Low	8.9	1.1	12.5	15
Medium	17.4	1.8	10.3	15
High	32.7	5.0	15.2	15

**16.2 TRUENESS: THE ASSAY TRUENESS HAS BEEN CHECKED BY THE LINEARITY TEST AND THE RECOVERY TEST.**

**Linearity (Parallelism)**

Serial Dilution Study of 3 Unknown Samples (values = nmol/L)

Sample Number	Undiluted	1/2	1/4
1	21.8	26.5	29.6
2	29.2	28.2	26.4
3	37.5	34.3	37.0

**Recovery**

Recovery Study (values = nmol/L). Percent recovery is calculated as measured value ÷ expected value x 100.

Background	Calibrator Added	Expected Value	Measured Value	Percent Recovery
Set No. 1				
15.7	1.2	16.9	16.9	100
15.7	3.4	19.1	18.3	96
15.7	7.7	23.4	24.6	106
Set No. 2				
20.7	1.2	21.9	17.7	81
20.7	3.4	24.1	21.3	88
20.7	7.7	28.4	28.5	101

**16.3 Analytical Sensitivity**

When defined as the apparent concentration at 3 standard deviations from the counts at maximum binding, the minimum detectable amount is <2.0 nmol/L.

**16.4 Analytical Specificity**

Comparison of the reactivity of IGF-I antibody at the 50% suppression B/B<sub>0</sub> was made with the following peptides:

Peptide	%Cross-Reactivity
IGF-II	<1.0
Human Growth Hormone	<1.0
Fibroblast Growth Factor	<1.0
Transforming Growth Factor	<1.0
Platelet-Derived Growth Factor	<1.0

**SEE LAST PAGE FOR REFERENCES**

### SCHEME OF THE ASSAY

1. Reconstitute the lyophilized reagents and allow any frozen specimens to thaw completely.
2. Extract the controls and unknown samples according to section 9 of the product insert.
3. Identify tubes in duplicate.
4. Dispense reagents according to the following scheme.

Tubes/Reagents	Total Counts	NSB	Cal 0-5	Controls and unknown samples
Calibrator 0	-	50 $\mu$ L	50 $\mu$ L	-
Calibrators (1-5)	-	-	50 $\mu$ L	-
Extracted Controls	-	-	-	50 $\mu$ L
Extracted Unknown Samples	-	-	-	50 $\mu$ L
Antiserum	-	-	200 $\mu$ L	200 $\mu$ L

5. Mix well; incubate for 2 hours (+/- 15 minutes) at 2-8°C.
6. Dispense 200  $\mu$ L of Tracer into all wells.
7. Mix well; incubate for 20 hours (+/- 1 hour) at 2-8°C.
8. Dispense 500  $\mu$ L Precipitating Complex into all wells, except the Total Count tubes.
9. Mix well, incubate for 2 hours (+/- 15 minutes) at 2-8°C.
10. Centrifuge using  $760 \times g^*$  for 20 minutes.
11. Decant the supernatants.
12. Count each tube in a gamma counter for 60 seconds or longer.

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
$$*g = (1118 \times 10^{-8}) (\text{radie i cm}) (\text{rpm})^2$$

#### REFERENCES/BIBLOGRAPHIE/LITERATUR/REFERÊNCIAS/REFERENSER

1. Bala, R.M., J. Lopatka, A. Leung, E. McCoy and R.G. McArthur, "Serum Immuno-reactive Somatomedin Levels in Normal Adults, Pregnant Women at Term, Children at Various Ages and Children with Constitutionally Delayed Growth," **Journal of Clinical Endocrinology and Metabolism**, 52(3):508, (1981).
2. Borges, J.L., M.C. Gelato, A.D. Rogol, M.L. Vance, R.M. Macleod, D.L. Loriaux, J. Rivier, R.M. Blizzard, R.F. Furlanetto, W.S. Evans, D.L. Kaiser, G.R. Merriam, J. Speiss and W. Vale, "Effects of Human Pancreatic Tumor Growth Hormone Releasing Factor on Growth Hormone and Somatomedin C Levels in Patients with Idiopathic Growth Hormone Deficiency," **The Lancet**, 199, July, (1983).
3. Hayek, A. and G.T. Peake, "Growth and Somatomedin-C Responses to Growth Hormone in Dwarfed Children," **The Journal of Pediatrics**, 99(6):868, (1981).
4. Klapper, D.V., M.E. Svoboda and J.J. Van Wyk, "Sequence Analysis of Somatomedin-C: Confirmation of Identity with Insulin-Like Growth Factor I," **Endocrinology**, 112(6):2215, (1983).
5. Rieu, M., F. Girard, H. Bricaire and M. Binoux, "The Importance of Insulin-Like Growth Factor (Somatomedin) Measurements in the Diagnosis and Surveillance of Acromegaly," **Journal of Clinical Endocrinology and Metabolism**, 55(1):147, (1982).
6. Rinderknecht, E. and R.E. Humbel, "The Amino Acid Sequence of Human Insulin-like Growth Factor I and Its Structural Homology with Proinsulin," **The Journal of Biological Chemistry**, 253(8):2769, (1978).
7. Rosenfeld, R.G., S.F. Kemp and R.L. Hintz, "Constancy of Somatomedin Response to Growth Hormone Treatment of Hypopituitary Dwarfism and Lack of Correlation with Growth Rates," **Journal of Clinical Endocrinology and Metabolism**, 53:611, (1981).
8. Underwood, L.E., A.J. D'Ercole and J.J. Van Wyk, "Somatomedin-C and the Assessment of Growth," **Pediatric Clinics of North America**, 27:771, (1980).
9. Wass, J.A.H., D.R. Clemmons, L.E. Underwood, I. Barrow, G.M. Besser and J.J. Van Wyk, "Changes in Circulating Somatomedin-C Levels in Bromocriptine-Treated Acromegaly," **Clinical Endocrinology**, 17:369, (1982).
10. Wilson, D.M., A.M. Luna, R.L. Hintz and R.G. Rosenfeld, "Somatomedins in Adolescence: The Effect of Puberty on Plasma Levels of Insulin-Like Factors I and II," **Clinical Research**, 31(1):61A, (1983).

### SYMBOLS USED WITH DEVICES

English	Français	Deutsch	Português	Svenska
European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformidade com as normas europeias	Europeisk överensstämmelse
Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Prazo de validade	Utgångsdatum
Manufacturer	Fabricant	Hersteller	Fabricante	Tilverkare
Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consulte as instruções de utilização	Se bruksanvisningen
In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnóstico in vitro.	Diagnostik in vitro.
Lot No.	No. de lot	Chargen-Nr.	N.º do lote	Batch-nummer
Maximum temperature	Température maximale	Höchsttemperatur	Temperatura máxima	Max. temperatur
Ab	Antisérum	Antiserum	Anti-soro	Antiserum
Ab PEG	Réactif précipitant	Fällungsreagenz	Precipitado no reagente	Utfällningsreagens
Ag <sup>125</sup> I	Traceur : antigène marqué à l'iode <sup>125</sup>	Tracer: <sup>125</sup> I-markiertes Antigen	Traçador: antígeno marcado com <sup>125</sup> I	Spårelement: antigen betecknad med <sup>125</sup> I
CAL	Étalon	Kalibrator	Calibrador	Kalibrator
CONTROL	Sérum de contrôle	Kontrollserum	Soro de controlo	Kontrollserum
RCNS X mL	Reconstituer avec X mL	Mit X ml auflösen	Reconstituir com X mL	Rekonstituera med X mL
Radioactive	Radioactif	Radioaktiv	Radioactivo	Radioaktiv
Harmful	Nocif	Gesundheitsschädlich	Nocivo	Skadligt

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DiaSorin S.A.  
11, rue Georges Besse  
92160 Antony, France

DiaSorin Inc.  
1951 Northwestern Avenue  
P.O. Box 285  
Stillwater, MN 55082-0285 U.S.A.

Phone: 651-439-9710  
FAX: 651-779-7847

For Customer Service in the U.S. and Canada Call Toll Free: 800-328-1482
In the United Kingdom Call: 44 118 9364200 FAX: 44 118 9792061
10417
27277 9/03

PRINTED IN U.S.A.