Fairway

Balloon catheter

Cordynamic

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RAPID EXCHANGE DILATATION CATHETER FOR PERCUTANEOUS TRANSLUMINAL CORONARY ANGIOPLASTY Cordynamic® FAIRWAY

INSTRUCTIONS FOR USE

Note: The instructions given in this document include all lengths and diameters of the balloon. Also include each diameter of the shaft.

Sterilised by ethylene oxide gas. Non pyrogenic. For only one use. Do not re-sterilise or reprocess. Check for the integrity of the package. Use prior the past due date. Store in a dry, cool, dark place.

1.- Description

The dilatation catheter for PTCA, is a catheter with a balloon at the distal end. The shaft has a combination of a single lumen at the proximal end, and double lumen at the distal end. At the distal end, the outer lumen is to inflate the balloon and the inner lumen permits to use of guidewires to facilitate advancement of the catheter to and through the stenosis to be dilated. The balloon is designed to provide a inflatable segment of known diameter and length at recommended pressures. The catheter has a tapered atraumatic tip.

The balloon has radiopaque markers to facilitate the positionament of the balloon at the stenosis and also it has been designed to provide an expandable segment to a known diameter and length under the application of a specific pressure.

In the shaft of the catheter there are two markers for calculate the catheter position according with the tip of the guiding catheter (the marker which is nearest to the adapter is for the guiding femoral catheter, and the furthest is for brachial guiding catheters).

2.- Indications

The Cordynamic® Balloon catheter FAIRWAY is indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion.

3.- Contraindications

- Unprotected left main coronary artery.
- Coronary artery spasm in the absent of a significant stenosis.

4.- Warning

- This device is intended for one time use only. DO NOT RESTERILIZE. Use the catheter prior to the: "Use Before" date specified on the package.
- Do not use organic solvents or oil-based (such as contrast mediums Ethiodol or Lipiodol) or highly viscous products.
- For avoiding the damage of the blood vessels, the diameter of the inflated balloon must be as near as possible to the diameter of the artery before and after the stenosis.
- When the catheter is exposed to the vascular system, it should be manipulated while under fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated. If there were any resistance during handling, determine the cause of the resistance before proceeding.
- Balloon pressure should not exceed the rated burst pressure. It is recommended the use of a pressure control device to avoid the overpressurization. At least the 99.9% of the balloons (with a 95% confidence) will not burst at or below their rated burst pressure.
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potentially injurious of life-threatening complication. The procedures must be performed by physicians with suitable training in the use of the device.
- Never use air or any gaseous medium to inflate the balloon.

5.- Precautions

- Prior to the use, choose the appropriate size of the catheter. Inspect the packaging and remove any products whose packaging is in bad condition.
- This type of devices must be used only by physicians trained in the performance of PTCA.
- Before the insertion, administer anticoagulant and coronary vasodilator therapy.
- To avoid drawing blood into the catheter, release vacuum immediately after full deflation of the balloon, or while manipulating the catheter and/or guidewire in the body.

6.- Potential adverse effects

Possible adverse effects can be:

- Death
- Acute myocardial infarction
- Unstable angina pectoris
- Embolism
- Spasm in the coronary artery
- Arrhythmia
- Total occlusion of the coronary artery
- Coronary vessel dissection, perforation, rupture or injure
- Restenosis
- Hematoma
- Haemorrhage
- Infection
- Arteriovenous fistula
- Changes in the arterial pressure
- Allergic reactions to contrast medium

7.- Form of administration

Sterile: This product is sterilised with ethylene oxide and is pyrogen-free. For single use. Do not use if the packaging is not intact.

Contents: A Cordynamic® Balloon catheter FAIRWAY.

A handbook with instructions for use

Storage: Keep at room temperature, in a dry, cool place, protected from light.

Instructions for disposal: After use, dispose of the device and its package following hospital, administrative or local government rules.

8.- Instructions for the operator

Before to start with the procedure, carefully examine all equipment to be used, including the dilatation catheter. Check about possible damage of the packaging and do not use those whose package have been damaged.

1. Prepare the inflation device with the recommended contrast medium according with the manufacturer instructions. Purge air from the system.
2. Choose the correct length and diameter of the balloon. The diameter of the balloon must not exceed the diameter of the artery which must be dilated before and after of the stenosis.
3. Dilatation catheter preparation:
 - 3.1. Leave the catheter in its protective tubing. The distal end of the protective tubing may be removed to visualise the balloon. Slide the balloon protector distally, so as not to damage it.
 - 3.2. Prepare the balloon catheter for purging. Fill a 20 cc luer-lock syringe with 4 cc with the appropriate contrast medium. Do not use air or any gaseous medium to inflate the balloon.
 - 3.3. Connect the syringe to the balloon inflation lumen, placing the catheter with the distal tip and the balloon pointing down vertically .
 - 3.4. Apply negative pressure and aspirate for 15 s. Slowly release the embolus.
 - 3.5. Remove the syringe from the balloon inflation lumen and extract all the air.
 - 3.6. Reconnect the syringe to the catheter and maintain negative pressure in the balloon until air stops returning to the device. Release the pressure from the syringe to the neutral point and disconnect it.
 - 3.7. Immerse the catheter in sterile saline solution to reactivate the hydrophilic coating.
 - 3.8. Connect the inflation device to the catheter inflation port without introducing air into the system.

4. Insertion technique:
 - 4.1. Release the catheter from its protective point
 - 4.2. Insert the introducer and the guiding catheter using standard procedures.
 - 4.3. Before the introduction of the guidewire in the distal end of the catheter, irrigate the lumen of the catheter with sterile saline solution. After that introduce the guidewire from the distal end ensuring that the guidewire comes out from the skive which is located 25-30 cm from the tip of the catheter. For doing this operation it may be required to straighten the catheter to facilitate the advance of the guidewire. Pull the guidewire proximally so that the distal end of the guide wire is inside the distal of the catheter lumen for protection. Place the guiding catheter hemostatic adapter (with swivel) over the balloon dilatation catheter and advance the adapter approximately 40 cm up to the catheter shaft. Advance the catheter and guidewire slowly, while the balloon is fully deflated. If resistance is encountered, do not advance the dilatation catheter through the adapter. Caution should be taken not to over tighten a Touhy-Borst type hemostatic adapter around the dilatation catheter shaft as lumen constriction may occur, affecting inflation/deflation of the balloon.
 - 4.4. Aspirate and flush the guiding catheter in preparation for introduction of the dilatation catheter.
 - 4.5. Connect the side port of the guiding catheter hemostatic adapter to the proximal pressure recording/infusion line or manifold assembly, which permit proximal pressure recording or infusion through the guiding catheter. Empty thoroughly.
 - 4.6. With the balloon deflated, insert the dilatation catheter into de luer fitting of the guiding catheter.
 - 4.7. Advance the guidewire through the dilatation catheter. The guidewire should extend approximately 1 cm beyond the tip of the dilatation catheter.
 - 4.8. Using test injection of contrast medium to confirm the right position of the guiding catheter.
 - 4.9. Advance the dilatation catheter and the guidewire to the tip of the guiding catheter. Be careful do not rotate or twist the balloon around the guide wire in order to avoid its damaging.
 - 4.10. Advance the guidewire through the coronary artery while carefully rotating until cross the stenosis. Perform a control angiography to ensure that the stenosis has been crossed.
 - 4.11. Position the balloon in the desired position. Inflate the balloon to the appropriate pressure, do not exceed the recommended rated burst pressure.
 - 4.12. After the procedure has been completed, deflate the balloon and withdraw it from the stenosis site. Maintain the guidewire across the dilated stenosis.
 - 4.13. Perform angiography using the guiding catheter as an angiographic catheter to confirm dilatation.
 - 4.14. Maintain the guidewire across the stenosis dilated for 10 minutes. Once the angiography has confirm the dilatation of the artery, withdraw the guidewire and the balloon catheter through the connector.
 - 4.15. Withdraw the guiding catheter through the introducer.
 - 4.16. Leave the introducer sheath in situ until hemodynamic profile becomes normal. Close the skin in routine manner.
5. Dilatation catheter exchange procedure:

In case of needing to exchange the balloon catheter, the rapid exchange balloon Iberhospitex, has been designed to allow a rapid exchange with only one operator. To carry out the interchange:

 - 5.1. Loosen the hemostatic valve.
 - 5.2. Hold the wire and hemostatic valve in one hand, while grasping the balloon shaft in the other hand.
 - 5.3. Maintain the guidewire's position in the coronary artery by holding the wire stationary, and begin pulling the balloon catheter out of the guiding catheter.
 - 5.4. Withdraw the catheter until opening in the guidewire lumen is reached (at approximately 25 cm from the tip of the balloon). Carefully inch the flexible, distal portion of the balloon catheter off the guidewire maintaining the guidewires position across the lesion. Close the hemostatic valve.
 - 5.5. Prepare the new balloon as previously described.
 - 5.6. Back load the new balloon onto the guidewire. The proximal portion of the guidewire will exit at the opening 25 cm from the distal balloon tip.
 - 5.7. Open the hemostatic valve and advance the balloon while holding the guidewire and maintaining its position in the coronary artery. Be careful do not twist or rotate the balloon catheter around the guidewire.
 - 5.8. Advance the catheter to the guiding catheter tip. Continue the PTCA procedure, using the previously described method.

9.- Warnings about re-use

Contents are supplied STERILIZED by ethylene oxide (EO). Do not use if the sterile barrier is damaged. If damage is detected, call an Iberhospitex agent. For single use in a single patient. Do not re-use, re-process or re-sterilise. Re-use, re-processing or re-sterilisation can affect the structural integrity of the device and/or cause its failure, which in turn could cause injuries to the patient, disease or death. Re-use, re-processing or re-sterilisation can also create a risk of device contamination and/or cause infection or cross infection to the patient, including but not limited to the transmission of infectious diseases from one patient to another. Contamination of the device can cause injuries, disease or the death of the patient.

10.- Warranty

Iberhospitex is not liable for damages caused to patients by misuse of the catheter; therefore, no claims will be accepted except for those stated in this document. Under no circumstances will Iberhospitex be liable for any damage caused, except for those specifically established by the law. The information provided in this package insert is a general description of the product and does not entail any express warranty for use. This information is based on current knowledge of the product and medical science.

Explanation of the symbols used on the package labels:

	Valid for a single use.		External diameter
	See Instructions of use before using this product.		Manufacturer
	Sterilised with ethylene oxide.		Store in a dry place
	Reference number		Store in a dark place
	Batch number		Do not resterilize
	Use before		

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