
Vacuum Regulator, Push-To-Set™ Continuous (PTS-CVR)

Two-Mode

ANSI

Analog

Adult (Standard) - (1225)
Adult (High) - (1330)
C.A.S.S. - (1227)
Neonatal - (1231)
Pediatric - (1234)

Digital

Adult (Standard) - (1325)
Adult (High) - (1430)
C.A.S.S. - (1327)
Neonatal - (1331)
Pediatric - (1334)


ISO

Analog

Adult (High) - (1230)
Neonatal - (1232)
Pediatric - (1235)

Digital

Adult (High) - (1337)
Neonatal - (1332)
Pediatric - (1335)

Three-Mode

ANSI

Analog

Adult (Standard) - (1224)
Adult (High) - (1226)

Digital

Adult (Standard) - (1324)
Adult (High) - (1326)

ISO

Analog

Adult (High) - (1229)

Digital

Adult (High) - (1329)

EN	Instructions for Use.....	2
NL	Gebruiksaanwijzing	11
FR	Mode d'emploi	20
DE	Gebrauchsanweisung	29
IT	Istruzioni per l'uso	38
ES	Instrucciones de uso	47
SV	Bruksanvisning	56

Table of Contents

Safety Instructions	2
Intended Use	2
Receiving/Inspection	2
User Responsibility	2
Definitions	3
Specifications	5
Regulator Identification	6
Operation	7
Pre-use Checkout Procedure	9
Patient Set-up	9
Cleaning and Disinfection	10
Troubleshooting	10

Safety Instructions

This manual provides you with important information about the Push-To-Set™ Continuous (PTS-CVR). To ensure the safe and proper use of this device, READ and UNDERSTAND all of the safety and operating instructions. IF YOU DO NOT UNDERSTAND THESE INSTRUCTIONS, OR HAVE ANY QUESTIONS, REFER TO SERVICE MANUAL, CONTACT YOUR SUPERVISOR, DEALER OR THE MANUFACTURER BEFORE ATTEMPTING TO USE THE DEVICE.

Intended Use:

The vacuum regulator is intended to be used in the medical facility as a means to evacuate media (i.e. fluids) from the body. **DO NOT** use this vacuum regulator for anything other than its intended use.

Receiving / Inspection:

Remove product from package and inspect for damage. If product is damaged, DO NOT USE and contact your dealer or equipment provider.

User Responsibility

This Product will perform as described in this operating manual and accompanying labels and/or inserts, when assembled, operated, maintained and repaired in accordance with the instructions provided. This Product must be checked periodically. A defective product should not be used. Parts that are broken, missing, worn, distorted or contaminated should be replaced immediately. Should such repair or replacement become necessary, see the Ohio Medical service manual for service or repairs to this product. For service advice, Ohio Medical recommends that a telephone request be made to the nearest Ohio Medical Regional Service Center. This product and any of its parts should only be repaired using written instructions provided by Ohio Medical or by Ohio Medical trained personnel. The Product must not be altered without the prior written approval of Ohio Medical's Quality Assurance Department. The user of this Product shall have the sole responsibility for any malfunction which results from improper use, faulty maintenance, improper repair, damage, or alteration by anyone other than Ohio Medical.

Failure to follow these instructions may cause damage to the patient, user and/or vacuum regulator.

AAAA 12345 This alpha character indicates the year of product manufacture and when the serial number
_____ was assigned; "L" = 2007, "M" = 2008, "N" = 2009, etc. "I" and "O" are not used.

Definitions

WARNING	= possible injury to patient or operator
CAUTION	= possible damage to equipment
Note	= Provides additional information to clarify a point in the text.
Important	= Similar to a note but of greater emphasis
 O O (INT)	= intermittent; cycles ON and OFF
 (CONT)	= continuous, ON
O (OFF)	= OFF

 = Consult Instructions for Use

 = Serial Number

 = Manufacturer

Abbreviations

in	Inch	N-m	Newton-Meter (N-m x 0.737 = ft-lb)
CVR	Continuous Vacuum Regulator	ft-lb	Foot-Pound Force (ft-lb x 1.356 = N-m)
kPa	Kilo pascals (kPa x 7.50 = 1 mmHg)	in-lb	Inch-Pound Force (ft-lb x 12 = in-lb)
mmHg	Millimeters of mercury (mmHg x 0.133 = kPa)	DISS	Diameter Index Safety System
mm	Millimeters	NPTF	National Pipe Thread Female (USA)
mL	Milliliters	NEO	Neonatal
°C	Degrees Celsius	PED	Pediatric
°F	Degrees Fahrenheit	PTS	Push-To-Set™
		PTFE	Polytetrafluoroethylene

WARNINGS

The Pre-Use Checkout Procedure must be performed before using the equipment on each patient. If the regulator fails any part of the Pre-Use Checkout Procedure, it must be removed from service and repaired by qualified service personnel.

Do not use this device in the presence of flammable anesthetics. Static charges may not dissipate and a possible explosion hazard exists in the presence of these agents.

This device is to be used only by persons who have been adequately instructed in its use.

To reduce service personnel exposure to hazardous contamination, clean and disinfect all suction equipment before disassembly.

To reduce transportation personnel and/or service personnel exposure to hazardous contamination, DO NOT ship any suction equipment that has been contaminated.

Ohio Medical will assume no responsibility for incidents which may occur if the product is not used in accordance with product labeling.

Connection to positive pressure sources such as oxygen and medical air, even momentarily, could injure the patient or operator.

After patient use, regulators may be contaminated. Handle in accordance with your hospital's infection control policy.

To help prevent aspirate from entering the device, wall outlet and pipeline equipment, a safety trap should be attached prior to its use. Aspirate in the device, wall outlet and pipeline equipment may impair its operations. The use of the safety trap and suction filter will help prevent this and extend the life of suction equipment.

Definitions

CAUTIONS

Connection to **positive** pressure sources **such as oxygen and medical air**, even momentarily, could damage the equipment.

Detailed information for more extensive repairs is included in the service manual for users having proper knowledge, tools and test equipment, and for service representatives trained by Ohio Medical.

The suction control knob must be completely pushed in to adjust the vacuum level. Failure to do so may damage the vacuum regulator.

Not for field or transport use: The categories of Field and Transport Use are specifically defined in ISO 10079-3 "Field" means accidents or emergencies outside the hospital. "Transport" means use in ambulances, cars and airplanes. These situations may expose the equipment to uneven support, dirt, water, mechanical shock and temperature extremes. Ohio Medical suction equipment has not been tested to comply with the specific requirements of these categories.

Do not use harsh chemical or cleaning solution. Do not spray cleaners directly onto suction regulators. Only use chemical recommended in this manual.

Important: BATTERY LOW INDICATOR: When a battery icon appears on the gauge it indicates that the battery is low. Please take the unit out of service immediately and contact an Ohio Medical Customer Service Representative for battery replacement.

If the low battery condition is not addressed and the battery becomes fully depleted, the gauge will not show any readout, including the low battery icon or gauge pressure. If the gauge were to go blank during suctioning, the unit will continue to suction and the intermittent feature will continue to operate. Once completing that procedure, it is important to immediately take the unit out of service and contact an Ohio Medical Customer Service Representative for battery replacement.

Specifications

Performance					
	Standard	High	Pediatric	Neonatal	C.A.S.S.
Flow rate ¹	80 LPM	80 LPM	80 LPM	80 LPM	80 LPM
Positive Pressure Safety Relief Valve:	Located in the vacuum supply line to prevent pressurization of the patient connection by failed injector vacuum (Venturi) units, or inadvertent cross connection to pressured gases.				
Gauge Accuracy	Gauge Range		Analog ²	Digital ³	
Standard	0-200 mmHg (0-26.7 KPA)		±5% Full-scale	±1% Full-scale	
High	0-760 mmHg (0-101.3 KPA)		±5% Full-scale	±1% Full-scale	
Pediatric	0-160 mmHg (0-21.3 KPA)		±5% Full-scale	±1% Full-scale	
Neonatal	0-160 mmHg (0-21.3 KPA) 0-100 mmHg (0-13.3 KPA)		±5% Full-scale ±5% Full-scale	±1% Full-scale ±1% Full-scale	
C.A.S.S.	0-160 mmHg (0-21.3 KPA)		±5% Full-scale	±1% Full-scale	
Physical					
Dimensions	6.7”H x 2.8”W x 4.8”D (17.02 cm x 7.1 cm x 12.2 cm)				
Weight	1 lb 4 oz (0.57 kg)				
Battery					
Lithium	Two 2/3 AA batteries, 3.6V, 1.6 Ah, lithium				

¹without fittings at full increase setting depending on supply vacuum and open air flow.

²full scale deflection

³full range at 22°C

ENVIRONMENTAL SPECIFICATIONS

Operating Temperature Range: 50 to 104°F (10 to 40°C)

Storage Temperature Range: -13 to 158°F (-25 to 70°C)



Operating and storage Relative Humidity: 5 to 95% RH (Non-condensing)

Ingress Protection Ratings: IP 20

Regulator Identification

Push-To-Set™ Continuous Vacuum Regulator (PTS-CVR)

- ① Vacuum Gauge/Analog
- ② Mode Selector Switch
- ③ Suction Control Knob
- ④ Fitting Port (inlet)
- ⑤ Probe/Adapter Port (outlet)
- ⑥ Vacuum Gauge/Digital

⑦ The pediatric and neonatal models are identified by the Baby Icon on the front right when looking at the unit.

Operation

Equipment Set-up

Insert the adapter/probe into the vacuum wall outlet. If the regulator is mounted elsewhere, connect a vacuum supply hose between the regulator's adapter/probe and the wall outlet.

Use hospital-supplied suction tubing between the vacuum regulator and the collection container, as well as between the patient port of the collection container and the patient (minimum inside diameter is 6 mm [0.25 in]).

An Ohio Medical high flow suction filter and/or overflow safety trap (OST) should be used between the collection container and regulator to prevent contamination of the regulator, wall outlet and pipeline system.

ISO 10079-3 (Section 5.1.2) states that "the usable volume of the collection container shall not be less than 500 mL."

High Flow Suction Filters

Description	Part Number		
Disposable Hydrophilic Filter	6730-0350-800 20 Pack	6730-0351-800 200 Pack	
Disposable Hydrophobic Filter Tubing x Tubing nipple	6730-0570-800 3 Pack	6730-0571-800 10 Pack	6730-0572-800 50 Pack
Disposable Hydrophobic Filter 1/8 NPT x Tubing nipple	6730-0580-800 3 Pack	6730-0581-800 10 Pack	6730-0582-800 50 Pack
Tubing	6700-0005-300		

Note: For proper installation of adapters/probe and fittings, see PTS-CVR Service Manual.


Attaching the Overflow Safety Trap (OST)

CAUTION To help prevent aspirate from entering the regulator, wall outlet and pipeline equipment, a safety trap should be attached prior to its use. Aspirate in the regulator, wall outlet and pipeline equipment may impair its operation. The use of a safety trap and suction filter will help prevent this and extend the life of suction equipment.

Locking Gland Fitting

1. Raise the sleeve and insert the safety trap into the regulator fitting.
2. Turn the trap clockwise about one and a half turns to engage the threads. The trap does not need to be screwed tight; an O-ring in the regulator fitting provides a vacuum seal. The trap should rotate freely to allow the desired nipple positioning.
3. Lower sleeve to lock trap in position.

- ① Regulator
- ② Sleeve
- ③ Safety Trap


Operation

DISS Fitting

1. Insert the safety trap into the regulator fitting. Situate the tubing in the desired position.
2. Turn the DISS wing nut clockwise to engage threads and tighten (there is no O-ring, so the vacuum seal depends on a tight connection).


- ① Regulator
- ② Wing nut
- ③ Safety Trap


Mode Selection

2-Mode Continuous

| (ON) - Suction can be adjusted with the suction control


O (OFF) - No suction supplied to the patient

3-Mode Continuous


MAX - Maximum, full-line vacuum is supplied to the patient



Note: Available on Three-Mode Vacuum Regulators ONLY

Setting the Suction Level

1. Move the mode selector switch to | (ON).

2. Push and rotate the suction control knob until the vacuum gauge indicates the required setting.

CAUTION: The suction control knob must be completely pushed in to adjust the vacuum level. Failure to do so may damage the vacuum regulator.

Pre-Use Checkout Procedure/Patient Set-up

Pre-Use Checkout Procedure

WARNING The Pre-Use Checkout Procedure must be performed before using the equipment on each patient. If the regulator fails any part of the Pre-Use Checkout Procedure, it must be removed from service and repaired by qualified service personnel.

Important: All tests must be performed with supply vacuum of 66.7 kPa (500 mmHg) minimum.


Note: Pediatric/Neonatal/C.A.S.S. models Only: Rotate suction control knob fully clockwise (increase) to verify that the suction level does not exceed 135 mmHg $\pm$ 5 mmHg (18.0 kPa $\pm$ 0.7 kPa) on the pediatric model or 100 mmHg $\pm$ 5 mmHg (13.3 kPa $\pm$ 0.7 kPa) on the neonatal model or 140 mmHg $\pm$ 5 mmHg (18.7 kPa $\pm$ 0.7 kPa) on the C.A.S.S. model. When maximum vacuum is reached, the safety relief valve will emit an audible pulse.

1. Move the mode selector switch to O (OFF).
2. Push and rotate the suction control knob one full turn clockwise (increase).
3. Release the suction control knob. The gauge needle should not move (reading should remain at zero on units with a digital gauge).
4. Move the mode selector switch to | (ON). The gauge should indicate vacuum.
5. Push and fully rotate the suction control knob counter-clockwise (decrease) until it stops.
6. Release the suction control knob. The gauge needle should move to zero and remain there. (Readout should be zero on the digital gauge.)
7. Set the unit to the appropriate vacuum level. Slowly release and push the suction control knob to create various flow rates through the regulator. Check that the suction level is maintained when the knob is fully pushed in.
 - For analog units, the gauge should return to $\pm$ 10 mmHg ($\pm$ 1.3 kPa) for the standard model, $\pm$ 8 mmHg ($\pm$ 1.1 kPa) for the PED/NEO/C.A.S.S. models, and $\pm$ 38 mmHg ($\pm$ 5.1 kPa) for the high model.
 - For digital units, the gauge should return to $\pm$ 2 mmHg ($\pm$ 0.3 kPa) for the standard/PED/C.A.S.S. models, $\pm$ 1 mmHg ($\pm$ 0.1 kPa) for the NEO model, and $\pm$ 8 mmHg ($\pm$ 1.0 kPa) for the high model.
8. Move the mode selector knob to **MAX** on three mode units.
9. Push the suction control knob and check the gauge to ensure maximum suction is applied.
10. Release the suction control knob.
11. Push the suction control knob and reduce the suction to zero. Set the mode selector to O (OFF).

Patient Set-up

1. Make sure the Pre-Use Checkout Procedure has been performed.
2. Move the mode selector switch to | (ON).
3. Set the prescribed suction level by pushing and turning the suction control knob.

CAUTION: The suction control knob must be completely pushed in to adjust the vacuum level. Failure to do so may damage the vacuum regulator.
4. Move the mode selector switch to O (OFF).
5. Attach tubing from the vacuum regulator to the vacuum port of the collection container.
6. Move the mode selector switch to | (ON).

Cleaning and Disinfection

Cleaning and Disinfection

Routine cleaning of the regulator is recommended as a standard procedure after each use.

Wipe all exterior surfaces with a solution of water and mild detergent and/or an approved cleaning solution.

Should misuse occur resulting in accidental flooding of the regulator, the regulator may be sterilized after cleaning using ethylene oxide (ETO). See Section 4.0 (Sterilization) of the PTS-CVR service manual. After sterilization follow the service checkout procedures in Section 7.0 (Service Checkout Procedure) of the PTS-CVR service manual.

Approved Cleaning Solutions

- Sodium Hypochlorite 0.5% (Bleach): Mixture of 13 fl. oz. of bleach to 1 gallon (128 fl. oz.) tap water
- Isopropyl Alcohol 70%
- Hydrogen Peroxide 3%
- Cavicide® Ready to use full strength

WARNING **After patient use, regulators may be contaminated. Handle in accordance with your hospital's infection control policy.**

Troubleshooting

If the regulator does not operate and you have performed the Pre-Use Checkout Procedure, remove the regulator from service and refer to PTS-CVR service manual or refer servicing to qualified service personnel.

Installation Procedure for Adapters/Probes and Fittings

For adapter/probe installation, refer to PTS-CVR service manual, section 8.

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Notes

**North America
United States**

**Customer Service, Distribution Center
Technical Support, Sales and Service
Equipment Service Center**

Ohio Medical Corporation
1111 Lakeside Drive
Gurnee, IL 60031 USA
P: 866 549 6446
P: +1 847 855 0800
F: +1 847 855 6218

**Ohio Medical Corporation Authorized
Representative**

OxygenCare Ltd.
2 Holfeld Business Park
Kilmacanogue Co Wicklow
Ireland
Tel: +353 1 2769700
Fax: +353 1 2764970

© 2014 Ohio Medical Corporation. All rights reserved
This document contains information that is proprietary and confidential to Ohio Medical Corporation. Use of this information is under license from Ohio Medical Corporation.
Any use other than that authorized by Ohio Medical Corporation is prohibited.

Ohio Medical, Ohio Medical Corporation, and the Ohio Medical Corporation Logo are registered trademarks and Push-To-Set is a trademark of Ohio Medical Corporation.

Cavicide is a registered trademark of Metrex Research, LLC.