## **Visual Field Clearing System**

### **User Manual**

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## Foreword

This manual and the equipment it describes are for use only by qualified medical professionals who have been trained to use electrosurgical equipment.

Equipment covered in this manual  
Ultravision™ visual field clearing system, comprising:

Item	Part number
Ultravision™ Generator	DAD-001-010
Ultravision™ Battery (x2)	DSD-001-034
Ultravision™ Battery Recharging Station	DAD-001-024
Ultravision™ Power Supply Unit	DSD-001-035
Ultravision™ Patient Return Adaptor (Solid)	DAD-001-006
Ultravision™ Patient Return Adaptor (Split)	DAD-001-007
Ultravision™ Ionwand™ Sterile Pack*	DAD-001-003

\*supplied separately

### Patents pending

Worldwide patents are pending for the Ultravision™ system, including PCT/GB2010/051196 and PCT/GB2012/052707.

### Trademarks

Ultravision™ and Ionwand™ are trademarks of Alesi Surgical Limited.

### Manufacturer

Alesi Surgical Limited  
Cardiff Medicentre  
Heath Park  
Cardiff, CF14 4UJ, U.K.

<http://www.alesi-surgical.com>

For information call +44 (0) 29 2029 1022.

USA representative:

Issue v6.0

MEDICAL – GENERAL MEDICAL EQUIPMENT  
AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY  
IN ACCORDANCE WITH ANSI/AAMI ES60601-1 (2005), CAN/CSA-C22.2 No. 60601-1 (2008)

© Alesi Surgical Limited, 2015. All rights reserved.

## Terms Used in this Manual

### **Warning**

Indicates a situation which, if not avoided, could result in death or serious injury.

### **Caution**

Indicates a situation which, if not avoided, may result in minor or moderate injury.

### **Be aware**

Indicates that the operator should be aware of a hazard that may damage the product.

### **Advice**

Indicates advice to the operator to optimise performance of the system.

## Statement on Electromagnetic Compatibility

The Ultravision generator and battery recharging station require special precautions regarding Electromagnetic Compatibility (EMC) and need to be installed and put in to service according to the EMC information provided in Chapter 9.

Portable and mobile RF communications equipment can affect the Ultravision generator and battery recharging station.

## Ultravision™ system warning

No modification of this equipment is allowed.

## Warranty

Please consult term and conditions of sale for details of the warranty offered with the Ultravision system and its components.

## Table of contents

Foreword .....	2
Terms Used in this Manual .....	3
Statement on Electromagnetic Compatibility .....	3
Ultravision™ system warning .....	3
Warranty .....	3
Table of contents .....	4
Chapter 1: Ultravision™ – overview and general features .....	5
Indications for use. ....	5
Contraindication.....	5
Parts shipped with the Ultravision™ system .....	5
The Ultravision™ system – principle of operation .....	6
Operating Parameters .....	7
Transport and Storage.....	7
Patient and Operating Room Safety.....	7
Maintenance .....	8
Chapter 2: Controls, Sockets and Indicators .....	9
Ultravision™ generator – Front panel view .....	9
Ultravision™ generator – Side panel view.....	10
Ultravision™ generator – Rear view.....	11
Ultravision™ Battery Recharging Station – Top view .....	11
Chapter 3: Preparing the Ultravision™ system for use .....	13
Periodic Inspection .....	13
Charging the battery using the Ultravision™ Battery Recharging Station.....	13
Inserting the battery into the Ultravision™ generator .....	14
Preparing Ultravision™ for use - Instruments requiring a patient return electrode.....	14
Preparing Ultravision™ for use – Instruments that do not require a patient return path .....	15
Chapter 4: Using the Ultravision™ system.....	16
Inserting the Catheter and Ionwand .....	16
Use during surgery .....	16
End of surgery .....	17
Chapter 5: After surgery .....	18
Cleaning the Ultravision™ Patient Return Adaptor.....	18
Cleaning the Ultravision™ Generator.....	18
Cleaning the Ultravision™ Battery Recharging station.....	18
Recharging the Ultravision™ battery.....	18
Chapter 6: Troubleshooting .....	19
Chapter 7: Summary System Specifications .....	21
Chapter 8: Accessories and replacement items .....	22
Chapter 9: Electromagnetic guidance .....	23
Symbols used on Ultravision .....	132

# Chapter 1: Ultravision™ – overview and general features

This section describes the indications for use for the Ultravision surgical smoke visual field clearing system and its features.

## Caution

Read all warnings, cautions, and instructions provided with the Ultravision system before use. Failure to do so could result in death or serious injury to the patient or operator.

## Indications for use.

The Ultravision system is indicated for the clearance of smoke and other particulate matter that is created during laparoscopic surgery by monopolar, bipolar and ultrasonic instruments.

## Warning

Do not allow the lonwand to come in to contact with the uninsulated portion of an electrosurgical instrument during operation. In such a situation the Ultravision generator issues an continual audio-visual proximity indication. Continued contact for in excess of 10 minutes may result in an electrical discharge from the electrosurgical instrument that may be sufficient to cause atrial fibrillation depending on proximity to the thoracic cavity.

## Contraindication

The Ultravision system is an electrical energy-based system that must not be used on patients fitted with a pacemaker or implantable cardioverter defibrillator (ICD).

## Parts shipped with the Ultravision™ system

The Ultravision system comprises both reusable and single-use, sterile disposable items. The single-use, sterile disposable items are supplied separately. Contact your local sales representative in order to establish a supply of these items. All other items required for Ultravision to function are shipped with the Ultravision system.

Before use, confirm that the shipping carton contains the following items:

1. Generator	x1
2. Battery	x2
3. Recharging station	x1
4. Power supply unit for recharging station	x1
5. Solid-plate patient return adaptor	x1
6. Split-plate patient return adaptor	x1
7. User manual	x1

If any of the items listed above are not in the shipping carton, please contact your local sales representative before attempting to use the system.

Prior to use, inspect all of the above parts for any damage that may have occurred during shipping. If you have any concerns, please contact your local sales representative before using the system.

The following reusable items will require replacement:

Item	Approximate number of uses before requiring replacement
Battery	200 charging cycles at 20°C (68°F) or 1 year from last recharge at 20°C (68°F).
Patient return adaptor (solid)	500 mating cycles with Ultravision generator and electrosurgical generator
Patient return adaptor (split)	500 mating cycles with Ultravision generator and electrosurgical generator
Battery recharging station	1000 mating cycles with Ultravision battery

Replacement reusable items are available from Alesi Surgical or its local distributors.

## The Ultravision™ system – principle of operation

The Ultravision system removes the particulate material that is created by energy-based instruments during laparoscopic surgery from the visual field. It achieves this by electrostatic precipitation, a process that transiently charges the particulate matter and causes it to be rapidly precipitated inside the peritoneal cavity as it is created. The system has been specifically designed to improve the clarity of the visual field and to minimise the release of particulate matter and gases into the operating theatre environment during laparoscopic surgery.

The Ultravision system requires a patient return pathway, similar to monopolar diathermy. The system has been designed such it can be configured in two ways, depending upon the type of cutting instrument that is used:



### With monopolar instruments that require a patient return electrode

The Ultravision system is provided with a Patient Return Adaptor that allows the electrosurgical unit and the Ultravision generator to share a single patient return electrode with the monopolar instrument (Figure 1).

### With bipolar and ultrasonic cutting instruments that do not require a patient return pathway

The Ultravision patient return socket in the generator allows a standard patient return electrode to be used and connected directly to the Ultravision generator (Figure 2).

The Ultravision system is not supplied with a patient return electrode but has been designed to be compatible with patient return electrodes that have an “international” style connector. Please contact Alesi Surgical or your local distributor if you experience any problems with your preferred patient return electrode.

The following components are required for Ultravision to function and are supplied with the system:

### 1. The generator unit

The generator unit is the source of the high voltage power. The output power cannot be adjusted manually. The unit has no connection to the mains power socket and is entirely powered by a removable battery.

### 2. The battery

The generator is powered by a battery that allows approximately 10 hours of continuous generator operation. The system is supplied with two batteries and a recharging station. As with all batteries of this type, the life will depend upon the number of recharges. It is recommended that batteries should be replaced every 12 months or when a fully charged battery no longer delivers six hours of generator power.

### 3. Battery recharging station

The Ultravision system is supplied with a twin-bay recharging station and a power supply unit. The power supply unit has four different plugs that allows it to be used in most countries and has a universal input allowing 90-264V AC.

### 4. Patient return adaptor – solid plate

The Ultravision system requires a patient return electrode. The patient return electrode can be shared with the electrosurgical generator if one is required by using the patient return adaptor. The “solid plate” patient return adaptor is used with solid plate patient return electrodes. The blue end of the adaptor plugs into the socket on the front of the electrosurgical generator. The orange connector at the other end of the cable plugs into the front of the Ultravision generator.

## 5. Patient return adaptor – split plate

This is the same as the solid plate patient return adaptor but is used with split plates. It has the same central pin that activates the sensing function of the electrosurgical generator (if present). The split plate adaptor end is red in colour.

The following single-use, sterile disposable components are required for Ultravision to function and are supplied as separate items:

## 6. Ionwand™

The Ionwand is the active cable that carries the current from the Ultravision generator to the patient. It terminates in an atraumatic stainless steel brush that is responsible for the electrostatic charging of the particles created during surgery.

**7. Catheter and trocar.** The catheter allows the Ionwand to be inserted into the patient's peritoneal cavity during surgery. It is introduced percutaneously using the trocar.

## Operating Parameters

The Ultravision system may be used within the following ranges

- Ambient temperature range 10°C to 40°C (50°F to 104°F)
- Relative humidity 20% to 75%, non-condensing

## Transport and Storage

- Ambient temperature range -10°C to 40°C (14°F to 104°F) except batteries
- Relative humidity 20% to 85%, non-condensing
- Battery ambient temperature range 0°C to 40°C (32°F to 104°F). Prolonged exposure to higher temperatures will reduce battery shelf life.
- Battery shelf life – 1 year from last recharge at 20°C (68°F)

## Patient and Operating Room Safety

### Warning

Only surgical staff trained to use electrosurgical equipment may operate the Ultravision system.

During initial set up and before each use, inspect the Ultravision system and battery recharging station for any damage. If damaged, do not use or attempt to repair. Call Alesi Surgical or your local distributor for assistance.

The Ultravision generator enclosure is plastic which completely insulates it from any internal electrical hazards, and as such requires no connection to an equipotential point.

To maximise patient safety, contact between the Ionwand and metallic instruments should be minimised. Otherwise, the metal instrument may be charged resulting in an electrostatic discharge to the patient or operator. This may be sufficient to cause atrial fibrillation depending upon proximity to the thoracic cavity.

If the Ultravision generator or battery recharging station becomes wet for any reason, switch off and air-dry thoroughly before reusing.

Do not remove any covers or panels exposing the internal components of the Ultravision system and battery recharging station. Removal of covers or panels will void any Warranty.

Ensure that the battery is correctly inserted into the Ultravision generator. Failure to do so may result in a fire hazard.

Do not operate this device in potentially explosive environments, such as in the presence of flammable anaesthetics.

The battery recharging station has only been certified for use outside the immediate operating theatre environment. It does not have sufficient redundant voltage or current protection to allow safe contact with the patient.

Use only with the power supply unit and cable supplied with the battery recharging station.

Plug the power cable directly into the power socket without any extension cords or adaptor plugs.

The power cable assembly should be checked periodically for damaged insulation or connectors. Do not use a damaged power cable.

Do not connect a wet power cable to a power socket.

Inspect any sterile component for damage to the packaging before use. If any damage is present, do not use.

Do not reuse or re-sterilise the lonwand, trocar or catheter. Doing so may result in cross contamination, infection or injury to the patient or medical staff and may damage the Ultravision system.

**Caution**

Read all warnings, cautions, and instructions provided with the Ultravision system before using.

Do not stack equipment on top of the Ultravision system.

To maximise patient safety, contact between the lonwand and patient tissue should be minimised. Otherwise, performance will be compromised and visibility will be impaired.

**Be Aware**

Using any unapproved battery with the Ultravision system may cause damage and will void the warranty.

Connect the power cable to a wall socket having the correct voltage. Otherwise, product damage to the recharger or battery may result.

Power supply voltage below 90V AC will significantly reduce the rate of recharging.

**Maintenance**

**Warning**

Inspect the Ultravision generator and battery recharging station before each use. If there is evidence of damage, do not use. Call Alesi Surgical or its local distributor for assistance.

Always turn off the Ultravision generator before cleaning.

Always unplug the battery recharging station before cleaning.

**Caution**

Dispose of the lonwand, trocar and catheter along with other operative waste materials, following the procedures for your institution.

**Be Aware**

Use only detergents that are already routinely used in your hospital to clean non-sterile electrosurgical generators and cables to clean the Ultravision generator, battery recharging station and patient return adaptors. Do not apply alcohol, caustic, corrosive, or abrasive cleaning or disinfectant compounds, solvents, or other materials that could cause damage.

Do not attempt to sterilise the Ultravision generator or patient return adaptors in any way.

Keep the Ultravision generator and battery recharging station away from liquid. These may damage internal components.


## Chapter 2: Controls, Sockets and Indicators

This section describes all of the controls, sockets and indicators on the front and side panels of the Ultravision generator and on the Ultravision battery recharging station.

### Warning

Read all instructions before operating this equipment. Failure to do so may result in damage to the generator and/or injury to the patient or operator.

### Ultravision™ generator – Front panel view

- 1 – Power On/Off switch and indicator LED
- 2 – Patient return socket
- 3 – Ionwand socket
- 4 – Proximity indicator LEDs
- 5 – Battery level indicator
- 6 – Generator fault indicator LED

### Warning

All third-party equipment used with Ultravision must be certified according to the relevant particular standards associated with IEC 60601. The user is responsible for verifying that any accessory equipment complies with the requirements of IEC 60601.

### Be Aware

Connect accessories to the proper socket type. Otherwise, the system will not function properly.

#### 1 – Power On/Off Switch and indicator LED

The power switch that controls power to the generator is located in the lower-left corner of the front panel (see Front Panel View, above). To switch the generator on, place the power switch in the ON position (I). The indicator light should illuminate green.

Turn the generator power off by placing the power switch in the OFF position (O). The indicator light should no longer be lit.

#### 2 – Patient return socket and proximity indicator LED

The solid or split plate patient return electrode or Ultravision Patient Return Adaptor is inserted into this socket in order to establish the return path.

### Be Aware

The Ultravision Patient Return Adaptor must be inserted the correct way up. Failure to do so may cause damage to the connector or Socket.

#### 3 – Ionwand socket

The Ionwand connector is inserted into the Ionwand socket. This delivers the power to the Ionwand and causes the electrostatic precipitation of the particulates.

**Be Aware**

The Ultravision lonwand must be inserted in the correct orientation. Failure to do so may cause damage to the connector or socket.

**4 – Proximity indicator LEDs**

When the generator is switched on, the proximity indicator LEDs illuminate green. In order to function correctly, contact between the lonwand and patient tissue or other surgical instruments should be minimised. Both proximity indicator LEDs will flash when the lonwand comes within approximately 2mm (1/12”) of patient tissue. There is also an audible indicator that will sound. When the contact is removed the LEDs will stop flashing and the audible alarm will stop sounding.

**Warning**

Do not allow the lonwand to come in to contact with the uninsulated portion of an electrosurgical instrument during operation. Continued contact for in excess of 10 minutes may result in an electrical discharge from the electrosurgical instrument that may be sufficient to cause atrial fibrillation depending on proximity to the thoracic cavity.

**5 – Battery level indicator**

With a new, fully-charged battery, the battery level indicator displays the approximate power remaining as follows:

Green LED	>8hr power remaining
Amber LED	2-8hr power remaining
Red LED	<2hr power remaining
Red flashing LED	<30min power remaining

The red flashing LED is also accompanied by an audible indicator.

**Caution**

The battery should not be used for a surgical procedure when the red flashing LED is present and should be re-charged. Failure to do so may result in loss of power during the procedure, leading to impaired visibility.

**Be Aware**

The battery may lose power during transportation, storage and delivery. It is recommended that the battery be fully charged before first use of the Ultravision system.


It is recommended that the battery be replaced once it is no longer able to deliver 6 hours of power when fully charged. Contact Alesi Surgical or its local distributor to purchase replacements.

**6 – Generator fault indicator LED**

This LED is not illuminated during use. If there is a problem within the generator, the LED will illuminate red and flash. This is also accompanied by an audible indicator.

**Warning**

Do not attempt to use the Ultravision generator when the generator fault indicator is lit. Contact Alesi Surgical or its local distributor to arrange for repair or replacement of the generator unit.

**Ultravision™ generator – Side panel view**

1 – Battery insertion mechanism


**1 – Battery insertion slot and release mechanism**

The battery is inserted into the slot on the right-hand side of the Ultravision generator. The battery is ejected from the Ultravision generator by pulling the ejection trigger toward the front of the generator.

**Caution**

The battery must be inserted in the correct orientation. If not, damage to the battery and generator may result and create a fire hazard.

**Ultravision™ generator – Rear view**

1 – Volume control

**1 – Volume control for audible indicators**

The Ultravision generator has four audible indicators that accompany the visual indicators:

Indicator	Tone	Description	Frequency	Visual indicator
Low battery	High (1047 Hz)	A single 1 second tone	2 seconds	Low battery LED
Proximity	Medium (523.3 Hz)	A single 0.3 second tone	6 seconds	Proximity LEDs
Generator fault	Low (261.6 Hz)	A single 1 second tone	2 seconds	Fault (!) LED
Power ON reminder	Low (261.6 Hz)	A single 1 second tone	15 minutes	None


The generator sounds the three different tones in sequence when switched on. The volume of these indicators may be adjusted by twisting the volume control on the rear of the generator. The volume is increased by turning the volume control clockwise, and reduced by turning anti-clockwise.

**Warning**

Ensure that the audible indicator can be heard within the operating theatre environment, and that the output LED indicators are visible to the operating personnel. Failure to observe indications so may result in:

- Capacitive charging of electrosurgical instruments and the risk of atrial fibrillation.
- Loss of function caused by power failure or a fault within the generator.

**Ultravision™ Battery Recharging Station – Top view**

- 1 – Power and fault indicator
- 2 – Battery docking bay 1 and bay 2
- 3 – Battery charging status indicators
- 4 – Power supply socket

## **1 – Power/fault indicator**

The power indicator lights green when power is being correctly supplied to the recharging station. If there is a fault within the recharging station it may still charge the batteries but it may be possible to over-charge them. This is indicated by the light changing colour from green to red.

## **2 – Battery docking bay 1 and 2**

The recharging station is able to recharge up to two batteries at the same time. Batteries should be inserted such that the battery terminals engage with the terminals on the recharging station.

## **3 – Battery charging status indicators**

When power is being supplied to the recharging station and a battery is inserted in the correct orientation the status indicator will illuminate. The indicator will be orange until the battery is fully charged. When the battery is fully charged the light will extinguish.

## **4 – Power supply socket**

The battery recharging station is supplied with a power supply unit. The power supply unit has different plugs that can be attached to the power supply unit. Select the correct plug for your local power supply.

### **Be aware**

If the power supply unit supplied with the recharging station does not have a suitable plug that works with your local power supply, please contact Alesi Surgical or your local distributor.

This charger includes temperature-compensated characteristics and has been tailored for optimal charging and life preservation of only the batteries supplied with your Ultravision system. The use of other recharging stations, even if mechanically compatible, is likely to dramatically shorten the life of your Ultravision system batteries.

Only use the power supply unit provided with the battery recharging station. The use of other power supply units may damage batteries and recharging station.

When the power/fault indicator is red there is a fault with the recharger that may damage the batteries and the recharger should not be used. Contact Alesi Surgical to purchase a replacement recharger.

## Chapter 3: Preparing the Ultravision™ system for use

This section describes how to prepare the Ultravision system for use, specifically:

- Periodic inspection
- Charging a battery using the battery charging station
- Inserting and removing a battery
- Preparing the Ultravision system for use:
  - Use with instruments requiring a patient return electrode
  - Use with instruments that do not require a patient return electrode

### Warning

Read all warnings, cautions and instructions provided with the system and its accessories before use. Failure to do so could result in death serious injury to the patient or operator.

### Periodic Inspection

The Ultravision system should be visually inspected at least once every year regardless of use. This inspection should include checks for:

- Damage to the generator and recharging station.
- Damage to the battery recharging station power cable and plug.
- Damage to the patient return adaptors cable or connectors.
- Proper mating between the patient return adaptors and the Ultravision generator and electrosurgical generator.

### Warning

Do not use the Ultravision system with any damaged component. Damaged cables must be replaced before use. If there is any damage to the Ultravision generator or battery recharging station contact Alesi Surgical or its local distributor.

### Charging the battery using the Ultravision™ Battery Recharging Station

A battery should be recharged when purchased and whenever the Ultravision generator battery level indicator is red or flashing red. To fully charge a battery from zero to full power takes approximately 12 hours. A fully charged, new battery should provide at least 10 hours use of the Ultravision system.

1. Inspect the battery recharging station power cable for any signs of damage.
2. If there is no evidence of damage, attach the power cable to the recharging station.
3. Attach the correct face plate to the power cable plug and insert into a grounded, hospital-grade wall socket and confirm that the power/fault indicator light is illuminated green.
4. Insert a battery into a docking bay, ensuring that the battery terminals engage fully with the recharging station terminals. Confirm that the battery charging indicator is illuminated orange.
5. Allow the battery to charge until the charging indicator is extinguished. Remove the battery from the docking bay.
6. Unplug the power cable from the wall socket.

Batteries should be stored in a dry, clean environment when not in use.

### Warning

Connect the battery charger power cable to a properly grounded, hospital-grade power socket of the correct voltage. Plug the power cable directly into the power socket without any extension cables and/or adaptors. The use of extension cables may increase the risk of a fire hazard.

The battery recharging station has been certified for use outside the immediate operating theatre environment. It does not have sufficient redundant voltage or current protection to allow safe contact with the patient.

### Caution

Batteries must be recycled and not incinerated or disposed of in general hospital waste disposal.

Do not attempt to directly connect the negative and positive terminals of a battery.

The battery contains acid and no attempt should be made to disassemble a battery.

### Be aware

If the power supply unit supplied with the recharging station does not have a suitable face plate that works with your local power supply, please contact Alesi Surgical or your local distributor.

The recharging station includes temperature-compensated characteristics and has been designed for optimal charging and life preservation of the batteries supplied with your Ultravision system. The use of other recharging stations even if mechanically compatible is likely to dramatically shorten the operating life of your Ultravision system batteries. Do not use batteries other than those supplied by Alesi Surgical with the recharging station.

Only use the power supply unit provided with the battery recharging station. The use of other power supply units may damage batteries and recharging station.

Do not mechanically shock or drop a battery. This may damage or reduce the life of the battery.

Do not use the battery to power equipment other than the Ultravision generator. This may damage the battery and equipment.

When the power indicator is red there is a fault with the recharger that may damage the batteries and the recharger should not be used.

#### **Advice**

To maximise battery life, batteries should be returned to the recharging station at the end of use, such as at the end of each day, and should be stored in as cool an environment as is possible. Leaving a battery uncharged over a weekend or holiday break is not advised.

Do not discharge battery completely. This will greatly shorten battery life.

## **Inserting the battery into the Ultravision™ generator**

To insert a battery into the generator:

1. Locate the battery slot on the right hand side of the Ultravision generator.
2. Ensure that the battery bay is empty by pulling the battery ejection trigger towards the front of the generator.
3. Ensure that the battery is oriented with the battery terminals closest to the generator, with the arrow facing upwards and pointing to the generator.
4. Insert the battery until it engages with the contact terminals in the generator.
5. Switch on power to the generator by pressing the on switch located on the front panel.
6. Switch off the device until set up is complete.

#### **Caution**

Ensure that battery is always inserted into the Ultravision generator and battery recharging station in the correct orientation. Failure to do so may damage the battery and equipment and may cause a fire hazard.

## **Preparing Ultravision™ for use - Instruments requiring a patient return electrode**

Before use, confirm that the following components are available:

- Ultravision generator with an adequately charged battery
- Patient return electrode with "international" style connector: solid or split type
- Correct Ultravision Patient Return Adaptor: solid or split
- Electrosurgical generator unit compatible with the patient return electrode

Set up the system as follows:

1. Connect the patient return electrode to the patient using standard surgical practice. Either solid or split variants may be used.
2. Connect the patient return electrode to the patient return adaptor. Select the solid patient return adaptor for use with a solid return electrode. Select a split plate adaptor for use with a split plate.
3. Connect the patient return adaptor end (BLUE – solid plate; RED – split plate) to the electrosurgical generator.
4. Switch on the electrosurgical generator and ensure that the patient return path is recognised by the electrosurgical generator.
5. Connect the other end of the patient return adaptor (ORANGE) to the Ultravision generator. The connector must be inserted in the correct orientation.

If the electro-surgical generator does not recognise the return electrode when using the patient return adaptor, follow the following procedure:

- Attach the patient return electrode directly into the electro-surgical generator. If this resolves the problem the patient return adaptor is worn and should be replaced.
- If this does not resolve the problem then replace the patient return electrode and repeat steps 1-4 above.
- If this does not resolve the problem then there is a problem with the electro-surgical generator.

#### **Warning**

Ensure that the selected patient return adaptor (solid or split) matches the patient return electrode that is attached to the patient. The patient return adaptor allows a solid patient return electrode to be inserted into a split patient return adaptor. Incorrect use may bypass any return pad sensing functionality of the electro-surgical generator.

#### **Advice**

If the surgical procedure may require both instrument types i.e. those that require a return pathway (e.g. monopolar) and those that do not (e.g. bipolar or ultrasonic) it is recommended that the system is set up in this way. This avoids the need to delay the procedure whilst reconfiguring Ultravision so that it can share the same patient return electrode.

The Ultravision system has been verified to have no effect on the split plate patient return contact monitoring systems employed in the monopolar generators listed on page 21.

### **Preparing Ultravision™ for use – Instruments that do not require a patient return path**

Before use, confirm that the following components are available:

- Ultravision generator with charged battery
- Patient return electrode with “international” style connector (solid or split type)

Set up the system as follows:

1. Connect the patient return electrode to the patient using standard surgical practice. Either solid or split variants may be used.
2. Connect the other end of the patient return electrode directly to the Ultravision generator by inserting it into the return socket on the face of the Ultravision generator. The connector must be inserted in the correct orientation.

#### **Advice**

Patient return adaptors should be made readily available in case the surgical procedure demands the use of an instrument that requires a patient return electrode.

## Chapter 4: Using the Ultravision™ system

This section describes the use of the Ultravision system to clear the visual field of surgical smoke. This step requires the following consumable items:

- Ultravision catheter and trocar
- Ionwand

These items are supplied separately.

### Inserting the Catheter and Ionwand

The catheter is inserted using the trocar. The catheter and trocar are supplied assembled and ready for use.

The peritoneal cavity must be insufflated using standard practice before attempting to insert the catheter. Once insufflation has been achieved:

1. Open the Ionwand sterile packaging.
2. Remove the catheter/trocar assembly and the Ionwand from the packaging and check the integrity of the components.
3. Introduce the catheter by piercing the abdominal wall with the catheter/trocar assembly:
  - a. Always insert after the abdomen has been insufflated.
  - b. Always insert using direct vision to avoid iatrogenic injury to abdominal structures.
4. Remove the trocar from the catheter by holding the catheter whilst rotating the trocar anti-clockwise and dispose of safely.
5. Carefully insert the Ionwand down the catheter until it protrudes from the tip and can be visualised internally within the patient with the laparoscopic camera.
6. Secure the Ionwand in place by holding the catheter and rotating the hub of the Ionwand clockwise until it tightens. Do not overtighten.
7. Adjust the height of the catheter so that the Ionwand is free from contact with patient tissue and surgical instruments.
8. Connect the Ionwand connector to the active socket on the face of the Ultravision generator.

#### Warning

If the sterile packaging is open or damaged, do not use and discard.

If there is evidence of damage to any of the components, do not use and discard.

Always insert the catheter and trocar under direct vision once the abdominal cavity has been insufflated. Failure to do so may result in damage to internal tissue, blood vessels or organs.

Avoid major arteries and veins when inserting the catheter.

Do not switch on the Ultravision generator before inserting the Ionwand into the catheter.

#### Caution

The sterile components have been design specifically for use with this system. Only use items supplied by Alesi Surgical.

#### Advice

Insert the catheter in a position where it will not impair the surgical procedure and where it is least likely to come in to contact with other surgical instruments. Smoke clearing is most efficient when placed between the laparoscope and the site of surgery.

Ensure that the tip of the catheter is at least 25mm (1") from patient tissue to reduce the risk of accidentally activating the proximity indicator.

Cover the lumen of the catheter with a finger prior to inserting Ionwand to minimise loss of pneumoperitoneum.

Ensure that the Ionwand brush is straight before insertion into the catheter.

## Use during surgery

Before starting to cut or coagulate tissue and creating smoke:

9. Ensure that the lonwand is free of surgical debris and not touching patient tissue or a surgical instrument.
10. Switch on the Ultravision generator. If the battery status light indicator is flashing red and the battery alarm is sounding, replace the battery with a fully charged unit.
11. Perform surgical procedure.
12. In the event that the proximity indicator sounds when the lonwand™ is not touching patient tissue or a surgical instrument, perform the following steps:
  - a. Switch off the Ultravision generator.
  - b. Remove the lonwand™ from the catheter by holding the catheter and twisting the lonwand anti-clockwise, ensuring the catheter is not displaced.
  - c. Clean the lonwand™ and brush filaments by enclosing it within a sterile swab moistened with sterile saline and gently wiping it, moving the swab from the insulated portion of the lonwand towards the metal brush. Repeat until all contamination has been removed and the brush is clean.
  - d. Re-introduce the lonwand into the catheter and secure by holding the catheter and twisting the lonwand in a clockwise direction.
  - e. Switch on the Ultravision generator and confirm that the indicator condition has resolved and that smoke clearing has been resumed.
  - f. If the alarm continues to sound, remove and replace the lonwand.

### Warning

To maximise patient safety, contact between the lonwand and metallic instruments should be minimised. Otherwise, the metal instrument may be charged resulting in an electrostatic discharge to the patient or operator. This may be sufficient to cause atrial fibrillation depending upon proximity to the thoracic cavity.

Accidental activation of a monopolar electrode in contact with the lonwand should be avoided. Provided the lonwand cable is connected to the Ultravision generator there is no risk of accidental burn or shock to the patient or surgical staff, but there is no proven redundancy in this protective insulation.

If the Ultravision catheter becomes disengaged from its position during use replace with a new device.

## End of surgery

When the surgical procedure has been completed or there is no longer any requirement for the Ultravision to be active:

13. Turn off the Ultravision generator.
14. Disconnect the lonwand cable from the generator.
15. Carefully remove the lonwand from the catheter.
16. Remove the catheter and apply dressing to the wound site.
17. Dispose of lonwand cable and catheter safely.

### Warning

Ensure that the Ultravision generator is switched off before removing the lonwand. Failure to do this may lead to electrical charging of conductive objects and create a risk of static shock to the operator.

Do not reuse or re-sterilise the lonwand, trocar or catheter. Doing so may result in cross contamination, infection or injury to the patient or medical staff, or damage the Ultravision system.

Dispose of the lonwand, trocar and catheter along with other operative waste materials following the procedures for your institution.

### Advice

Removing the catheter as the last step in the procedure allows a surgeon to close the other surgical port sites without rapid loss of peritoneal pressure.

## Chapter 5: After surgery

This chapter describes how to clean the reusable components of the Ultravision system and how to recharge the battery when required.

The reusable components should be cleaned with a mild detergent that is approved and routinely used in your hospital for cleaning surgical generators and non-sterile medical cables and connectors.

### Cleaning the Ultravision™ Patient Return Adaptor

The patient return adaptors may require cleaning from time to time.

The cables and surfaces of the connectors should be wiped with a mild detergent.

Cables and connectors should be allowed to air dry thoroughly before being reused.

### Cleaning the Ultravision™ Generator

The generator may require cleaning from time to time:

1. Ensure that the generator is switched off ("O").
2. Eject the battery and remove.
3. Ensure that all cables are detached from the generator.
4. Wipe all external surfaces thoroughly with a soft cloth that has been dampened with a mild detergent.
5. Ensure that the generator is dry before reusing.

### Cleaning the Ultravision™ Battery Recharging station

The recharging station may require cleaning from time to time:

1. Ensure that the recharging station is unplugged from the power supply.
2. Ensure that the battery charging bays are empty.
3. Wipe the external surfaces thoroughly with a soft cloth that has been dampened with a mild detergent.
4. Ensure that the recharging station is dry before reconnecting the power supply.

#### Warning

Always disconnect the patient return adaptor from the electrosurgical generator and Ultravision generator before cleaning.

Always turn off and unplug the Ultravision generator before cleaning.

Allow the patient return adaptor and Ultravision generator to air dry thoroughly before reusing.

#### Be Aware

Do not attempt to clean the patient return adaptor or Ultravision generator with alcohol or caustic, corrosive or abrasive cleaning or disinfectant compounds, solvents, or any other material that could scratch or cause damage.

Do not attempt to sterilise the patient return adaptor or Ultravision generator.

Keep the patient return adaptor and Ultravision generator away from liquids. Liquids that enter the generator may damage internal components. Do not spray liquid cleaner into the Sockets of the Ultravision generator. Only apply appropriate cleaning liquids when applied to a cloth.

### Recharging the Ultravision™ battery

Follow the procedure described in Chapter 3 to recharge the battery.

## Chapter 6: Troubleshooting

This Chapter considers situations that may arise and suggests possible causes and recommended actions.

Situation	Possible cause(s)	Recommended action
<b>Generator does not power up when switched on</b>	No battery power remaining	Replace with fully charged battery
	Faulty generator	Contact manufacturer or local distributor
<b>Fully charged battery does not support generator operation for expected duration.</b>	Battery at end of life through use or damage.	Replace with new battery
<b>Battery not charging</b>	Power cable not supplying power to recharging station	Confirm battery recharger power indicator is green.
	Battery not inserted correctly	Confirm that battery is inserted correctly with the terminals in contact with the terminals on the battery recharging station.
	Battery exhausted	Purchase replacement battery.
<b>Battery recharger indicators (x3) not illuminated</b>	No connection to wall adaptor	Check connection
	Wall adaptor not connected to an acceptable mains voltage supply	Check mains supply acceptable
	Battery recharger or wall adaptor damaged	Contact manufacturer or local distributor
<b>Battery recharging lights fail to illuminate orange when battery still discharged</b>	Battery not properly seated in battery recharger cradle	Reposition battery until recharging light changes to orange
	Battery recharger damaged	Contact manufacturer or local distributor
<b>Battery recharging light(s) continuously illuminated orange after 12 hours of recharging</b>	Battery damaged and continuously dissipating charging energy	Replace with new battery
	Battery recharger or wall adaptor damaged	Replace with new recharger
<b>Electrosurgical generator does not recognise patient return electrode when using patient return adaptor <i>OR</i> No smoke clearing and electrosurgical instrument not working</b>	Incorrect patient return adaptor	Ensure chosen solid or split adaptor matches the patient return electrode being used in the procedure.
	Defective patient return electrode	Connect return electrode directly to electrosurgical generator – if the return electrode is recognised it is <b>NOT</b> defective.
	Defective patient return adaptor	Perform previous two actions. If these do not resolve the issue, replace the patient return adaptor.
<b>No smoke clearing but electrosurgical instrument working</b>	Ultravision generator switched off or has no power	Check generator is switched on and charged battery correctly inserted.
	Ultravision system incorrectly configured	Ensure system configured correctly for the electrosurgical instrument being used.


	Defective lonwand	Confirm patient return path is correctly configured. Deliberately touch patient tissue with the lonwand. If the proximity alarm <b>DOES NOT</b> sound, replace the lonwand and repeat.
	Defective patient return path	If above action does not resolve the problem replace the patient return adaptor and patient return electrode.
<b>Proximity alarm sounding frequently</b>	lonwand is too close to patient tissue	Locate lonwand with laparoscope and confirm it is not too close to patient tissue. If too close, reposition the lonwand.
	lonwand is too close to instrument or port	Locate lonwand with laparoscope and confirm it is not too close to port or trocar. If too close, reposition the lonwand.
	lonwand partially contaminated with fluid or smoke particles	Remove the lonwand and wipe with sterile surgical wipe. Refer to "Use during surgery" in Chapter 4: Using the Ultravision™ system
	lonwand is defective	If above actions do not resolve the situation, replace the lonwand.
<b>Proximity alarm sounding constantly</b>	Rate of smoke production is too high to allow a smoke-precipitating electric field to be maintained.	Reduce the rate of smoke production by cutting more slowly.
	lonwand is heavily contaminated with smoke particles	Remove the lonwand and wipe with sterile surgical wipe. Refer to "Use during surgery" in Chapter 4: Using the Ultravision™ system.
	Catheter is contaminated with smoke particles and is providing a return path from the lonwand	Insert the catheter further into the abdominal cavity.
	lonwand is defective	If above actions do not resolve the situation, replace the lonwand.

## Chapter 7: Summary System Specifications

Item	Specification
Storage requirements	
Generator	-10°C to 40°C at maximum 85% relative humidity
Battery recharging station	-10°C to 40°C at maximum 85% relative humidity
Battery	-10°C to 40°C at maximum 85% relative humidity
Patient return adaptor	-10°C to 40°C at maximum 85% relative humidity
Sterile consumables	-10°C to 40°C at maximum 85% relative humidity
Operating requirements	
Generator	+10°C to 40°C at 20-75% relative humidity
Battery recharging station	+10°C to 40°C at 20-75% relative humidity
Battery	+10°C to 40°C at 20-75% relative humidity
Patient return adaptor	+10°C to 40°C at 20-75% relative humidity
Sterile consumables	+10°C to 40°C at 20-75% relative humidity
Compatible patient return electrodes	Valleylab Adult Non-REM Polyhesive (Catalogue No. E7506) Valleylab Adult REM Polyhesive II (Catalogue No. E7507) Erbe MONOPlate™ Single Pad (Catalogue Number 20193-076) Erbe NESSY (Catalogue No. 20193-083 (-084 in USA))
Compatible electrosurgical generators	Valleylab Force Fx Erbe VIO 300D Karl Storz Autocon II 400 Olympus ESG 100
Battery capacity and operating life	12 Volt, 2.0 Amp hours Estimated shelf life of 12 months from purchase or 200 charging cycles at 20C
Recharger input voltage range	90-264V AC
Maximum battery recharger current consumption	0.5A at 230V AC
Battery run time (hours)	Approximately 10 hours
Number of battery recharges before replacement	200 at 20°C
Generator and battery recharging station EMC compliance standards	EN 60601-1-2: 2007 (See Chapter 9)
Transit packaging rating	ISTA Class 2A
Ingress protection rating for generator	1litre saline poured onto top surface of generator (IEC 60601-2-2: 2009 5 <sup>th</sup> Edition, clause 201.11.6.3).
Patient Return Adaptor rating	Comply with IEC60601-2-2:2009 ed. 5th clause 201.15.101.2 and 201.15.101.3.
Generator maximum audible alarm (dBA)	60dBA
Insulation type/shock prevention	Class 2
Weight/Dimensions: Generator	2.43kg; 395mm x 350mm x 100mm
Weight/Dimensions: Battery recharger	175g; 140mm x 100mm x 50mm

## Generator output characteristics

The graph below illustrates the relationship between voltage and current at the end of the lonwand cable.

The operating current is influenced by the distance of the lonwand from the patient return path (patient tissue). As the lonwand approaches the patient tissue the impedance (resistance) falls. This causes the current to increase and the voltage to fall. The generator monitors the current and an active feedback mechanism collapses the current as it approaches 10µA, which is the maximum DC current that can safely be applied to a patient. As a result the voltage falls below a level that is sufficient to cause electrostatic precipitation. This is accompanied the proximity indicator, which informs the operator that the lonwand is touching or very close to the patient tissue and requires repositioning.

## Chapter 8: Accessories and replacement items

The following replacement items are available from Alesi Surgical or its local distributor:

Item	Number per pack	Product code
Ultravision™ Battery Recharging Station (recharging station, power supply unit and power cable)	1	DAD-001-024
Ultravision™ Power Supply Unit (power supply unit and power cable)	1	DSD-001-035
Ultravision™ Battery	2	DSD-001-034
Ultravision™ Patient Return Adaptor (solid)	1	DAD-001-006
Ultravision™ Patient Return Adaptor (split)	1	DAD-001-007
Ultravision™ lonwand™ Sterile Pack	10	DAD-001-003

## Chapter 9: Electromagnetic guidance

### Electromagnetic emissions – Ultravision Generator

The Ultravision generator is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultravision generator should assure that it is used in such an environment.

Emissions test Compliance Electromagnetic environment – guidance	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Ultravision generator uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. RF emissions
RF emissions CISPR 11	Class A	The Ultravision generator is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

### Electromagnetic immunity – Ultravision Generator

The Ultravision generator is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultravision generator should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Not applicable		
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable		
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Electromagnetic immunity – Ultravision Generator

The Ultravision generator is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultravision generator should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment – guidance
<b>Conducted RF</b> <b>IEC 61000-4-6</b>  <b>Radiated RF</b> <b>IEC 61000-4-3</b>	3 V <sub>rms</sub> 150 kHz to 80 MHz  3 V/m 80 MHz to 2.5 GHz	3V  3V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Ultravision generator, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = 1.17\sqrt{P}$ $d = 1.17\sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = 2.33\sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Ultravision generator is used exceeds the applicable RF compliance level above, the Ultravision generator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Ultravision generator.			
<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

## Recommended separation distances between portable and mobile RF communications equipment and the Ultravision Generator

The Ultravision generator is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Ultravision generator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Ultravision generator as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.68	3.68	7.37
100	11.66	11.66	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### Electromagnetic emissions – Ultravision Battery Recharging Station

The Ultravision battery recharging station is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultravision battery recharging station should assure that it is used in such an environment.

Emissions test Compliance Electromagnetic environment – guidance	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Ultravision battery recharging station uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. RF emissions
RF emissions CISPR 11	Class A	The Ultravision battery recharging station is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

### Electromagnetic immunity – Ultravision Battery Recharging Station

The Ultravision battery recharging station is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultravision battery recharging station should assure that it is used in such an environment.

IMMUNITY test	IEC 60601test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 s	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Ultravision battery recharging station requires continued operation during power mains interruptions, it is recommended that the Ultravision battery recharging station be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Electromagnetic immunity – Ultravision Battery Recharging Station

The Ultravision battery recharging station is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultravision battery recharging station should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	3V	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Ultravision battery recharging station, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = 1.17\sqrt{P}$ $d = 1.17\sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = 2.33\sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup> should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Ultravision battery recharging station is used exceeds the applicable RF compliance level above, the Ultravision battery recharging station should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Ultravision battery recharging station.

<sup>b</sup>Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

## Recommended separation distances between portable and mobile RF communications equipment and the Ultravision Battery Recharging Station

The Ultravision battery recharging station is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Ultravision battery recharging station can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Ultravision battery recharging station as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.68	3.68	7.37

100	11.66	11.66	23.33
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance <math>d</math> in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

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