LANGE 

The Lange logo symbol is a red square containing a white stylized target or bullseye design.

DOC022.98.90248

sensION™ + PH 3

Manual
Bedienungsanleitung
Manuel d'utilisateur

07/2010, Edition 1

07/2010, Heft 1

07/2010, Edition 1

ENGLISH

DEUTSCH

FRANÇAIS

CONTACT

Index

1. Technical Specifications	2
2. General information	3
2.1. Safety information	3
2.1.1. Use of hazard information	3
2.1.2. Precautionary labels	3
2.2. General product information	3
2.3. Instrument composition	3
3. Installation	4
3.1. Assemblage	4
3.2. Connections	5
4. Operation	5
4.1. Description	5
4.2. Start-up	6
4.3. General diagram	6
4.4. pH and mV calibration	7
4.5. pH and mV measurements	8
4.6. Other options	8
4.7. Recognized buffers	9
5. Maintenance	9
6.1. Clean the instrument	9
6.2. Sensor cleaning	9
6. Warning messages	10
7. Accessories and spares	10
8. Warranty, liability and complaints	11

Contact Information

1. Technical Specifications

Measuring ranges:

pH -2.00 to 19.99

mV ± 2000

Temp -20 to 150 °C (-4 to 302 °F)

Resolution

pH 0.01

mV 1

Temp 0.1 °C (0.18 °F)

Measuring error (± 1 digit)

pH ≤ 0.01

mV ≤ 1

Temp ≤ 0.2 °C (≤ 0.4 °F)

Reproducibility (± 1 digit)

pH ± 0.01

mV ± 1

Temp ± 0.1 °C (± 0.1 °F)

Automatic temperature compensation

Manual entrance.

With Pt 1000 temperature probe (A.T.C.).

With NTC 10 K Ω probe.

pH calibration

With 1, 2 or 3 buffers selectable inside the range.

Automatic recognition of technical buffers pH 2.00, 4.01, 7.00, 9.21, 10.00, values at 25 °C (77 °F).

Programmable calibration validity between 0 h and 7 days.

Automatic recalibration warning.

Automatic rejection of electrodes in poor conditions.

Theoretical calibration.

pH calibration criteria

Asymmetry potential ± 58 mV.

Slope 41 to 62 mV/pH, at 25 °C (77 °F)

(sensitivity 70 to 105 %).

mV calibration

Automatic recognition of redox standard 220 mV at 25 °C (77 °F).

Measuring mode

By stability.

Continuous.

Languages

English, German, Spanish, French, Italian, and Portuguese.

Display

Graphic, backlit liquid crystal, 128 x 64 dots.

Inputs and outputs

Combined or indicator electrode, BNC connector (Imp. $>10^{12}\Omega$).

Reference electrode, banana connector.

A.T.C. type Pt 1000, banana or telephonic connector.

Stirrer, RCA connector.

Ambient conditions

Working temperature 5 to 40 °C (41 to 104 °F).

Storage temperature -15 to 65 °C (5 to 149 °F).

Relative humidity $< 80\%$ (not condensed).

Power supply

External plug-in power supply 100-240 VAC 0.4 A 47-63 Hz.

Materials

Enclosure, ABS and PC.

Keypad, PET with protective treatment.

Physical parameters

Weight: 1100 g.

Size: 350 x 200 x 110 mm.

2. General information

As a result of constant improvements to our products sometimes differences may exist between this manual and the instructions supplied with the instrument.

2.1. Safety information

Please, read carefully this information before installing and using the instrument !

Pay attention to all danger and caution statements.

2.1.1. Use of hazard information

DANGER

Indicates a potentially or imminently hazardous situation that, if not avoided, will result in death or serious injury.

WARNING

Indicates a potentially or imminently hazardous situation that, if not avoided, may result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation that, if not avoided, may result in minor or moderate injury.

Important note: Indicates a situation that, if not avoided, could lead to damage to the instrument. Important information that requires special emphasis.

Note: Information that supplements points in the main text.

2.1.2. Precautionary labels

Read carefully all labels and tags attached to the instrument.

Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of life equipment to the Producer for disposal at no charge to the user.

Note: For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.

2.2. Description

The **sensION™ + PH3** is an instrument to measure pH, mV or temperature.

2.3. Instrument composition

Version	Electrode	Accessories	Manual
LPV2000.98.0002	—	✓	✓
LPV2010T.98.002	LZW5010T.97.002	✓	✓
LPV2011T.98.002	LZW5011T.97.002	✓	✓
LPV2014T.98.002	LZW5014T.97.002	✓	✓
LPV2021T.98.002	LZW5021T.97.002	✓	✓

Accessories:

- Buffer solutions, 135 ml flasks.
- Magnetic stirrer.
- Calibration flasks, with magnetic bar inside.
- Electrolyte for the electrode, KCl 3M, 50 ml flask.
- Plug-in power supply.
- Electrode's stick and support.

3. Installation

⚠ DANGER

Only qualified personnel, following the local security regulations, should conduct the tasks described in this section of the manual.

3.1. Assemblage

1. Unpack the instrument.

2. Verify that the shift is in the same position than in the picture.

6. Place the electrode on the electrode's support. Insert the electrode's cable through the groove.

3. Drive the stick into the hole.

Note: Stirrer replacement

Before replacing the stirrer, disconnect it from the rear panel and release the cable from the bottom of the instrument.

4. Turn the shift in a counterclockwise until the stick is well fixed.

Place the o-ring at 135-140 mm from the base in order to avoid the shock of the sensor against the bottom of the glass.


1. Turn the stirrer to release or fix from base.

5. Place the electrode's support on the stick by softly pressing the button.

Holder for 3 sensors.
PN= LZW9321.99

2. Pull up the stirrer and replace it. Place the cable of the new stirrer through the bottom groove of the instrument.


3. Installation

3.2. Connections

DANGER

Make sure the power adapter matches your local power supply.

1. Magnetic stirrer, RCA connector.
2. Combined pH electrode (or indicator), BNC connector.
3. Reference electrode (measuring with separated electrodes), banana connector.
4. Temperature probe for the electrodes version **T**, banana connector.
5. Separated temperature probe, telephone connector.
6. Plug-in power supply (12 V).

4. Operation

4.1. Description

Keys

Switch on / Switch off.

pH measurement.

mV measurement.

Calibration and calibration data.

Select between temperature (°C, °F) and stirrer (St) to modify the temperature or the stirring speed.

Selection of numeric values.

Access to:

- Language.
- Date and hour.
- Display contrast.
- Calibration frequency.

4. Operation


4.2. Start-up

Connect to the rear panel of the instrument:


- the pH electrode
- the temperature probe (A.T.C.) built-in the pH electrode or as separated sensor
- the magnetic stirrer
- the plug-in power supply

Follow the next steps:

For the first time...


Successive times...

4.3. General diagram

This is a general view over the way to access to all possibilities of the instrument:

- DISPLAY CONTRAST
- LANGUAGE
- DATE AND TIME
- CHANGE OF TEMPERATURE UNITS
- pH AND mV CALIBRATION
- pH AND mV MEASUREMENT

4. Operation

4.4. pH and mV calibration

pH Calibration

To measure pH correctly, the electrode and the instrument set needs to be calibrated regularly with buffer solutions. This compensates deviations of asymmetry potential and slope of the electrode over time i.e. the calibration corrects variations in the electrode response.

The instrument allows one-, two- and three- point calibration. The calibration parameters remain in memory until a new calibration is performed.

One-point calibration

This type of calibration is acceptable when measuring pH values similar to those of the buffer solution.

It corrects the asymmetry potential of the electrode.

Two-point calibration

This is the most common calibration.

It is recommended to start with pH 7 as first buffer. pH 4 or pH 9 buffers can be used as a second buffer, depending on the measured zone (acid or alkaline).

The two-point calibration compensates the asymmetry potential and the electrode's sensitivity (slope).

Three-point calibration

This type of calibration is recommended when measurements are made across the whole pH scale.

As a first point pH 7 is recommended.

For the second and third points, two of the remaining values must be chosen (2.00, 4.01, 9.21, 10.00 at 25°C (77 °F)).

The three-point calibration compensates the asymmetry potential and electrode's sensitivity in the acid (a) and alkaline zone (b).

mV calibration

The metallic electrodes do not suffer significant variation in their potentials, for this reason normally they are not calibrated. Nevertheless, some potential deviations can be produced due to the continuous use of the electrode.

The **sensION™ + PH3** allows calibration of the metallic electrodes, with platinum or gold indicator metal, using redox standard of 220 mV (value at 25°C (77 °F)).

Calibration data

They can be viewed at any moment. See diagram.

Calibration frequency

Select between 0 h and 7 days. The instrument by default is programmed for daily calibration and indicates the remaining time to the new calibration. Selecting 0 h the automatic re-calibration warning is deactivated.

The calibration frequency depends on the required accuracy and the effect of the samples on the electrode. Each recommends daily calibration, but users should decide on the right time based on their own experience.

Theoretical calibration

Selecting "theoretical calibration" the calibration data of the electrode in use are replaced by the theoretical values at 25°C (77 °F).

This step is prior to the obtention of a calibration certificate.


Preparation

Fill the flasks with the corresponding buffer solution (40 ml.)

These buffers can be used to perform between 5 and 10 calibrations.

Every calibration flask is supplied with a magnetic bar.

Diagram for calibration

4. Operation

4.5. pH and mV measurements

Measurement by stability

The value is fixed on the display when during 6 seconds the reading variation is less than 0.01 pH (0.5 mV).


Continuous measurement

The instrument displays the value obtained by the electrode all the time.

Warning:


The temperature resolution is 1°C (1 °F) when no temperature sensor is connected. The resolution is 0.1°C (0.1 °F) when there is a temperature sensor connected.

During mV measurement if temperature sensor is not connected the temperature is not shown.


4.6. Other options

- DATE AND TIME
- LANGUAGE
- DISPLAY CONTRAST

- CHANGE OF TEMPERATURE UNITS, °C to °F.

4. Operation

4.7. Recognized buffers

Buffer solutions

The **sensION™ + PH3** has in memory the following data (pH values as a function of the temperature).

Technical buffers (DIN 19267)

°C	°F	pH					mV
0	32	2.01	4.01	7.12	9.52	10.30	--
10	50	2.01	4.00	7.06	9.38	10.17	245
20	68	2.00	4.00	7.02	9.26	10.06	228
25	77	2.00	4.01	7.00	9.21	10.01	220
30	86	2.00	4.01	6.99	9.16	9.96	212
40	104	2.00	4.03	6.97	9.06	9.88	195
50	122	2.00	4.06	6.97	8.99	9.82	178
60	140	2.00	4.10	6.98	8.93	9.76	160
70	158	2.01	4.16	7.00	8.88	--	142
80	176	2.01	4.22	7.04	8.83	--	--
90	194	2.01	4.30	7.09	8.79	--	--

6. Maintenance

6.1. Clean the instrument

Important note: Never use cleaning such as turpentine, acetone or similar products to clean the instrument including the display and accessories.

Only clean the housing and accessories using a soft, damp cloth. Mild soap solution may also be used. Dry the cleaned parts carefully with a soft cotton cloth.

CAUTION

When using chemicals or solvents, comply with the instructions of the producer and all local safety regulations.

6.2. Sensor cleaning

See the manual of the sensor in use.

6. Warning messages

During calibration

CALIBRATION 1st buffer
 ATENTION !
 Asymmetry > 58 mV
 Repeat with 

CALIBRATION 1st buffer
 ATENTION !
 Buffer or electrode
 in poor conditions.

CALIBRATION 2nd buffer
 ATENTION !
 Sens. (a) < 70 %
 Repeat with 

CALIBRATION X buffer
 UNSTABLE READING
 Time t> 100 s
 Repeat with 

CALIBRATION X buffer
 ATENTION !
 Check temperature.
 Check the electrode.

CALIBRATION X buffer
 BUFFER NOT RECOGNISED
 Repeat with 

CALIBRATION X buffer
 ATENTION !
 Electrode in poor
 conditions.

CALIBRATION 2/3 buffer
 SAME BUFFERS
 Repeat with 

CALIBRATION 2/3 buffer
 ATENTION !
 Check the electrode.

During measurement

MEASURE STABLE
 pH 12.78 19°C
 Calibrate again

MEASURING by stability
 pH out of range
 Check the electrode.

MEASURING by stability
 Temp out of range °C
 Check temp. sensor

MEASURING by stability
 UNSTABLE READING
 Time > 150 s
 Check temperature.
 Check the electrode.

7. Accessories and spares

Code	Description
LZW9110.98	Set of 3 x 50 ml printed flasks for benchtop pH-meters calibration
LZW9322.99	Magnetic stirrer
LZW9400.99	Standard redox solution of 220 mV, 250 ml flask.
LZW9463.99	pH buffer solution 4.01, 250 ml flask.
LZW9464.98	pH buffer solution 7.00, 250 ml flask.
LZW9465.99	pH buffer solution 9.21, 250 ml flask.
LZW9500.99	Electrolytic solution, KCl 3M, 250 ml flask.
LZW9321.99	Three-sensor holder.

For electrodes please ask for specific brochure.

8. Warranty, liability and complaints

The manufacturer warrants that the product supplied is free of material and manufacturing defects and undertakes the obligation to repair or replace any defective parts at zero cost.

The warranty period for instruments is 24 months. If a service contract is taken out within 6 months of purchase, the warranty period is extended to 60 months.

With the exclusion of the further claims, the supplier is liable for defects including the lack of assured properties as follows: all those parts that can be demonstrated to have become unusable or that can only be used with significant limitations due to a situation present prior to the transfer of risk, in particular due to incorrect design, poor materials or inadequate finish will be improved or replaced, at the supplier's discretion. The identification of such defects must be notified to the supplier in writing without delay, however at the latest 7 days after the identification of the fault. If the customer fails to notify the supplier, the product is considered approved despite the defect. Further liability for any direct or indirect damages is not accepted.

If instrument-specific maintenance and servicing work defined by the supplier is to be performed within the warranty period by the customer (maintenance) or by the supplier (servicing) and these requirements are not met, claims for damages due to the failure to comply with the requirements are rendered void.

Any further claims, in particular claims for consequential damages cannot be made.

Consumables and damage caused by improper handling, poor installation or incorrect use are excluded from this clause.

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DR. LANGE NEDERLAND B.V.