www.baylismedical.com

Instructions for Use

Toronto Transseptal Catheter

English..... 1

Français..... 4

Deutsch 7

Italiano 11

Español..... 14

Português.....17

Baylis Medical Company Inc.
 5959 Trans-Canada Highway
 Montreal, Quebec, Canada, H4T 1A1
 Tel: (514) 488-9801/ (800) 850-9801 Fax: (514) 488-7209
 www.baylismedical.com

EU Authorized Representative:
 Quality First International
 20 Eversley Road, Bexhill-on-Sea, East Sussex
 TN40 1HE, United Kingdom
 Tel: +44-(20)-8-522-1937 Fax: +44-(20)-8-522-1937

© Copyright Baylis Medical Company Inc., 2007-2015

Baylis Medical Company logo is a trademark and/or registered trademark of Baylis Medical Company Inc. in the United States of America and/or other countries. Patents pending and/or issued.

English

Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these instructions. Failure to do so may result in patient complications.

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

I. DEVICE DESCRIPTION

The Toronto Transseptal Catheter delivers radiofrequency (RF) power in a monopolar mode between its distal electrode and a commercially available external Disposable Indifferent (Dispersive) Patch (DIP) Electrode, which is in compliance with IEC 60601-2-2. The Toronto Transseptal Catheter is loaded through a Transseptal Sheath/Dilator set, and is connected at its proximal end to the BMC Radiofrequency Puncture Generator via the BMC Connector Cable and to an external pressure monitoring system via a luer connection. Detailed information concerning the BMC Radiofrequency Puncture Generator is contained in a separate manual that accompanies the Generator (entitled "BMC Radiofrequency Puncture Generator Instructions for Use").

The dimensions for the Toronto Transseptal Catheter can be found on the device label. The distal end of the catheter contains holes to facilitate injection of contrast solution and the monitoring of cardiac pressures. As well, the distal portion of the catheter is flexible and the active tip is specially shaped to be atraumatic to the cardiac tissue unless RF energy is applied. When used with the Baylis Medical TorFlex Transseptal Guiding Sheath, the markings on the proximal end of the catheter indicate the location of the active tip in relation to the dilator tip.

II. INDICATIONS FOR USE

The Toronto Transseptal Catheter is used to create an atrial septal defect in the heart. Secondary indications include monitoring intracardiac pressures, sampling blood, and infusing solutions.

III. CONTRAINDICATIONS

The Toronto Transseptal Catheter is not recommended for use with any conditions that do not require cutting or coagulation of soft tissue.

IV. WARNINGS

- Only physicians with a thorough understanding of angiography and percutaneous interventional procedures should use this device.
- The Toronto Transseptal Catheter is supplied STERILE using an ethylene oxide process. Do not use if the package is damaged.
- Laboratory staff and patients can undergo significant x-ray exposure during radiofrequency puncture procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure.
- The Toronto Transseptal Catheter is intended for single patient use only. Do not attempt to sterilize and reuse the catheter. Reuse can cause the patient injury and/or the communication of infectious disease(s) from one patient to another. Failure to do so may result in patient complications.
- The Toronto Transseptal Catheter must be used with the BMC Connector Cable. Attempts to use it with other connector wires can result in electrocution of the patient and/or operator.
- For RFP-100: Do not attempt to puncture with an initial power setting of greater than 10 Watts. The initial attempt should be made with a setting of 10 Watts. In subsequent punctures, the power setting can be increased, if necessary.
- The active tip of the Toronto Transseptal Catheter is fragile. Be careful not to damage the tip while handling the catheter. If the tip becomes damaged discard the catheter immediately.
- The pressure transducer system used with the Toronto Transseptal Catheter must comply with the electrical safety requirements of IEC 60601. Failure to use compliant pressure transducers may result in patient or operator injury.
- If the active tip of the catheter becomes bent at any time during its use, dispose of the catheter immediately. Do not attempt to straighten the active tip.

V. PRECAUTIONS

- Do not attempt to use the Toronto Transseptal Catheter or ancillary equipment before thoroughly reading the accompanying Instructions for Use.
- Radiofrequency puncture procedures should be performed only by physicians thoroughly trained in the techniques of radiofrequency powered puncture in a fully equipped catheterization laboratory.

- The sterile packaging should be visually inspected prior to use to detect any compromise. Ensure that the packaging has not been damaged. Do not use the equipment if the packaging has been compromised.
- Visually inspect the catheter prior to use. Do not use the catheter if there is any damage.
- Do not use the Toronto Transseptal Catheter after the "Use By" date indicated on the label.
- The Toronto Transseptal Catheter is intended for use with only those devices listed in section VII "Equipment Required."
- Read and follow the manufacturer's instructions for use of the Disposable Indifferent (Dispersive) Patch (DIP) electrode. Always use DIP electrodes that meet or exceed IEC 60601-2-2 requirements.
- Placement of the dispersive electrode on the thigh could be associated with higher impedance.
- In order to prevent the risk of ignition make sure that flammable material is not present in the room during RF power application.
- Take precautions to limit the effects that the electromagnetic interference (EMI) produced by the Generator may have on the performance of other equipment. Check the compatibility and safety of combinations of other physiological monitoring and electrical apparatus to be used on the patient in addition to the Generator.
- Adequate filtering must be used to allow continuous monitoring of the surface electrocardiogram (ECG) during radiofrequency power applications.
- Careful catheter manipulation must be performed to avoid cardiac damage, or tamponade. Catheter advancement should be done under fluoroscopic guidance. If resistance is encountered, DO NOT use excessive force to advance or withdraw the catheter.
- Do not attempt to puncture until firm position of the active tip has been achieved against the atrial septum.
- It is not recommended to exceed five (5) radiofrequency power applications per Toronto Transseptal Catheter.
- Do not bend the Toronto Transseptal Catheter. Excessive bending or kinking of the catheter shaft may damage the integrity of the catheter and may cause patient injury. Care must be taken when handling the catheter.
- Do not exceed the maximum recommended pressure of 100 psi (0.69 MPa)
- The Generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the catheter and DIP electrode, particularly when operating the device.
- During power delivery, the patient should not be allowed to come in contact with ground metal surfaces.
- Apparent low power output or failure of the equipment to function properly at normal settings may indicate faulty application of the DIP electrode, failure to an electrical lead, or poor tissue contact at the active tip. Check for obvious equipment defects or misapplication. Attempt to better position the tip of the catheter against the atrial septum. Only increase the power if low power output persists.
- Baylis Medical Company relies on the physician to determine, assess and communicate to each individual patient all foreseeable risks of the Baylis Medical Radiofrequency Puncture System.
- Do not use the device for the treatment of children or treatment of pregnant or nursing women.

VI. ADVERSE EVENTS

Adverse events that may occur while using the Baylis Medical Radiofrequency Puncture System include:

Tamponade	Vascular thrombosis
Vessel perforation	Allergic reaction to contrast medium
Vessel spasm	Arteriovenous fistula
Hemorrhage	Thromboembolic episodes
Hematoma	Myocardial Infarction
Pain and Tenderness	Atrial Flutter
Sepsis/Infection	Perforation of the myocardium
Atrial Fibrillation	Ventricular Tachycardia
Sustained arrhythmias	

VII. EQUIPMENT REQUIRED

Intracardiac puncture procedures should be performed in a specialized clinical setting equipped with a fluoroscopy unit, radiographic table, physiologic recorder, emergency equipment and instrumentation for gaining vascular access. Ancillary materials required to perform cardiac Puncture include:

- BMC Radiofrequency Puncture Generator
- Baylis Connector Cable (RFP-102 for use with RFP-100 Generator, or RFX-BAY-TS-12-SU for use with RFP-100A Generator).
- Transseptal Sheath/Dilator set, such as the Baylis Medical Company TorFlex Transseptal Guiding Sheath.
- Disposable Indifferent (dispersive) Patch (DIP) electrode must meet or exceed IEC 60601-2-2 requirements for electrosurgical electrodes.

VIII. INSPECTION PRIOR TO USE

Prior to use of the Baylis Medical Radiofrequency Puncture System, the individual components including the BMC Radiofrequency Puncture Generator, Toronto Transseptal Catheter, and the BMC Connector Cable should be carefully examined for damage or defects, as should all equipment used in the procedure. Do not use defective equipment.

IX. DIRECTIONS FOR USE

- All instructions for equipment required should be carefully read, understood, and followed. Failure to do so may result in complications.
- The Toronto Transseptal Catheter is supplied sterile. Use aseptic technique when opening the package and handling the product in the sterile field.
- A Transseptal Sheath and Dilator is usually inserted through the right femoral vein and is then advanced over a guidewire to be positioned into the superior vena cava (SVC) under fluoroscopic guidance. The Baylis Medical TorFlex Transseptal Guiding Sheath is recommended for this purpose.
- Insert the Toronto Transseptal Catheter through the sheath/dilator set.
- With the Toronto Transseptal Catheter inside the dilator, position the catheter/sheath/dilator assembly in the right atrium against the fossa ovalis under fluoroscopic guidance using standard technique.
- Connect the Toronto Transseptal Catheter to the pressure monitoring system, if using, by joining its clear luer connector to a luer lock and turning the connector to ensure a secure connection.
- Connect the Toronto Transseptal Catheter to the BMC Connector Cable ensuring that the catheter is inserted properly into the connector. Make sure that the Connector Cable is plugged into the appropriate port on the BMC Radiofrequency Puncture Generator. Be sure to carefully follow the Instructions for Use provided with the Generator and Cable.
- Advance the Toronto Transseptal Catheter so that the active tip is engaging the septum at the fossa ovalis, but still within the dilator. Once appropriate positioning has been achieved, radiofrequency power can be delivered via the BMC Radiofrequency Puncture Generator to the distal tip. This results in the puncture of the targeted cardiac tissue. Please refer to the Generator Instructions for Use before using the Generator.
- **NOTE: If the active tip of the catheter becomes bent at any time during its use, dispose of the catheter immediately. Do not attempt to straighten the active tip.**
- For RFP-100 Generator: A power setting of 10 Watts has been experimentally determined to be sufficient for successful puncture.
- **NOTE: Use the lowest appropriate power setting to achieve the desired puncture. Do not use an initial power setting greater than 10 Watts.**
- For RFP-100A Generator: An initial RF setting between one (1) second on "PULSE" mode to two (2) seconds on "CONSTANT" mode has been shown to be sufficient for successful puncture.
- Apply firm pressure to the Toronto Transseptal Catheter during the application of radiofrequency energy to successfully advance it through the tissue.
- Radiofrequency power delivery can be terminated by pressing the RF ON/OFF button on the Generator if the timer has not expired.
- Entry into the left atrium can be confirmed by monitoring the Toronto Transseptal Catheter under fluoroscopy, by injecting contrast through the catheter, and/or by obtaining hemodynamic pressure waveforms.
- If septal puncture is not successful after five (5) radiofrequency power applications, it is advised that the user proceed with an alternate method for the procedure.
- Once the initial puncture is successfully completed, the Toronto Transseptal Catheter should be mechanically advanced without any radiofrequency power.

- The transeptal dilator can then be used to enlarge the puncture. Connections (pg. Error! Bookmark not defined.)

X. CLEANING AND STERILIZATION INSTRUCTIONS

The Toronto Transeptal Catheter is intended for single use only. Do not clean or re-sterilize the Toronto Transeptal Catheter.

XI. TROUBLESHOOTING

The following table is provided to assist the user in diagnosing potential problems.

PROBLEM	COMMENTS	TROUBLESHOOTING
Generator Error Messages	In order to successfully puncture tissue using RF energy, the entire system must be connected and all devices must be in good working order.	Ensure that all connections are made: - catheter to connector cable - connector cable to generator - generator to power outlet - generator to grounding pad Visually inspect the catheter or cable for damage. Immediately discard any damaged equipment. If the problem persists discontinue use. For error messages encountered while attempting radiofrequency puncture, refer to the operator's manual that accompanies the Generator.
Inaccurate Pressure Readings	In order to accurately monitor pressure, the entire system must be properly connected and all devices must be in good working order.	Ensure that the following connections are made: -catheter to pressure transducer -pressure transducer to monitoring system Ensure that the transducer is zeroed. Ensure that the transducer is leveled with the phlebostatic axis Perform a 'fast-flush test' to determine the dynamic response. Visually inspect the catheter for any damage. If there are any breaks or kinks, discard immediately
Catheter breaks or kinks.	Breaks and kinks in the catheter are a potential cause of patient injury.	Discard immediately

XII. CUSTOMER SERVICE AND PRODUCT RETURN INFORMATION


If you have any problems with or questions about Baylis Medical Equipment contact our technical support personnel.

NOTES:

- In order to return products you must have a return authorization number before shipping the products back to Baylis Medical Company.
- Baylis Medical will not accept any piece of used equipment without a sterilization certificate. Ensure that any product being returned to Baylis Medical has been cleaned, decontaminated and sterilized as indicated in the user instructions before returning it for warranted service.

XIII. LABELING AND SYMBOLS

	Manufacturer
---	--------------

	EU Authorized Representative
	Sterile using ethylene oxide
	Use By
	Caution
	Consult Instructions for Use
	Model number
Rx ONLY	Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
	Single Use – Do not reuse
	Lot Number
	Only for EU member states: Use of this symbol indicates that the product must be disposed of in a way that complies with local and national regulations. For questions regarding recycling of this device please contact your distributor
	Do Not Use if Packaging is Damaged
	Keep Away from Sunlight
	Contains Phthalates: Bis(2-ethylhexyl)phthalate (DEHP)
	Non-Pyrogenic

XIV. LIMITED WARRANTY – DISPOSABLES AND ACCESSORIES

Baylis Medical Company Inc. (BMC) warrants its Disposable and Accessory products against defects in materials and workmanship. BMC warrants that sterile products will remain sterile for a period of time as shown on the label as long as the original package remains intact. Under this Limited Warranty, if any covered product is proved to be defective in materials or workmanship, BMC will replace or repair, in its absolute and sole discretion, any such product, less any charges to BMC for transportation and labor costs incidental to inspection, removal or restocking of product. The length of the warranty is: (i) for the Disposable products, the shelf life of the product, and (ii) for the Accessory products, 90 days from shipment date.

This limited warranty applies only to new original factory delivered products that have been used for their normal and intended uses. BMC's Limited Warranty shall not apply to BMC products which have been resterilized, repaired, altered, or modified in any way and shall not apply to BMC products which have been improperly stored or improperly cleaned, installed, operated or maintained contrary to BMC's instructions.

DISCLAIMER AND LIMITATION OF LIABILITY

THE LIMITED WARRANTY ABOVE IS THE SOLE WARRANTY PROVIDED BY SELLER. SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE OR PURPOSE.

THE REMEDY SET FORTH HEREIN SHALL BE THE EXCLUSIVE REMEDY FOR ANY WARRANTY CLAIM, AND ADDITIONAL DAMAGES, INCLUDING CONSEQUENTIAL DAMAGES OR DAMAGES FOR BUSINESS INTERRUPTION OR LOSS OF PROFIT, REVENUE, MATERIALS,

ANTICIPATED SAVINGS, DATA, CONTRACT, GOODWILL OR THE LIKE (WHETHER DIRECT OR INDIRECT IN NATURE) OR FOR ANY OTHER FORM OF INCIDENTAL, OR INDIRECT DAMAGES OF ANY KIND, SHALL NOT BE AVAILABLE. SELLER'S MAXIMUM CUMULATIVE LIABILITY RELATIVE TO ALL OTHER CLAIMS AND LIABILITIES, INCLUDING OBLIGATIONS UNDER ANY INDEMNITY, WHETHER OR NOT INSURED, WILL NOT EXCEED THE COST OF THE PRODUCT(S) GIVING RISE TO THE CLAIM OR LIABILITY. SELLER DISCLAIMS ALL LIABILITY RELATIVE TO GRATUITOUS INFORMATION OR ASSISTANCE PROVIDED BY, BUT NOT REQUIRED OF SELLER HEREUNDER. ANY ACTION AGAINST SELLER MUST BE BROUGHT WITHIN EIGHTEEN (18) MONTHS AFTER THE CAUSE OF ACTION ACCRUES. THESE DISCLAIMERS AND LIMITATIONS OF LIABILITY WILL APPLY REGARDLESS OF ANY OTHER CONTRARY PROVISION HEREOF AND REGARDLESS OF THE FORM OF ACTION, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND FURTHER WILL EXTEND TO THE BENEFIT OF SELLER'S VENDORS, APPOINTED DISTRIBUTORS AND OTHER AUTHORIZED RESELLERS AS THIRD-PARTY BENEFICIARIES. EACH PROVISION HEREOF WHICH PROVIDES FOR A LIMITATION OF LIABILITY, DISCLAIMER OF WARRANTY OR CONDITION OR EXCLUSION OF DAMAGES IS SEVERABLE AND INDEPENDENT OF ANY OTHER PROVISION AND IS TO BE ENFORCED AS SUCH. IN ANY CLAIM OR LAWSUIT FOR DAMAGES ARISING FROM ALLEGED BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, PRODUCT LIABILITY OR ANY OTHER LEGAL OR EQUITABLE THEORY, THE BUYER SPECIFICALLY AGREES THAT BMC SHALL NOT BE LIABLE FOR DAMAGES OR FOR LOSS OF PROFITS, WHETHER FROM BUYER OR BUYER'S CUSTOMERS. BMC'S LIABILITY SHALL BE LIMITED TO THE PURCHASE COST TO BUYER OF THE SPECIFIED GOODS SOLD BY BMC TO BUYER WHICH GIVE RISE TO THE CLAIM FOR LIABILITY.

No agent, employee or representative of Baylis Medical has the authority to bind the Company to any other warranty, affirmation or representation concerning the product.

This warranty is valid only to the original purchaser of Baylis Medical products directly from a Baylis Medical authorized agent. The original purchaser cannot transfer the warranty.

Use of any BMC product shall be deemed acceptance of the terms and conditions herein.

The warranty periods for Baylis Medical products are as follows:

Disposable Products	The shelf life of the product
Accessory Products	90 days from the shipment date

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