LANGE 

DOC22.98.90251

sensION™ + EC 7

Manual
Bedienungsanleitung
Manuel d'utilisateur

07/2010, Edition 1

07/2010, Heft 1

07/2010, Edition 1

ENGLISH

DEUTSCH

FRANÇAIS

CONTACT

Index

1. Technical Specifications	2
2. General information	3
2.1. Safety information	3
2.1.1. Use of hazard information	3
2.1.2. Precautionary labels	3
2.2. General product information	3
2.3. Instrument composition	3
3. Installation	4
3.1. Assemblage	4
3.2. Connections	5
4. Operation	5
4.1. Description	5
4.2. Start-up	6
4.3. General diagram	6
4.4. Calibration	7
4.5. EC and NaCl measurements	8
4.6. Other options	8
4.7. Recognized standards	9
5. Maintenance	10
6.1. Clean the instrument	10
6.2. Sensor cleaning	10
6. Warning messages	10
7. Accessories and spares	11
8. Warranty, liability and complaints	11

Contact Information

1. Technical Specifications

Measuring ranges (resolution depending on scale)

Conductivity	0.01 μ S to 500 mS/cm
Resistivity	0 to 100 M Ω
Salinity	5.85 mg/l to 311.1 g/l NaCl
Temp	-20.0 to 150.0 °C (-4 to 302 °F)

Measuring error (± 1 digit)

Conductivity	$\leq 0.5\%$
Resistivity	$\leq 0.5\%$
Salinity	$\leq 0.5\%$
Temp	$\leq 0.2\text{ °C}$ (0.4 °F)

Reproducibility (± 1 digit)

Conductivity	$\pm 0.1\%$
Resistivity	$\pm 0.1\%$
Salinity	$\pm 0.1\%$
Temp.	$\pm 0.1\text{ °C}$ (0.1 °F)

Automatic temperature compensation

Manual entrance.

With Pt 1000 temperature probe (A.T.C.)

With NTC 10 K Ω probe.

Temperature coefficient, TC

Lineal, TC = 0.00 to 9.99 % / °C

Temperature reference, TRef

20 °C or 25 °C.

Conductivity calibration (EC)

With 1, 2 or 3 standards selectable inside the range.

Manual introduction of cell constant.

Programmable calibration validity between 0 and 99 days.

Automatic recalibration warning

Accepted cell constants

With calibration: values between 0.05 to 50 cm⁻¹.

Manual introduction: values between 0.05 to 15 cm⁻¹.

Measuring mode

By stability or in Continuous.

Languages

English, German, Spanish, French, Italian, and Portuguese.

Display

Graphic, backlit liquid crystal, 128x64 dots.

Inputs and outputs

Cell with built-in A.T.C. type Pt 1000, telephonic connector.

Stirrer, RCA connector.

Ambient conditions

Working temperature 5 to 40 °C (41 to 104 °F).

Storage temperature -15 to 65 °C (5 to 149 °F).

Relative humidity < 80% (not condensed).

Power supply

External plug-in power supply 100-240 VAC 0.4 A 47-63 Hz.

Materials

Enclosure, ABS and PC.

Keypad, PET with protective treatment.

Physical parameters

Weight: 1100 g.

Size: 350 x 200 x 110 mm.

2. General information

As a result of constant improvements to our products sometimes differences may exist between this manual and the instructions supplied with the instrument.

2.1. Safety information

Please, read carefully this information before installing and using the instrument !

Pay attention to all danger and caution statements.

2.1.1. Use of hazard information

DANGER

Indicates a potentially or imminently hazardous situation that, if not avoided, will result in death or serious injury.

WARNING

Indicates a potentially or imminently hazardous situation that, if not avoided, may result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation that, if not avoided, may result in minor or moderate injury.

Important note: Indicates a situation that, if not avoided, could lead to damage to the instrument. Important information that requires special emphasis.

Note: Information that supplements points in the main text.

2.1.2. Precautionary labels

Read carefully all labels and tags attached to the instrument.

Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of life equipment to the Producer for disposal at no charge to the user.

Note: For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.

2.2. Description

The **sensION™+ EC7** is an instrument to measure conductivity, salinity and temperature.

2.3. Instrument composition

Version	Cell	Accessories	Manual
LPV3010.98.0002	—	✓	✓
LPV3070.98.0002	LZW5070.97.0002	✓	✓

Accessories:

- Standard solutions 147 µS/cm , 1413 µS/cm and 12.88 mS/cm, 135 ml flasks.
- Magnetic stirrer.
- Calibration flasks, with magnetic bar inside.
- Plug-in power supply.
- Electrode's stick and support.

3. Installation

DANGER

Only qualified personnel, following the local security regulations, should conduct the tasks described in this section of the manual.

3.1. Assemblage

1. Unpack the instrument.

2. Verify that the shift is in the same position than in the picture.

6. Place the electrode on the electrode's support. Insert the electrode's cable through the groove.

3. Drive the stick into the hole.

4. Turn the shift in a counterclockwise until the stick is well fixed.

Place the o-ring at 135-140 mm from the base in order to avoid the shock of the sensor against the bottom of the glass.

Note: Stirrer replacement

Before replacing the stirrer, disconnect it from the rear panel and release the cable from the bottom of the instrument.


1. Turn the stirrer to release or fix from base.

5. Place the electrode's support on the stick by softly pressing the button.

Holder for 3 sensors.
PN LZW93 21.

2. Pull up the stirrer and replace it. Place the cable of the new stirrer through the bottom groove of the instrument.

3. Installation

3.2. Connections

DANGER

Make sure the power adapter matches your local power supply.

1. Magnetic stirrer, RCA connector.
2. Conductivity cell, with built-in ATC, telephonic connector.
3. Plug-in power supply (12 V).

4. Operation

4.1. Description

Keys

Switch on / Switch off.

Electric conductivity measurement.

Salinity measurement.

Calibration and calibration data.

Select between temperature (°C, °F) and stirrer (St) to modify the temperature or the stirring speed.
Access to TC and TRef.

Selection of numeric values.

Access to:
- Language.
- Date and hour.
- Display contrast.
- Calibration frequency.

4. Operation


4.2. Start-up

Connect to the rear panel of the instrument:


- the conductivity cell
- the magnetic stirrer
- the plug-in power supply

Follow the next steps:

For the first time...


Successive times...

4.3. General diagram

This is a general view over the way to access to all possibilities of the instrument:

- DISPLAY CONTRAST
- LANGUAGE SELECTION
- DATE AND TIME
- CHANGE OF TEMPERATURE UNITS
- EC and NaCl CALIBRATION
- EC, RESISTIVITY and NaCl MEASUREMENT

4. Operation

4.4. Calibration

This process involves adjusting the values read by a conductivity measuring system (instrument and cell) according to the values of standard solutions.

Calibration is very important for highly accurate readings.

The equipment is capable of one-, two- and three-point calibration using 147 $\mu\text{S}/\text{cm}$, 1413 $\mu\text{S}/\text{cm}$, 12.88 mS/cm and 111.8 mS/cm standards.

The calibration parameters are stored in the memory until a new calibration is performed.

It is recommended to use the standards with values closest to the expected one in the samples.

One-point calibration

This form of calibration is used when measuring conductivity values around that of the standard used.

This is the most common type of calibration. In this type, 1413 $\mu\text{S}/\text{cm}$ is the most commonly-used standard.

Two-point calibration

This type of calibration is recommended for work in high or low conductivity regions where accuracy is required.

If you are calibrating with more than one standard, it is necessary to start with the one with the lowest conductivity to avoid contamination problems.

Three-point calibration

This type of calibration is recommended when the conductivity of the samples to be measured ranges from low to high.

If you are calibrating with more than one standard, it is necessary to start with the one with the lowest conductivity to avoid contamination problems.

Calibration data

They can be viewed at any moment. See diagram.

Calibration frequency

Select between 0 and 99 days. The instrument by default is programmed for calibration every 15 days, and indicates the remaining time to the new calibration.

Selecting 0 h the automatic re-calibration warning is deactivated

Theoretical calibration

Selecting "theoretical calibration" the calibration data of the cell in use are replaced by $C=1,000\text{ cm}^{-1}$

This step is prior to the obtention of a calibration certificate.

C introduction

Manual introduction of the cell constant. Useful when the cell was previously calibrated, for example in calibration laboratory.

Press **EC** at any point to exit.


Preparation

Fill the flasks with the corresponding standard solution (50 ml.)

The times of use of these standards portions depends on the cell cleaning. Crison recommends replacing the solution in the flasks after 2-3 times of usage.

Every calibration flask is supplied with a magnetic bar.

Diagram for calibration

4. Operation

4.5. EC and NaCl measurements


Measurement by stability

The value is fixed on the display when during 6 seconds the reading variation is less than ± 1 digit.

During measurement, the measured value appears on the display from time to time to orient the user.

Continuous measurement

The instrument displays the value obtained by the electrode all the time.

4.6. Other options

TRef and TC selection, change of temp. units

By default the instrument is programmed with TRef 25°C and TC 2.00 %/°C. To modify these values, follow the scheme.

5. Maintenance

6.1. Clean the instrument

Important note: Never use cleaning such as turpentine, acetone or similar products to clean the instrument including the display and accessories.

Only clean the housing and accessories using a soft, damp cloth. Mild soap solution may also be used. Dry the cleaned parts carefully with a soft cotton cloth.

CAUTION

When using chemicals or solvents, comply with the instructions of the producer and all local safety regulations.

6.2. Sensor cleaning

See the manual of the sensor in use.

6. Warning messages

During calibration

EC CALIBRATION
UNSTABLE READING
Time t> 100 s
Repeat with

EC CALIBRATION
ATTENTION !
Check temperature.
Check the cell.

EC CALIBRATION
ATTENTION !
Difference C > 30 %
Repeat with

EC CALIBRATION
ATTENTION !
Check the cell.
Check the standards.

EC CALIBRATION
SAME BUFFERS
Repeat with

EC CALIBRATION
ATTENTION !
Check the cell.

EC CALIBRATION
STANDARD TEMPERATURE
< 15 °C or > 30 °C
Repeat with

EC CALIBRATION
ATTENTION !
C<0.05 or C>50.00
Repeat with

EC CALIBRATION
ATTENTION !
Buffer or cell in
poor conditions.

During measurement

MEASURE STABLE
728 23.2 °C
uS/cm @ 25 °C
Calibrate again

MEASURING by stability
EC out of range
Check the cell.

MEASURING by stability
out of range °C
Check temp. sensor

MEASURING by stability
ATTENTION ! TC = 0
Impossible to measure
salinity
Exit with

MEASURING by stability
UNSTABLE READING
Time > 150 s
Check temperature.
Check the cell.

7. Accessories and spares

Code	Description
LZW9111.99	Set of 3 x 50 ml printed flasks for benchtop EC-meters calibration
LZW9322.99	Magnetic stirrer
LZW9700.99	Standard 147 $\mu\text{S}/\text{cm}$ (at 25 °C, 77 °C), 250 ml flask.
LZW9710.99	Standard 1413 $\mu\text{S}/\text{cm}$ (a 25 °C, 77 °C), 250 ml flask.
LZW9720.99	Standard 12.88 mS/cm (a 25 °C, 77 °C), 250 ml flask.
LZW9740.99	Standard 111.8 mS/cm a (25 °C, 77 °C), 250 ml flask.

For cells please ask for specific brochure.

8. Warranty, liability and complaints

The manufacturer warrants that the product supplied is free of material and manufacturing defects and undertakes the obligation to repair or replace any defective parts at zero cost.

The warranty period for instruments is 24 months. If a service contract is taken out within 6 months of purchase, the warranty period is extended to 60 months.

With the exclusion of the further claims, the supplier is liable for defects including the lack of assured properties as follows: all those parts that can be demonstrated to have become unusable or that can only be used with significant limitations due to a situation present prior to the transfer of risk, in particular due to incorrect design, poor materials or inadequate finish will be improved or replaced, at the supplier's discretion. The identification of such defects must be notified to the supplier in writing without delay, however at the latest 7 days after the identification of the fault. If the customer fails to notify the supplier, the product is considered approved despite the defect. Further liability for any direct or indirect damages is not accepted.

If instrument-specific maintenance and servicing work defined by the supplier is to be performed within the warranty period by the customer (maintenance) or by the supplier (servicing) and these requirements are not met, claims for damages due to the failure to comply with the requirements are rendered void.

Any further claims, in particular claims for consequential damages cannot be made.

Consumables and damage caused by improper handling, poor installation or incorrect use are excluded from this clause.

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DR. LANGE NEDERLAND B.V.