Bio-Rad
Medical Diagnostics GmbH

Industriestr. 1, D-63303 Dreieich,
Germany
Telefon : +49 (0) 6103 3130-0
Telefax :+49 (0) 6103 3130-724

Bio-Rad Medical Diagnostics GmbH · P.O.Box 10 21 50 · 63267 Dreieich

To whom it may concern

Tel.: +49 (0) 6103 3130-807
Fax: +49 (0) 6103 3130-724
daniel_czachurski@bio-rad.com

11. Juni 2015

Urgent Field Safety Notice
Solidscreen II Anti-D Blend 806530, lot 2429090-00 and 2416111-00

Addressor:

Bio-Rad Medical Diagnostics GmbH, Dr. Daniel Czachurski, Industriestraße 1, 63303 Dreieich, Germany, Tel.: +49 (0) 6103 3130-807, Fax: +49 (0) 6103 3130-724, E-Mail: daniel_czachurski@bio-rad.com

Details of affected devices:

Product	Ref.-N°	Lot N°	Date of Production	Expiry
Solidscreen II Anti-D Blend	806530	2429090-00	2014-07-16	2016-07-16
Solidscreen II Anti-D Blend	806530	2416111-00	2014-04-16	2016-04-16

Description of the problem:

The above mentioned product lots contain a weak anti-Wr^a. In cells carrying a Wr^a antigen, the product may give a weak false-positive reaction. Usually, due to the weak reactivity, this reaction will not be interpreted as a weak D antigen. However, in case of stronger reactions, Wr(a+) samples might be misinterpreted as weak D. For patients there is no risk as these will be transfused with D-negative blood. However, there might be a risk for pregnant women for not receiving a Rhesus prophylaxis in case of misinterpretation of weak D.

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Advise on action to be taken by the user:

Please do not use mentioned product lots any longer. Affected products will be replaced free of charge. Please contact your local Bio-Rad representative.

Transmission of this Field Safety Notice:

This notice needs to be passed on to all those who need to be aware within your organization or to any organization where the potentially affected devices have been transferred.

Please maintain awareness of this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.

We apologize for any inconvenience and would like to thank you in advance for your cooperation. If you have any additional questions please contact us.

Contact reference person:

Manuela Hornstein
Bio-Rad Medical Diagnostics GmbH
Industriestraße 1
63303 Dreieich
Germany
E-Mail: contact_bmd@bio-rad.com

Best regards,

Bio-Rad Medical Diagnostics GmbH

Dr. Rolf Vornhagen
Safety Officer

Dr. Daniel Czachurski
International Product Manager

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International Product Manager

Solidscreen™ II Anti-D (RH1) Blend

For the determination of the red cell antigen weak D and partial D of RhD negative blood samples

Solid phase antiglobulin test Solidscreen™ II Strip

(REF) 806 530 5 ml

(PRES) Preservative: < 0.1 % NaN₃

(IVD) Diagnostic reagent for in vitro use only

To be used by trained laboratory personnel only. Symbols are used according to regulation EN 980. Symbols not contained in regulation EN 980 are described.

Test purpose

Solidscreen™ II Anti-D (RH1) Blend is used with Solidscreen™ II for the detection of weak D and partial D antigens (DVI and DVII). It is used to test samples which have been tested RhD negative with a direct agglutination method.

Test principle

Solidscreen™ II Anti-D (RH1) Blend reacts with weak D and partial D positive red cells in solidphase method with Solidscreen™ II Strip. The active ingredients Anti-Human Globulin Solidscreen™ II-Donor bind to the cavity wall activated with protein A. The bound antiglobulin molecules bind to sensitized test red cells thus forming a link between the coated test red cells and the cavity wall. A positive reaction shows a layer of cells on the plates' cavities. Uncoated red cells react with sedimentation after the last centrifugation on the cavity's bottom¹.

Reagent

As reactive component Solidscreen™ II Anti-D (RH1) Blend contains human monoclonal antibodies of the immunoglobulin class IgG (clone BS221 and H41 11B7) and therefore is suited for testing with Solidscreen™ II. The antibodies are derived from cell culture supernatant and demonstrate the consistent specificity and reproducibility characteristic for monoclonal antibodies. Due to the selection of appropriate clones, weak D and partial D antigens (DVI and DVII) can be detected.

Solidscreen™ II Anti-D (RH1) Blend clones BS221 and H41 11B7 (IgG)

Material needed but not supplied

- Anti-Human Globulin Solidscreen™ II Donor (REF) 806 516)

For further materials please refer to the instructions for use Solidscreen™ II Strip (REF) 806 521).

Sample material

Blood samples should be prepared following general blood sampling guidelines. Fresh, non-hemolytic samples (EDTA blood) should be used. If testing is delayed the EDTA anticoagulated samples should be stored at 2...8°C. EDTA anticoagulated whole blood samples may be tested for up to 7 days following collection. Donor blood stored in citrate anticoagulant may be stored at 2...8°C until the expiration date of the donor unit. Samples to be tested must be centrifuged prior to testing (analogous to the TANGO® optimo and TANGO infinity® user manual). TANGO® optimo and TANGO infinity® systems prepare the red blood cells suspension automatically. Blood sample segments of the blood unit can be used, but have to be prepared according to the TANGO® optimo and TANGO infinity® user manual. ◀

Test procedure

For each test 50 µl Solidscreen™ II Anti-D (RH1) Blend are pipette to the respective cavity of Solidscreen™ II Strip.

The exact test procedure is described in the instructions for use of Solidscreen™ II Strip (REF) 806 521).

Performance characteristics and limitations of the method

- Solidscreen™ II Anti-D (RH1) Blend is only for testing blood samples which have been previously typed Rh D negative with a direct agglutination technique. RhD positive blood samples with normal D antigen expression can be interpreted as false negative with Solidscreen™ II Anti-D (RH1) Blend due to agglutination of the red cells instead of solidphase binding.
- Solidscreen™ II Anti-D (RH1) Blend might not react with cells in the combination C in trans to D cat VI Type 1 (Cde/cD^{cat VI} e, extremely rare), due to the reduced number of antigens.
- Rh33 cannot be detected with Solidscreen™ II Anti-D (RH1) Blend.
- Insufficient or inappropriate washing can lead to false negative reactions. Small amounts of residual sera/plasma can neutralize the Anti-Human Globulin.
- Strongly lipemic, icteric or microbiologically contaminated samples may lead to false results.
- In case of ambiguous results it is recommended to wash red blood cells at least 2 times.
- Conditions that may cause false positive result are:
 - Contamination of sample or reagents
 - Autoantibodies
 - Improper storage or preparation of red blood cells

- Antibodies to antibiotics or other reagents
- Cold antibodies
- Positive Direct Anti-Globulin Test

Please refer to the instructions for use of Solidscreen™ II Strip (REF) 806 521). In case of questionable results of unknown origin our Bio-Rad Service will assist you (phone: +49-6103-3130-611).

Shelf life

After opening the vial the product can be stored under proper storage conditions (2...8°C) until the expiry date. The expiry date is printed on the label. Do not use damaged vials.

Warning and precautions

- Turbidity or other visible changes may indicate a bacterial contamination. In this case the reagent must be discarded. The cause for turbidity must be examined by the manufacturer.
- All techniques are to be performed according to the manufacturer's instructions. Each deviation from these instructions is the sole responsibility of the user and needs to be validated by the user.
- Solidscreen™ II Anti-D (RH1) Blend contains human monoclonal antibodies. Because of biotechnological manufacturing of human monoclonal antibodies, a contamination with infectious diseases can be considered impossible.
- All materials of human origin used for this product were tested for HbsAg, Anti-HCV, and Anti-HIV-1/2 and proved to be non-reactive. Nevertheless, all products of human origin must be regarded as potential transmitters of hepatitis, HIV or other infectious diseases.
- All test reagents of biological origin must be regarded as potential transmitters of infectious agents. Appropriate safety precautions are recommended.
- Used test material must be discarded as hazardous material. Waste management information are in the safety data sheet.
- The bovine albumin used for the production of this reagent is purchased from BSE-free sources.

Internal quality controls according to national guidelines are recommended at regular intervals.

Please also refer to the instructions for use of Solidscreen™ II Strip (REF) 806521).

References

- (1) AABB Technical Manual, 17th edition.

Key: Underline = Addition or significant change ◀ = Deletion of text

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