Reflotron® Amylase

REF 11200658

15 tests

English

Intended use

Test for the quantitative determination of the total activity of α -amylase (EC 3.2.1.1) in blood, serum or plasma with Reflotron systems.


Summary

The α -amylase (1,4- α -D-glucanohydrolase, EC 3.2.1.1) catalyzes the hydrolytic breakdown of polymeric carbohydrates like amylose, amylopectin and glycogen by cleavage of 1,4- α -glucosidic bonds. In poly- and oligosaccharides, several glucosidic bonds are simultaneously hydrolyzed. Being the smallest unit, maltotriose is cleaved into maltose and glucose, at a much slower rate.

Two types of α -amylase are distinguished: the pancreatic type (P-type) and the salivary gland type (S-type). The P-type comes almost exclusively from the pancreas, and is thus organ-specific; the S-type can have various sources. Apart from in the salivary glands, it can also occur in tears, sweat, breast milk, amniotic fluid, lungs, testes and the epithelium of the oviducts. On account of the rather non-specific clinical symptoms of pancreatic diseases, enzyme determinations are of particular value in the diagnosis of these disorders. They are performed particularly for the diagnosis and monitoring of acute pancreatitis. However, hyperamylasemia does not only occur in acute pancreatitis, or in the inflammatory phase of chronic pancreatitis but also in renal failure as a result of reduced glomerular filtration, in tumours of the lung, or the ovaries, in disorders of the salivary glands, diabetic ketoacidosis, cerebral trauma, after surgery or in macroamylasemia or in pneumonia. To confirm the pancreatic origin, additional determination of a further pancreas-specific enzyme like pancreatic α -amylase is recommended.

Test principle
After application to the test strip, the sample flows into the reaction zone, where, in the case of blood samples, the separation of the erythrocytes from the plasma occurs. The α -amylase contained in the sample and the enzyme α -glucosidase contained in the test cleave the substrate indolyl- α -D-maltoheptaose yielding indoxyl and glucose.

The indoxyl is coupled with a diazonium derivative to form a purple dye. The amount of dye formed per unit of time is directly proportional to the α -amylase activity:

The enzyme activity is measured kinetically at a wavelength of 567 nm and 37 °C, and is displayed after about 170 seconds in U/L or $\mu\text{katal/L}$.

Reagents
Components per test: indolyl- α -D-maltoheptaose 81 μg ; α -Glucosidase (yeast rec.) ≥ 3.1 U; 2-methoxy-4-morpholinophenyl diazoniumtetrachlorozincate 6.8 μg ; buffer.

Precautions and warnings
For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request.

Disposal of all waste material should be in accordance with local guidelines.

Avoid any contact to the application zone of a test strip (e.g., during pipetting of sample).

Reagent handling

Test strips are ready for use.

Storage and stability

Store at 2–30 °C. Do not use the test strip after the specified expiry date.

Specimen collection and preparation

Capillary blood; whole blood collected in standard sample collection tubes; serum; heparinized blood, or heparinized plasma.

Use fresh capillary or venous blood immediately after collection.

Use only heparin (preferably Li-heparin) as anticoagulant. Do not use any other anticoagulants or additives.

Heparinized blood kept in closed containers must be used within 8 hours. After precipitation of the cellular components, the supernatant plasma can be used.

If heparinized single-use containers or capillary pipettes are used, please observe the stability date given by the manufacturer.

Serum and heparinized plasma collected in standard sample collection tubes is stable for 7 days at 20–25 °C, or for 1 month at 2–8 °C.

Sample volume: 30 μl .

Materials provided
1 container with 15 test strips

Materials required (but not provided)

- Reflotron instrument
- Reflotron pipette
- Reflotron pipette tips
- Reflotron capillary tubes
- Reflotron Precinorm U, Reflotron Precipath U, or Reflotron Clean + Check
- General laboratory equipment
- 0.9 % NaCl

Assay procedure

For optimum performance of the assay, follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

- Remove a test strip from the container. **Tightly recap the container immediately after removing a test strip.**
- Peel off the aluminum protective foil, taking care not to bend the test strip.
- All Reflotron tests require a sample volume of 30 μl .
- Apply the required volume of sample onto the **centre** of the reagent application zone using a pipette (e.g., Reflotron pipette) – being careful not to touch the application zone. Avoid air-bubbles.
- Open the flap or sliding cover. Within 15 seconds of applying the sample, place the test strip onto the guide, and slide it forward horizontally until it locks into place. Close the sliding cover or flap.
- The test parameter abbreviation is shown on the display, if the test strip has been correctly inserted and the magnetic code has been read. The result is displayed depending on the setting of the instrument.

Calibration

The function curve for the Reflotron Amylase assay for converting reflectance values into activities is defined for each lot using the α -amylase liquid method from Roche

Diagnostics. The parameters of the curve are automatically transferred to the instrument via the magnetic strip during testing.

Quality control

For quality control, use Reflotron Precinorm U, Reflotron Precipath U or Reflotron Clean + Check. The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Follow the applicable government regulations and local guidelines for quality control.

Calculation

The amylase activity is calculated automatically from the measurements taken, and function and conversion factors read from the magnetic strip on the lower face of each test strip. The enzyme activity is shown for 25 °C, 30 °C or 37 °C in U/L or $\mu\text{katal/L}$, depending on the way the instrument has been set.

Limitations - interferences^{3,4}

The following had no influence on the results in the concentration ranges tested (criterion: recovery ± 10 % of baseline): haemocrit values up to 55 %, hemoglobin up to 0.37 mmol/L (6 g/L), bilirubin up to 205 $\mu\text{mol/L}$ (12 mg/dL), triglycerides up to 22.8 mmol/L (1995 mg/dL).

Of 28 drugs tested, only paracetamol in high concentrations led to a reduction in the α -amylase activities measured.

Ascorbic acid concentrations above 10 mg/dL lead to increased recovery. Saliva and sweat contain α -amylase. Therefore, take care not to touch the application zone or test area.

High maltose concentrations of > 300 mg/dL lead to falsely low test results. As a result, patients undergoing peritoneal dialysis cannot be tested with the Reflotron Amylase or Reflotron Pancreatic Amylase assays because of their high maltose levels.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Measuring range

29–860 U/L or 0.48–14.33 $\mu\text{katal/L}$. If the α -amylase activity measured is above the measuring range for the Reflotron Amylase assay (indicated by an asterisk next to the displayed value or the message DILUTE AMYL), the sample may be diluted 1 + 2 with physiological saline solution. Multiply the result by a factor of 3.

Expected values⁵

Adults: Blood, serum, plasma: < 100 U/L or < 1.65 $\mu\text{katal/L}$ (37 °C, 30 °C, 25 °C). For technical reasons, the value is the same at 37 °C, 30 °C and 25 °C.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

The data for the Reflotron Amylase assay were determined in evaluation studies. The majority of the test results were within the given ranges.

Precision

Repeatability (within-run precision): CV (coefficient of variation) 2.9 % in the normal range, 3.2 % in the pathological range; sample material: serum.

Intermediate precision (between-day precision): CV 2.9 % in the normal range, 3.6 % in the pathological range; sample material: control sera.

Method comparison

A comparison of the Reflotron Amylase assay (y) with the α -amylase PNP/ α -amylase liquid method (x) using either serum, heparinized plasma or heparinized blood gave the following correlations: $y = 0.981x + 9.7$ and $y = 1.03x - 6.2$; $n = 50$; $r = 0.996$

References

1. Rothe A et al. Clin Chem 1987; 33:996.
2. Tietz NW. Clinical Guide to Laboratory Tests. 3rd edition. Philadelphia, PA: WB Saunders Company; 1995: 46.
3. Bais R et al. Clin Chem 1989; 35: 317-320.
4. Koller PU et al. Lab med 1989; 13: 399-402.
5. Heil W et al. Reference Ranges for Adults and Children, Roche Diagnostics GmbH, 1999.

Last update: 05/2010

Symbol table

REF	Catalogue number
LOT	Batch code/Lot number
IVD	In vitro diagnostic medical device
	Manufacturer
	Use by
	Temperature limitation (Store at)
	Consult instructions for use
	This product fulfills the requirements of the European Directive 98/79/EC for in vitro diagnostic medical devices

Significant additions or changes are indicated by a change bar in the margin.

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