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
ISO HYDROCYCLONE AND ITAS APPLICATIONS

INDEX

	<i>Page</i>
— REFERENCES	13
— ISO HYDROCYCLONE	
— TECHNICAL FEATURES	14
— INTRODUCTION	14
— SIGNALS	14
— RECOMMENDATION	15
— OPERATION	15
— WARNING AND AMALGAM LEVEL ELECTRONIC CONTROL SYSTEMS DESCRIPTION	15
— AMALGAM LEVEL VISUAL CONTROL	15
— AMALGAM CONTAINER REPLACEMENT	16
— INSTALLATION	16
— MAINTENANCE AND PERIODIC CHECKS RESERVED TO AUTHORIZED ENGINEERS	17
— NOTES	17
— GUARANTEE	17
— TRANSPORT OF SECONDHAND APPLIANCES	17
— ISO HYDROCYCLONE APPLICATIONS	
— ISO HYDROSEPARATOR	17
— SHORT INSTRUCTIONS	18
— TURBO-JET modular AND ISO HYDROCYCLONE	18
— PAL 22-50 VE PUMPS AND ISO HYDROCYCLONE	18
— LARGE SYSTEMS WITH SEPARATING TANKS AND ISO HYDROCYCLONE	18
— SINGLE-SURGERY ISO HYDROCYCLONE	19
— MAINTENANCE AND PERIODICAL CHECKS RESERVED TO AUTHORISED ENGINEERS	19
— GENERAL RECOMMENDATIONS AND BIOLOGIC HAZARD	19
— MAINTENANCE	19
— IMPORTANT RECOMMENDATION	20

REFERENCES (*Tables pag. 52-53-54*)

- 1 – Debris inlet
- 2 – Cyclone inlet holes
- 3 – Cone bottom
- 4 – Debris outlet
- 5 – Amalgam container
- 6 – Power-on signal
- 7 – Amalgam level at 95%
- 8 – Amalgam level at 100%
- 9 – Not in use
- 10 – Buzzer reset
- 11 – Amalgam max. level
- 12 – Liquid draining tap
- 13 – Amalgam container handle
- 14 – Amalgam container cover
- 15 – Amalgam container cover valve
- 16 – Screws
- 17 – Fixing holes
- 18 – AC 55 Power printed board
- 19 – AC 54 Logic printed board
- 20 – AC 37 Control panel
- 21 – AC 56 Optical reader
- 22 – Photodiodes

ISO HYDROCYCLONE

Technical features

Classification according to ISO 11143

Type 2

Rated tension

24 V ~ AC

Frequency

50/60 Hz

Max flow according to ISO 11143

5.5 l/min

Separation rate according to ISO 11143

95.6%

Power

8 W

Fuse

T 1.6 A

Protection class

I

	Alternating current	IEC 417-5032
	Protective earth (ground)	IEC 417-5019

Introduction

The following presentation is aimed at making operators and installers familiar with the equipment and facilities, to advise about maintenance and to inform engineers and users about hazards and necessary precautions for preventing accidents.

Signals

- **Electrical shock risk: even 230 V power can be lethal**
- **Biological danger, danger of infections from epidemic diseases.**
- **High temperature**
- **Generic danger sign**
- **Compulsory direction of flow and rotation sense**

Not all hazards and compulsory instructions can be indicated by signals, so the user must read all notices and consider them carefully. Failing to observe a signal or a notice can be harmful for the operator or damage the appliance.

Protection devices must not be removed. Structural or functioning features of the appliance or of the system must not be modified.

Recommendations

The retailer or the installation engineer will take care to train the surgery staff with trials on a brand new, non-contaminated appliance.

The installation of the amalgam separator is reserved to dental engineers authorized by the manufacturer. Aspirated debris is always contaminated and infected: for this reason, we stress that the greatest care must be used to prevent contamination of operators or the environment. Contamination may also result from an appliance in bad working order so we recommend to contact only dental engineers whose teaching and training is duly certified by the manufacturer. Any modification of the appliance must be agreed upon with the manufacturer.

Operation (Fig. 1, page 52)

Fluid passing through the inlet (1) **gets to the diffusion holes** (2); they distribute uniformly the fluid inside the cone and make it whirl downwards. When the fluid gets to **the bottom of the cone** (3), **it is subjected to a reversal of rotation and to a backflow**: at the point of the rotation reversal the speed is null and **the amalgam particles fall into the collecting container whilst the liquid rises upwards in the middle of the cone to be drained off** from the outlet (4).

Warning and amalgam level electronic control systems description (Fig. 2-3-4, page 53-54)

Next to the amalgam container (5) **two photodiodes (22) gauge the amalgam level in the container**: a yellow warning light (7) and a buzzer indicate that the container has been filled up to **95%** of its capacity. To silence the buzzer push the button (10). The Hydrocyclone will keep on running as usual whilst the warning light (7) will stay on. The warning light (8) and a second beep mean that **the amalgam container is filled at 100%**. The buzzer can be reset again by pushing button (10) and the Hydrocyclone will work until it is switched off. **This time it cannot be started again unless the amalgam container is replaced**. The warning system (see fig. 4) consists of a power printed board AC 55 and of a logic one AC 54 (superposed to form a unit), of a control panel AC 37, and of a max. of three level sensors AC 56 (connected in series) as the AC 54 Printed Board can control up to three separators; however it is enough that one of them gets to a critical level of 95% or of 100% to alarm the whole system. The AC 55 printed board is fitted with four relay connecting points; the first one for the 95% level, the second one for the 100% level, the third one for the buzzer, and the fourth for the contact (normally closed) stopping the aspiration when the max. amalgam level has been reached as described above. The same AC 55 P.B. is also fitted with a contact for the remote reset. To facilitate the checking, it is possible to install a remote control panel close to the operative area by using the relay connecting points also in case of a centralized separation system. **All the printed boards must be fed with low tension 24 V ~ AC**. They can be supplied without protecting cover or fitted with their IP20 protection -degree boxes, see fig. 4.

Amalgam level visual control (Fig. 2-3 , page 53)

The visual control of the amalgam level must be carried out very often and methodically.

As the demolition of a reconstructed molar tooth can raise the amalgam level by one millimetre, we deem necessary to check the actual status of the container before every operation requiring the use of the amalgam separator. **The red mark (11) shows the max. level compatible with the separator:** it must never be exceeded. You could wait and replace the container when the red mark has been reached, but we recommend to replace it when the level is 5 mm below the red mark (11) and to **check the amalgam level every evening** after the usual cleaning of the aspiration system.

Amalgam container replacement (Fig. 1-2, page 52-53)

Wear disposable gloves, an operating mask, and glasses as **protections against any accidental spurt before starting to replace the amalgam container. Disconnect the appliance** where the separator is installed on **from the power supply**. Open the tap (12) to drain the liquids left in the separator and in the container. Once the amalgam separator has been closed, suck up the drained fluid to separate amalgam particles that may be present in the fluid. Hold firmly the container with one hand to prevent it falling, lift the handle (13) upright, slide off the container and fill it with undiluted Puli-Jet + or with other suitable disinfectant. Open the valve (15) on the container lid (14), place the lid on the container and close the valve. The container is manufactured in shock-proof material and it is sealing closed. The container is full of toxic, harmful, and epidemic diseases contaminating material: it must be handled with care. To obtain the list of firms or institutes **authorized for disposal** or further information about the disposal, contact to the professional associations and unions or to local authorities. **In compliance with local regulations**, the container must be given to the care of firms or institutes specifically authorized for transport, disposal or recycling. A new container must be slid in as the previous one and the handle turned down up to touch the container. To reset the appliance, close the tap (12) and connect the appliance to the power supply again. Once the amalgam level has reached 95% of container capacity, in case of electronic control fitted appliances, the operator can wait for the max. level but the aspiration may be completely stopped, as the switching off of the appliance is enough to block it. For this reason we recommend to wash carefully the dental unit aspiration system according to the manufacturer directions once the amalgam has reached the 95% container capacity or is 5mm below the red mark and then replace the container without waiting for the maximum level. This procedure reduces noticeably the risks of contamination and of unwanted forced interruptions.

Installation (Fig. 1-2 , page 52-53)

The Hydrocyclone must be installed downstream an air-water separator and in an upright position. The appliance is fitted with four fixing holes (17); suitable brackets are available on request. **Debris inlet and outlet are marked by raised arrows: they must not be inverted.** The separator must be easy to reach from any side. An open space of at least 3 cm is needed below the amalgam container to slide it off and there must be enough space above to dismantle the separator. **Inside the Hydrocyclone, the max. pressure acceptable is 80 kPa.** A higher pressure would cause a flow exceeding the allowed 5.5 l/min and would impair the separation and jeopardize the integrity of the separator. If a pump is used to force the fluid into the separator, its max. pressure must be checked, to prevent it from exceeding the mentioned value. To avoid the clogging of the holes (2) that could only reduce the flow, the fluid must be filtered through a filter-net not larger than 0.8 mm before entering the separator. The AC 56 P.B. (21) must be fixed by means of the holes on the rear part of the cone.

Maintenance and periodic checks reserved to authorized engineers (Fig. 1-2-5 , page 52-53-55)

Besides the amalgam container replacement, **carry out a functioning test of the photodiodes** (22) by placing a metal tool between the two functioning photodiodes. If necessary, **clean the liquid inlet holes** (2) and the cyclone cone at least once a year. **After disconnecting from the power supply the appliance** where the separator is installed on empty the separator from the liquids, replace the amalgam container. The holes can be reached by unscrewing the two screws (16) and dismantling the separator as per split-up drawing, see fig. 5: the holes can be seen and inspected then. Gauge the diameter of the holes (it must not exceed 2.8 mm) with a caliper. Check the cleaning and integrity of the cone; replace seals, lubricating new ones with Lubri-Jet before assembling them.

Notes

On request, we can supply kits for the installation of the separator - outside the main appliances present on the market - and others - certified by dental unit manufacturers - for the installation inside the dental unit water-systems. We would like to remind that in some countries, besides the compliance to the ISO 11143 regulation on the amalgam separator, the installation of an electronic amalgam level control is compulsory. The manufacturer, dealers, and agents are at disposal for further details.

Guarantee

The appliance is guaranteed for a year from the date of sale, provided that the guarantee card addressed to the manufacturer is returned with date of sale, retailer's stamp and customer's name.

The guarantee and manufacturer's liability cease in case the appliance and the plants have been tampered with operations of any kind, carried out by unqualified persons and unauthorized by the manufacturer.

Transport of secondhand appliances

- Prior to packing, clean and sanitize with Puli-Jet+ or with other suitable disinfectant.
- Place the appliance into a polyethylene bag, seal it and pack in a 3-layer corrugated board.

ISO HYDROCYCLONE APPLICATIONS

ISO Hydroseparator (Mini-Separator and ISO Hydrocyclone)

The ISO Hydroseparator is manufactured in two versions: single-surgery and multi-surgery with and without electronic control of the amalgam separator.

Short instructions (fig. 6-7-8 pag. 56-59)

The appliance is to be assembled inside the dental unit (or inside a floor cabinet, next to the dental unit) and the aspirated fluid (filtered through a 0,8 mm filter) is conveyed to the Mini-Separator, by means of a slightly sloping pipe (fig. 6 part. A) to enhance its operation. Air exits from the passage (B), when the liquid level has reached the two long probes (part. C) the drainage pump (D) sends the liquid to the ISO Hydrocyclone (E), the amalgam separates from the liquid and collects in the container (F) whilst the liquid is sent to sewers (outlet 4). When the liquid reaches the shorter probe (H), the suction stops and starts up again when the liquid retracts from the two longer probes.

The spittoon drainage unit (I) starts the aspiration when liquid is detected. The spittoon rinsing is to be timed at 20 seconds. Filter (L) must be cleaned everyday. Through the electronic control of the amalgam level, over-filling of the Mini-Separator and amalgam container will be detected. AC 35 electrical diagram is shown on page 59 (fig. 13) together with electrical connections of Mini-Separator, draining pump, spittoon draining unit and amalgam level control.

Turbo-Jet modular and ISO Hydrocyclone

(fig. 8-9-10-14-16 pages 57-58-60-62)

Turbo-Jet 1 and 2 modular can be completed with the ISO Hydrocyclone, (fig. 8). Turbo-Jet 3 modular is combined with two ISO Hydrocyclones (fig. 9). The liquid coming out of the centrifugal separator creates pressure into the ISO Hydrocyclone, which separates amalgam. A relief valve detail A, fig. 10, opens and recycles the excess liquid, detail B fig. 10. The ISO Hydrocyclone can be equipped with an electronic control of the amalgam level, relevant electric diagrams and connections of the level control to the electric circuits AC 54/55 and AC 48 of the Turbo-Jet 1 and 2 modular, are illustrated on page 60, to circuits AC 54/55 of Turbo-Jet 3 modular, page 62.

Detailed instructions shall be supplied with the equipment.

***(General recommendations and biologic hazard)**

Pal 22 and 50 VE Pumps and ISO Hydrocyclone (fig. 11-15 page 58-61)

The centrifugal separator receives the fluid that is conveyed to the Hydrocyclone ISO B through pipe A. Detailed instructions shall be supplied with the equipment.

Detailed instructions shall be supplied with the equipment.

***(General recommendations and biologic hazard)**

Large systems with separating tanks and ISO Hydrocyclone

(fig. 12-16 pages 58-62)

The draining pump (F fig. 12), started by the presence of liquid in the Separator-Tank (D), raises the pressure in the ISO Hydrocyclones (C) that separate amalgam and drain the liquid.

The P.B.AC 54/55 (page 62) shows the connections between the amalgam level electronic control (if present) and the electric circuit controlling the Separator-Tank.

Detailed instructions shall be supplied with the equipment.

***(General recommendations and biologic hazard)**

Single-surgery ISO Hydrocyclone (fig. 17-18 page 63)

The single-surgery ISO Hydrocyclone is an amalgam separator for air systems with air-water separator in the dental unit. The single-surgery ISO Hydrocyclone unit shall be placed at the bottom and on the left of the chair.

The spittoon shall be timed for a maximum of 20 seconds. The liquid arriving from the spittoon and from the air-water separator enters the filter (A), possible air bubbles are exhausted through the pipe (B) and the solenoid valve (C), the pump (D) pushes the liquid into the ISO Hydrocyclone (I) that separates the amalgam and drains the liquid (E). To replace the amalgam container (F), ("General recommendations and biologic hazard") open the tap (G), remove the cover (H) to drain the liquid contained in the ISO Hydrocyclone (I) and in the amalgam container (F), then follow the instructions given on page 15.

Detailed instructions shall be supplied with the equipment.

***General recommendations and biologic hazard)**

Maintenance and periodical checks reserved to authorised engineers

Test photodiodes operation (detail 22 page. 53) by inserting a metal object between the two photodiodes. Ensure lights are working.

Short-circuit sensors and check operation of:

- aspiration system in the single-surgery units
- closure of electro-pneumatic valve in multi-surgery units
- draining pump for both units

Remove the lid (detail D page 56) pull out the **draining pump, ensure it is working** and clean the area before replacing it. **Carry out operation tests of the printed board, tracks, relays and any part that is subject to wear and tear.**

For every replacement use original spare parts. **Check the condition of pipes** carrying aspirated liquids, in particular the part between the draining pump (or the separator) and the Hydrocyclone (detail 3 page 56). In case of cracks, however small, pipes should be replaced.

Make sure that staff in charge of cleaning and sanitising the aspirator remember their tasks and that they use the recommended products and procedures. **It is the responsibility of the technician** in charge of general maintenance **to instruct surgery staff as to the daily maintenance routines required.**

***General recommendations and biologic hazard**

Before servicing any used equipment, carry out some washing operations with Puli jet+, or other disinfectant appropriate for this purpose. **Disconnect electric supply** and padlock mains switch if provided. **Before servicing the equipment wear disposable gloves, goggles, mask and disposable overall.** **Aspirated debris is infected and contaminated.** The draining pump (and the separator) put the aspirated liquid under pressure and push it into the Hydrocyclone (detail E page 56). **Should the pipe or the equipment under pressure break the biologic hazard increases.**

Maintenance

After every surgical and particularly long operations, rinse the system thoroughly by aspirating clean water. At the end of every working day, cleanse and sanitise the equipment by aspirating a Puli-Jet + and hot water (50°) solution, according to the

instructions, and replace dirty filters. In the morning, before starting work, make sure that the filter contains sanitising anti-foam, make sure also that liquid anti-foaming is available in the technical room.

The use of products not recommended by the manufacturer may jeopardise the correct operation of the equipment and cause damage.

Important Recommendations

- *The manufacturer, dealers, agents and authorised technicians are available for advice, instructions and the supply of documentation, spare parts and whatever may be of help.*
- *Technical requirements, product improvements, legislative and operational problems, difficulty in finding finished or semi-manufactured products may lead the manufacturer to introduce changes to the production without previous notice.*

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CATTANI S.p.A.
ESAM S.p.A.