**FLEX**  
**Monoclonal Mouse**  
**Anti-Human**  
**Cytokeratin, High Molecular Weight**  
 Clone 34 $\beta$ E12  
**Ready-to-Use**  
 (Dako Autostainer/Autostainer Plus)

**Code IS051**

**Intended use**

For in vitro diagnostic use.

FLEX Monoclonal Mouse Anti-Human Cytokeratin, High Molecular Weight (HMW), Clone 34 $\beta$ E12 (1), Ready-to-Use (Dako Autostainer/Autostainer Plus), is intended for use in immunohistochemistry together with Dako Autostainer/Autostainer Plus instruments. This antibody is useful for the identification of basal cells and squamous epithelium in various tissues. Positive results aid in the differentiation of benign prostatic glands from prostate adenocarcinoma and for detecting squamous differentiation in poorly differentiated tumors. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

**Summary and explanation**

Cytokeratins are intermediate filament cytoskeletal proteins essential to development and differentiation of epithelial cells. Approximately twenty different cytokeratins have been identified and are classified and numbered according to molecular weight and isoelectric points (2). In general, most low molecular weight cytokeratins (40–54 kDa) are distributed in non-squamous epithelium, Moll's catalog numbers 7-8 and/or 17-20 (3). High molecular weight cytokeratins (48–67 kDa) are found in the squamous epithelium, Moll's catalog numbers 1-6 and/or 9-16 (3).

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures.

**Reagent provided**

Ready-to-use monoclonal mouse antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L sodium azide.

Clone: 34 $\beta$ E12(1) Isotype: IgG<sub>1</sub>, kappa.

**Immunogen**

Solubilized immunogen keratin extracted from human stratum corneum (1).

**Specificity**

Anti-Cytokeratin, High Molecular Weight (Anti-CK HMW), Clone 34 $\beta$ E12 has been shown to react with the 66, 57, 51 and 49 kDa proteins in Western blotting corresponding to cytokeratins 1, 5, 10 and 14 of the Moll catalog (1, 2, 4).

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused solution should be disposed of according to local, State and Federal regulations.

**Storage**

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

**Quick guide**

Step		Comments
Fixation	Formalin	
Pre-treatment	EnVision FLEX™, High pH (Code K8004)	20 min HIER, 3-in-1 using PT Link and PT Link Rinse Station
Dilution Range	Ready-to-use	20 min incubation
Dilution buffer	Pre-diluted	
Negative Control	FLEX Negative Control, Mouse (Code IS750)	20 min incubation
Visualization	EnVision™ FLEX, High pH (Code K8010)	20 min incubation, 2x5 min DAB+ incubation
Counterstain	EnVision™ FLEX Hematoxylin (Code K8018)	5 min incubation
Control Tissue	Tonsil	Cytoplasmic staining

Slides/Mounting	FLEX IHC Microscope Slides (Code K8020)	Recommended for greater adherence. Permanent mounting required.
Instrumentation	Dako Autostainer/Autostainer Plus	Use instrument-specific vials (S3425)

*\*The user must always read the package insert for detailed instructions of the staining procedure and handling of the product.*

### Specimen preparation

**Paraffin sections:** The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm.

**Pre-treatment:** Pre-treatment of formalin-fixed, paraffin-embedded tissue sections with heat-induced epitope retrieval (HIER) is required. Optimal results are obtained by pretreating tissues with HIER using diluted EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8004). Deparaffinization, rehydration and epitope retrieval can be performed in Dako PT Link (Code PT100/PT101). For details, please refer to PT Link User Guide. The following parameters should be used for PT Link: Pre-heat temperature: 65 °C; epitope retrieval temperature and time: 97 °C for 20 (±1) minutes; cool down to 65 °C. Remove slide rack from PT tank and immediately dip slides in jar/tank (e.g., PT Link Rinse Station, Code PT109) containing diluted room temperature EnVision™ FLEX Wash Buffer (20x) (Code K8007). Leave slides in Wash Buffer for 1–5 minutes.

The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended. After staining the sections must be dehydrated, cleared and mounted using permanent mounting medium.

### Staining procedure

**Visualization:** The recommended visualization system is EnVision™ FLEX, High pH (Dako Autostainer/Autostainer Plus) (Code K8010).

**Program:** Use the protocol given below. Reagents must be appropriately diluted before use. All steps should be performed at room temperature (20–25 °C).

1. Incubate tissue section with 200 µL EnVision™ FLEX Peroxidase-Blocking Reagent (Code DM821) for 5 (±1) minutes
2. Rinse in EnVision™ FLEX Wash Buffer (Code K8007) for 1–5 minutes
3. Incubate tissue section with 200 µL Primary Antibody (Code IS051) for 20 (±1) minutes
4. Rinse in EnVision™ FLEX Wash Buffer (Code K8007) for 1–5 minutes
5. Incubate tissue section with 200 µL EnVision™ FLEX/HRP (Code DM822) for 20 (±1) minutes
6. Rinse in EnVision™ FLEX Wash Buffer (Code K8007) for 1–5 minutes
7. Incubate tissue section with 200 µL EnVision™ FLEX Wash Buffer (Code K8007) for 5 (±1) minutes (Auxiliary step)
8. Rinse in EnVision™ FLEX Wash Buffer (Code K8007) for 1–5 minutes
9. Incubate tissue section with 200 µL EnVision™ FLEX Substrate Working Solution (Codes DM823 and DM827) for 5 (±1) minutes
10. Incubate tissue section with 200 µL EnVision™ FLEX Substrate Working Solution (Codes DM823 and DM827) for 5 (±1) minutes
11. Rinse in EnVision™ FLEX Wash Buffer (Code K8007) for 1–5 minutes
12. Counterstain slide with EnVision™ FLEX Hematoxylin (Code K8018) for 5 (±1) minutes
13. Rinse in deionized water for 1–5 minutes
14. Incubate tissue section with 200 µL EnVision™ FLEX Wash Buffer (Code K8007) for 5 (±1) minutes
15. Rinse in deionized water for 1–5 minutes
16. Mount slides with permanent mounting medium

Programming grid for recommended assay:

The Auxiliary step should be set to “rinse buffer” in staining runs with ≤10 slides. For staining runs with >10 slides the Auxiliary step should be set to “none”. This ascertains comparable wash times.

For details, please refer to the Operator’s Manual for the dedicated instrument. If the protocols are not available on the Dako Autostainer instrument, please contact Dako Technical Services.

**Counterstaining:** The recommended counterstain is EnVision™ FLEX Hematoxylin (Dako Autostainer/Autostainer Plus) (Code K8018). Non-aqueous, permanent mounting medium is recommended.

**Controls:** Positive and negative controls should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include squamous epithelial cells and the cells/structures should display reaction patterns as described for this tissue in “Performance characteristics” in all positive specimens. The recommended negative control reagent is FLEX Negative Control, Mouse (Dako Autostainer/Autostainer Plus) (Code IS750).

**Staining interpretation** The cellular staining pattern is cytoplasmic.

**Performance characteristics**

Normal tissues (5):

Tissue Type (# tested)	Positive Tissue Elements	Tissue Type (# tested)	Positive Tissue Elements
Adrenal (3)	0/3	Ovary (3)	3/3 Surface epithelium (80–100%), cytoplasmic
Bone marrow (3)	0/3	Pancreas (3)	3/3 Pancreatic duct epithelium (50–80%), cytoplasmic
Breast (2)	2/2 Gland and duct epithelium (80%), cytoplasmic	Parathyroid (3)	0/3
Cerebellum (3)	0/3	Pituitary (3)	0/3
Cerebrum (3)	0/3	Prostate (3)	3/3 Basal and ductal epithelium (90%), cytoplasmic
Cervix (2)	2/2 Surface squamous epithelium (100%), cytoplasmic	Salivary gland (3)	3/3 Excretory duct epithelium (100%), cytoplasmic
Colon (3)	3/3 Surface epithelium (10%), cytoplasmic	Skin (3)	3/3 Squamous epithelium (100%), cytoplasmic
Esophagus (3)	3/3 Squamous epithelium (100%), cytoplasmic		3/3 Sweat duct and gland epithelium (80–100%), cytoplasmic
Heart (3/3)	0/3		1/3 Sebaceous gland epithelium (100%) cytoplasmic
Kidney (3)	3/3 Renal tubule epithelium (10–30%), cytoplasmic	Small intestine (3)	3/3 Epithelium (20–50%), cytoplasmic
Liver (3)	3/3 Bile duct epithelium (50%), cytoplasmic	Spleen (3)	0/3
Lung (3)	3/3 Bronchial epithelium (10%), cytoplasmic	Stomach (3)	3/3 Gastric crypt gland epithelium (10–30%), cytoplasmic
	3/3 Alveolar epithelium (50%) cytoplasmic	Testis (3)	0/3

Mesothelial cells (2)	2/2 Mesothelium (100%), cytoplasmic	Thymus (3)	3/3 Squamous epithelium (100%), cytoplasmic
Muscle, cardiac (3)	0/3	Thyroid (3)	1/3 Follicular cells (5%), cytoplasmic
Muscle, skeletal (3)	0/3	Tonsil (3)	3/3 Squamous epithelium (100%), cytoplasmic
Nerve, peripheral (3)	0/3	Uterus (2)	2/2 Endometrial epithelium (5–80%), cytoplasmic

Abnormal tissues: Positive immunoreactivity with anti-CK HMW, 34 $\beta$ E12 antibody has been reported in squamous cell, ductal and transitional cell carcinomas including: squamous cell carcinoma of the skin, lung and nasopharynx; ductal carcinoma of the breast, pancreas, bile duct and salivary gland; transitional cell carcinomas of the bladder and nasopharynx and thymomas (6-11). Anti-CK HMW, 34 $\beta$ E12 has also been shown to positively identify epithelial mesotheliomas, but was unreactive with sarcomatoid or desmoplastic mesotheliomas (12). Variable positivity with anti-CK HMW 34 $\beta$ E12 has been found in adenocarcinomas of ovary, gastrointestinal tract and thyroid. Tumors reported to be largely unreactive with clone 34 $\beta$ E12 include adenomas of endocrine organs, carcinomas of liver (hepatocellular carcinoma), endometrium, kidney and neuroendocrine tumors (6, 7, 11). In prostate, clone 34 $\beta$ E12 positively labels the basal cells of benign lesions including atrophy, atypical adenomatous hyperplasia and post-sclerotic hyperplasia. The absence of staining of the basal layer with anti-CK HMW, 34 $\beta$ E12 is suggestive of the diagnosis of prostate cancer, but not diagnostic of malignancy and requires other morphologic data (3, 7, 13-15). Anti-CK HMW clone 34 $\beta$ E12 has been reported to label tumor cells in a subset acinar adenocarcinomas(3,15). Spurious immunostaining of astrocytomas has also been reported with clone 34 $\beta$ E12 (7,16).

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**References  
Bibliographie  
Literaturangaben**

- Gown AM, Vogel AM. Monoclonal antibodies to intermediate filament proteins of human cells: Unique and cross-reacting antibodies. J Cell Biol 1982; 95:414-24.
- Moll R, Franke WW, Schiller DL. The catalog of human cytokeratins: Patterns of expression in normal epithelia, tumors and cultured cells. Cell 1982; 31:11-24.
- Miettinen M. Keratin immunohistochemistry: update of applications and pitfalls. Pathol Ann 1993; 28:113-43.
- Gown AM, Vogel AM. Monoclonal antibodies to human intermediate filament proteins II. Distribution of filaments in normal human tissues. Amer J Pathol 1984; 114:309-21.
- IR051(IHC003-D09346).
- Gown AM, Vogel AM. Monoclonal antibodies to human intermediate filament proteins. III. Analysis of tumors. Am J Clin Pathol.1985 Oct;84(4):413-24.
- Diagnostic immunohistochemistry: theranostic and genomic applications. Dabbs, David J. 3rd ed. Philadelphia, PA : Saunders /Elsevier, 2010.
- Hurlimann J, Gardiol D. Immunohistochemistry in the differential diagnosis of liver carcinomas. Amer J Surg Pathol 1991; 15:280-8.
- Dairkee SH, Puett L, Hackett AJ. Expression of basal and luminal epithelium-specific keratins in normal, benign and malignant breast tissue. J Nat Can Inst 1988; 80:691-5.
- Thike AA, Cheok PY, Jara-Lazaro AR, Tan B, Tan P, Tan PH. Triple-negative breast cancer: clinicopathological characteristics and relationship with basal-like breast cancer. Mod Pathol 2010;23(1):123-33.
- Gown AM, Vogel AM. Anti-intermediate filament monoclonal antibodies; tissue-specific tools in tumor diagnosis. Surv Synth Pathol Res. 1984;3:369-385.
- Bolen JW, Hammar SP, McNutt MA. Reactive and neoplastic serosal tissue: a light-microscopic, ultrastructural, and immunocytochemical study. Amer J Surg Pathol 1986; 10:34-47.
- Varma M, Amin MB, Linden MD, Zarbo RJ. Discriminant staining pattern of small glandular and preneoplastic lesions of the prostate using high molecular weight cytokeratin antibody—A study of 301 consecutive needle biopsies. Mod Pathol 1997; 10:93A.
- O'Malley FP, Grignon DJ, Shum DT. Usefulness of immunoperoxidase staining with high-molecular-weight cytokeratin in the differential diagnosis of small-acinar lesions of the prostate gland. Virch Arch Pathol Anat 1990; 417:191-6.
- Brimo F, Epstein JI. Immunohistochemical pitfalls in prostate pathology. Hum Pathol 2012;43:313-24.
- Bacchi CA, Zarbo RJ, Jiang JJ, Gown AM. Do glioma cells express cytokeratin? Appl Immunohistochem 1995; 3:45-53.

Catalogue number Référence du catalogue Bestellnummer	Temperature limitation Limites de température Zulässiger Temperaturbereich	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	 Dako North America, Inc. 6392 Via Real Carpinteria, California 93013 USA  Tel 805 566 6655 Fax 805 566 6688 Technical Support 800 424 0021 Customer Service 800 235 5763	 Dako Denmark A/S Produktionsvej 42 DK-2600 Glostrup Denmark  Tel +45 4485 9500 Fax +45 4485 9595  www.dako.com
Manufacturer Fabricant Hersteller	Batch code Code du lot Chargenbezeichnung	Contains sufficient for <N> tests Contenu suffisant pour <n> tests Inhalt ausreichend für „n“ Ansätze		
Use by Utiliser jusque Verwendbar bis	Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU		

PT0020/ Rev C

Edition 07/12

(115639-002)

307118EFG\_002\_IS051/2012.07 12/12