USB2.0 microscope "Pen Scope"

Product-No:93180 Manual no:93180-a www.delock.com

Description

This handy USB2.0 microscope can be applied in different fields, e.g. science, industry or home use. Equipped with an aspherical high-tech lense and an built-in LED illumination X3 you can zoom in objects from 50x to 200x, take snap-shots and time-lapsed videos. Due to its small size and light weight, the microscope is very flexible in its application and can also display positions which are difficult to reach.

Specification

- USB2.0 specification
- Built-in LED illumination X3
- Aspherical high-tech lense
- Features snap-shots and time-lapsed video with image capture and image saving
- Zoom: 50x to 200x
- Pixel: 645(H) x 485(V) (300k Pixel)
- CMOS Sensor
- Power consumption: 65~130mA
- No extra power supply required
- Size: 165 x 38 mm (HxR))
- White balance: AWB (autom.)
- Automatic Iris Aperture
- For Windows 98/ME/2000/XP/Vista

System requirements

- · PC or Notebook with a free USB Port
- Windows 98/ME/2000/XP/Vista

Package content*

- USB2.0 microscope
- Stand
- Software
- User manual

*Make sure that the product package contains all items before operation. If any item is missing or damaged, then please contact your dealer immediately.

Note:

Follow the general instructions for installation of the USB device. If necessary, read the manual of your computer, or ask a specialist for help.

Driver Installation Windows Vista

- 1. Switch on your computer.
- 2. Connect the USB cable with an available USB port of your computer.
- 3. Put the driver CD in the CD ROM drive.
- 4. Afterwards, Windows displays the Hardware wizard.
- 5. Choose "Don't show this message again for this device".
- 6. Open Windows Explorer, select: CD-drive\PenScope\Setup.exe.
- 7. Choose your language and click "Next".
- 8. Click on "Next".
- 9. Choose "Install the driver software anyway".
- 10. Click on "Finish", to restart the computer.
- 11. After restart you can start the Software "Toprocap", and you can use the device.

- 1. Switch on your computer.
- 2. Connect the USB cable with an available USB port of your computer.
- 3. Put the driver CD in the CD ROM drive.
- 4. Afterwards, Windows displays the Hardware wizard. Click on "Cancel".
- 5. Open Windows Explorer, select: CD-drive\PenScope\Setup.exe.
- 6. Choose your language and click "Next".
- 7. Click on "Next".
- 8. When driver is located, click "Continue installation".
- 9. Click on "Finish", to restart the computer.
- 10. After restart you can start the Software "Toprocap", and you can use the device.

Support Delock

If you have further questions, please contact our customer support support@delock.de / www.delock.com or phone +49 30 81789324. You can contact our service hotline during the following times: Mon – Fri: 9:00 – 17:30 pm charge 0,12 € / min*

* A connection fee applies according to fee schedule of Deutsche Telekom AG or your local telephone provider.

You can find current product information on our homepage: www.delock.com

Final clause

Information and data contained in this manual are subject to change without notice in advance. Errors and misprints excepted.

Copyright

No part of this user manual may be reproduced, or transmitted for any purpose, regardless in which way or by which means, electronically or mechanically, without explicit written approval of Delock.

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