DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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Thanks for purchasing this equipment. We hope you enjoy the precision balance Serie 5153. We advise you to care the equipment following the instructions of this manual.

Nahita develops its products according to directives CE and with special emphasis on ergonomics and user's safety.

The quality of used materials and a right procedure will assure a long life of the equipment.

An inappropriate used may cause accidents, electric shocks, short circuits, fires, injuries, etc. Please read Maintenance section.

TO GET THE BEST RESULTS AND A HIGHER LIFE OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING.

Please consider the following points:

- ◆ This manual is an inseparable part of Nahita digital balance series 5153, so it should be available for all users.
- ◆ Handle with care to avoid sudden movements, knocks, falls of heavy or sharp objects; avoid spilling liquids inside.
- ◆ Never dismantle the equipment to repair it by yourself, you will lose the warranty and it might cause inadequate operation as well as injuries to the people who manipulate it.
- ◆ To prevent fire or electric shocks, avoid dry and dusty environments. If this happens, disconnect immediately the equipment from the main current.
- ◆ If you have any doubt please contact with your wholesaler (installations, setting up, working). You can also send your doubts to Nahita Technical Service (asistencia@auxilab.es) or calling 807117040 (0.35 euros/min).
- ◆ This equipment is protected under Warranty and consumer goods regulations (10/2003).
- ◆ Warranty does not cover overhauls.
- ◆ Manipulation by unauthorized personal means total loss of the warranty.
- ◆ Warranty does not cover fuses, accessories or spare parts with normal wear.
- ◆ Please keep the purchasing invoice to have the right to claim or asking for warranty coverage. If you send the equipment to the Technical Service, attach the original or a copy of the invoice.
- ◆ Fill in and send the warranty document before 15 days after purchasing.
- ◆ Manufacturer reserves the right to introduce modifications or improvements.

Cleanliness

- ◆ To clean metallic parts, stainless steel, aluminum, paintings, etc, do not use scourers or products that may scratch.
- ◆ We recommend use a free fuzzy cloth wet with soapy water without abrasive products.

ATTENTION! AUXILAB DOES NOT ADMIT REPARATIONS IF THE EQUIPMENT IS NOT CLEAN AND DISINFECTED.

NOTE: In accordance with the legislation currently in force, Nahita balances series 5153 cannot be used for:

- .- Make commercial transactions.*
- Calculate taxes, duties, remunerations, compensations and other similar canons.*
- Judicial inspections.*
- Prepare medicines to order, as well as analysis made at medical and chemist laboratories.*
- Determinate retail price or total amount of prepared precanning.*

INSTRUCTIONS ABOUT PROTECTION OF THE ENVIRONMENT

Do not deposit this equipment in the ordinary garbage; deposit it in a recycle collect point. The balance does not contain dangerous or toxic elements for humans but it is detrimental for environment. Materials may be recycled, with this act you contribute to protect the environment. If you need more information about recycle collect points, please contact with the administration of your community.

| | | |
|---------------------------------|--------------------------------------------------------------------|----------------------------------------------------------------------------------|
| Value changes constantly | 1 Balance suffers vibrations | 1 Put the balance in a place without vibrations. |
| | 2 Door of the balance is closed. | 2 Close all doors. |
| | 3 There is an stranger object and the balance does not stabilize. | 3 Check internal part of the balance is empty and there is no stranger objects. |

| | | |
|--------------------|-----------------------------------------------|----------------------------------------------------------------------|
| Wrong value | 1 Balance is not calibrated | 1 Calibrate the balance |
| | 2 Balance has not been tared before weighing | 2 Tare the balance before each weighing |
| | 3 Balance is not leveled | 3 Level the balance until the bubble placed in the back is centered |

9. MAINTENANCE AND CLEANNESS

Note: all these utilization rules have no value if you do not maintenance the balance continually.

For a good working is required follow some advices:

- ◆ Follow the instructions and advices given in this manual.
- ◆ This manual should be available for all users.
- ◆ Avoid sudden movements or knocks as well as direct sun light or draughts. Handle the balance carefully, like a precision instrument.
- ◆ The balance has an adapter; it has to be connected to a grounded current, it should be visible to unplugged it in case of emergency.
- ◆ Unplugged the adapter pulling its base, never the wire.
- ◆ Leave space enough around the balance for a good ventilation. Never use the balance squashed, for example a shelf.
- ◆ Do not use sharp objects as pens, etc. to touch the buttons of the front panel; only use your fingers.
- ◆ Do not place inside the balance an object heavier than the indicated on the range of the balance, the sensor may be damaged.
- ◆ Do not immerse the balance and do not spill liquids on it.
- ◆ If you work with batteries, remove them when you are not going to use the balance for a long time.
- ◆ If some liquid come into contact with electrical components, turn off and unplugged it. Send the balance to the technical service as soon as possible.
- ◆ Use original components and spare parts. Maybe you can find similar devices, but use them may damage the equipment.

ATTENTION! AUXILAB DOES NOT ADMIT REPARATIONS IF THE EQUIPMENT IS NOT CLEAN AND DISINFECTED.

INDEX OF LANGUAGES

| | |
|---------------|-------|
| Spanish | 2-11 |
| English | 12-21 |
| French | 22-31 |

INDEX OF CONTENTS

| | |
|--------------------------------------|----|
| 1. USES OF THE EQUIPMENT | 13 |
| 2. DESCRIPTION | 14 |
| 3. TECHNICAL SPECIFICATIONS | 15 |
| 4. INSTALLATION / SETTING UP | 15 |
| 5. CALIBRATION | 17 |
| 6. FUNCTIONS OF THE BALANCE | 18 |
| 7. SOFTWARE | 19 |
| 8. TROUBLESHOOTING | 19 |
| 9. MAINTENANCE AND CLEANLINESS | 20 |
| ANEX I: CE CERTIFICATE | 32 |


1. USES OF THE EQUIPMENT

Precision electronic balances are basic instruments at any laboratory when weighing tasks are needed to be carried out. Nahita balances are developed with the most advanced technology, they are reliable and durable instruments, with solid, clean and ergonomic design suitable for educational, industrial or research use.


2. DESCRIPTION

- 1.1 Front panel
- 1.2 Camera
- 1.3 Pan
- 1.4 Outlet RS232 (back part)
- 1.5 Connection adapter (back part)

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to weigh from the original container to other container. When the decanting is finished, its weight (with negative sign) will be show on the screen of the balance.

◆ If you want to compare the weights of two different substances, put one of them on the balance and press "TARE". When you see the value "0.000 g" removes it and the balance will show with negative sign the value of the weight. Put the second substance on the pan and directly you will see the difference between both weighings.

6.2 Counting pieces (Mode "PCS")

In mode "PCS" the balance calculate the number of pieces put on the pan. Please follow these steps:

- ◆ Press button "MODE" (1.5), select mode "PCS".
- ◆ Screen shows "-COU-" and finally the flashing value "10". This value is showed by default, it indicates the number of pieces that have to be placed on the pan and they will be the reference for next counting.

NOTE: If you want to modify this value, press button "CAL" successively and select one of these values 10, 20, 50, 100, 150, 200, 250 and 500 pieces.

- ◆ Place 10 pieces or selected quantity on the pan.
- ◆ Press "MODE" and the value is fixed as a reference.
- ◆ Remove the pieces.
- ◆ Screen shows "0 pcs".
- ◆ Place on the pan the pieces you want to count and the quantity will be show.

7. SOFTWARE

Balance is equipped with universal output RS232 and it can be connected to computers or printers.

The format is:

- ◆ Bauds 9600
- ◆ 8 bits (1 start bit (0), 8 data bits (code ASCII) and 1 stop bit (1)).

8. TROUBLESHOOTING

| PROBLEM | CAUSE | SOLUTION |
|------------------------|-----------------------------------------------------------------------|-----------------------------------------|
| Digits are not visible | 1 Check that the balance receive current. The adapter can be damaged | 1 Adjust the pins. Change the adapter. |

Note: If the screen does not show "0.000 g", press button TARE (1.4) until the mentioned value is fixed and stabilized led is illuminated (0).

- ◆ Press button "CAL" (1.6) for five seconds.
- ◆ Screen shows "–CAL–" after this, it indicates the calibration weight required (flashing).
- ◆ Put into the balance the calibration weight.
- ◆ Screen shows "-----" for several seconds and finally the respective weight of the calibration weight (Ex.: "100.000"), the reference is registered automatically.
- ◆ Remove the calibration weight from the balance.
- ◆ When screen shows "0.000 g", the calibration has finished correctly.

It is required to make the verification of the calibration:

When screen shows "0.000 g", put again the same weight on the pan and check that the value on the screen coincides with the value of the weight.

6. FUNCTIONS OF THE BALANCE

6.1 Weighing

Take into account the temperature of the balance; if you have placed it in another place, wait at least one hour to stabilize the temperature.

Calibrate the balance if you think it is required, follow "Calibration" section.

- ◆ Switch on the balance pressing button "ON/OFF" (1.2), the pan should be empty.
- ◆ When screen shows the value "0.000 g" with stabilized led (0), you can make the reading (if you want to use other measure unit, you can select it pressing the button UNITSW (1.3)).
- ◆ If you want to tare the balance to make the weighing in a container, put the empty container into the balance and press button "TARE" (1.4). Wait until you see the value "0.000 g" stabilized and then you can make the measures.

Note:

1. *The value of the subtractive tare is recorded in the memory of the balance until you press button "TARE" (1.4) again.*

2. *If you want to weigh more than one substance in the same container and you need independent weights, add the first substance and when you obtain the desired weight press "TARE" again (1.4). Screen will show "0.000 g" and you can weigh another substance.*

- ◆ If you want to weigh a substance contained in other container, put the full container in the balance. When its weight is stabilized, press "TARE" (1.4) and screen will show the value "0.000 g". With a spatula or something like that, decant the substance that you want

3. TECHNICAL SPECIFICATIONS

| | | | |
|----------------------------|----------------|----------|----------|
| <i>Code</i> | 55153210 | 55153220 | 55153230 |
| <i>Capacity</i> | 100 g | 200 g | 300 g |
| <i>Readability</i> | 0.001 g | | |
| <i>Repeatability</i> | ± 0.002 g | | |
| <i>Linearity</i> | ± 0.002 g | | |
| <i>Plate diameter</i> | 80 mm | | |
| <i>Internal dimensions</i> | 180x180x175 mm | | |
| <i>External dimensions</i> | 270x190x265 mm | | |
| <i>Weight</i> | ~2.6 kg | | |

Among its main characteristics we emphasize the following:

- ◆ Quickly setting up and stabilization.
- ◆ Easy use, clear reading on screen LCD.
- ◆ Weight and counting functions.
- ◆ Selection of weight units:
- ◆ Grams (g), carats (ct), pounds (lb), ounces (oz), Americans ounces (ozt), grains (GN), pennyweights (dwt).
- ◆ Subtractive tare.
- ◆ Calibration through calibration weights (included).
- ◆ Connection to main supply through adapter (included).

4. INSTALLATION / SETTING UP

Preliminary check

Unpack the Balance, remove the plastic and the polispan protection.

Make sure that it does not present any damages because of the shipment. If the balance presents damages, communicate them immediately to your transport agent.

If you want to return the equipment, please do it within 15 days after purchasing with original packing, all accessories and documents included

Check the accessories you have to receive with the equipment:

- Pan
- Source adapter
- Calibration weight (included)
- User's manual
- Warranty certificate

Installation

Before operating with this equipment, it is advisable to check components and basic

foundations, as well as functions of its controls.

TO GET THE BEST RESULTS AND A HIGHER LIFE OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING.

- ◆ Put the balance on an horizontal and stable table, leave space enough to work.
- ◆ Do not place the equipment close to heating sources (burners, gas welding torches...) or magnetic material and do not expose directly to sun light, etc.
- ◆ Level the balance through screw legs, center the level bubble.
- ◆ Nahita digital balances series 5153 can work with batteries or connect them to the power supply through the AC adapter.
- ◆ Place the AC adapter
- ◆ Take into account that polarity of the adapter should be to avoid  to avoid damages in the electrical system.
- ◆ Insert the adapter connector in the hole that you find on the back of the balance, and plug into grounded main source.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installations or damages to people suffered by improper use of the electric connection. Tension should be 220 V, 50 Hz +/- 10%.

- ◆ Battery installation
 - Open the lid of the battery compartment placed at the lower.
 - Put eight alkaline batteries 1.5 V type AA and put them correctly according to its polarity.
 - Replace the lid of the compartment.
- ◆ Put the stainless steel pan. Both pieces have to be fitted carefully in the right order, first the hoop and after the pan.
 - Balance is ready for its setting up.

Setting up

NOTE: Before starting to use the balance you have to calibrate it following "Calibration" section of this manual. When you observe important variations or weight errors you should calibrate the balance again.

- ◆ Be sure the pan is empty and press ON/OFF (1.2) to switch on the balance.
- ◆ Screen (1.1) is on and make the autotest, at the beginning, it shows the capacity of the balance and finally "0.000 g".

- ◆ Select by MODE button (1.5) the function you want to work with; weight or counting pieces.
- ◆ Press button UNITSW (1.3) you can select the unit you want to weigh with: grams (g), carats (ct), pounds (lb), ounces (oz), troy ounces (ozt), grains (GN), pennyweights (dwt).
- ◆ The balance is ready for working.
- ◆ Switch off the balance pressing ON/OFF button when you had finished the work.
 - The following table shows unit conversion.

Table 1: "UNIT CONVERSION"

| <i>Abbreviation</i> | <i>Name</i> | <i>Conversion</i> |
|---------------------|--------------|-------------------|
| g | Gram | 1 g |
| ct | Caract | 0.1999694 g |
| lb | Pound | 453.592 g |
| oz | Ounce | 28.3495231 g |
| ozt | Troy ounce | 31.1034 g |
| GN | Grain | 0.064799 g |
| dwt | Pennyweights | 1.55517 g |

NOTE: To prevent oscillations in weighing, it is important to keep close the balance as well as avoid vibrations on the working table while the measuring.

5. CALIBRACIÓN

When you observe important variations in the weighing or wrong weighing you have to calibrate the balance. Factors such as gravity variation depending on the latitude, the change of working place as well as an inappropriate or sudden use lead to the necessity of calibrate the equipment again.

Important!: do not move the balance during the calibration

Here you have the weights to calibrate your balance according to the model

Table 2: "CALIBRATION WEIGHT"

| | |
|----------|-------|
| 55153210 | 50 g |
| 55153220 | 100 g |
| 55153230 | 200 g |

Proceed like this:

- ◆ Switch on the balance pressing ON/OFF (1.2), the pan of the balance should be empty.
- ◆ The balance makes the autotest and it stabilizes at "0.000 g".