Pacific Hemostasis® Thrombin Time Reagent (For Thrombin Time Determinations)

- I. Intended Use**
The Thrombin Time (TT) is a simple test to screen for conditions that can interfere with the conversion of fibrinogen to fibrin.
- II. Summary and Principle**
Fibrinogen (Factor I) is a soluble plasma protein that is instrumental in the normal coagulation process. Following trauma or injury, fibrinogen is converted to an insoluble fibrin clot by a two-stage process. In stage one, thrombin cleaves fibrinogen to form a fibrin monomer. In stage two, these fibrin monomers aggregate to form the insoluble fibrin polymers that are recognized as the end point in thrombin clotting assays.
A low potency thrombin is added to undiluted plasma and clot formation is timed. The TT may be prolonged when any of the following conditions exist:
- Decreased fibrinogen levels
 - Dysfunctional fibrinogen molecules (dysfibrinogenemia)
 - Heparin therapy
 - Increased Fibrinogen Degradation Products (FDP)
 - Presence of abnormal serum globulins or increased immunoglobulins ^{1,3}
- III. Reagent**
For *in vitro* diagnostic use.
Composition
Thrombin Time Reagent is prepared from *bovine* source material with buffers. Store between 2-8°C.
Reconstitute with 1.0 mL distilled water. According to procedure below agitate gently until solution is complete. The reconstituted material is stable for 4 hours when stored in original glass container at 2-8°C. Reconstituted material is stable for 8 hours when stored in plastic at 2-8°C. ⁴
- IV. Specimen Collection**
Anticoagulant: 3.2% (0.109M) sodium citrate is recommended for coagulation assays. A ratio of 9 parts blood to 1 part citrate is appropriate for most samples. However for the occasional patient with a hematocrit of less than 20% or greater than 55%, this ratio should be adjusted to ensure valid results. The following formula can be used to calculate the amount of anticoagulant needed for a specific hematocrit and blood volume:
 $0.00185 \times \text{blood (mL)} \times (100 - \text{hematocrit}) = \text{anticoagulant (mL)}$ ¹
Sample Collection: Blood samples can be obtained either by venipuncture or from an indwelling catheter. Venipuncture should be non-traumatic, avoiding hemolysis and contamination by tissue fluids. If blood is drawn from an indwelling catheter, the line should be flushed with saline and the first 5 mL of blood must **not** be used for coagulation assays. Samples that have less than 90% of the expected fill volume in the collection tube should be rejected. ²
Specimen Processing: Blood must be gently mixed with anticoagulant immediately after collection and transported to the lab as quickly as possible. Centrifuge with the minimum force necessary to yield platelet poor plasma (PPP) or platelet free plasma (PFP). Contamination of the plasma with excess platelets could cause falsely elevated levels of some factors, including fibrinogen (Factor I). Remove plasma without disturbing the buffy coat or red cells unless testing will be performed immediately. Cover samples to prevent pH changes that may affect test results. Specimens maintained at 22-24°C must be tested within 2 hours, or within 4 hours if held at 2-4°C. For longer periods of time, samples can be frozen at -20°C for 2 weeks, or at -70°C for 6 months. Thaw samples rapidly at 37°C and test immediately. ¹
- V. Test Procedure**
Materials Provided
Thrombin Time Reagent, 10 x 1 mL
Incubate 0.2 mL of undiluted plasma at 37°C for 3 minutes. Add 0.1 mL of Thrombin Time Reagent and time clot formation.
- VI. Results**
As with any test, normal ranges will vary according to the laboratory and the instrument used. However, when tested with Thrombin Time Reagent following the procedure described in Section V, normal individuals will generally exhibit a Thrombin Time of approximately 8-14 seconds.

- VII. Limitation**
As with any test normal ranges will vary according to the laboratory and the instrument used. Always follow instrument and reagent manufacturer's guidelines.
- VIII. Expected Values**
Normal individuals will generally exhibit a Thrombin Time of approximately 8-14 seconds.
- IX. References**
1. Bick, R.L., et al. *Hematology: Clinical and Laboratory Practice*. Mosby: 1993.
 2. NCCLS. 1991. *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays; Approved Guideline 3rd Edition*. NCCLS document H21-A3. NCCLS, Wayne, PA, 1998.
 3. Powers, L.W. *Diagnostic Hematology: Clinical and Technical Principles*. Mosby: 1989.
 4. Stability data found in DHF.

ORDERING INFORMATION		
Cat. No.	Description	Contents
100011	Thrombin Time Reagent	10 x 1 mL

FISHER DIAGNOSTICS® LIMITED WARRANTY

Fisher Diagnostics (FD) warrants to the purchaser only that FD products will perform as described on their labeling and product literature. Purchaser must determine the suitability of FD products for their specific applications. FD's sole obligation will be, at its option, to either replace a non-conforming or defective product, or return the purchase price. FD DISCLAIMS ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR ANY PARTICULAR PURPOSE. Neither FD nor its affiliates shall, in any event, be liable for incidental or consequential loss or damage.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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IX. Bibliographie

- Bick, R.L., et al. *Hematology: Clinical and Laboratory Practice*. Mosby: 1993.
- NCCLS. 1991. *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays; Approved Guideline 3rd Edition*. NCCLS document H21-A3. NCCLS, Wayne, PA, 1998.
- Powers, L.W. *Diagnostic Hematology: Clinical and Technical Principles*. Mosby: 1989.
- Données de stabilité provenant de DHF.

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- Bick, R.L., et al. *Hematology: Clinical and Laboratory Practice*. Mosby: 1993.
- NCCLS. 1991. *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays; Approved Guideline 3rd Edition*. NCCLS document H21-A3. NCCLS, Wayne, PA, 1998.
- Powers, L.W. *Diagnostic Hematology: Clinical and Technical Principles*. Mosby: 1989.
- Los datos de estabilidad proceden de DHF.

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IX. Riferimenti

1. Bick, R.L., et al. *Hematology: Clinical and Laboratory Practice*. Mosby: 1993.
2. NCCLS. 1991. *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays; Approved Guideline 3rd Edition*. NCCLS document H21-A3. NCCLS, Wayne, PA, 1998.
3. Powers, L.W. *Diagnostic Hematology: Clinical and Technical Principles*. Mosby: 1989.
4. Dati relativi alla stabilità del prodotto rinvenuti nell'archivio storico.

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EC REP
WMDE
Bergenweg 18
6085 AT Horn
The Netherlands

Fisher Diagnostics
a division of Fisher Scientific Company, LLC
a part of Thermo Fisher Scientific Inc.
Middletown, VA 22645-1905 USA
Phone: (800) 528-0494
(540) 869-3200
Fax: (540) 869-8132

JL840807 (R0)

Thermo
S C I E N T I F I C

Symbols Key	
	Manufacturer Hersteller Fabricant Fabricante Fabbricante
IVD	<i>In Vitro</i> Diagnostic Medical Device Medizinprodukt für die <i>in-vitro</i> -Diagnostik Matériel médical pour utilisation diagnostique <i>in vitro</i> Dispositivo médico para diagnóstico <i>in vitro</i> Dispositivo medico per diagnosi <i>in vitro</i>
LOT	Lot Number Chargennummer Numéro de lot Número de lote Numero di lotto
	Use By Verfallsdatum Utiliser jusque Fecha de caducidad Da utilizzare entro
	Temperature Limitation Temperatureinschränkungen Limite de température Limite de temperatura Limiti di temperatura
CE	CE Mark CE-Markierung Marquage CE Marca CE Marchio CE
REF	Catalogue Number Katalognummer Référence catalogue Número de catálogo Numero di catalogo
	Consult Instructions for Use Bedienungsanleitung lesen Consulter le manuel d'utilisation Consultar las instrucciones de uso Consultare le istruzioni per l'uso
EC REP	Authorized Representative in the European Community Autorisierte Vertretung in der Europäische Gemeinschaft Représentant agréé pour la Communauté européenne Representante autorizado en la Comunidad Europea Rappresentante autorizzato nella Comunità europea
	Biological Risks Biologische Risiken Risques biologiques Riesgos biológicos Rischi biologici