ENGLISH

INTENDED USE

The PreVu® Non-Invasive Skin Cholesterol Point of Care (POC) Test (PreVu POC Test) is an *in vitro* diagnostic test for the quantification of skin cholesterol. Like serum cholesterol, elevated skin cholesterol is associated with an increased risk of coronary artery disease. Thus, skin cholesterol can be used as part of risk assessment for coronary artery disease (CAD).

SUMMARY AND EXPLANATION OF THE TEST

The assessment of skin cholesterol allows for the quantification of cholesterol present in the tissues. Skin contains approximately 11%, by weight, of all body cholesterol and has been cited as “mirroring” vascular changes associated with age and atherosclerosis¹. A number of studies have demonstrated a relationship between skin cholesterol levels and CAD (reviewed in Lopukhin¹) and more recently studies showing a relationship between skin cholesterol levels and treadmill stress test outcome², and skin cholesterol correlated with angiographically proven CAD³ have been reported.

The PreVu POC Test affords a practical and effective method of determining skin cholesterol. Skin cholesterol as measured by the PreVu POC Test can be used in conjunction with clinical evaluation and other lipoprotein tests as part of risk assessment for coronary artery disease.

PRINCIPLE OF THE PROCEDURE

The PreVu POC Test is a rapid two-step test for the detection of skin cholesterol. The test is performed on the palmar surface of the hand and requires no prior fasting. The principle of the procedure is as follows:

- A solution (Detector) containing digitonin and horseradish peroxidase enzyme (HRP) conjugated to a copolymer is first placed on the palm of the hand for one minute. During this time, digitonin forms complexes with cholesterol in the skin. The greater the amount of cholesterol present in the epidermis, the greater the amount of Detector will attach.
- Unbound Detector is blotted off the palm with a Blotting Stick.
- A second solution, Indicator, is added to the test sites for two minutes. The Indicator is a substrate that reacts with HRP in the Detector to produce a blue-green colour. The greater the amount of Detector bound to skin cholesterol, the darker the colour will become.
- A hand-held spectrophotometer is used to determine an epidermal skin cholesterol value by measuring the hue² of the colour developed in the test well.

Negative and positive controls applied to the palm in parallel with the Detector solution serve as reagent and procedural controls to monitor test performance.

MATERIALS

Materials Supplied in the PreVu POC Test Kit
(Contains 40 tests)

Alcohol Swab

Forty four (44) isopropyl alcohol (70% v/v) swabs

Foam Pad

Forty four (44) foam applicator pads with medical-grade adhesive protected by a paper liner.

Wet-Nap® moist towelette

Forty four (44) pre-moistened hand cleaning wipes.

Blotting Stick

Forty four (44) blotting sticks (cellulose fibre rods).

NOTE: Additional supply volumes provided for practice/unintentional wastage.

Detector (green cap)

One dropper bottle with 2.5 mL of digitonin horseradish peroxidase conjugate in StabilZyme SELECT® Stabilizer⁴ as a buffer solution. Store refrigerated. Stable at 2-8°C (36-46°F) for duration of labelled expiry date.

PC - Positive Control (red cap)

One dropper bottle with 2.5 mL of digitonin horseradish peroxidase conjugate in StabilZyme SELECT® Stabilizer⁴ as a buffer solution. Store refrigerated. Stable at 2-8°C (36-46°F) for duration of labelled expiry date. The PC solution is 10X more concentrated than the Detector solution.

Indicator (blue cap)

One dropper bottle with 6.0 mL of 3, 3’, 5, 5’–tetramethylbenzidine (TMB) – hydrogen peroxide solution. Store refrigerated. Stable at 2-8°C (36-46°F) for duration of labelled expiry date.

Materials Required But Not Provided

PreVu POC Handheld Spectrophotometer (PVU-SPECT-101-UNI)

WARNINGS

Methylisothiazolone warning:

Methylisothiazolone in combination with methylchloroisothiazolinone (MI/MCI) is routinely used in the cosmetics industry as a biocide. MI/MCI has been reported to be slightly irritating to the skin at some concentrations. The PreVu POC Test reagents contain only 100ppm of MI and are not expected to cause skin irritation. However, use of the PreVu POC Test should be discontinued in individuals who experience skin irritation.

NOTE: For External Use Only!

PRECAUTIONS

- Dispense hanging drop of each reagent. In order to prevent contamination, do not touch dropper bottle nipple assembly to surface of skin.
- Do not use kit components beyond expiry date and do not mix reagents from kits with different lot numbers.
- Avoid exposure of reagents to excessive heat or light.
- Observe good laboratory practices, i.e., do not eat, drink, smoke or apply cosmetics in areas where biological reagents are used.
- Bring reagents to room temperature (18-25°C or 64-77°F) prior to testing. Return all reagents to refrigerator after use.

STORAGE INSTRUCTIONS

Store reagents between 2-8°C (36-46°F) for duration of labelled expiry date.

TEST PROCEDURE

Preparation for the Assay

- Bring all reagents to room temperature (18-25°C or 64-77°F).
- Activate the PreVu POC Handheld Spectrophotometer and consult User Manual.
- Finish assay without interruption.

Assay Procedure

Follow test directions using the PreVu POC Handheld Spectrophotometer.

- Hands must be clean prior to the test. Wash hands with soap and water or vigorously clean outside portion of palm of one hand with Wet-Nap® moist towelette provided. Dry the palm using a clean paper towel or other absorbent, lint-free tissue. Then wipe the outside portion of the palm of the hand vigorously with an Alcohol Swab and wait about 30 seconds for the cleaned area to air dry. Discard Wet-Nap® moist towelette but save Alcohol Swab for use after test.
NOTE: Do not use other brands of moist towelette to clean hands.
- Air-dry palm. Palm must be completely dry before proceeding with the next step.
- Remove adhesive backing from Foam Pad. Firmly place Foam Pad on alcohol-wiped portion of hand with the diamond well towards fingers.
- Rest subject’s hand on flat surface.
- Add one drop of **Detector** (green cap) to round **test well**.
- Add one drop of **PC** (red cap) to diamond well.
- Do not add any liquid to square well.
- Wait 1 minute.
- Slowly insert one end of the double-ended Blotting Stick into the round, centre well of the Foam Pad to absorb the Detector reagent. After insertion press firmly and hold for 3 seconds to ensure that all Detector reagent has been absorbed, then invert Stick and insert the second (clean) end into the same well (centre well) of the Foam Pad. Press to make firm contact and gently twist the Blotting Stick approximately one quarter of a turn forward and then back while still making contact with the palm. Once the centre well has been blotted dry, remove the Blotting Stick and slowly insert the second end into the diamond-shaped well to absorb the PC reagent. Discard the Blotting Stick. **NOTE:** Visually ensure that Foam Pad and test wells are completely dry.
- Add one drop of **Indicator** (blue cap) to all wells, including the square well.
- Wait two minutes.
- Immediately read **test well** with the PreVu POC Handheld Spectrophotometer.

Blot excess liquid and remove Foam Pad. Clean palm with Alcohol Swab from Step 1.

Test Validity

- Liquid in square well is colourless.
- Liquid in diamond well is blue-green.

Repeat test on other hand if assay is invalid. Consult Technical Service at 1-888-928-7205 (Canada) if assay fails the repeat test.

QUALITY CONTROL PROCEDURES

Instrument Calibration

When using the PreVu POC Handheld Spectrophotometer, the instrument must be calibrated prior to each test. The instrument will prompt the operator to open the alignment shoe completely and leave it open until a beep is heard. If you have a tethered version of the spectrophotometer please consult the instructions provided with that device for the proper calibration procedure. The instrument will then indicate whether calibration was successful. If calibration was not successful, the instrument will indicate to repeat the calibration procedure. If the calibration error persists, please refer to the User Manual.

LIMITATIONS OF THE PROCEDURE

- The wells in the Foam Pad must not be placed over prominent creases on the palm as this may allow the liquid to leak from the wells during incubation.
- Patients using topical medications and individuals with skin disease (e.g. eczema, psoriasis) on the hands should not be tested.
- Patients who have recently used a dermatological preparation, such as a hand lotion or barrier cream, may produce aberrant skin cholesterol values.

- In the clinical studies supporting the PreVu POC Test safety and effectiveness only a small number of African American subjects were evaluated in relation to other patient subgroups. Consequently, use of the PreVu POC Test for risk assessment of CAD in African American subjects must be carefully considered with all other risk factors before any follow-up patient management decisions are made.
- Consult the PreVu POC Handheld Spectrophotometer User Manual prior to use.

NOTE: It is recommended to test the hypothenar eminence of the palm.

INTERPRETATION OF RESULTS

Skin cholesterol values less than 81 are below the detection limit of the PreVu POC Test. Individuals with skin cholesterol measurements equal to or below 134 are considered to be at normal risk of developing CAD, whereas those with skin cholesterol values greater than 134 are at increased risk (Table 1).

Table 1. Skin Cholesterol Ranges and Interpretation		
Skin Cholesterol Range	Interpretation	CAD Risk Category
<81	Below limit of detection	Normal risk
81-134	Normal range of skin cholesterol	Normal risk
>134	Elevated skin cholesterol	Increased risk

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

Table 2 shows the reproducibility of the PreVu POC Test within-day, between days, and from batch to batch.

Table 2. Reproducibility of PreVu POC: Mean CV (range)			
	Within- Day	Day-to- Day	Batch-to- Batch ^a
Low skin cholesterol	11% (5-19%)	8% (4-17%)	10% (1-17%)
High skin cholesterol	7% (2-13%)	3% (1-5%)	
ⁱ 3 lots of PreVu POC were compared			

CLINICAL DATA

Skin cholesterol levels were measured using the PreVu POC Test, in 750 individuals at three angiography clinics (649 cases) and three out-patient facilities (101 controls). Skin cholesterol levels were significantly lower in control subjects compared to the case population (105 vs. 123, p<0.001). The distribution of controls and cases in the different skin cholesterol ranges is shown below (Table 3).

Table 3. Distribution of Controls and Angiography Clinic Patients by Skin Cholesterol Range		
Skin Cholesterol Range	Controls (n=101)	Angiography Patients (n=649)
<81	12% (n=12)	10% (n=62)
81-134	79% (n=80)	56% (n=363)
>134	9% (n=9)	34% (n=224)

Table 4 shows the relationship between skin cholesterol and CAD as determined by coronary angiography. Skin cholesterol was correlated with CAD even after adjusting for Framingham risk factors. The strongest relationship was with multi-vessel disease. This shows that skin cholesterol provides new information about CAD risk that is independent of Framingham risk factors: age, gender, serum LDL, serum HDL, hypertension status, diabetes status, and smoking status.

Table 4. Relationship Between Skin Cholesterol and CAD Determined by Coronary Angiography (odds ratio, 95% confidence interval, p-value)			
Parameter	Insignifi- cant CAD ^a	Any CAD ^b	Multi- Vessel CAD ^c
Skin cholesterol (10 unit increase)	1.06 (0.99-1.12) 0.06	1.07 (1.01-1.13) 0.02	1.08 (1.03-1.14) 0.004
Skin cholesterol Adjusted ^d (10 unit increase)	1.05 (0.99-1.12) 0.13	1.06 (1.01-1.13) 0.03	1.08 (1.02-1.15) 0.009
^a ≤ 50% stenosis in any coronary arteries ^b >50% stenosis in at least one coronary artery ^c >50% stenosis in more than one coronary artery ^d adjusted for Framingham Global Risk and trial site			

REFERENCES

- Lopukhin YuM. The skin and atherosclerosis (a three-drop test). *Physiochemical Aspects of Medicine Reviews*. 1992; 3:1-124.
- Zawyidiwski R, Sprecher DL, Eveleigh MJ, *et al*. Cholesterol 1,2,3TM: A novel test for the determination of skin cholesterol levels. *Clinical Chemistry* 2001; 47:1302-1304.

- Sprecher DL, et al. Skin Tissue Cholesterol (Skin Tc) is related to angiographically-defined cardiovascular disease. *Atherosclerosis* 2003; 171:255-258.

There are no conditions, warranties or representations, expressed or implied, statutory or otherwise, made by Miraculins Inc., its officers, agents or distributors, regarding the performance of this kit when:

- the validated test procedure described in this package insert is not followed,
- the kit components are used beyond the indicated expiry date, or
- components of different lot numbers are mixed, whether such components are purchased individually or as part of a kit.

Patent Information
The PreVu Non-Invasive Skin Cholesterol test technology is protected by patents and patents pending in North America and Internationally.

Questions?
1-888-928-7205 (Canada)
www.prevu.com

Manufactured By: Miraculins Inc. 6-1250 Waverley Street Winnipeg, Manitoba Canada R3T 6C6	EC Rep: Emergo Europe Molenstraat 15 2513 BH The Hague The Netherlands
---	--

Manufactured for Miraculins by Fisher Diagnostics, part of Thermo Fisher Scientific.

PreVu® is a registered trademark of Miraculins Inc. All rights reserved. 2012.

† All colours can be characterized by hue (the dominant shade), saturation (how much colour of any hue is present) and lightness (the degree of lightness or darkness of a particular colour) in a three-dimensional space using cylindrical coordinates. Lightness is the centre vertical axis and saturation is the horizontal axis that extends from the lightness axis. Hue is the angle at which the saturation axis extends from the lightness axis.

‡ StabilZyme SELECT® is an aqueous solution containing a purified bovine protein and other non-toxic stabilizing chemicals in MOPS (3-(N-Morpholino) propanesulfonic acid), pH 6.2-6.7, buffer. The solution is preserved with the following components: 0.02% methylisothiazolone, 0.02% bromonitrodioxane and 20 ppm Proclin 300 (Rohm and Haas Company).

KIT CONTENTS

(Contains 40 tests)

- 44 Alcohol Swabs – for external use only!
- 44 Foam Pads
- 44 Wet-Nap® moist towelettes
- 44 Blotting Sticks
- Detector Solution (2.5 mL)
- PC (This is the positive control solution) (2.5 mL)
- Indicator Solution (6.0 mL)

NOTE: Additional supply volumes provided for practice/unintentional wastage.

SYMBOL LEGEND

	: Contains N Tests		: Store at
	: Expiration Date		: Lot #
	: Manufacturer		: <i>In Vitro</i> Diagnostic
	: See instructions		: European Conformity

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RÉFÉRENCES

- Lopukhin YuM. The skin and atherosclerosis (a three-drop test). Physiochemical Aspects of Medicine Reviews. 1992; 3:1-124.
- Zawdydiwski R, Sprecher DL, Evelegh MJ, et coll. Cholesterol 1,2,3™: A novel test for the determination of skin cholesterol levels. Clinical Chemistry 2001; 47:1302-1304.
- Sprecher DL, et coll. Skin Tissue Cholesterol (Skin Tc) is related to angiographically-defined cardiovascular disease. Atherosclerosis 2003; 171:255-258.

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PreVu® est une marque déposée de Miraculins Inc. Tous droits réservés. 2012.