Dade® Ci-Trol® Coagulation Control Level 2

Intended Use

For use in citrated coagulation systems

Summary and Principle

The use of controls in the coagulation laboratory is an established procedure. Dade® Ci-Trol® Coagulation Control Level 2, is recommended as an abnormal control with patient citrated plasma for prothrombin time and activated partial thromboplastin time determinations.

This control provides a means of assuring the accuracy of prothrombin time test results in those ranges usually sought during anticoagulant therapy. In activated partial thromboplastin time testing, the control represents valuable and reliable data in the lower abnormal ranges as seen in first and severe second stage coagulation disorders.

Reagent

Materials provided

Dade® Ci-Trol® Coagulation Control Level 2

Pack for 20 x 1,0 ml, Code No. B4244-20

Composition

Dade® Ci-Trol® Coagulation Control Level 2: a lyophilized preparation of human plasma, stabilizers and buffers.

Preparation of the Control

Add exactly 1 mL of distilled water, gently tilt the vial and allow 15 minutes for reconstitution of the plasma in the **closed** vial. Before use, gently mix once more. **Do not shake.**

Note: Do not use distilled water containing preservative.

Storage and Stability

Store unopened vials at +2 to +8°C and use by date given on the label.

Stability after reconstitution:

6 hours at +2 to +8°C (closed vial).

Indication of deterioration: lack of vacuum upon opening vial or inability to obtain reproducible values.

Warnings and Precautions

1. For *in vitro* diagnostic use only.
2. Human source material. Treat as potentially infectious. **Each donor unit used in the preparation of this product has been tested by an FDA-approved method and found non-reactive for the presence of HBsAg, antibodies to HIV and HCV. Because no known test method can offer complete assurance that hepatitis B virus, hepatitis C virus, Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all human blood-based products should be handled in accordance with good laboratory practices using appropriate precautions.¹**

Procedure

After reconstitution, use Ci-Trol® Level 2, in the same manner as freshly drawn patient citrated plasma using the same test procedure. Ci-Trol® Level 2 can be used in the manual tilt tube method as well as most semi-automated and automated instrument techniques. See instrument operator's manual for detailed instructions and limitations of procedures.

Controls such as Dade® Ci-Trol® Coagulation Control Level 2, should be tested at the initiation of testing, upon reagent changes, and at least once each 8-hour shift. The control material should be run in the same manner as the test sample. Each laboratory should establish a range for control values. If control values are outside of determined range check controls, reagents, and instrument. It is recommended before reporting any patient data to document any actions taken to identify and correct the problem. New control ranges should be established for each lot of reagent or control material.

Limitations of the Procedure

Inability to obtain proper control values may be a sign of product deterioration. However, when proper control values are not obtained, a study of each component of the test system (i.e., reagents, control and technical conditions) is indicated to ascertain that all other components are functioning properly.

Expected Values

Prothrombin Time (PT)

Ci-Trol® Level 2 is specifically designed for use with Dade Behring PT reagents: Dade® Innovin®, Thromboplastin C Plus and Thromborel® S. Using Dade® Innovin®, the following performance-guidelines apply for Ci-Trol® Level 2:

Level 2

BCS® (US-setting)	23 – 31 seconds
Sysmex® CA-1500	23 – 31 seconds

Activated Partial Thromboplastin Time (APTT)

Ci-Trol® Level 2 is specifically designed for use with Dade Behring APTT reagents: Actin®, Actin® FS, Actin® FSL and Pathromtin® SL. With these reagents, the APTT values expected with Ci-Trol® Level 2 are in the intermediately elevated (abnormal) range.

Specific Performance Characteristics

Studies of Ci-Trol® Level 2 in normal clinical laboratory usage show an intralaboratory variation resulting in a total CV of approximately 3% for prothrombin times and approximately 4% for activated partial thromboplastin times.

Since laboratory control materials are used for the effective monitoring of the performance of a coagulation test, each laboratory should establish its own level of performance to monitor quality assurance.

Bibliography

1. U.S. Department of Health and Human Services CDC, Biosafety in Microbiological and Biomedical Laboratories, HHS Publication (CDC) 93-8395; 1999; Section II; 8-16

This product is warranted to perform as described in its labeling and in the product literature, and **Dade Behring disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event shall Dade Behring be liable for any consequential damages arising out of the aforesaid express warranty.**

Dade, Ci-Trol, Innovin and Actin are registered trademarks of Dade Behring, Inc. in the USA, Germany and other countries.

Thromborel and BCS are registered trademarks of Dade Behring Marburg GmbH in the USA, in Germany and other countries.

Sysmex is a registered trademark of SYSMEX CORPORATION in the USA, in Germany and other countries.

* Pathromtin is a registered trademark of Dade Behring Marburg GmbH in Germany and other countries.

Dade Behring Marburg GmbH
Emil-von-Behring-Str. 76
D-35041 Marburg
www.dadebehring.com

USA Distributor: Dade Behring Inc.
Newark, DE 19714 U.S.A.

Edition March 2003

Dade® Ci-Trol® Coagulation Control Level 2

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