Medtronic

BIO-PUMP[®]

BP-50, CBBP-50

Centrifugal Blood Pump

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Medtronic[®], Bio-Pump[®], and Bio-Console[®] are all trademarks of Medtronic, Inc.

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1 DESCRIPTION

The Medtronic® Bio-Pump® centrifugal blood pump with or without Carmeda®^{1, 2} BioActive surface has been designed to utilize the constrained forced-vortex pumping principle. In action, blood is allowed to pass through a vortex created by the rotation of a series of smooth-surfaced rotating cones. Energy is transferred from the cones to the blood in the form of pressure and velocity as the blood is gently accelerated toward the outlet of the pump. The Bio-Pump blood pump is fabricated from sturdy, thromboresistant, biocompatible acrylics which are mounted on precision bearings for dependable operation. The non-occlusive hemodynamic design of the pump promotes laminar flow, allowing for improved blood handling capabilities and decreasing trauma which may be associated with extracorporeal circulatory support during cardiopulmonary bypass.

The Bio-Pump blood pump couples to a safety-locked magnetic drive unit called the Medtronic® Bio-Console® pump speed controller, which serves as the hardware component of the Medtronic extracorporeal blood delivery system.

All Bio-Pump blood pumps are sterilized prior to shipment. Depending on the method of sterilization a slight amber cast may be present. This discoloration poses no risks to the patient.

The Bio-Pump blood pump with Carmeda® BioActive surface (CBBP-50) is coated with a non-leaching bioactive surface (heparin) which provides thromboresistant blood contact surfaces.

¹ Carmeda® is a registered trademark of Carmeda A.B.

² This product is coated with the Carmeda® BioActive Surface under license from Carmeda A.B. and is licensed for use only as part of an extracorporeal blood circulation system or circuit which includes an oxygenator or blood pump.

2 INDICATIONS FOR USE

The Medtronic Bio-Pump centrifugal blood pump is indicated for use only with the Medtronic Bio-Console pump speed controller to pump blood through the extracorporeal bypass circuit for extracorporeal circulatory support for periods appropriate to cardiopulmonary bypass procedures (up to six hours).

3 CONTRAINDICATIONS

The Medtronic blood pumping system is contraindicated as a cardiotomy suction device.

This device used for any other purposes than the intended use is the responsibility of the user.

4 WARNINGS

- The pump has not been qualified through *in-vitro*, *in-vivo*, or clinical studies for long-term use (longer than six hours) as a bridge-to-transplant or for pending recovery of the natural heart.
- Read all warnings, precautions and instructions for use carefully prior to use. Failure to read and follow instructions or failure to observe the stated warnings could cause serious injury or death to the patient.
- Possible side effects include, but are not limited to, infections, mechanical failure, hemolysis, and thromboembolic phenomena. These are potential side effects with all extracorporeal blood systems.
- Ensure that the pump and circuit have been debubbled and primed properly prior to beginning bypass to minimize the risk of air reaching the patient. The use of an arterial filter is recommended.
- Massive air entry into the pump will cause the pump to deprime and blood flow to stop. Stop the pump and remove air prior to resuming circulation.

- Do not use alcohol or alcohol based solutions on any surface of the pump.
- To prevent backflow of the patient's blood when the pump outlet tubing is open, establish and maintain a minimum pump speed that overcomes line and patient resistance. Failure to do this could allow retrograde flow and exsanguinate the patient.
- Arterial lines must be clamped when the pump is at rest to prevent retrograde flow.
- The Bio-Pump blood pump is designed to be operated only with the Bio-Console pump speed controller. There is no safety or performance data known to Medtronic which establishes compatibility of any other manufacturer's devices or components to the Medtronic system. Run the Bio-Pump blood pump only on a properly maintained Bio-Console pump speed controller.
- Potential risk to the patient should be evaluated prior to changing the pump.
- Frequent patient and device monitoring is recommended; do not leave device unattended while in use. Monitor carefully for signs of occlusion throughout the circuit.
- Do not operate the pump for more than thirty (30) seconds in the absence of flow. The temperature within the pump will rise and increased cellular damage may result.
- The pump must be handled and used in an aseptic manner.
- Do not operate the Bio-Pump blood pump with its inlet clamped as a negative pressure would be generated in the pump and air bubbles may be formed in the blood.

5 PRECAUTIONS

USA Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

- This device should be used only by persons thoroughly trained in extracorporeal circulation procedures.

- Do not hit or strike the Bio-Pump blood pump with hands or instruments. Shock may cause damage to the device, which may cause device malfunction.
- The Bio-Pump blood pump is sterile and nonpyrogenic in an unopened and undamaged unit package. Inspect device and package carefully prior to use. Do not use if the unit package or the product has been damaged or soiled, or if caps are not in place.
- Each pump is intended for single use only. Do not resterilize. Resterilization by any means may cause severe damage to the device or its components. Dispose of safely after single use to avoid risk of infection.
- Attach tubing in such a manner as to prevent kinks or any restrictions that may alter flow.
- Do not run the pump unprimed as damage to the internal components may occur.
- Do not use pump if dropped. Dropping or other severe shock may cause damage which could lead to device malfunction.
- Do not use excessive force to install tubing on pump as damage to pump may occur.
- Take care to prevent damage to connectors when repositioning or attaching tubing.
- Always have a spare Bio-Pump blood pump and backup equipment available with appropriate protocol for change out.
- Do not place the pump by items adversely affected by magnetic fields.
- A strict anticoagulation protocol should be followed and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician.

6 RECOMMENDED PRIMING PROCEDURES

1. With accepted aseptic technique, attach the appropriate tubing to the inlet and outlet of the Bio-Pump blood pump. (Applying a sterile saline solution to the inlet/outlet connectors may aid tubing attachment.)
2. Where practical, flush the circuit and the pump with CO₂ (from filter or reservoir bag).
3. Fill the Bio-Pump blood pump with prime, by gravity, to a point beyond the flow probe and clamp. Eliminate air from the device by “walking” the air out the outlet.

Caution: Do not hit or strike the Bio-Pump blood pump with hands or instruments. Shock may cause damage to the device, which may cause device malfunction.

Warning: Ensure that the pump and circuit has been debubbled and primed properly prior to beginning bypass to minimize the risk of air reaching the patient. The use of an arterial filter is recommended.

Warning: Massive air entry into the pump will cause the pump to deprime and blood flow to stop. Stop the pump and remove air prior to resuming circulation.

4. With the outlet clamped, turn the Bio-Console pump speed controller power “ON.”
5. Zero (gain and balance) the flow probe according to the Bio-Console Operator and Reference Manual.
6. With the outlet still clamped, turn the RPMs of the Bio-Console pump speed controller to maximum RPM for 30 seconds. This delivers maximum pressure within the pump which aids in assuring its integrity. Observe the pump for leakage or other anomalies.

Warning: Do not operate the pump for more than thirty (30) seconds in the absence of flow. The temperature within the pump will rise and increased cellular damage may result.

7. Bring the RPMs to zero and recheck the integrity of the pump.

Warning: If leaks or other anomalies are found, remove the Bio-Pump blood pump and replace with a new, sterile pump, repeating the above steps to prime.

8. If no anomalies are noted, continue the priming of your circuit in the usual fashion.
9. Verify all connections and the integrity and flow of your circuit prior to use.

Warning: Do not operate the Bio-Pump blood pump with its inlet clamped as a negative pressure would be generated in the pump and air bubbles may be formed in the blood.

7 PRE-BYPASS CHECKLIST

Information from this blood pump bypass checklist can be used to revise, upgrade, or expand existing checklists appropriately.

1. EQUIPMENT ASSEMBLY

- Mount motor drive correctly, relative to venous reservoir.
- Check that all electrical connections are secure.
- Test control module power and display.
- Check date and integrity of sterile blood pump (and disposable probe) package(s).
- Check that flow transducer / sensor (and disposable probe) are sized properly.
- Assemble perfusion circuit in a sterile manner.
- Allow sufficient tubing length for standby pumping unit.
- Connect flow transducer / sensor (and disposable probe) to circuit in correct location and flow direction.

2. PRIME PUMP AND CIRCUIT

- CO₂ flush pump and circuit, if indicated; turn off CO₂.
- Gravity-prime and debubble pump and perfusion circuit.
- Check pump for leaks, irregular motion, and noise.
- Check circuit for visible air.

- Check that all tubing connections are secure.
- Clamp pump outlet line completely.
- Clamp venous return line completely.

3. OPERATING PARAMETERS

- Calibrate transducers / sensors according to the manufacturer's instructions.
- Set low / high flow alarms.
- Verify flow alarms.

4. EMERGENCY BACKUP EQUIPMENT

- Backup power available.
- Hand crank available.
- Pump (and disposable probe) available.
- Control module or roller pump available.

5. PERFUSION

- Achieve minimum pump speed prior to unclamping lines.
- Monitor control module for messages and alarms.
- Monitor perfusion circuit for visible air.
- Maintain minimum pump speed prior to clamping lines.

6. CLEANUP

- Turn off power.
- Properly discard disposable components.
- Cover magnetic motor drive.
- Clean console and transducer / sensor.

7. CHECK EQUIPMENT

- Inspect and verify that equipment is operational.
- Maintain indicated preventative maintenance schedule.
- Recharge batteries to full capacity (if necessary).

8 INSTRUCTIONS FOR PRODUCTS WITH CARMEDA® BIOACTIVE SURFACE

The blood-contact surface of the product (CBBP-50) is coated with Carmeda® BioActive Surface (heparin) that provides thromboresistant blood-contact surfaces.

Warning: A product coated with a Carmeda® BioActive Surface is intended for single use only. Reesterilization may adversely affect the Carmeda® BioActive Surface.

A strict anticoagulation protocol should be followed and anticoagulation should be routinely monitored during all procedures. The Carmeda® BioActive Surface provides a nonleaching, nonthrombogenic coating on all blood-contact surfaces. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician.

Caution: Do not store a Carmeda®-coated product above 40°C.

9 SPECIFICATIONS

Priming Volume (approximate)	48 mL	Maximum flow rate	1.5 L/min.
Inlet/Outlet I.D.	6.35 mm (1/4 inch)	Maximum pump speed	4,500 rpm
Maximum Operating Pressure	900 mm Hg (120 kPa)	Compatible consoles	All Medtronic Bio-Consoles

9.1 STORAGE ENVIRONMENT

Storage temperature:

- BP-50: -40° to 50°C (-40° to 122°F)
- CBBP-50: -40° to 40°C (-40° to 105°F)

Storage humidity: 15% to 95% noncondensing.

10 LIMITED WARRANTY*

(For Single Use Product)

The following LIMITED WARRANTY Applies to United States Customers Only:

- A. This LIMITED WARRANTY provides the following assurance to the patient who receives the Medtronic® Bio-Pump® Centrifugal Blood Pump, Models BP-50 and CBBP-50, hereafter referred to as the “Product”:
- (1) Should the Product fail to function within normal tolerances due to a defect in materials or workmanship prior to its “Use By” date, Medtronic will at its option:
 - (a) issue a credit to the purchaser equal to the Purchase Price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.
 - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.
- B. To qualify for the LIMITED WARRANTY, these conditions must be met:
- (1) The Product must be used prior to its “Use By” date.
 - (2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.
 - (3) The Product must not have been altered or subjected to misuse, abuse or accident.
 - (4) The Product must be used in accordance with the labeling and instructions for use provided with the Product.

* This LIMITED WARRANTY is provided by Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432-5604. It applies only in the United States. Areas outside of the United States should contact their local Medtronic representative for exact terms of the LIMITED WARRANTY.

- C. This LIMITED WARRANTY is limited to its express terms. In particular:
- (1) Except as expressly provided by this LIMITED WARRANTY, MEDTRONIC IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.
 - (2) This LIMITED WARRANTY is made only to the patient in whom the product was used. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOM OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.
 - (3) The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of The LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid. This LIMITED WARRANTY gives the patient specific legal rights. The patient may also have other rights which vary from state to state.

- (4) No person has any authority to bind Medtronic to any representation, condition or warranty except this LIMITED WARRANTY.

11 LIMITED WARRANTY*

(For Single Use Product / Outside the United States)

The following LIMITED WARRANTY Applies to Customers Outside the United States:

- A. This LIMITED WARRANTY provides assurance for the purchaser who receives a Medtronic® Bio-Pump® Centrifugal Blood Pump, Models BP-50 and CBBP-50, hereafter referred to as “Product,” that should the product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient. THE WARNINGS CONTAINED IN THE PRODUCT LABELING ARE CONSIDERED AN INTEGRAL PART OF THIS LIMITED WARRANTY. CONTACT YOUR LOCAL MEDTRONIC REPRESENTATIVE TO OBTAIN INFORMATION ON HOW TO PROCESS A CLAIM UNDER THIS LIMITED WARRANTY.
- B. To qualify for the LIMITED WARRANTY, these conditions must be met:
- (1) The Product must be used prior to its “Use By” date.
 - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
 - (3) The Product may not have been used for any other patient.

* This LIMITED WARRANTY is provided by Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432-5604. It applies only outside the United States.

- C. This LIMITED WARRANTY is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation or material alteration of the replaced Product.
 - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

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All Rights Reserved
© Medtronic, Inc. 2006