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© 2008 Handicare

All rights reserved.

The information provided herein may not be reproduced and/or published in any form, by print, photoprint, microfilm or any other means whatsoever (electronically or mechanically) without the prior written authorization of Handicare.

The information provided is based on general data concerning the constructions known at the time of the publication of this manual. Handicare executes a policy of continuous improvement and reserves the right to changes and modifications.

The information provided is valid for the product in its standard version. Handicare cannot be held liable for possible damage resulting from specifications of the product deviating from the standard configuration.

The available information has been prepared with all possible diligence, but Handicare cannot be held liable for possible errors in the information or the consequences thereof. Handicare accepts no liability for loss resulting from work executed by third parties.

Names, trade names, etc. used by Handicare may not, as per the legislation concerning the protection of trade names, be considered as being available.

2008-01

Introduction	24
This manual	24
Product reference	24
Service and technical support	24
Identification plate	24
Approval	25
CE declaration of conformity	25
Warranty	25
Liability	27
Used wheelchairs and the environment	27
1 Safety	28
1.1 General safety instructions	28
1.2 Decals and instructions on the wheelchair	28
2 General description	29
2.1 Product description Ibis X-Series	29
2.2 Main components	29
3 Operating instructions	30
3.1 Inspection before use	30
3.2 Getting in and out	30
3.3 Parking	30
3.4 Driving	31
3.5 Obstacles	31
3.6 Slopes	32
3.7 Reducing the wheelchair	32
3.8 Transport	33
3.9 Adjustment options (dealer)	34
3.10 Other adjustments (user)	36
4 Maintenance	38
4.1 Maintenance table	38
4.2 Checking and inflating the tyres	38
4.3 Cleaning the wheelchair	38
5 Technical specifications	39


Introduction

This manual

With this manual you can safely use and maintain (clean) the wheelchair. This manual should be used in combination with Handicare's Sedeo seating system manual (supplement). In case of doubt always contact your dealer.

Strictly follow the procedures with a caution symbol! These procedures can if they are not executed properly - result in human injury or damage to the product or the environment.


Product reference

- A Seat tilt adjustment
- B Backrest recline with gas spring (option)
- C Transport brackets
- D Identification plate
- E Tyre pressure rear wheels

Service and technical support

For information concerning specific settings, maintenance or repair work, please contact your dealer. Make sure you always mention the model, year of manufacture and identification number. This information is provided on the identification plate of the product. A service manual is available for the dealer.

Identification plate

- A Model
- B Year of manufacture
- C Identification number
- D User area indoors or outdoors
- E Maximum load in kg

Approval

The product meets the following standards:

1. EN12183 Hand-operated wheelchairs.
2. ISO7176-8 Requirements for impact, static and fatigue strength.
3. ISO7176-9 Climate tests for wheelchairs and scooters.
4. ISO7176-16 Requirements of resistance to fire.
5. ISO 7176-19 Crash test.

CE declaration of conformity

The product is in conformity with the provisions of the Medical Devices Directive and thus has CE marking.

Warranty

In the following warranty and liability stipulations the terms and definitions as summed up hereafter must be explained as follows:

- **Product:** The hand-operated or electric wheelchair or scooter manufactured and delivered by Handicare.
- **Customer:** The person who directly obtains a Product from Handicare.
- **Dealer:** The person who delivers a Product obtained from Handicare to customers or third parties.
- **User:** The person who uses a Product manufactured by Handicare.

Irrespective of what is determined concerning warranty conditions applicable to the Product, in any case the following applies with regard to the warranty:

1. Except insofar as described otherwise Handicare guarantees the Product for its suitability for the purpose for which the Product is

intended - all of these points as described in this manual - and for the quality of the material of which the Product is made and the manner in which the Product is manufactured.

2. Repairs or replacement of parts of the Product that may be necessary as a result of faults that are based on qualitatively faulty material or manufacturing errors will be executed free of charge, as long as such faults occurred within one (1) year after the date of delivery of the Product to the Customer. The parts to be replaced must be shipped post-paid to Handicare. Disassembly or assembly of these parts is for the expense of the Customer. Therefore the following cases are not eligible for free repair or replacement as meant above:

- The repair or replacement that is necessary in connection with faults that arise after one (1) year after the date of delivery of the Product to the Customer;
- The repair or replacement that is required in connection with faults due to improper or careless use of the Product or that are based on the Product being used for another purpose than the one for which it is intended, in which regard it applies that if the Customer is a Dealer, this Dealer will indemnify Handicare against possible claims from Users or other third parties for faults based on an incorrect or careless use of the Product;
- Parts that are subject to wear, and the need for repair or replacement of the parts are the actual consequence of normal wear.

3. Irrespective of that stipulated under 2, it applies as far as an electric Product is concerned that with regard to the battery that forms part of the Product warranty is only given in case of faults or non-functioning of the battery that are demonstrably the direct consequence of material or manufacturing errors. A fault or non-functioning of the battery as a result of normal wear is not covered by the warranty as meant in these warranty stipulations. Similarly not covered by the warranty are faults or non-functioning that are the consequence of improper or unprofessional use of the Product or the battery that is part of the Product, including the incorrect charging of the battery and the failure to carry out timely and good maintenance, in which context it also applies that if the Customer is a Dealer, this Dealer will indemnify Handicare against possible claims from Users or other third parties that are based on the above mentioned improper or unprofessional use of the Product or the battery that is part of the Product.
4. The warranty conditions as mentioned in above stipulations in any case become null and void if:
 - The guidelines of Handicare for the maintenance of the Product are not, or insufficiently, followed;
 - A necessary repair or replacement of parts is based on neglect, damage or abuse of the Product or a use of the Product for another purpose than the one for which it was intended;
 - Parts of the Product are replaced by parts of another origin than those which Handicare uses and/or parts of the Product are replaced without the permission of Handicare.
1. The warranties as mentioned in stipulations 1 through 3 also become null and void if what is involved is re-used by a new user within the warranty period and that such re-use necessitates adaptations to the product and those adaptations are carried out without the instructions and/or at the order of Handicare.
2. To retain rights under the above delineated warranties the Customer must, in case of damage or other calamities, as rapidly as possible contact Handicare and inform them as fully as possible. The possibility of taking recourse to the above-mentioned warranty conditions becomes null and void for the Customer in any case after 20 workdays after the claim situation or the calamity arises that was the reason for the recourse to the guarantee.
3. The replacement of a part or the repair or the reconditioning of the Product within a current warranty period does not extend the warranty period.
4. Handicare gives no warranty on repair to or reconditioning of the Product carried out other than under order of and/or at the instructions of Handicare. If repairs and/or reconditioning are executed by or on behalf of a Customer, the Customer indemnifies Handicare with respect to the claims of third parties who result, in the broadest sense of the word, from such repairs or reconditioning.

Liability

Irrespective of what is determined regarding liability in the general conditions applicable to the Product, with regard to liability in any case the following applies:

1. Taking into consideration the following stipulations, Handicare only accepts liability for loss due to death or physical injury that is the result of a defect in the Product for which Handicare is responsible and for damage to another object that is the private property of the user of the Product, as long as said loss is the direct result of a fault in the Product.
2. Handicare accepts no other or further liability than delineated under 1. In particular Handicare accepts no liability for consequential damage, in any form whatsoever.

Used wheelchairs and the environment

If your wheelchair has become superfluous or needs to be replaced, it can usually be taken back by your dealer. If this is not possible, please contact your local authorities for the possibilities of recycling or an environmentally friendly way of disposing of the used materials. For the production of the wheelchair several plastics and metals have been used.

1 Safety

Handicare accepts no liability whatsoever for damage or injury caused by the failure to (strictly) adhere to the safety guidelines and instructions or else due to carelessness during the use and cleaning of the wheelchair and possible accessories. Depending on the specific work circumstances or the accessories used, additional safety instructions may be required. Please contact your dealer immediately if you observe a potential danger during the use of the product.

1.1 General safety instructions

- The user of the wheelchair is at all times fully responsible for the fulfilment of the locally applicable safety regulations and guidelines.
- The technical specifications must not be altered.
- Modifications to (parts of) the wheelchair are not permitted.

1.2 Decals and instructions on the wheelchair

- Signs, symbols and instructions placed on this wheelchair comprise part of its safety facilities. They must never be covered or removed and must remain present and clearly legible throughout the entire lifespan of the wheelchair.
- Immediately replace or repair illegible or damaged signs, symbols and instructions. Contact your dealer with regard to this.

2 General description

2.1 Product description Ibis X-Series

The Ibis X-Series is a range of wheelchairs designed for users who need excellent seating posture and maximum comfort. Every model in the Ibis X-Series range has tilt-in-space adjustment, allowing the user to rest in the wheelchair.

The Ibis X-Series provides a range that extends from a push wheelchair (Ibis XC) and manual or foot-propelling wheelchair (Ibis XA) to an electric model (Ibis XP) for use in and around the house. As a result, a wheelchair offering maximum comfort and safety is available for every situation.

Depending on the model, the user will be able to push the wheelchair or an attendant will be required to push and operate the wheelchair.

The frame height of the wheelchair can easily be adjusted, ensuring a ideal seat height for everyone. Additionally, the seating system can be fully customized to satisfy the specific wishes and requirements of the individual.


The Ibis X-Series has been designed for:

- Transport of persons to a maximum weight of 160 kg.
- Use on paved roads, sidewalks, foot- and bicycle paths.
- Use in and around the house.

If you use the wheelchair for other purposes than those for which it is intended, Handicare accepts no liability whatsoever for damage or injury resulting from such use other than for which the wheelchair was developed and designed.


Please study this manual carefully before you use the wheelchair.

2.2 Main components

- A Seat
- B Backrest
- C Seat tilt adjustment
- D Armrest
- E Legrest
- F Push bracket
- G Brake
- H Castor wheels
- I Rear wheels
- J Anti-tip brackets
- K Headrest (option)
- L Lap belt, to be provided by the dealer (option)

3 Operating instructions

3.1 Inspection before use

Check if the tyres are sufficiently inflated before you start using the wheelchair.

3.2 Getting in and out

Transfer forwards

Apply the wheelchair brake.
Fold up the foot support (A)
Seat yourself in the chair.

Transfer sideways

1. Apply the wheelchair brake.
2. Loosen the cap (B) and remove the armrest (C).
3. Seat yourself in the chair.

3.3 Parking

Ibis XC




Apply the brake to the wheelchair by pushing the footbrake (A) downward or pulling the handle (B) backward.

Ibis XA

Apply the brake to the wheelchair by pushing the handle (B) downward.

- Make sure the brake handle is pushed all the way backward while driving the wheelchair, so that the brake cannot be applied accidentally.
- If the wheelchair is stopped on a slope, the brake should be applied to prevent the wheelchair from moving down the slope.

3.4 Driving

- As a wheelchair user, you could be vulnerable in traffic. When travelling in or across a road, keep in mind that other traffic participants may not always notice you. Keep to the applicable traffic regulations.
- Avoid out-of-the way routes, so that help can be on its way quickly when needed.

Driving the Ibis XC

1. Adjust the height of the push handles.
2. Release the brake from the rear wheels.
3. The wheelchair can now be pushed.

Driving the Ibis XA

1. Release the brake from the rear wheels.
2. Take the hoops of the rear wheels and push the hoops forward or backward with your hands.
3. The wheelchair can also be moved by foot-propelling.

Braking with the Ibis XA

Tightly grasp both hoops.

By braking this way, friction is caused between your hands and the hoops, which will create heat. If you brake too fast, your hands may become very hot.

Miscellaneous

- Drive carefully on slippery roads, resulting from rain, ice or snow!
- Prevent the wheelchair from coming into contact with seawater: seawater is caustic and may damage the wheelchair.
- Prevent the scooter from coming into contact with sand: sand may affect the moving parts of the wheelchair, resulting in unnecessarily rapid wear.
- Do not put the wheelchair in direct sunlight: the upholstery may become too hot, which may cause burning or an allergic reaction of the skin.
- Do not drive the wheelchair when you are under the influence of drugs, alcohol or medication that could affect your ability to drive.
- You must have sufficient eye sight to be able to safely drive the wheelchair.
- Be sure that no items of clothing hang loose. These could be caught between the wheels.

3.5 Obstacles

The method described below is to be used when the wheelchair is pushed by an attendant and will make the user feel safe.

Driving up obstacles

1. Push the kick-tube downward and pull the wheelchair slightly backward, so that the castor wheels are released from the ground.
2. Push the wheelchair onto the elevation and release the kick-tube.

Operating instructions

3. Push the wheelchair further and lift the rear wheels over the obstacle using the push bracket.

- Always descend from a sidewalk back to front, so that the wheelchair cannot tip over and make the user fall out of the chair.
- Do not drive down stairs.

3.6 Slopes

Driving up slopes


- Slopes with an angle of more than 5° must not be taken.
- If you move the wheelchair yourself, always make sure that you hold on to the hoops, so the wheelchair cannot start to move by itself.
- The wheelchair is less stable on a slope.
- Do not drive on slopes when the seat is tilted.
- The attendant pushing the wheelchair must continue to push to prevent the wheelchair from moving downhill on its own accord.
- Apply the brake of the wheelchair if it is not pushed or not pushed hard enough.

Driving down slopes

- The wheelchair will start to move by itself, caused by gravity.
- If you push the wheelchair yourself, use both hands to control the speed with the hoops.


- The attendant pushing the chair must hold on to the push handles to apply some brake force.
- Prevent the wheelchair from gaining too much speed.

3.7 Reducing the wheelchair

Removing the seat

1. Pull up the seat (A) in the front.
2. When replacing the seat, the clips (B) must click around the frame (C) of the chair.


010407010

Check if the seat is clicked in place securely.

Folding down the backrest

1. Flip the lock clip (A) off the adjusting bar.
2. Pull the lock clip from the hinged joint.
3. Fold down the backrest.

Removing the armrest

Loosen the cap and remove the armrest.

Folding up the footplate


Fold up the footplate.

3.8 Transport

The wheelchair can be placed in a car. The best way to do this is by reducing the wheelchair or using steel ramps.

Handicare recommends wheelchair users to travel in a regular car seat whenever possible. A wheelchair cannot offer the level of safety that is offered by standard car seats, no

matter how well the wheelchair is fixed to the vehicle concerned. As it is not possible for every wheelchair user to make a transfer, the wheelchair should be secured as described below.

- The Ibis X-Series has been crash tested and has met all the pass requirements of ISO 7176-19.
- Your safety during transport largely depends on the person securing the wheelchair. It is your responsibility to ensure the wheelchair is secured as safely as possible.
- Always ask for confirmation of the transporter that the vehicle is suitably designed, insured and equipped to transport a person in a wheelchair.
- Parts of the wheelchair that can be easily detached should be removed and stored elsewhere during transportation.
- The wheelchair should always be transported in a forward facing direction.
- Handicare recommends an Unwin tie-down system that is tested and approved according to ISO10542.

Operating instructions

- The wheelchair should be secured only at the tie-down points on the frame, not onto any accessories (armrest, anti-tip brackets, etc).
- Apply the wheelchair brake.
- Handicare recommends that the wheelchair be supplied with a headrest and that this should be used during transportation.
- During transportation it is essential that you are secured by a three point lap belt, which is anchored to the vehicle wall and floor rails. A lap belt in your wheelchair must not be considered as an alternative to the occupant restraint system.


3.9 Adjustment options (dealer)

The seat tilt range can vary depending on the way the wheelchair is adjusted.

Seat height


Before adjusting the seat height, remove as many heavy seat parts as possible (seat, backrest, etc.) to limit the weight to be lifted.


1. Loosen the 4 bolts (A) using a 13 mm spanner, and remove them.
2. Lift the seat post off the bracket (B) and place it back in the desired position on the bracket.
3. Tighten the bolts (A) securely.

Seat angle


1. Tilt the seat with the seat tilt adjustment, so that the adjusting plate becomes free.
2. Loosen the bolt (A) with a 4 mm Allen key, and remove it.
3. Loosen the bolt (B) a few turns.
4. Turn the plate © in the desired position (angle).
5. Tighten both bolts (A and B) securely.

Frame width Ibis XA

1. Loosen the bolts (A) with a 13 mm spanner and remove them.
2. Put the frame tube (B) in the desired position, so that the holes are level with the frame.
3. Tighten the bolts (A) securely.

Position of the wheels Ibis XC

1. Put the frame on a solid base, so that the wheels are off the ground.
2. Remove the manual control (A).
3. Remove the foot bracket (B).
4. Loosen the nut (C) with a 19mm spanner.
5. Move the wheel into the desired position.


6. Tighten the nut (C) securely.

- Changing the position of the wheels may affect the stability of the wheelchair. Check the stability of the chair and make sure the anti-tip brackets are well adjusted.
- Check the brake.

Position of the wheels Ibis XA

1. Put the frame on a solid base, so that the wheels are off the ground.
2. Loosen the bolts (A) with a 19 mm spanner and remove them.
3. Move the wheel plate (B) into the desired position.
4. Tighten the bolts (A) securely.


- Changing the position of the wheels may affect the stability of the wheelchair. Check the stability of the chair and make sure the anti-tip brackets are well adjusted.
- Check the brake.

Anti-tip brackets

1. Loosen the bolts (B).
2. Put the anti-tip brackets in the desired position.
3. Tighten the bolts (B) securely.

Check if the anti-tip brackets work properly.

Wheel size Ibis XA

For the various wheel sizes, specific holes are available.


Readjust the brake, so that the brake power is sufficient.

Brake



1. Loosen the Allen bolt (D) a few turns with a 15 mm Allen key.
2. Put the brake in the desired position.
3. Tighten the Allen bolt (D) securely.

Check the brake.

3.10 Other adjustments (user)

Seat tilt adjustment

1. Apply the wheelchair brake.
2. Grasp the push handles of the push bracket and squeeze the handle (A).
3. Push the push handles downward until the desired angle is reached.

Backrest recline with gas spring

1. Apply the wheelchair brake.
2. Grasp both handles of the push bracket and squeeze the handle (B).
3. Push the push handles downward until the desired angle is reached.

Push bracket height

1. Loosen both wing nuts ©.
2. Push the push bracket (D) in the desired position.

The seat tilt adjustment cable should not be too tight.

4 Maintenance

4.1 Maintenance table

It is recommended to have your wheelchair serviced by your dealer once a year or, in case of intensive use, every six months. The table below shows the maintenance activities you can perform yourself.

Frequency	Description
Weekly	Check tyre pressure and inflate the tyres if necessary. See 4.2
Monthly	Clean the wheelchair. See 4.3.
	Lubricate the turning system of the seat.

4.2 Checking and inflating the tyres

- **Make sure the tyres are kept at the correct pressure.**
- **The tyre pressure may never exceed the maximum value indicated on the tyre.**

Inflate the tyres with a foot pump or a bicycle pump, in combination with the reducing nipple included in the delivery. You can also have your tyres inflated at your local petrol station.

- Tyre pressure 12" rear wheel: 2.5 bar

4.3 Cleaning the wheelchair

Removing dry dirt

Upholstery, metal parts and frame parts can usually easily be cleaned with a dry soft cloth.

Removing mud and other wet dirt

The dirty parts can best be cleaned by first wiping them with a wet sponge and then wiping them dry with a dry soft cloth.

- **Never use abrasive or aggressive cleansers. They can scratch the wheelchair.**
- **Do not use organic solvents such as thinner, benzene or white spirit.**

5 Technical specifications

Description	
Model	Ibis XC / XA
Max. user weight	160 kg
Total length	1040 mm
Total width Ibis XC	625 mm
Total width Ibis XA	660-700-740 mm
Total weight	40 kg
Total weight of the heaviest part Ibis XC	26 kg
Total weight of the heaviest part Ibis XA	28 kg
Static stability downward	15°
Static stability upward	15°
Static stability sideways	15°
Seat tilt adjustment	30°
Seat angle	0 - 6°
Effective seat depth	440 - 520 mm
Effective seat width Ibis XC	380 - 555 mm
Effective seat width Ibis XA	380 - 555 mm
Seat height front	390 - 540 mm
Backrest angle	89 - 120°
Backrest height	520 - 570 mm
Lower leg length	390 - 500 mm
Armrest height	180 - 280 mm
Front armrest to backrest	370 - 450 mm
Minimum turning radius	810 mm
Ground clearance	70 mm
Push bracket height	830 - 1200 mm

Wheels	
Diameter front wheel	7"
Diameter rear wheels Ibis XC	12" / 16"
Diameter rear wheels Ibis XA	20" / 22" / 24"
Tyre pressure rear wheels	12": 2.5 bar 20" / 22" / 24": 7.5 bar

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Handicare B.V.