Automated Impella® Controller

**USER MANUAL
GEBRAUCHSANWEISUNG
MANUEL D'UTILISATION
MANUAL DEL USUARIO
MANUALE D'USO
ANVÄNDARHANDBOK
BRUGERVEJLEDNING**

ABIOMED®
Recovering hearts. Saving lives.™

TABLE OF CONTENTS

English..... 1

Deutsch 21

Français..... 41

Español..... 61

Italiano 81

Svenska..... 101

Dansk 121

INTENDED AUDIENCE

This manual is intended for use by cardiologists, surgeons, nurses, perfusionists, and cardiac catheterization lab technicians who have been trained in the use of the Impella® Circulatory Support System.

WARNINGS

The Impella® System is intended for use only by personnel trained in accordance with the Abiomed® Training Program.

During defibrillation, do **NOT** touch the Impella® Catheter, cables, or Impella® Console.

Power the Automated Impella® Controller using its internal battery if the integrity of the protective earth conductor is questionable.

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the electromagnetic compatibility (EMC) information provided in this manual.

During transport, the Automated Impella® Controller may be exposed to stronger electromagnetic disturbance than during in-hospital use. Strong electromagnetic disturbance may cause the Automated Impella® Controller to display soft button menu selections that were not selected by the user. Operators should be aware that, under these conditions, the operating parameters are not affected. No user intervention is required. Monitor Impella® Catheter flow and patient hemodynamics to confirm normal operation. The condition will resolve itself once the Automated Impella® Controller is no longer exposed to the disturbance.

Portable and mobile RF communications equipment can affect medical electrical equipment.

The equipment or system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

Use of cables, other than those sold by Abiomed, may result in increased emissions or decreased immunity of the Automated Impella® Controller.

The Automated Impella® Controller uses RFID (radio frequency identification) to identify and communicate with the purge cassette. Other equipment may interfere with the Automated Impella® Controller even if that other equipment complies with CISPR emission requirements.

CAUTIONS

To prevent malfunction of the Automated Impella® Controller, avoid long-term exposure to direct sunlight and excessive heat (40 °C).

To prevent overheating and improper operation, do NOT block the cooling vents of the Automated Impella® Controller while it is operating.

The Li-Ion batteries must be charged for 5 hours prior to system operation in order to meet the runtime requirement of 1 hour. Failure to do so will yield a shorter runtime. After being unplugged, the Automated Impella® Controller will operate for at least 60 minutes after the batteries have been fully charged.

Minimize exposure of Impella® System components to sources of electromagnetic interference (EMI). Exposure to sources of EMI, such as cell phones and two-way radios, may cause operational interference. To clear interference, either increase the distance between system components and the EMI source or turn off the EMI source.

Operation of Impella® System components may interfere with the operation of other devices. If interference occurs, increase the distance between the device and system components.

Have a backup Impella® Controller, purge cassette, connector cable, and Impella® Catheter available in the unlikely event of a device failure.

Do **NOT** use the bed mount as a handle.

INTENDED USE (EU AND CANADA)

The Automated Impella Controller is intended for exclusive use of Impella catheters and accessories. The controller may only be connected with devices which are listed in these instructions.

CONTRAINDICATIONS (EU AND CANADA)

There are no contraindications for use of the Automated Impella Controller, which operates the Impella pumps. Please see the respective catheter user manuals for indications, contraindications and possible complications.

OVERVIEW

REUSABLE SYSTEM COMPONENTS

The Impella® System consists of the following reusable components:

- Automated Impella® Controller—The controller provides an interface for monitoring and controlling the function of the Impella® Catheters. The controller provides a fluid purge to the Impella® Catheters. The controller provides backup power when the Impella® System is operated away from AC power. The controller weighs 26 lbs (11.8 kg) and can operate on its internal battery for at least 60 minutes when fully charged.
- Automated Impella® Controller Cart—The Automated Impella® Controller cart holds the controller. The cart has wheels for easy transport of the controller and a storage basket.


SINGLE-USE SYSTEM COMPONENTS

The Impella® System also includes the following single-use components:

- Impella® Catheter—microaxial blood pump which is operated by the Automated Impella Controller. There are different types of left ventricular Impella pumps. Due to a programmable memory chip integrated in the red plug, an Automated Impella Controller can recognize the type of connected Impella catheter and can operate it as intended.
- Purge Cassette—The purge cassette delivers rinsing fluid to the Impella® Catheters. The purge fluid (typically 20% dextrose solution) flows from the purge cassette through the catheter to the microaxial blood pump to prevent blood from entering the motor.
- Connector cable—The white connector cable connects the Impella® Catheter to the Automated Impella® Controller. Clips on the cable are used to secure the purge tubing to the cable.

AUTOMATED IMPELLA® CONTROLLER FEATURES



IMPORTANT NOTE: The underside of the Automated Impella® Controller has a battery switch to turn on the batteries. This switch is turned off for shipping purposes. Before operating the Automated Impella® Controller for the first time, make sure you turn this switch on. If the battery switch is not turned on, the Automated Impella® Controller will not be able to operate on battery power.

Feature	Description
Display screen	Displays user information, including the labels for the soft buttons. (Display screen elements described in detail later in this section.)
Soft buttons	Display, open, and close menus. The function for each soft button is defined by labels adjacent to the button on the display screen; function changes depending on the screen. When the Impella® Catheter is running, the default soft button labels are as follows: • MUTE ALARM • P-LEVEL • DISPLAY • PURGE SYSTEM • MENU
Power indicator	LED light above the selector knob; indicates the power status of the Automated Impella® Controller. • Green light—controller is on and plugged into AC power or running on battery power • Amber light—controller is off but plugged into AC power • No light—controller is off and not plugged into AC power
Selector knob	Rotating push button; turn clockwise and counterclockwise to navigate through menu items; push to make a selection.
Purge pressure transmitter	A flexible diaphragm on the purge cassette tubing that applies pressure to the sensor in the controller so that purge pressure can be measured.
Catheter plug	Connection point on the controller for the connector cable that connects to the Impella® Catheter.
Purge cassette	Contains the components for delivering the purge fluid; maintains the pressure barrier between the blood and the motor to prevent blood from entering the motor. (The purge cassette and its components are described in section 3 of this manual.)
Purge cassette door	Spring-loaded door that opens to provide access to the purge cassette.

Selector Knob Function


Rotate the selector knob on the controller to navigate through menu items. Push the selector knob to confirm your selection.

Feature	Description
Bed mount	Metal bracket on the back of the controller; attaches controller to the cart or bed
Purge cassette door release	Button located on the left side of the controller; press to open the purge cassette door
RS-232 service jack	Interface for data transfer by Abiomed maintenance or service personnel
USB connector	Connection for data downloading by Abiomed maintenance or service personnel
AC fuses	Electrical safety device in the event of current overload
AC plug	Connection point on the controller for the AC power cord
Power switch	Button that turns the controller on or off • ON: Press and hold the power switch for 3 seconds • OFF: (1) Disconnect the Impella® Catheter from the Automated Impella® Controller (2) Press and hold the power switch for 3 seconds (3) A pop-up confirmation box will appear (4) Press OK using the selector knob to confirm that the controller should be turned off NOTE: Holding down the power switch for longer than 30 seconds during operation will cause the controller to initiate an emergency shutdown
Equipotential ground stud	Used to ground the Automated Impella® Controller according to hospital procedures
Ethernet jack	Connection for downloading data or software upgrades

AUTOMATED IMPELLA® CONTROLLER HOME SCREEN

The home screen displays operating parameters and information for the entire Impella® System.

Display Element	Description
Alarm window	The alarm window displays up to 3 alarms simultaneously, in order of priority from top to bottom. For each alarm, the alarm window displays: - Alarm header – displayed in the left column; window is color-coded red for critical alarms, yellow for serious alarms, white for advisory notifications, gray for resolved alarms - Alarm subhead (if applicable) – further describes the alarm condition - Detailed text – up to 3 lines of instructions for resolving the alarm condition are displayed in the right column of the alarm window next to the alarm header and subhead information
Catheter serial number	Displayed in the upper left of the display screen if a catheter is connected to the controller.
System date and time	The current date (DD-MM-YYYY) and time (24-hour format; HH:MM) are displayed in the upper center of the screen display. (In this example it is December 21, 2010 at 5:30am.)
Mute alarm indicator	Displayed in place of the words "MUTE ALARM" when an alarm is silenced. - Yellow bell with red X displayed when an alarm is muted - Not displayed when an alarm is active (but not muted) or when there are no active alarms
Soft button labels	The soft buttons on the Automated Impella® Controller have corresponding labels adjacent to them on the display screen. These labels change depending on the type of screen displayed. MUTE ALARM - Mutes (silences) active alarms P-LEVEL (or NEXT) P-LEVEL – Allows you to set the performance level for the Impella® Catheter NEXT – Advances to the next screen DISPLAY (or BACK) DISPLAY – Brings up the Display menu for viewing waveforms and navigating to other screen displays BACK – Returns to the previous screen PURGE SYSTEM (or EXIT) PURGE SYSTEM – Brings up the Purge System menu for changing the purge fluid, purge cassette, or purge system; de-airing the purge system; or transferring to the standard configuration EXIT – Exits the current procedure MENU (or Exit Repositioning Guide) MENU – Brings up a menu of options related to controller settings, alarm history, repositioning, and starting a case Exit Repositioning Guide – Exits the repositioning guide

Display Element	Description
System power area	System power information is displayed to the right of the purge system information on the bottom of the display screen. Battery status – bar within battery symbol indicates the overall remaining capacity of the batteries • Full green bar for fully charged battery • Partial green bar for battery that is at least 50% charged • Partial yellow bar for battery that is between 16% and 50% charged • Partial red bar for battery that is less than or equal to 15% charged • Moving gray bar for battery that is in charging mode • Numeric percentage of battery power remaining displayed below the battery icon AC plug indicator • Green plug indicates that the controller is running on AC power • Gray plug with a red X indicates no AC power detected and the controller is running on battery power
Purge system area	Information about the purge system is displayed to the right of the flow area at the bottom of the display screen. Purge system marquee – scrolls from left to right when purge system is operating • Slow scrolling represents normal purge flow rate • Fast scrolling represents bolus flow rate Y connector icon • Appears above the purge system marquee when the Impella® System is configured using the Y connector in the set-up configuration Purge flow • Current purge flow displayed in mL/hr below the purge system marquee if the purge flow is known • Not displayed when the purge system is stabilizing, when there is no purge cassette, or when the procedure has not yet started Time remaining in reservoir • Displays the remaining runtime based on the volume of the purge fluid bag at the start of the case and the cumulative time and flow rate of delivery of purge fluid • Not displayed if there is no purge cassette or the procedure has not yet started
Flow area	Information about Impella® Catheter flow is displayed in the lower left corner of the display screen. Current flow rate • Mean catheter flow displayed in liters per minute (L/min) – the numbers appear in <i>white</i> if the catheter position is correct; <i>yellow</i> if the catheter position is incorrect or unknown • If the system is unable to calculate flow, a yellow triangular caution icon is displayed with the message “Flow Calculation Disabled” Catheter operation icon • The circular catheter operation icon rotates when the Impella® Catheter is running Low flow alarm limit icon • Indicates whether the alarm is enabled or disabled • Enabled: bell symbol and numeric alarm limit are displayed • Disabled: bell symbol is crossed out
Central display area	On the home screen, the central display area displays a heart pictogram and Impella® Catheter position indicator message. Heart pictogram appears in the center of the home screen display. • Provides a visual representation of the current Impella® Catheter position • Overlaid with a translucent yellow “?” when the controller detects an incorrect catheter position or cannot determine catheter position Impella® Catheter position indicator message displayed to the left of the heart icon. • Displays “Impella Position OK” in green when catheter position is correct • Displays “Impella Position Unknown” in yellow when catheter position is unknown • Displays “Impella Position in Ventricle” in yellow when catheter is in the ventricle • Displays “Impella Position Wrong” in yellow when catheter position is incorrect • Displays “Placement Monitoring Suspended” in yellow when there is a fault in the sensor • Displays “Placement Monitoring Disabled” in yellow when you turn off placement monitoring through the menu

To select a performance level:

1. Press the **P-LEVEL** soft button to open the performance level icon.
2. Turn the selector knob to increase or decrease the performance level.
3. Press the selector knob to select the new performance level.

AUTOMATED IMPELLA® CONTROLLER WAVEFORM SCREEN

The waveform screen displays real-time operating data for the system. The screen displays the placement signal and motor current waveforms as well as the maximum/minimum and average values for each waveform in the central display area of the screen.

AUTOMATED IMPELLA® CONTROLLER INFUSION SCREEN

The infusion screen displays purge system data. In the central display area of the screen, the purge flow rate and purge pressure are plotted as a function of time. To the right of the plots, the current purge flow rate and purge pressure are displayed.

MOBILE OPERATION

The Li-Ion batteries must be charged for 5 hours prior to system operation in order to meet the runtime requirement of 1 hour. Failure to do so will yield a shorter runtime. After being unplugged, the Automated Impella® Controller will operate for at least 60 minutes after the batteries have been fully charged.

The Automated Impella® Controller can be operated on internal battery power when it is not connected to AC power.

1. Disconnect the Automated Impella® Controller from AC power.
2. The controller beeps once every 5 minutes to alert you that it is running on battery power and a white advisory notification appears in the alarms area on the screen. The AC power icon turns gray with an X through it.
3. When the controller is connected back to AC power, the white advisory notification turns gray and the AC power icon turns green.

TRANSPORT WITHIN THE HOSPITAL

Patients supported with the Impella® System may require transport within the hospital for various reasons. Transport can be safe and simple for patients supported with Impella® Catheters.

Considerations for transport within the hospital:

- The Automated Impella® Controller and Impella® Catheters are designed to operate on battery power for at least 1 hour.
- Confirm that the battery capacity displayed on the controller is 100%.
- If transport time might be longer than 1 hour, bring an extension cord or confirm that you will be able to connect the controller to AC power once you arrive at your destination.
- Use care when rolling the Automated Impella® Controller cart, and pay close attention when going over thresholds and through elevator doors.
- Do not stress the connector cable from the controller to the Impella® Catheter.

TRANSPORT BETWEEN HOSPITALS

During transport, the Automated Impella® Controller may be exposed to stronger electromagnetic disturbance than during in-hospital use. Strong electromagnetic disturbance may cause the Automated Impella® Controller to display soft button menu selections that were not selected by the user. Operators should be aware that, under these conditions, the operating parameters are not affected. No user intervention is required. Monitor Impella® Catheter flow and patient hemodynamics to confirm normal operation. The condition will resolve itself once the Automated Impella® Controller is no longer exposed to the disturbance.

The Automated Impella® Controller is qualified for safe use by healthcare professionals to facilitate the transport of patients supported by the Impella® System from one medical facility to another. The Impella® System performs life-sustaining functions. To use the system during transport you must understand and follow the instructions in this manual. The transport team should include person(s) fully trained in the use of the Automated Impella® Controller and the Impella® Catheter.

Battery Power Note

If the Automated Impella® Controller is allowed to discharge completely and the system shuts down due to low battery, the controller will need to charge for an extended period of time before it will turn on again.

Patients are placed on the Automated Impella® Controller for circulatory support for periods compliant with the intended use duration of the Impella catheters. If, during this time, a patient requires additional resources and specialized teams at another location (eg, a larger facility such as a transplant center), the patient may be transferred safely to such a location using the Automated Impella® Controller.

GUIDELINES FOR PATIENT TRANSPORT

Maintaining optimal patient hemodynamic status and correct Impella® Catheter position are two key factors in managing patients supported with the Impella® System during transport. Steps should be taken to eliminate or minimize any aspect of the transport that might adversely affect these factors.

The Automated Impella® Controller is designed to operate for 60 minutes on battery power. Transport teams should take this into consideration when planning the transport. If the total transport time is expected to include more than 60 minutes during which the system will be disconnected from AC power, arrangements should be made to use a vehicle with a built-in DC to AC power inverter.

IMPORTANT TRANSPORT CONSIDERATIONS

1. Planning is critical to success. Abiomed representatives can help with planning for transport. They can be contacted 24 hours a day at the Abiomed 24hr hotline (see back cover).
2. The Automated Impella® Controller should be fully charged prior to transport. Keep the Automated Impella® Controller connected to AC power (or an AC inverter) whenever possible.
3. Do not stress the connector cable from the controller to the Impella® Catheter. Such tension could move the catheter out of correct position and compromise patient circulatory support.
4. Carefully monitor purge pressures during changes in altitude.
5. The controller should be positioned to allow easy access to the display screen and soft buttons to view alarms and make any necessary changes.
6. Maintain ACTs between 160 and 180 or at the level recommended by the physician responsible for the patient.

DEFIBRILLATION

During defibrillation, do NOT touch the Impella® Catheter, cables, or Automated Impella® Controller.

If emergency defibrillation is required during support with the Impella® System, universal safety precautions should be followed. To avoid the risk of electrical shock, caregivers should use caution during defibrillation to ensure they do not touch the Impella® Catheter, cables, or controller.

ECG INTERFERENCE

Operating the Automated Impella® Controller may cause interference with electrocardiogram (ECG) signals. Please check the electrode pads and leads for good fixation and contact. If interference persists, activate the 50/100 Hz band-elimination filter or the 60/120 Hz band-elimination filter (also known as notch filter) on your ECG device. The filter frequency will be based on the AC power frequency for the country in which you are operating the equipment. If your ECG device does not have the appropriate filters, disconnect the Automated Impella® Controller temporarily from AC power to obtain an undisturbed signal. Please observe the battery status while running the Automated Impella® Controller on battery power.

LATEX

The Automated Impella® Controller is 100% latex free.

DATA SNAP SHOT RECORDING

The Automated Impella® Controller can hold up to 24 hours of real-time data. Once memory is full, the controller starts overwriting the old data. The Data Snap Shot feature allows you to permanently save real-time operating data for later analysis. Data Snap Shot is automatically turned on during certain alarm conditions to capture data for analysis. You can also manually turn on the feature at any time to capture data for later analysis.

To manually access the Data Snap Shot feature:

1. Press **MENU** and scroll to "Data Snap Shot." Press the selector knob.
2. The controller records data for a predefined period of 10 minutes.

INFUSION HISTORY

The Automated Impella® Controller has an Infusion History screen that displays the infusion volume as well as the amount of heparin and dextrose infused each hour. The current time period is displayed at the top of the list. The calculations begin when the case start procedure is completed and Impella® Catheter performance level is greater than P0.

HOW TO CHANGE TO A BACKUP CONTROLLER

A backup Automated Impella® Controller should be available at all times when a patient is on support. In the event that the controller fails, follow the steps below to transition the Impella® Catheter to the backup controller.

1. Confirm that the backup controller is powered on and ready.
2. Disconnect the yellow luer connector from the Impella® Catheter to release the pressure and then immediately reconnect it.
3. Transfer the purge cassette and purge solution from the original controller to the backup controller.
4. Remove the white connector cable from the original controller and plug it into the catheter plug on the front of the backup controller.
5. Once the Impella® Catheter is connected to the backup controller, a message will appear on the screen asking you to confirm re-starting the Impella® Catheter at the previously set performance level.
6. Press **OK** within 10 seconds to confirm restarting the Impella® Catheter at the previously set performance level.

EMERGENCY SHUTDOWN PROCEDURE

In the unlikely event that the Automated Impella® Controller software stops responding, follow the procedure below to restart the controller without stopping the Impella® Catheter.

1. Press and hold the power switch for 30 seconds.
2. An "Emergency Shutdown Imminent" alarm will sound at 15 seconds.
3. The controller will shut down after 30 seconds.
4. Restart the controller.

ALARMS OVERVIEW

The Automated Impella® Controller monitors various functions to determine whether specific operational parameters are within expected limits. When a parameter goes outside of its specified limits, the controller sounds an alarm tone and displays an alarm message that can be viewed on the display screen on the front of the controller. The alarm tone indicates the severity of the alarm. The alarm message on the display screen is color-coded for severity and provides details on the cause of the alarm and how to resolve the alarm.

Change Purge Fluid to Obtain Accurate Purge Values

To get accurate purge values after changing to a backup controller, perform the Change Purge Fluid procedure (described in section 5 of this manual) and replace the purge fluid bag.

ALARM LEVELS

Alarms are divided into three levels of severity:

- Advisory (white)
- Serious (yellow)
- Critical (red)

Category	Description	Audible Indicator*	Visual Indicator
Advisory	Notification	1 beep every 5 minutes	Alarm header on white background
Serious	May become harmful or life-threatening if not addressed immediately	3 beeps every 15 seconds	Alarm header on yellow background
Critical	Immediately harmful or life-threatening	10 beeps every 6.7 seconds	Alarm header on red background
* Sound pressure of audible alarm indicators is >80 dBA			

MUTE ALARM FUNCTION

Pressing the **MUTE ALARM** button on the upper right of the Automated Impella® Controller display screen will silence the audible alarm indicator for 2 minutes (for red or yellow alarms) or 5 minutes (for white advisory alarms). When an alarm is silenced, the words "MUTE ALARM" next to the button are replaced by the mute alarm indicator, a crossed-out bell icon.

ALARM HISTORY SCREEN

The alarm history screen may be accessed through the **MENU**. This screen contains a log of the alarms that occurred during the case. This log is not maintained when the Automated Impella® Controller is powered down or after a power failure. The controller does, however, maintain a long-term log that is saved after the Automated Impella® Controller is powered down or after a power failure and this information may be downloaded by Abiomed personnel.

Alarms That Resolve On Their Own

The audible indicator will shut off if an alarm condition is resolved before you press **MUTE ALARM**. The visual message, however, will continue to be displayed, with the alarm header on a gray background, for 20 minutes or until you press **MUTE ALARM**. This allows you to identify the alarm that occurred.


ALARM MESSAGE SUMMARY

Severity	Alarm Header	Action	Cause
Critical Alarms	Impella Stopped	1. Replace controller.	There is a problem with the controller electronics.
	Controller Failure	2. Restart Impella.	
	Controller Failure	Switch to backup controller.	There is a problem with the controller electronics.
	Impella Stopped	1. Restart Impella. 2. Replace Impella after 3rd unsuccessful restart attempt.	There may be a mechanical or electrical problem in the Impella® Catheter.
	Impella Disconnected	1. Check cable connection to console. 2. Check Impella connection to cable.	Running Impella® Catheter disconnected.
	Emergency Shutdown Imminent	Release ON/OFF push button.	Power switch pressed for 15 seconds while Impella® Catheter still connected.
	Battery Failure	1. Plug controller into AC power. 2. Press switch located on the underside of the controller. 3. Switch to backup controller.	A battery switch is turned off or there is a malfunction of the switch.
	Battery Temperature High	Switch to backup controller.	Battery temperature is greater than 60 °C.
	Battery Critically Low	Plug controller into AC power.	Battery power has 15% remaining capacity.
	Battery Failure	Plug controller into AC power.	One of the batteries has failed.
	Air in Purge System	The purge system has stopped. Follow the instructions in the De-air Tool to remove the air from the system.	There is air in the purge tubing.
	Purge System Failure	Switch to backup controller and replace purge cassette.	There is a problem with the purger unit driver.
	Impella Stopped Reverse Flow	Restart Impella, or remove Impella from ventricle.	Impella® Catheter is not running; possible reverse flow through Impella® Catheter.
	Impella Position Wrong	1. Confirm Impella position with imaging. 2. Follow repositioning guide if needed.	Controller has detected that Impella® Catheter is in the wrong position.
	Impella Failure	Replace Impella.	There is a problem with the Impella® Catheter motor.
	Impella Motor Current High	Replace Impella.	There is a problem with the Impella® Catheter motor.
	Purge System Open	1. Check the purge system tubing for open connections or leaks. 2. Replace purge cassette.	Purge pressure has dropped below 100 mmHg for 20 seconds or longer.
	Purge Pressure Low	1. Check purge system tubing for leaks. 2. Increase concentration of dextrose in the purge solution. 3. Replace purge cassette.	Purge pressure has dropped below 300 mmHg with the purge flow ≥30 mL/hr for 30 seconds or longer.
	Purge System Blocked	1. Check all purge system tubing for kinks or blockages. 2. Decrease concentration of dextrose in the purge solution.	Purge flow has dropped below 1 mL/hr. Kinked or blocked purge connecting tube. Kinked or blocked purge lumen in Impella® Catheter.
	Purge Flow Low	1. Check purge system tubing for kinks. 2. Decrease concentration of dextrose in the purge solution.	Purge pressure is ≥1100 mmHg with the purge flow <2 mL/hr.
	Purge Line Click-On Not Detected	Check the purge line click-on and make sure it is fully inserted.	The controller is not detecting that the purge pressure transmitter is clicked into the front of the controller.

Severity	Alarm Header	Action	Cause
Serious Alarms	Controller Error	Switch to backup controller.	There is a problem with the controller electronics.
	Impella Defective	Do not use Impella. Replace Impella.	There is a problem with the Impella® Catheter electronics.
	Impella Position Wrong	1. Confirm Impella position with imaging. 2. Pull Impella back 2 cm. 3. Follow repositioning guide if needed.	Controller has detected that the Impella® Catheter is in the wrong position, with the outlet area too close to the aortic valve.
	Suction	1. Check left side filling and volume status. 2. Check Impella position. 3. Reduce flow or performance level.	Suction is detected.
	Battery Temperature High	1. Check controller for blocked air vents. 2. Switch to backup controller.	Battery temperature is greater than 50 °C and less than or equal to 60 °C.
	Impella Flow Low	1. Adjust low flow limit value. 2. Check for suction. 3. Check for high afterload pressure.	Actual flow is below the target P-level.
	Impella Sensor Failure	Placement monitoring is suspended. 1. Monitor patient hemodynamics. 2. Monitor Impella position with imaging.	There is a problem with the Impella® Catheter sensor signal.
	Battery Level Low	Plug controller into AC power.	Battery power has 50% remaining capacity.
	Battery Comm. Failure	Plug controller into AC power.	Loss of communication to the battery.
	Purge Cassette Failure	Replace purge cassette.	There is a problem with the purge cassette software.
Advisory Alarms	Purge Volume Critically Low	1. Open the PURGE SYSTEM menu and select Change Purge Fluid. 2. Follow the instructions to change the purge fluid.	There are 15 mL (in addition to 5% of the starting bag volume) or fewer remaining in the purge fluid bag.
	Purge Volume Low	1. Open PURGE SYSTEM menu and select Change Purge Fluid. 2. Follow the instructions to change the purge fluid.	There are 30 mL (in addition to 5% of the starting bag volume) or fewer remaining in the purge fluid bag.
	Purge Flow Increased	The purge flow has increased by 2.5 mL/hr or more. This is a notification only; no action is required.	Purge flow has increased by ≥2.5 mL/hr.
	Purge Flow Decreased	The purge flow has decreased by 2.5 mL/hr or more. This is a notification only; no action is required.	Purge flow has decreased by ≥2.5 mL/hr.
	Impella Position Unknown	Impella Catheter position unknown due to low pulsatility. Assess cardiac function.	Impella Catheter position unknown due to low pulsatility; or catheter position unknown detected by algorithm.
	AC Power Disconnected	Controller is running on battery power.	AC power was disconnected.
	Transfer to Standard Configuration	Follow instructions under Purge System to transfer to Standard Configuration.	Follow instructions or press MUTE ALARM to clear the alarm for 30 minutes.
	Complete Procedure	1. Follow the steps on the screen, or 2. Exit the procedure.	User has not responded to a particular screen for more than 1 minute (in de-air and purger procedures) or 5 minutes (in transfer to standard configuration procedure).

TERMINOLOGY, ABBREVIATIONS, AND SYMBOLS

SYMBOLS

	Caution; consult instructions for use
	Defibrillator-proof type CF equipment
	Keep dry
	Storage temperature (eg, 10 °C to 30 °C)
	Declares conformity with directive 93/42/EEC for medical devices
	Date of manufacture (eg, February 2009)
	Protect from sunlight
	Symbol for lot designation; the manufacturer's lot designation must be stated after the LOT symbol
REF 123456	Abiomed part number (eg, part number 123456)
SN 123456	Manufacturer's serial number (eg, serial number 123456)
Non Sterile!	The product is not sterile
	Use-by date (eg, use before June 2012)
	Do not reuse
	Sterilized using ethylene oxide
	Electric scrap; must be disposed of separately. Must not be disposed of as domestic waste.
	Protective Earth
	ON / OFF
	Alternating current (AC) only
	Equipotentiality
	Fuse
	Non-ionizing electromagnetic radiation
	USB port
	CAT 5 Port (Ethernet)

AUTOMATED IMPELLA® CONTROLLER MECHANICAL SPECIFICATIONS

Parameter	Specification	
Temperature	Operating:	10 °C to 40 °C (50 °F to 104 °F)
	Storage:	–15 °C to 50 °C (5 °F to 122 °F)
Relative Humidity	Operating:	95%
	Storage:	95%
Atmospheric Pressure	Operating:	8000 ft (750 hPa) to –1000 ft (1050 hPa)
	Storage:	18,000 ft (500 hPa) to –1000 ft (1050 hPa)
Dimensions	Height:	351 mm (13.8 in)
	Width:	443 mm (17.4 in)
	Depth:	236 mm (9.3 in)
Dimensions – Packaged	Height:	508 mm (20.0 in)
	Width:	559 mm (22.0 in)
	Depth:	406 mm (15.0 in)
Weight	Maximum:	11.8 kg (26.1 lbs)
Weight – Packaged	Maximum:	13.6 kg (30 lbs)
Maintenance and repair interval	12 months (Work must be performed by technicians authorized by Abiomed)	

AUTOMATED IMPELLA® CONTROLLER ELECTRICAL SPECIFICATIONS

AC operation	100-230 V AC (nominal); 47-63 Hz; 1.1 A
Internal battery operation	14.4 V DC (nominal); lithium ion
Characteristic values	
Max. power consumption under load	120 VA
9.7 fuses	2 Amp. 250 V. 5 mm x 20 mm, slow-blow fuses
Running time without AC power with fully charged batteries	At least 60 minutes (charging duration of at least 5 hours)
Electrical system	Installation in accordance with pertinent regulations is required for use in medical facilities (eg, VDE 0100, VDE 0107, or IEC stipulations). Observe country-specific regulations and national deviations.

EQUIPMENT DESIGN

The Automated Impella® Controller conforms to the applicable requirements of the following standards:

- UL 60601-1 (2003), 1st Edition *Medical Electrical Equipment, Part 1: General Requirements for Safety*
- CAN/CSA C22.2 No 601.1-M90 (1990; Reaffirmed 2005), 2nd Edition *Medical Electrical Equipment, Part 1: General Requirements for Safety*
- EN 60601-1 (1990), 2nd Edition *Medical Electrical Equipment, Part 1: General Requirements for Safety + A1(93) + A2(95) + A1.3(96)*
- IEC 60601-1 (1988), 2nd Edition *Medical Electrical Equipment, Part 1: General Requirements for Safety + A1(91) + A2(95)*
- IEC 60601-1-2 (2001), *Medical Electrical Equipment, Part 1-2: General Requirements for Safety – Collateral Standard: Electromagnetic Compatibility – Requirements and Tests*
- IEC 60601-1-4 (2000), Edition 1.1 Consolidated Edition *Medical Electrical Equipment, Part 1-4: General Requirements for Safety – Collateral Standard: Programmable Electrical Medical Systems*
- IEC 60601-1-1, (2000/12/01), 2nd Edition *Medical Electrical Equipment, Part 1-1: General Requirements for Safety – Collateral Standard: Safety Requirements for Medical Electrical Systems*
- IEC 60601-1-8, (2003/08/01), Edition 1 *Medical Electrical Equipment, Part 1-8: General Requirements for Safety – Collateral Standard: General Requirements, Tests and Guidance for Alarm Systems in Medical Electrical Equipment and Medical Electrical Systems*
- IEC 60601-1-6, (2004/06/01), Edition 1 *Medical Electrical Equipment, Part 1-6: General Requirements for Safety – Collateral Standard: Usability*

EQUIPMENT CLASSIFICATIONS

Type of protection against electric shock	IEC 60601-1: Class I degree of protection: CF defibrillation-proof and internally powered. Relies not only on basic insulation against shock but also includes additional protection. Accomplished by providing means for connecting the equipment to the protective earth conductor of the fixed wiring of the installation in a way that prevents accessible metal parts from becoming live if basic insulation fails.
Degree of protection against electric shock for Automated Impella® Controller	Class I Equipment
Mode of operation	Continuous
Degree of protection against explosion hazard	Not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. Also not suitable for use in an oxygen-enriched atmosphere.
Degree of protection against harmful ingress of water	IEC 60529: IPX1 protected against dripping water.

FEDERAL COMMUNICATIONS COMMISSION (FCC) NOTICE

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by Abiomed, Inc. could void the user's authority to operate this device.

ELECTROMAGNETIC COMPATIBILITY

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the electromagnetic compatibility (EMC) information provided in this document.

Portable and mobile RF communications equipment can affect medical electrical equipment.

The equipment or system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

Use of cables, other than those sold by Abiomed, may result in increased emissions or decreased immunity of the Automated Impella® Controller.

The Automated Impella® Controller uses RFID (radio frequency identification) to identify and communicate with the purge cassette. Other equipment may interfere with the Automated Impella® Controller even if that other equipment complies with CISPR emission requirements.

During transport, the Automated Impella® Controller may be exposed to stronger electromagnetic disturbance than during in-hospital use. Strong electromagnetic disturbance may cause the Automated Impella® Controller to display soft button menu selections that were not selected by the user. Operators should be aware that, under these conditions, the operating parameters are not affected. No user intervention is required. Monitor Impella® Catheter flow and patient hemodynamics to confirm normal operation. The condition will resolve itself once the Automated Impella® Controller is no longer exposed to the disturbance.

NOTE: The EMC tables and other guidelines that are included in this manual provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use permit the equipment or system to perform to its intended use without disturbing other equipment and systems or non-medical electrical equipment.

TABLE 201

Guidance and Manufacturer's Declaration – Emissions, All Equipment and Systems

The Automated Impella® Controller is intended for use in the electromagnetic environment specified below. The customer or user of the Automated Impella® Controller should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Enforcement – Guidance
RF Emissions CISPR 11	Group 1 Class A	The Automated Impella® Controller uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Harmonics IEC 61000-3-2	Class A	The Automated Impella® Controller is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Flicker IEC 61000-3-3	Complies	

TABLE 202
Guidance and Manufacturer's Declaration – Immunity

The Automated Impella® Controller is intended for use in the electromagnetic environment specified below.
The customer or user of the Automated Impella® Controller should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are synthetic, the relative humidity should be at least 30%.
Electrical Fast Transient/burst IEC 61000-4-4	±2 kV Mains ±1 kV for input/output lines	±2 kV Mains ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV Differential ±2 kV Common	±1 kV Differential ±2 kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	> 95% dip for 0.5 cycle 60% dip for 5 cycles 30% dip for 25 cycles > 95% dip for 5 seconds	> 95% dip for 0.5 cycle 60% dip for 5 cycles 30% dip for 25 cycles > 95% dip for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Automated Impella® Controller requires continued operation during power mains interruptions, it is recommended that the Automated Impella® Controller be powered from an uninterruptible power supply or battery.
Power Frequency 50/60 Hz Magnetic Field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be that of a typical location in a typical commercial or hospital environment.

TABLE 203
Guidance and Manufacturer's Declaration – Emissions, Equipment and Systems that are Life-Supporting

The Automated Impella® Controller is intended for use in the electromagnetic environment specified below.
The customer or user of the Automated Impella® Controller should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
			Portable and mobile RF communications equipment should be separated from the Automated Impella® Controller by no less than the recommended separation distances calculated/listed below:
Conducted RF IEC 61000-4-6	10 Vrms 150 kHz to 80 MHz	10 Vrms	$d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	20 V/m	$d = 0.6\sqrt{P}$ 80 MHz to 800 MHz $d = 1.2\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum power rating in watts and d is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey ^(a) , should be less than the compliance level in each frequency range. ^(b) Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

^(a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Automated Impella® Controller is used exceeds the applicable RF compliance level above, the Automated Impella® Controller should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Automated Impella® Controller.

^(b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

TABLE 205
Recommended Separation Distances Between Portable and Mobile
RF Communications Equipment and the Automated Impella® Controller,
Equipment and Systems that are Life-Supporting

The Automated Impella® Controller is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the Automated Impella® Controller can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the Automated Impella® Controller as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output Output Power of Transmitter (Watts)	Recommended Separation Distances for the Automated Impella® Controller (m)		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 0.6\sqrt{P}$	800 MHz to 2.5 GHz $d = 1.2\sqrt{P}$
0.01	0.12	0.06	0.12
0.1	0.38	0.19	0.38
1	1.2	0.6	1.2
10	3.8	1.9	3.8
100	12	6.0	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

RFID Transmitter / Receiver Specifications

Frequency	13.56 MHz
Receiver bandwidth	14 kHz
Effective radiated power	30 nW
Modulation	ASK

WHITE CONNECTOR CABLE

Length	2.5 m
Service life	Single use only

CLEANING

- Clean the Automated Impella® Controller keypad and display with either 70% isopropyl alcohol or soap and water. (NOTE: Be aware that soft buttons may be activated when you spray or wipe the display.)
- Clean the Automated Impella® Controller housing with mild detergent.
- Do not allow any fluids to enter the connector sockets.
- Clean the connector cable with 70% isopropyl alcohol.

STORING THE AUTOMATED IMPELLA® CONTROLLER

The Li-Ion batteries must be charged for 5 hours prior to system operation in order to meet the runtime requirement of 1 hour. Failure to do so will yield a shorter runtime. After being unplugged, the Automated Impella® Controller will operate for at least 60 minutes after the batteries have been fully charged.

- Place the Automated Impella® Controller on a horizontal surface to prevent falling.
- Connect the AC power cord to an AC outlet.
- The battery may be destroyed if the Automated Impella® Controller is stored with a depleted battery.

Storing the Controller

To keep the Automated Impella® Controller battery charged, the controller should be plugged into an AC outlet. When plugged into an AC outlet, the controller battery will charge whether the controller is on or off.

IMPELLA® SYSTEM LIMITED SERVICE WARRANTY

IMPELLA® SYSTEM LIMITED SERVICE WARRANTY (EUROPE)

Please contact your EU representative for warranty information.

TECHNICAL SAFETY INSPECTIONS, MAINTENANCE, AND REPAIR

TECHNICAL SAFETY INSPECTIONS

In accordance with the specification issued by Abiomed as per §6 MPBetreibV (German regulation), the Automated Impella® Controller must be subject to yearly technical safety inspections. Technical safety inspections may only be carried out by authorized technicians in accordance with the technical safety inspection requirements presented below and must be documented in the medical product logbook in accordance with §7 MPBetreibV.

A sticker on the device indicates the date for the next required inspection. This example shows a sticker indicating that inspection is required in May 2012. However, the stipulations in the medical product logbook are binding in any case.

			6	7	8	
4			00			9
3			00			10
2						11
1	EXPIRE DATE					12
07	08	09	10	11		13

The following technical safety inspections are required for the Automated Impella® Controller:

- Inspection of labeling and instructions for use
- Visual inspection of the device and its accessories for any signs of damage
- Testing for electrical safety as per DIN VDE 751 or DIN EN 60 601
- Leakage current test
- Dielectric strength test
- Functional testing of all switches, keys, rotary knobs, sockets, and control lights on the device
- Checking battery operation

If defects become apparent during the technical safety inspections that could endanger patients, employees, or third parties, then the device must not be operated until the defects have been remedied by proper technical servicing.

MAINTENANCE AND REPAIR

The Automated Impella® Controller is subject to 1 year maintenance and repair intervals. The corresponding work must be performed by technicians authorized by Abiomed, and must be documented in the medical product logbook in accordance with §7 MPBetreibV (German regulation).

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Abiomed Europe GmbH