
Distribution

P&R
P&R Medical
High Tech Science Park 1/13
University Hasselt 3590 B-Diepenbeek
Belgium
Tel: 0032(0) 11 23 46 81
Fax: 0032(0)11 23 46 82

Manufacturer

MHMedical Tec GmbH
D-63739 Aschaffenburg

AFTER INJECTION. Disconnect the patient line and leave the proximal part of the Transflux™ connected to the Transset™ as a safety cap. Close the clamp. Remove this part before installing the new Transflux™ for a new patient.

4. Contraindications

Not elsewhere specified.

5. Warnings symbols

Transflux™ is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Not for reuse
Any reuse of this product may cause a change in its mechanical and biological characteristics and could lead to product failure, allergic reactions or bacterial infections.

Do not use when package is damaged!
Patient or operator injury may result if package is opened or damaged or if damaged components are used. For such reason, it is very important that you visually inspect contents and package before each use.

Do not re-sterilize the product!


2. Description of the Transflux™ system

Transflux™ is a patient line with non-return valves, safety zones, Clip, 3-LL, Reflux control system.


The anti-reflux valves offer the security of an anti contamination barrier between the patient and the transfer system.

3. Application


These systems are designed only for single use (one-way) to be used with a power injector up to a maximum pressure of 300 psi (20,615 bar or 2068,5 kPa) in radiology.

1. Open the peel package, take out the Transflux™ and remove the protective caps.
2. Keep Transflux™ vertically and fill the saline emulsion into the Transflux™ either manually with the syringe or automatically with the injector in order to remove any air bubbles.

3. Insert the needle into the vein or into a port-a-cath and make sure that the luer fittings are firmly connected.
4. Connect Transflux™ with the needle.

5. Avoid subcutaneous injection or extravasations.
Briefly disconnect the connection at the one-way valve proximal to the patient to check the blood reflux.
6. Connect Transflux™ with the transfer set and the injector. Verify and make sure that the connecting tubing is not kinked or obstructed. Prime the connector tube and make certain that the air is expelled.

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4. Precautions

Please take care that all connections are properly installed and assure they are all firm and secure. Please do not overtighten! This will help minimize leaks, disconnections and system damage.

5. Adverse effects

Not elsewhere specified

8. Storage

The product should be stored in its original package at a cool and dry place under normal temperature. Avoid direct sunlight.


9. Disclaimer

The manufacturer will not be liable for any damage caused by reuse of the product, abuse, improper handling and non-compliance with warning and handling instructions.

Identification of the graphic symbols

- Caution ! Consult accompanying documents
- Consult instructions for use !
- Use by YYYY-MM-DD
- Batch code
- Order number
- Do not use when package is damaged
- Do not reuse !
- Sterilized using ethylene oxide
- Do not re-sterilize the product

Devices are worldwide patented and have Trade Market.
*In certain jurisdictions based upon PCT/EP2006/062554, PCT/EP2009/066383 and PCT/EP2011/062179

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