**Manuale di Istruzione Uso e Manutenzione**  
**User and Maintenance manual**  
**Manuel d'Utilisation et d'Entretien**  
**Gebrauchs- und Wartungsanleitung**  
**Manual de Instrucciones de Uso y Mantenimiento**

**Micromotore da laboratorio**  
**Laboratory micromotor**  
**Micromoteur de laboratoire**  
**Labor-Mikromotor**  
**Micromotor de taller**

**carlo de giorgi** s.r.l.

---

## Dati modello

### Model data / Donnés du modèle / Modelldaten / Datos del modelo

| <b>Tipo</b><br><b>Type</b><br><b>Type</b><br><b>Typ</b><br><b>Tipo</b> | <b>Modello</b><br><b>Model</b><br><b>Modèle</b><br><b>Modell</b><br><b>Modelo</b> | <b>Alimentazione</b><br><b>Power supply</b><br><b>Alimentation</b><br><b>Betriebs-spannung</b><br><b>Alimentación</b> | <b>Anno di costruzione</b><br><b>Year of manufacture</b><br><b>Année de fabrication</b><br><b>Baujahr</b><br><b>Año de fabricación</b> |
|------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| 2900                                                                   | Da tavolo                                                                         | <input type="checkbox"/> 230 V 50 Hz                                                                                  | <input type="checkbox"/> 2005                                                                                                          |
| Green Star                                                             | Da terra (R)                                                                      | <input type="checkbox"/> 120 V 60 Hz                                                                                  | <input type="checkbox"/> 2006                                                                                                          |
| Green Star Gold                                                        |                                                                                   |                                                                                                                       | <input type="checkbox"/> 2007                                                                                                          |
| 4000 light                                                             |                                                                                   |                                                                                                                       | <input type="checkbox"/> 2008                                                                                                          |
| 3000                                                                   |                                                                                   |                                                                                                                       | <input type="checkbox"/> 2009                                                                                                          |
| 5000                                                                   |                                                                                   |                                                                                                                       |                                                                                                                                        |

**carlo de giorgi** s.r.l.

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Tel. +39.02.356.15.43 r.a. – Fax +39.02.356.18.08

## CONTENTS

|                                                             |    |
|-------------------------------------------------------------|----|
| 1 – STANDARDS AND GENERAL PRECAUTIONS .....                 | 4  |
| 1.1 - INTRODUCTION.....                                     | 4  |
| 1.2 – DESCRIPTION OF THE APPLIANCE.....                     | 4  |
| 1.3 – REFERENCE STANDARDS .....                             | 4  |
| 1.3.1 – Obligatory compliance .....                         | 4  |
| 1.3.2 – Voluntary compliance.....                           | 4  |
| 1.4 – RESPONSIBILITIES OF THE CUSTOMER.....                 | 4  |
| 1.5 – REQUESTS FOR SERVICE AND SPARE PARTS .....            | 5  |
| 2- TECHNICAL FEATURES.....                                  | 5  |
| 2.1 – TECHNICAL FEATURES AND PERFORMANCE.....               | 5  |
| 3 - INSTALLATION.....                                       | 5  |
| 3.1 - POSITIONING.....                                      | 5  |
| 3.2 – CONNECTIONS AND PLACING INTO SERVICE .....            | 6  |
| 4 – OPERATION AND USE .....                                 | 7  |
| 4.1 - OPERATION.....                                        | 7  |
| 4.2 - OPERATOR.....                                         | 7  |
| 4.3 – INTENDED USE.....                                     | 7  |
| 4.4 – SAFETY WARNINGS .....                                 | 7  |
| 5 – OPERATOR INSTRUCTIONS .....                             | 8  |
| 5.1 –CONTROLS AND POWER SUPPLY UNIT .....                   | 8  |
| 5.2 – INSTRUCTIONS FOR USE.....                             | 8  |
| 5.2.1 - Startup .....                                       | 8  |
| 5.2.2 – Normal stopping and emergency stop procedures ..... | 9  |
| 5.2.3 - Warnings .....                                      | 9  |
| 6 - MAINTENANCE .....                                       | 10 |
| 6.1 – CLEANING AND LUBRICATION .....                        | 10 |
| 6.2 – SPECIAL MAINTENANCE .....                             | 10 |
| 6.2.1 – Collet fails to open .....                          | 10 |
| 6.2.2 – Replacing the handpiece .....                       | 11 |
| 6.2.3 – Replacing the motor.....                            | 11 |
| 7 - DIAGNOSTICS .....                                       | 12 |
| 8 - ACCESSORIES .....                                       | 13 |
| 8.1 – REDUCER BUSHES.....                                   | 13 |
| 8.2 – SCREWDRIVER FOR COLLET REPLACEMENT .....              | 13 |
| 8.3 - 3 mm Ø COLLET .....                                   | 13 |
| 8.4 – HANDPIECE SERVICING KIT .....                         | 13 |
| 8.5 – JEWELLERS HANDPIECE.....                              | 13 |
| 8.6 – HAMMER HANDPIECE.....                                 | 14 |
| 8.7 – INDUSTRIAL HANDPIECE .....                            | 15 |

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declares, under its own responsibility, that the apparatus:

Type: **LABORATORY MICROMOTOR**

described as:

laboratory micromotor for roughing, polishing and finishing materials utilized in orthodontic laboratories, jewellers workshops and tool shops,

conforms to the Laws implementing the requirements of the "Machinery Directive" 98/37/CE, the "Low Voltage Directive" 73/23/EEC and subsequent amendments, and the "Electromagnetic Compatibility" Directive 89/336/EEC.

Signed: De Giorgi Ariberto  
Managing Director

## **1 – STANDARDS AND GENERAL PRECAUTIONS**

### **1.1 - INTRODUCTION**

This manual is the property of CARLO DE GIORGI S.r.l. The contents of this manual cannot be reproduced or disclosed to other parties. All rights reserved.

### **1.2 – DESCRIPTION OF THE APPLIANCE**

The appliance in question is a laboratory micromotor designed for roughing, polishing and finishing of the materials utilized in orthodontic laboratories, jewellers workshops and tool shops.

### **1.3 – REFERENCE STANDARDS**

#### **1.3.1 – Obligatory compliance**

- EEC directive n.89/392 – Machinery Directive - (D.P.R. n° 459/1996)
- EEC directive n.73/23 – Low voltage Directive (Italian legislation: Law n° 791/1977, Decree Law n° 626/1996, Decree Law n° 277/97)
- EEC directive n.89/336 relative to Electromagnetic Compatibility (EMC) - (Italian legislation: Decree Law n° 615/1996)
- DPR 27.4.1955 n.547, "Norme per la prevenzione degli infortuni sul lavoro" (work accident prevention regulations)
- DPR 27.4.1956 n.303, "Norme generali per l'igiene del lavoro" (general work hygiene regulations)
- Italian Decree Law 12 August 1991 n.277 implementing EEC directives n.80/605; n.83/477; 86/188; 88/642 concerning the protection of workers against risks deriving from exposure to chemical, physical and biological agents in the workplace, as defined by art. 7 of Law 30 July 1990, n.212

#### **1.3.2 – Voluntary compliance**

- EN 292 (1992) Machine safety – Basic concepts; general design principles – Part 1 – Basic terminology and methods (EN 292-1) /Part 2 – Specifications and technical principles (EN 292-2).
- EN 60204-1 Machine safety – Machine electrical equipment - Part 1: General requisites

### **1.4 – RESPONSIBILITIES OF THE CUSTOMER**

The user must install the laboratory micromotor in a suitable place equipped with an electrical system that complies with established legislation.

The appliance must be installed in a dry environment that is illuminated in accordance with the requirements of established legislation.

NOTE: The expression "established legislation" refers to the laws in force in the user's country.

## 1.5 – REQUESTS FOR SERVICE AND SPARE PARTS

For all electrical maintenance requirements, contact CARLO DE GIORGI di Baranzate di Bollate, Via Tonale, 1, Milan, Italy.

In the case of faults, notify CARLO DE GIORGI technical personnel of the exact nature of the problem, specifying also the lot number shown on the label affixed to the rear of the appliance.

## 2- TECHNICAL FEATURES

### 2.1 – TECHNICAL FEATURES AND PERFORMANCE

Performance: intermittent duty

|                                         |                  | 2900                                                                              | 3000    | 4000<br>light | 5000    | Green<br>Star |
|-----------------------------------------|------------------|-----------------------------------------------------------------------------------|---------|---------------|---------|---------------|
| Voltage                                 | V                | 230/120                                                                           | 230/120 | 230/120       | 230/120 | 230/120       |
| Frequency                               | Hz               | 50/60                                                                             | 50/60   | 50/60         | 50/60   | 50/60         |
| Installed power                         | W                | 70                                                                                | 120     | 100           | 300     | 90            |
| V <sub>out</sub> to motor               | V                | 59                                                                                | 24      | 30            | 40      | 30            |
| Max. tool speed                         | rpm              | 25,000                                                                            | 30,000  | 30,000        | 30,000  | 30,000        |
| Delayed acting fuse                     | A                | 1                                                                                 | 1.25    | 1.25          | 2       | 1.25          |
| Mass                                    | kg               | 2.7                                                                               | 3.2     | 3.3           | 4.3     | 3.3           |
| Transmitted vibration<br>(acceleration) | m/s <sup>2</sup> | < 2.5                                                                             | < 2.5   | < 2.5         | < 2.5   | < 2.5         |
| Class II appliance                      |                  |  |         |               |         |               |

## 3 - INSTALLATION

### 3.1 - POSITIONING

The appliance is supplied assembled and ready for use after connecting the control devices and the micromotor complete with handpiece to the power supply. The package contains:


- power supply unit, complete with mains cable;
- micromotor complete with handpiece, equipped with 2.35 mm Ø collet (excluding Gold version);
- footpedal (bench top versions only);
- rubber mat;
- user instructions.

The appliance can be moved by hand without requiring lifting devices.

### 3.2 – CONNECTIONS AND PLACING INTO SERVICE

With reference to Figure 1, for bench top models, and Figure 2 for floor types, after positioning the appliance and checking that the main power switch (1) is set to OFF, make the following connections:

- footpedal (2) to connector (3) on the rear of the appliance (bench top models only)
- micromotor complete with handpiece (4) to connector (5) on the front of the appliance (on the rear in the case of floor versions – see Fig. 2)
- mains power cable, designed for connection to a mains socket outlet with overcurrent protection devices.

**Fig. 1 (bench top model)**

DO NOT START THE APPLIANCE UNLESS THE FOOTPEDAL AND HANDPIECE ARE CONNECTED. DISCONNECT FOOTPEDAL AND HANDPIECE ONLY AFTER THE APPLIANCE HAS BEEN DISCONNECTED FROM THE MAINS POWER SUPPLY.

**N.B.** DURING THE INITIAL PERIOD OF USE THE HANDPIECE MAY TEND TO BECOME HOT DURING OPERATION. THIS PHENOMENON IS INDUCED BY INITIAL ADJUSTMENT OF THE TRANSMISSION ELEMENTS AND WILL DISAPPEAR WITHIN THE SPACE OF ABOUT ONE WEEK OF NORMAL USE. IF THE PROBLEM SHOULD PERSIST, CONSULT CARLO DE GIORGI.

## 4 – OPERATION AND USE

### 4.1 - OPERATION

Bench top models should be placed on a work table in front of the operator; floor models must be positioned in such a way that they cannot be inadvertently displaced by the pressure exerted by the operator's foot – the unit should ideally be located against a contrast such as a pillar.

Check that the appliance is working by setting the main power switch (1) to ON and then pressing footpedal (2) (or by operating lever (9) in the case of floor models); the motor should begin to rotate. If the motor fails to start set the main power switch to OFF, unplug the appliance and then check that the various connectors are in the correct positions.

**Always switch OFF the appliance after use.**

### 4.2 - OPERATOR

The operator starts the appliance by means of a voluntary command (footpedal or floor lever) thereby assuming direct control of its operation.

### 4.3 – INTENDED USE

The laboratory micromotor is designed for use in orthodontic laboratories, jewellers workshops and tool shops, for machining various materials including resins, ceramics, gold, chromium cobalt, etc..

The appliance is not designed for use in applications that fall within the definitions of the Medical Equipment Directive 93/42/EEC, or any uses other than those specified by the manufacturer.

### 4.4 – SAFETY WARNINGS

- The appliance must be used exclusively by skilled personnel who are aware of the potential hazards;
- do not use the appliance until you have connected the footpedal (2) to connector (3) located on the rear of the unit (bench top models only);
- do not use the appliance in the vicinity of inflammable materials;
- do not apply oil to the handpiece from the collet side (part "B", Fig. 3);
- do not open the power supply case under any circumstances;
- do not attempt to change the tool while the motor is running;
- always check that the diameter of the tool shank is compatible with the supplied collets.

## 5 – OPERATOR INSTRUCTIONS

### 5.1 –CONTROLS AND POWER SUPPLY UNIT

The laboratory micromotor is equipped with the following command and control devices located on the power supply unit:

\* for bench top models refer to Fig. 1 – for floor models refer to Fig. 2\*

- a) main power switch with green indicator (1);
- b) black tool rotation direction selector (reverse) (6). When reverse tool rotation is selected (normal=clockwise, reverse=counterclockwise) amber LED (10) illuminates;
- c) tool rotation footpedal (2) (bench top models only);
- d) tool speed slider control (7) (bench top models only);
- e) motor rpm display (R.P.M. x 100) (8) (bench top models 3000 and 5000 only);
- f) footpedal type speed control (9) (floor models only);
- g) power circuit fuse (11).


Fig. 2 (floor model)

## 5.2 – INSTRUCTIONS FOR USE

### 5.2.1 - Startup

With the appliance powered down for safety, fit the required tool into the handpiece. This is achieved by turning lever (A - Fig. 3) through 90° either towards the right or towards the left, fitting the tool inside collet (B). Now return lever (A) to its original position to lock the tool.

The appliance can now be powered up by means of main power switch (1).


Start the handpiece motor by means of footpedal (2) for bench top models or pedal control (9) for floor models.

When the tool is stationary the direction of rotation can be inverted by means of reverse selector (6); when reverse mode is selected amber LED (10) will illuminate.

The tool speed can be set in relation to the type of work by means of slide control (7) on bench top models, or pedal control (9) on floor models.

Furthermore, bench top models 3000 and 5000 are equipped with a display (8) showing motor speed in rpm (x 100).

**Fig. 3**

### **5.2.2 – Normal stopping and emergency stop procedures**

The motor can be stopped as follows:

- by releasing footpedal (2) (or control (9) for floor models);
- by setting main power switch (1) to OFF, thereby disconnecting electrical power to the power supply unit;
- by unplugging the power cable from the electrical socket outlet.

### **5.2.3 - Warnings**

In relation to the material being machined, the user must adopt personal protection devices (facemask, gloves, etc.) and/or an efficient exhaustor system.


## 6 - MAINTENANCE

All work on electrical components accessible only after removing the casing must be performed exclusively by properly authorized and trained personnel. Certain parts of the unit carrying dangerously high voltages are accessible when the casing is removed, representing a source of danger during maintenance work if performed when the appliance is powered up.

### 6.1 – CLEANING AND LUBRICATION

The appliance requires the following maintenance work, especially when subject to continuous use (every day); before cleaning the appliance disconnect it from the power supply.

- clean felt (F) located inside head (E) as shown in Fig. 4;

**Fig. 4**

- blow compressed air through the collet in point (B) with the collet open (locking lever (A) turned through 90°); (See Fig. 3)
- lubricate the handpiece once a month by removing screw (D) and applying a drop of oil through the threaded hole. Never attempt to oil the handpiece by way of the collet (part B, Fig. 3).

### 6.2 – SPECIAL MAINTENANCE

#### 6.2.1 – Collet fails to open

With reference to Fig. 3:

if collet (B) should fail to open after lever (A) has been turned through 90°, or if the tool is insufficiently clamped when the lever is in the "locked" position, adjust the lever as follows:

- turn lever (A) through 90° to cause the collet to extend as far as possible;

- immobilize handpiece spindle (C) with a flat tool (2 mm);
- insert the special screwdriver (available on request – see chapter 8) into collet (B), which is secured to the end of spindle (C) by means of a screw;
- loosen the screw by one half turn in a counterclockwise direction. This operation will release collet (B);
- to adjust clamping force, turn the collet clockwise with your fingers to increase or counterclockwise to reduce;
- use the special screwdriver to tighten the collet screw;
- remove the tool used to immobilize the spindle.

### 6.2.2 – Replacing the handpiece

When replacing handpiece (13) by screwing it onto micromotor (14), take care that the raised features on handpiece spindle (C) (see fig. 3) enter their locations without having to force them; otherwise the motor coupling (15, Fig. 5) may be damaged.

As soon as you detect excess resistance, unscrew handpiece (14) and retighten it gradually; repeat this operation as many times as necessary until you find the correct position.

### 6.2.3 – Replacing the motor

After approximately 2000 hours of duty it may prove necessary to change the motor in order to restore the original performance characteristics of the new unit.

Proceed as follows:

- unscrew rear head (G) shown in Fig. 6 after releasing the cable to prevent it becoming twisted;
- unscrew the three screws (H) shown in Fig. 5;



Fig. 5

- extract the motor from its location after desoldering the wires in models 3000 and 5000, or after detaching connector (L) as shown in Fig. 6.
- in the case of soldered connections take care to reconnect the new motor observing the correct polarity; if you are in doubt, connect the motor temporarily and check the direction of rotation.

**Fig. 6**

When refitting the rear cap after installing a new motor take care not to twist the connecting wires.

## **7 - DIAGNOSTICS**

If the appliance is not working first check the condition of the fuse in the rear.

If necessary, fit a replacement fuse rated as specified in table 2.1 (technical data).

Before checking the fuse, set the main switch to 0 (OFF) and disconnect the appliance from the electrical power supply; now check the fuse and replace if necessary.

After the screwcap has been refitted the power cable can be reconnected.

**If malfunctions and faulty operation is persistent, set the main power switch to OFF, unplug the appliance and consult CARLO DE GIORGI.**

**When communicating a fault, ensure you explain the exact nature of the problem to our technical personnel.**

**WARNING!**

ALL WORK ON THE APPLIANCE MUST BE PERFORMED WITH  
THE POWER SUPPLY DISCONNECTED

**8 - ACCESSORIES****8.1 – REDUCER BUSHES**

The reducer bushes make it possible to use tools with different diameter shanks.

To use the bushes follow the instructions in heading 5.2.1 regarding tool insertion; in this case first insert the reducer bush and then insert the tool.

| Article | Diameter (mm) |
|---------|---------------|
| 500/11  | 1             |
| 500/09  | 1.6           |
| 500/10  | 2.35          |

**8.2 – SCREWDRIVER FOR COLLET REPLACEMENT**

This tool is the same one that is required to adjust the collet if it fails to open, as described in heading 6.2.1.

**8.3 - 3 mm Ø COLLET****8.4 – HANDPIECE SERVICING KIT**

The kit consists in a set of wrenches covering all user service interventions regarding the micromotor with handpiece.


**8.5 – JEWELLERS HANDPIECE**

Designed for use in jewellery workshops, this handpiece can be fitted exclusively to the Green Star Gold model micromotor.

As shown in Fig. 7, the handpiece comes complete with four collets (1) for different tool shank diameters (0 - 1 – 2.35 - 3 mm)

To insert the tool:

- a) press button (2) on the handpiece;
  - b) unscrew head (3);
  - c) fit the tool and retighten the head while keeping button (2) depressed;
- To fit a new collet perform the above procedure to open the head and then replace the collet (4) as shown in the figure on the right.


**Fig. 7**

## 8.6 – HAMMER HANDPIECE

This unit is compatible with all models of laboratory micromotors.

Tool (2) is inserted simply by screwing it into the handpiece head (1) (see Fig. 8).

**N.B.** The supplied bits (2) must be profiled in accordance with the specific requirements of the user.

**Fig. 8**

## 8.7 – INDUSTRIAL HANDPIECE


Designed for jewellers, toolmakers and die machining, this unit is compatible exclusively with model 5000.

As shown in Fig. 9, the handpiece comes complete with four collets (1) for different tool diameters (2.35 - 3 - 4 - 6 mm)

To fit the tool:

- insert the supplied pin (4) in hole (2);
- unscrew head (3);
- fit the tool, retighten the head and then remove the pin;

To fit a new collet perform the above procedure to remove the head and then replace the collet (5) as shown in the figure on the right.

**Fig. 9**

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