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
PNEUMATIC ASPIRATOR  
ASPIRATEUR PNEUMATIQUE  
LUFTDRUCKABSAUGER  
ASPIRADOR NEUMATICO

- MANUALE ISTRUZIONI
- OPERATOR'S HANDBOOK
- MANUEL D'UTILISATION
- GEBRAUCHSANWEISUNG
- MANUAL DE INSTRUCCIONES

CATTANI S.p.A.

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## **INDICE**

## **Pagina**

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# ***PNEUMATIC ASPIRATOR***

## ***INDEX***

	<b><i>Page</i></b>
— GENERAL RUNNING DATA .....	8
— LEGEND OF COMPONENTS .....	9
— GENERAL FEATURES .....	10
— INSTALLATION .....	10
— FUNCTIONING AND USE .....	10
— NOTICE .....	11
— MAINTENANCE .....	11
— MAINTENANCE OPERATIONS FOR TECHNICAL STAFF .....	11
— TRANSPORT AND STORAGE .....	11
— TRANSPORT OF SECOND-HAND APPLIANCES .....	11

## **General running data**

**Air-depression aspirator**  
**(testing carried out with compressor at 5 bar)**

		<b>SMALL TERMINAL</b>	
<b>PRESSURE bar</b>	<b>AIR CONSUMPTION l/min</b>	<b>ASPIRATED AIR l/min</b>	<b>DELIVERY HEAD mm H<sub>2</sub>O</b>
1.5	3.4	62.5	275
2	4.2	75	326
2.5	9.1	87.5	428
3	10.2	108.3	591
3.5	11.9	120.8	673
4	13.1	125	734
4.5	14.2	133.3	826
5	15.6	133.3	826

**The conversion factor NL / AIR liters at 5 bar is 5.84**

**Manufactured by CATTANI S.P.A. - PARMA - ITALY**

## **LEGEND OF COMPONENTS** (See page 30)

<b>040880</b>	<b><i>Outside pipe Ø 11</i></b>
<b>040634</b>	<b><i>Small terminal</i></b>
<b>041050</b>	<b><i>Secretion filter</i></b>
<b>201109</b>	<b><i>Small terminal seat</i></b>
<b>183010</b>	<b><i>Seat microswitch</i></b>
<b>201002</b>	<b><i>Cover</i></b>
<b>202450</b>	<b><i>Plug</i></b>
<b>201110</b>	<b><i>Pneumatic aspirator body</i></b>
<b>200443</b>	<b><i>Tube holder for exhausted air outlet</i></b>
<b>210909</b>	<b><i>O-ring</i></b>
<b>240500</b>	<b><i>Tube holder for compressed air inlet</i></b>
<b>040009</b>	<b><i>Pressure reducer</i></b>
<b>040007</b>	<b><i>Manometer</i></b>
<b>182020</b>	<b><i>Electrovalve 24 V~</i></b>
<b>201029</b>	<b><i>Secretion filter cover</i></b>
<b>201130</b>	<b><i>Secretion filter body</i></b>
<b>210900</b>	<b><i>O-ring</i></b>
<b>210300</b>	<b><i>O-ring</i></b>
<b>202520</b>	<b><i>Tube holder for outside pipe Ø 11</i></b>

## **General features**

*The pneumatic aspirator is a low-cost extremely reliable appliance, small in size and easy to install; it requires no maintenance. Its performance, low operating head and flow, does not allow for use where traditional high-speed suction is needed; on the contrary, it works very well as an emergency aspirator or else in the preliminary examination room where a weaker suction is sufficient.*

**See table at page 8 for air consumption value.**

## **Installation (illus. 1 page 30)**

*We suggest mounting the body of the pneumatic aspirator inside the dental unit (201110 with screws A to the respective holes, D), so that the bag container 201130 is outside the unit and filter 041050 can be easily inspected. Tube holder C is to be connected to the waste outlet creating a siphon.*

*Seat 201109 is to be fitted within the reach of the surgeon and his assistant with screws F. The aspirator-filter set 201110-201130 must be at a lower level compared to the seat 201109. A Rilsan tube Ø 6x4 will feed air to the pneumatic aspirator through electrovalve 182020 and tube holder 240500. Aspirated air will be exhausted through outlet E, when the unit is on. There are two options: you can allow the air to disperse in the room or convey it outside. In the first instance, in order to reduce noise amplification caused by metal plate or other objects, we suggest you to use a conveyer lined with foam rubber of 1-2 cm depth; it should reduce the noise considerably. We suggest you to direct the outlet tube upwards, so that condensed water (due to rarefaction) is free to move back to the aspirator and can be drained out by tube C. In case it should not be possible to direct the tube with the wanted slope, nor to create a siphon, it will be necessary to poke a hole in the lowest part of the tube in order to collect condensed water which will eventually evaporate.*

## **Functioning and use (illus. 1 pag. 30 - table page 8)**

*We suggest to rate the pressure level in the compressor tank at 5-6 bar, while the pressure reducing valve, 040009, must be rated so that the gauge 040007 is at 2.0-2.5 bar when the unit is in use and 2.5-3.0 when the unit is off. This way, we will obtain about 75/87,5 l/min aspirated air on the small terminal 040634, with the compressor using up 4,2/9,1 l/min.*

**Table a pag. 8 shows other possible working conditions of the aspirator.**

*Item 201109, single seat is available only for the small terminal, the most suitable for this kind of suction. By lifting the terminal from its seat 040634, microswitch 183010 opens electrovalve 182020 and aspiration starts. Aspirated liquid goes through filter 041050 and enters the unit, then it is drained to outlet C.*

## **Notice**

***Prior to starting any servicing operation on appliances that have been used, clean with Puli-Jet a few times as explained in the maintenance section. Use disposable GLOVES, GOGGLES, MASK and OVERALL.***

## **Maintenance (illus. 1 page 30)**

*In order to clean the pneumatic aspirator inside and avoid the formation of bacterial flora, every evening you must aspirate a Puli-Jet solution in hot water, as indicated on directions for Puli-Jet. A glass of Puli-Jet solution is enough for a thorough cleaning.*

## **Maintenance operations meant for technical staff**

*The manufacturer is at your disposal to supply spare parts, literature, instructions and anything that might be useful. Our Concessionaries, Agents, authorized Retailers and Service Engineers are provided with split up drawings, electrical diagrams, instructions and updatings about maintenance and service.*

*The appliance is guaranteed for a year from the date of sale provided that the guarantee card addressed to the manufacturer is returned with date of sale, retailer's stamp and customer's name.*

*The guarantee and manufacturer's liability cease in case the appliances and/or the plants have been tampered with operations of any kind done by persons unfit and so unauthorized by the manufacturer.*

*For any usage not contemplated or specified in this manual, please contact the manufacturer.*

## **Transport and storage**

*Packed equipment can be transported and stored at a temperature range of -10 °C +60 °C.*

*Packages must be kept away from water and splashing and cannot tolerate humidity >70%.*

*Packages with the same weight can be stored in piles of three only.*

## **Transport of second-hand appliances**

*Prior to packing cleanse and sanitize with Puli-Jet (see maintenance and cleaning) place unit into a polyethylene bag, seal and pack in 3-layer corrugated board.*

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ESAM S.p.A.