# HbA1c II

Tina-quant®  Hemoglobin A1c II

Standardized according to IFCC; transferable to DCCT/NGSP

● Indicates Roche/Hitachi analyzer(s) on which kit(s) can be used

Cat. No.	Bottle	Contents	717	747	902	904	911 912	917	MOD P	MOD D
11822039	1	REAGENT 4 x 17 mL								
	2	REAGENT 4 x 4 mL								
	3a-d	CALIBRATOR → 4 x 2 mL	●	●	●	●	●	●	●	
	4	REAGENT 4 x 17 mL								
11822080	1	REAGENT 2 x 66 mL								
	2	REAGENT 2 x 16 mL								
	3a-d	CALIBRATOR → 4 x 2 mL		●				●	●	
	4	REAGENT 2 x 66 mL								
11822098	1	REAGENT 2 x 71 mL								
	2	REAGENT 2 x 17 mL								
	3a-d	CALIBRATOR → 4 x 2 mL	●	●			●			
	4	REAGENT 2 x 71 mL								
11822497	1	REAGENT 40 x 66 mL								
	2	REAGENT 40 x 16 mL								
	3a-d	CALIBRATOR → 32 x 2 mL		●				●	●	
	4	REAGENT 40 x 66 mL								

Some analyzers and kits shown may not be available in all countries. For additional system applications, contact your local Roche Diagnostics representative.

## English

### System information

For Roche/Hitachi 904/911/912: ACN 027, ACN 125, ACN 162, ACN 165.

For Roche/Hitachi 917/MODULAR P: ACN 027, ACN 125, ACN 802

(US code 812), ACN 803 (US code 813).

### Intended use

Immunological assay for the in vitro quantitative determination of hemoglobin A1c in whole blood on Roche automated clinical chemistry analyzers. HbA1c determinations are utilized in the long-term monitoring of glycemia.

### Summary<sup>1-9</sup>

Hemoglobin (Hb) consists of four protein chains with four heme portions, and is the red-pigmented protein located in the erythrocytes. Its main function is the transport of oxygen and carbon dioxide in blood. Each Hb molecule is able to bind four oxygen molecules. Hb consists of a variety of subfractions and derivatives. Among this heterogeneous group of hemoglobins HbA1c is one of the glycosylated hemoglobins, a subfraction formed by the attachment of various sugars to the Hb molecule. HbA1c is formed in two steps by the non-enzymatic reaction of glucose with the N-terminal amino group of the β-chain of normal adult Hb (HbA). The first step is reversible and yields labile HbA1c. This slowly rearranges in the second reaction step to yield stable HbA1c.

In the erythrocytes, the relative amount of HbA converted to stable HbA1c increases with the average concentration of glucose in the blood. The conversion to stable HbA1c is limited by the erythrocyte's life span of approximately 100 to 120 days. As a result, HbA1c reflects the average blood glucose level during the preceding 2 to 3 months. HbA1c is thus suitable to monitor long-term blood glucose control in individuals with diabetes mellitus. More recent glucose levels have a greater influence on the HbA1c level.<sup>1</sup> The approximate relationship between HbA1c and mean blood glucose value during the preceding 2 to 3 months has been analyzed by several studies. A recent study obtained the following correlation:

### DCCT / NGSP Standardization<sup>8</sup>

$$\text{Mean plasma glucose [mmol/L]} = 1.98 \times \text{HbA1c (\%)} - 4.29$$

$$\text{Mean plasma glucose [mg/dL]} = 35.6 \times \text{HbA1c (\%)} - 77.3$$

### IFCC Standardization (recalculated acc. to Ref. 8)

$$\text{Mean plasma glucose [mmol/L]} = 1.73 \times \text{HbA1c (\%)} + 0.20$$

$$\text{Mean plasma glucose [mg/dL]} = 31.2 \times \text{HbA1c (\%)} + 3.51$$

The risk of diabetic complications, such as diabetic nephropathy and retinopathy increases with poor metabolic control. In accordance with its function as an indicator for the mean blood glucose level, HbA1c predicts the risk for the development of diabetic complications in diabetes patients.<sup>3,5</sup>

For routine clinical use, testing every 3 to 4 months is generally sufficient. In certain clinical situations, such as pregnancy diabetes or after a major change in therapy, it may be useful to measure HbA1c in 2 to 4 week intervals.<sup>7</sup>

### Test principle<sup>10</sup>

This method uses TTAB\* as the detergent in the hemolyzing reagent to eliminate interference from leukocytes (TTAB does not lyse leukocytes). Sample pretreatment to remove labile HbA1c is not necessary.

All hemoglobin variants which are glycosylated at the β-chain N-terminus and which have antibody-recognizable regions identical to that of HbA1c are determined by this assay. Consequently, in contrast to chromatographic methods, the metabolic state of diabetic patients having hemoglobinopathies or uremia can be determined using this assay.

### Hemoglobin A1c

The HbA1c determination is based on the turbidimetric inhibition immunoassay (TINIA) for hemolyzed whole blood.

- Sample and addition of R1 (buffer/antibody)  
Glycohemoglobin (HbA1c) in the sample reacts with anti-HbA1c antibody to form soluble antigen-antibody complexes. Since the specific HbA1c antibody site is present only once on the HbA1c molecule, complex formation does not take place.
- Addition of R2 (buffer/polyhapten) and start of reaction:  
The polyhapten react with excess anti-HbA1c antibodies to form an insoluble antibody-polyhapten complex which can be determined turbidimetrically.

Hemoglobin  
The hemoglobin concentration is determined in a second channel. Liberated hemoglobin in the hemolyzed sample is converted to a derivative having a characteristic absorption spectrum which is measured bichromatically.

### Reagents – working solutions

#### HbA1c

R1 MES buffer<sup>\*\*</sup>: 0.025 mol/L; TRIS buffer<sup>\*\*\*</sup>: 0.015 mol/L, pH 6.2; HbA1c antibody (ovine serum): ≥ 0.5 mg/mL; stabilizers

R2 MES buffer<sup>\*\*</sup>: 0.025 mol/L; TRIS buffer<sup>\*\*\*</sup>: 0.015 mol/L, pH 6.2; HbA1c polyhapten: ≥ 8 μg/mL; stabilizers

3a-d Hemolysate derived from human blood and ovine blood; TTAB<sup>\*</sup>: 9 g/L; stabilizer

#### Hemoglobin

R1 Phosphate buffer: 0.02 mol/L, pH 7.4; stabilizers

\* TTAB = Tetradecyltrimethylammonium bromide

\*\* MES = 2-morpholinoethane sulfonic acid

\*\*\* TRIS = Tris(hydroxymethyl)-aminomethane

### Precautions and warnings

For in vitro diagnostic use.

# HBA1C II

Tina-quant®  Hemoglobin A1c II

Exercise the normal precautions required for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines. All human material should be considered potentially infectious. All products derived from blood are prepared exclusively from the blood of donors tested individually and shown by FDA-approved methods to be free from HBsAg and antibodies to HCV and HIV.

However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be treated just as carefully as a patient specimen. In the event of exposure the directives of the responsible health authorities should be followed.<sup>11,12</sup>

This kit contains components classified as follows according to the European directive 88/379/EEC.

Xi - Irritant (tetradecyltrimethylammonium bromide, TTAB), R36/38; S24/25.

Irritating to eyes and skin. Avoid contact with skin and eyes.

Contact 'phone: all countries: +49-621-7590, USA: +1-800-428-2336

Safety data sheet available for professional user on request.

## Reagent handling

HbA1c

R1: Ready for use

R2: Ready for use

Hemoglobin

R1 (bottle 4): Ready for use

Calibrators 3a-d:

Carefully open each bottle, avoiding the loss of lyophilizate. Pipette exactly 2.0 mL of distilled or deionized water into each bottle. Mix gently, avoiding the formation of foam. Let stand 30 minutes before use.

## Storage and stability

Unopened kit components: up to the expiration date at 2-8°C

HbA1c

R1: 28 days opened and refrigerated on the analyzer

R2: 28 days opened and refrigerated on the analyzer

Hemoglobin

R1: 28 days opened and refrigerated on the analyzer

Calibrators: 8 hours at 15-25°C

2 days at 2-8°C

3 months at (-15) - (-25)°C

Freeze the calibrators immediately after reconstitution if they are to be used again. Freeze in discrete aliquots if necessary.

## Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

Capillary, EDTA- or heparinized blood

Stability:<sup>13</sup> 3 days at 15-25°C

7 days at 2-8°C

6 months at (-15) - (-25)°C (freeze once only)

## Sample preparation

Capillary blood

Obtain blood sample with a disposable capillary tube, and place the filled tube in a test tube with 100-times as much hemolyzing reagent (e.g. 20 µL capillary tube in 2000 µL hemolyzing reagent). Close the test tube and wash the blood out of the capillary tube by shaking vigorously. Then mix the solution by either swirling gently or by using a vibration mixer. The hemolysate can be used after the foam has subsided and the solution has changed color from red to brown-green (approx. 1-2 min).

Stability of the hemolysate: 24 hours at 2-8°C

6 months at (-15) - (-25)°C

Bring the hemolyzing reagent to room temperature prior to use.

EDTA- or heparinized blood

Immediately prior to pipetting the sample, mix the sample thoroughly to obtain uniform distribution of the erythrocytes. Avoid the formation of foam.

Pipette into a test tube or sample cup:

Hemolyzing reagent	1000 µL
EDTA- or heparinized blood	10 µL

Bring the hemolyzing reagent to room temperature prior to use.

Mix using a vibration mixer or by gentle swirling. Avoid the formation of foam. The hemolysate can be used after the solution has changed color from red to brownish-green (approx. 1-2 min).

Stability of the hemolysate: 4 hours at 15-25°C

24 hours at 2-8°C

6 months at (-15) - (-25)°C

## Materials provided

• See "Reagents-working solutions" section for reagents and calibrators.

## Materials required (but not provided)

- Controls: Precinorm HbA1c, Cat. No. 11488422 (USA # 1488422); Precipath HbA1c, Cat. No. 11488449 (USA # 1488449).
- Hemolyzing reagent for HbA1c, Cat. No. 11488457 (USA # 1488457)
- 0.9% NaCl
- General laboratory equipment

## Assay

For optimal performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

## Calibration

Traceability:

This method has been standardized against the approved IFCC reference method for the measurement of HbA1c in human blood<sup>14,15</sup> and can be transferred to results traceable to DCCT/NGSP by calculation.

HbA1c

S1: 0.9% NaCl

S2-5: Calibrators 3a-d

Enter the lot-specific calibrator values to two decimal places.

Hemoglobin

S1: 0.9% NaCl

S2: Calibrator 3a

Enter the lot-specific calibrator values to one decimal place.

The calibrators do not require pretreatment with the hemolyzing reagent for HbA1c.

Calibration frequency

Full calibration is recommended

- after reagent lot change
- after cuvette change
- as required following quality control procedures

Calibration verification: Not necessary.

## Quality control

For quality control, use control materials as listed above. Other suitable control material can be used in addition.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Precinorm HbA1c and Precipath HbA1c do not require pretreatment with the hemolyzing reagent for HbA1c.

## Calculation

Protocol 1 according to IFCC

Calculation of the hemoglobin A1c concentration in percent is made via the HbA1c/Hb quotient, e.g.

$$\% \text{ HbA1c} = \frac{\text{HbA1c [g/dL]}}{\text{Hb [g/dL]}} \times 100$$

Protocol 2 according to DCCT/NGSP

Calculation of the hemoglobin A1c concentration in percent is made via the correction formula:

$$\% \text{ HbA1c} = (87.6 \times \frac{\text{HbA1c [g/dL]}}{\text{Hb [g/dL]}}) + 2.27$$

# HbA1c II

Tina-quant®  Hemoglobin A1c II

For automatic calculation, this equation must be entered under PROGRAM "Calculated TEST" on the Roche/Hitachi 717, "Test Calculation" on the Roche/Hitachi 902, "TEST CALCULATION 04-03" on the Roche/Hitachi 911, "Main/Utility, Calcu Test, Calc" on the Roche/Hitachi 912/917 and "Utility, Calc Test" on the Roche/Hitachi MOD P.

## Limitations – interference<sup>16–22</sup>

- For diagnostic purposes, %HbA1c values should be used in conjunction with information from other diagnostic procedures and clinical evaluations.
- The test is designed only for accurate and precise measurement of %HbA1c. The individual results for total Hb and HbA1c concentration should not be reported.
- The test is not intended for the diagnosis of diabetes mellitus or for judging day-to-day glucose control and should not be used to replace daily home testing of urine or blood glucose.
- Any cause of shortened erythrocyte survival will reduce exposure of erythrocytes to glucose with a consequent decrease in %HbA1c values even though the time-averaged blood glucose level may be elevated. Causes of shortened erythrocyte lifetime might be hemolytic anemia or other hemolytic diseases, homozygous sickle cell trait, pregnancy, recent significant or chronic blood loss, etc. Caution should be used when interpreting the HbA1c results from patients with these conditions.

Criterion: Recovery within  $\pm 10\%$  of initial value.

Icterus: No significant interference up to 50 mg/dL or 855  $\mu\text{mol/L}$  (conjugated and unconjugated bilirubin).

Lipemia (Intralipid): No significant interference up to a triglycerides concentration of 800 mg/dL or 9.12 mmol/L. There is poor correlation between turbidity and triglycerides concentration.

Rheumatoid factor < 750 IU/mL and ascorbic acid < 50 mg/dL (2.84 mmol/L) do not interfere.

In very rare cases gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

## Measuring range<sup>16</sup>

The measuring range for HbA1c lies between 0.3 g/dL (0.124 mmol/L) and the concentration of the highest standard. This corresponds to a measuring range of 2–16% hemoglobin A1c at a typical hemoglobin concentration of 15 g/dL or 9.3 mmol/L (IFCC values; corresponding NGSP/DCCT values 4–16%). As the concentration of the highest standard is lot-specific, this should - where appropriate - be taken into account in the instrument settings for the upper limit of the measuring range.

If the HbA1c concentration is above that of the highest standard, dilute the hemolysate 1 + 1 (or the original sample 1 + 200) with hemolyzing reagent and repeat the two determinations (HbA1c and Hb). No conversion factor is required for the % HbA1c result.

If the HbA1c concentration is below 0.3 g/dL, then the original sample must be diluted 1 + 50 with hemolyzing reagent and the two determinations (HbA1c and Hb) repeated. No conversion factor is required for the %HbA1c result.

Measuring range for hemoglobin: 6–30 g/dL (3.73–18.6 mmol/L).

## Expected values

Metabolically healthy patients.

Protocol 1 (according to IFCC):  
2.9–4.2% HbA1c<sup>23</sup>

Protocol 2 (according to DCCT/NGSP): 4.8–5.9 % HbA1c.<sup>23</sup>

HbA1c levels above the established reference range are an indication of hyperglycemia during the preceding 2 to 3 months or longer. HbA1c levels may reach 20% or higher in poorly controlled diabetes. Therapeutic action is suggested at levels above 8%. Diabetes patients with HbA1c levels below 7% meet the goal of the American Diabetes Association.<sup>19,20</sup>

HbA1c levels below the established reference range may indicate recent episodes of hypoglycemia, the presence of Hb variants, or shortened lifetime of erythrocytes.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range.

## Specific performance data<sup>16</sup>

The data determined using a Roche/Hitachi system are given below. Results obtained in individual laboratories may differ.

## Precision

Reproducibility was determined using controls and pooled whole blood in an internal protocol (n = 21). The following results were obtained.

Hb	Within run			Between run		
	Mean	CV		Mean	CV	
Sample	g/dL	mmol/L	%	g/dL	mmol/L	%
Precinorm HbA1c	13.5	8.38	0.4	14.0	8.69	1.4
Precipath HbA1c	13.6	8.45	0.5	14.1	8.76	1.4
Whole blood	18.7	11.6	0.6	10.6	6.58	0.9

HbA1c	Within run			Between run		
	Mean	CV		Mean	CV	
Sample	g/dL	mmol/L	%	g/dL	mmol/L	%
Precinorm HbA1c	0.56	0.35	2.3	0.61	0.38	4.3
Precipath HbA1c	1.42	0.88	1.5	1.53	0.95	3.5
Whole blood	1.71	1.06	2.1	0.89	0.55	3.1

%HbA1c	Within run			Between run		
	Mean	CV		Mean	CV	
Sample	g/dL	mmol/L	%	g/dL	mmol/L	%
Precinorm HbA1c	4.1	2.55	2.2	4.4	2.73	4.4
Precipath HbA1c	10.5	6.52	1.8	10.9	6.77	3.0
Whole blood	9.2	5.71	2.3	8.4	5.22	3.2

## Analytical sensitivity (lower detection limit)

HbA1c: 0.2 g/dL (0.124 mmol/L)

Hemoglobin: 0.3 g/dL (0.186 mmol/L)

The detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying three standard deviations above that of the lowest standard (standard 1 + 3 SD, within-run precision, n = 21).

## Method comparison

Comparison of the % hemoglobin A1c determination using Roche Tina-quant  HbA1c II with the methods listed below gave the following correlations:

(x) = IFCC reference method

(y) = Tina-quant  HbA1c II, Roche/Hitachi 917, IFCC values

Passing/Bablok<sup>24</sup>

$y = 1.009x - 0.03$

$\tau = 0.950$

SD (md 95) = 0.393

Linear regression

$y = 0.979x + 0.14$

$r = 0.955$

$Sy.x = 0.190$

Number of samples measured: 64

The sample concentrations were between 3.3 and 13.5%.

(x) = Tina-quant  HbA1c II, Roche/Hitachi 917, IFCC values

(y) = Integra 700 hemolysate application, IFCC values

Passing/Bablok<sup>24</sup>

$y = 1.000x + 0.08$

$\tau = 0.951$

SD (md 95) = 0.349

Linear regression

$y = 0.984x + 0.20$

$r = 0.996$

$Sy.x = 0.18$

Number of samples measured: 62

The sample concentrations were between 3.2 and 13.9%.

# HbA1c II

Tina-quant®  Hemoglobin A1c II

(x) = Tina-quant  HbA1c II, Roche/Hitachi 917, DCCT/NGSP values  
(y) = Integra 700 hemolysate application, DCCT/NGSP values

Passing/Bablok <sup>24</sup>	Linear regression
y = 1.000 x + 0.07	y = 0.984 x + 0.21
τ = 0.951	r = 0.996
SD (md 95) = 0.306	Sy,x = 0.16

Number of samples measured: 62  
The sample concentrations were between 5.0 and 14.5%.

## Specificity<sup>16</sup>

No cross-reactivity with HbA0, HbA1a, HbA1b, acetylated hemoglobin, carbamylated hemoglobin, glycated albumin, labile HbA1c and HbA1d was observed for the anti-HbA1c antibodies used in this kit.

## References

- Goldstein DE, Little RR, Lorenz RA, Malone JL, Nathan D, Peterson CM. Tests of glycemia in diabetes. *Diabetes Care*. 1995;18:896-909.
- Goldstein DE, Little RR. More than you ever wanted to know (but need to know) about glycohemoglobin testing. *Diabetes Care*. 1994;17:938-939.
- The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl J Med*. 1993; 329:977-986.
- Santiago JV. Lessons from the diabetes control and complications trial. *Diabetes*. 1993;42:1549-1554.
- UK Prospective Diabetes Study (UKPDS) group. Intensive blood glucose control with sulfonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). *Lancet* 1998, 352, 837-853.
- Flückiger R, Mortensen HB. Review: glycated haemoglobins. *J Chromatogr*. 1988;429:279-292.
- Goldstein DE, Little RR, Wiedmeyer HM, England JD, McKenzie EM. Glycated hemoglobin: methodologies and clinical applications. *Clin Chem*. 1986;32:B64-B70.
- Rohlfing CL, Wiedmeyer HM, Little RR, England JD, Tennill A, Goldstein DE. Defining the relationship between plasma glucose and HbA1c. *Diabetes Care* 2002, 25, 275-278.
- Bunn HF, Gabbay KH, Gallop PM. The glycosylation of hemoglobin: relevance to diabetes mellitus. *Science*. 1978; 200:21-27.
- Karl J et al. Development and Standardization of a New Immunturbidimetric HbA1c Assay. *Klin Lab* 1993;39:991-996.
- Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR 1910.1030). *Federal Register*. July 1, 2001;17:260-273.
- Directive 2000/54/EC). *Official Journal of the European Communities* No. L262, from October 17, 2000.
- Guder WG, Narayanan S, Wisser H, Zawta B. The Quality of Diagnostic Samples. Brochure in: *Samples: From the Patient to the Laboratory*, 2<sup>nd</sup> edition. Darmstadt: GIT Verlag 2001.
- Kobold U, Jeppsson JO, Duellfer T, Finke A, Hoelzel W, Miedema K. Candidate reference methods for hemoglobin A1c based on peptide mapping. *Clin Chem* 1997, 43, 1944-1951.
- Jeppsson JO, Kobold U, Finke A, Hoelzel W, Hoshino T, Miedema K, Mosca A, Mauri P, Paroni R, Thienpont L, Umemoto M, Weykamp C. Approved IFCC reference method for the measurement of HbA1c in human blood. *Clin Chem Lab Med* 2002, 40, 78-89.
- Data on file at Roche Diagnostics.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986, 32, 470-474.
- Greiling H, Gressner AM, eds. *Lehrbuch der Klinischen Chemie und Pathobiochemie*, 3<sup>rd</sup> ed. Stuttgart/New York: Schattauer, 1995:279.
- American Diabetes Association. Standards of Medical Care for patients with diabetes mellitus. *Diabetes Care [Suppl]* 18/1, 1995;8-15.
- Sacks BW, Bruns DE, Goldstein DE, Maclaren NK, McDonald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. *Clin Chem* 2002,48, 436-472.
- Martina WV, Martijn EG, van der Molen M, Schermer JG, Muskiet FAJ. β-N-terminal glycohemoglobins in subjects with common

hemoglobinopathies: relation with fructosamine and mean erythrocyte age. *Clin Chem*. 1993;39:2259-2265.

- Weykamp CW, Penders TJ, Muskiet FAJ, van der Slik W. Influence of hemoglobin variants and derivatives on glycohemoglobin determinations, as investigated by 102 laboratories using 16 methods. *Clin Chem*. 1993;39:1717-1723.
- Junge W, Wilke B, Halabi A et al. Determination of reference intervals in adults for hemoglobin A1c (HbA1c). Poster presentation 18th International Diabetes Federation Congress, Paris, 2003.
- Bablok W et al. A General Regression Procedure for Method Transformation. *J Clin Chem Clin Biochem* 1988;26:783-790.

## Instrument settings

If assays are not performed in batch mode (test-oriented), the determination of Hb must be performed immediately before the determination of HbA1c to avoid possible interference by Creatinine-Jaffé. This can be achieved on the Roche/Hitachi 717/747/902 analyzers by assigning these tests to two consecutive channels (e.g. Hb in channel 14 and HbA1c in channel 15). When running large batches of HbA1c it may be necessary to implement a special washing procedure. Please contact your local Roche representative for further details.

## US users

Refer to application sheet for additional operating information.

## Roche/Hitachi 717

If the HbA1c II assay is run in batch mode (test-oriented), rinse the cuvette once weekly (cell-wash) with 0.1 N sodium hydroxide solution.

## Roche/Hitachi 904, 911, 912, 917 and Modular users

Enter the application parameters via the application diskette, settings sheet or barcode sheet, as appropriate.

## Roche/Hitachi 717

Temperature: 37°C

PROGRAM 2 CHEMISTRY PARAMETERS		
	Hemoglobin A1c	Hemoglobin
TEST	[HBA1C]	[HB]
ASSAY CODE	[2(2POINT)]- [24]-[50]	[1(1POINT)]- [24]-[0]
SAMPLE VOLUME (µl)	[10] - [10]	[20] - [20]
R1 VOLUME	[250] - [50] - [NO]	[230] - [50] - [NO]
R1 VOLUME	[50] - [20] - [NO]	[0] - [20] - [NO]
WAVELENGTH	[700] - [340]	[660] - [570]
CALIB. METHOD	[3] - [1] - [5]	[1] - [0] - [0]
STD (1) CONC RACK POS	[0.00] - [ ]	[0.0] - [ ]
STD (2) CONC RACK POS	[ ] - [ ]	[ ] - [ ]
STD (3) CONC RACK POS	[ ] - [ ]	[0] - [0]
STD (4) CONC RACK POS	[ ] - [ ]	[0] - [0]
STD (5) CONC RACK POS	[ ] - [ ]	[0] - [0]
STD (6) CONC RACK POS	[0] - [0]	[0] - [0]
SD LIMIT	[200]	[0.1]
DUPLICATE LIMIT	[350]	[50]
SENSITIVITY LIMIT	[2000]	[250]
ABS. LIMIT (INC/DEC)	[0] - [INC]	[0] - [INC]
PROZONE LIMIT	[32000] - [UPPER]	[0] - [UPPER]
EXPECTED VALUES	[ ] - [ ]	[ ] - [ ]
PANIC VALUE	[ ] - [ ]	[ ] - [ ]
INSTRUMENT FACTOR	[1.00]	[1.00]

— Data entered by the operator

# HBA1C II

Tina-quant®  Hemoglobin A1c II

Roche/Hitachi 747

Temperature: 37°C

Hemoglobin A1c

Hemoglobin

PROGRAM 4.2 CHEMISTRY PARAMETERS			
TEST	[HBA1C]	-	
ASSAY CODE	[2(2POINT)]	-[22]-[50]	
WAVELENGTH	[700(SUB)]	-[340 (MAIN)]	
	<b>SERUM</b>	<b>URINE</b>	
SAMPLE VOLUME (µl)	[10]	-[5]	[ ] -[ ]
EXPECTED VALUE	[ ]	-[ ]	[ ] -[ ]
PANIC VALUE	[ ]	-[ ]	[ ] -[ ]
ABS LIMIT (INC/DEC)	[0]	-[INCREASE]	[ ] -[ ]
PROZONE LIMIT	[32000]	-[UPPER]	[ ] -[ ]
	<b>R1</b>	<b>R2</b>	
R1/R2 VOLUME (µl)	[250]	[50]	
R1/R2 DUMMY INTERVAL	[0]	[0]	
DILUTION VOLUME (µl)	[0]		
CALIB. METHOD	[NONLINEAR]		
POINTS	[5]		
STD 1 CONC RACK POS	[ ]	-[ ]-[ ]	
STD 2 CONC RACK POS	[ ]	-[ ]-[ ]	
STD 3 CONC RACK POS	[ ]	-[ ]-[ ]	
STD 4 CONC RACK POS	[ ]	-[ ]-[ ]	
STD 5 CONC RACK POS	[ ]	-[ ]-[ ]	
STD 6 CONC RACK POS	[ ]	-[ ]-[ ]	
SD LIMIT	[200]		
DUPLICATE LIMIT	[350]		
SENSITIVITY LIMIT	[2000]		
STD 1 ABS. LEVEL	[ ]	-[ ]	
INSTRUMENT FACTOR	[1.00]		

\_\_\_ Data entered by the operator


PROGRAM 4.2 CHEMISTRY PARAMETERS			
TEST	[HB]	-	
ASSAY CODE	[1(1POINT)]	-[22]-[0]	
WAVELENGTH	[660(SUB)]	-[570 (MAIN)]	
	<b>SERUM</b>	<b>URINE</b>	
SAMPLE VOLUME (µl)	[20]	-[20]	[ ] -[ ]
EXPECTED VALUE	[ ]	-[ ]	[ ] -[ ]
PANIC VALUE	[ ]	-[ ]	[ ] -[ ]
ABS LIMIT (INC/DEC)	[0]	-[INCREASE]	[ ] -[ ]
PROZONE LIMIT	[32000]	-[UPPER]	[ ] -[ ]
	<b>R1</b>	<b>R2</b>	
R1/R2 VOLUME (µl)	[250]	[0]	
R1/R2 DUMMY INTERVAL	[0]	[0]	
DILUTION VOLUME (µl)	[0]		
CALIB. METHOD	[LINEAR]		
POINTS	[0]		
STD 1 CONC RACK POS	[ ]	-[ ]-[ ]	
STD 2 CONC RACK POS	[ ]	-[ ]-[ ]	
STD 3 CONC RACK POS	[ ]	-[ ]-[ ]	
STD 4 CONC RACK POS	[ ]	-[ ]-[ ]	
STD 5 CONC RACK POS	[ ]	-[ ]-[ ]	
STD 6 CONC RACK POS	[ ]	-[ ]-[ ]	
SD LIMIT	[0.1]		
DUPLICATE LIMIT	[50]		
SENSITIVITY LIMIT	[250]		
STD 1 ABS. LEVEL	[ ]	-[ ]	
INSTRUMENT FACTOR	[1.00]		

# HBA1C II

Tina-quant® Hemoglobin A1c II

Roche/Hitachi 902

Tina-quant, Precinorm and Precipath are trademarks of a member of the Roche Group.  
Intralipid is a trademark of KabiPharmacia, Inc.  
Significant additions or changes (where present) are indicated by a change bar in the margin.  
© 2004 Roche Diagnostics

 Roche Diagnostics GmbH, D-68298 Mannheim  
for USA: US Distributor:  
Roche Diagnostics, Indianapolis, IN  
US Customer Technical Support 1-800-428-2336

No.	<Chemistry>	HbA1c	Hb
1	Test Name	HbA1c	Hb
2	Assay Code (Mthd)	2 Point End	1 Point
3	Assay Code (2. Test)	0	0
4	Reaction Time	10	10
5	Assay Point 1	17	17
6	Assay Point 2	35	0
7	Assay Point 3	0	0
8	Assay Point 4	0	0
9	Wavelength (SUB)	700	660
10	Wavelength (MAIN)	340	570
11	Sample Volume	10.0	20.0
12	R1 Volume	250	230
13	R1 Pos.	.....	.....
14	R1 Bottle Size	Large	Large
15	R2 Volume	0	0
16	R2 Pos.	0	0
17	R2 Bottle Size	Small	Small
18	R3 Volume	50	0
19	R3 Pos.	.....	0
20	R3 Bottle Size	Small	Small
21	Calib. Type (Type)	Logit-Log(4P)	Linear
22	Calib. Type (Wght)	0	0
23	Calib. Conc. 1	0.00	0.00
24	Calib. Pos. 1	.....	.....
25	Calib. Conc. 2	.....	.....
26	Calib. Pos. 2	.....	.....
27	Calib. Conc. 3	.....	0
28	Calib. Pos. 3	.....	0
29	Calib. Conc. 4	.....	0
30	Calib. Pos. 4	.....	0
31	Calib. Conc. 5	.....	0
32	Calib. Pos. 5	.....	0
33	Calib. Conc. 6	.....	0
34	Calib. Pos. 6	.....	0
35	S1 ABS	0	0
36	K Factor	10000	10000
37	K2 Factor	10000	10000
38	K3 Factor	10000	10000
39	K4 Factor	10000	10000
40	K5 Factor	10000	10000
41	A Factor	0	0
42	B Factor	0	0
43	C Factor	0	0
44	SD Limit	200	0.1
45	Duplicate Limit	350	50
46	Sens. Limit	3000	250
47	S1 Abs. Limit (L)	-32000	-32000
48	S1 Abs. Limit (H)	32000	32000
49	Abs. Limit	0	0
50	Abs. Limit (D/I)	Increase	Increase
51	Prozone Limit	32000	32000
52	Proz. Limit (Upp/Low)	Upper	Upper
53	Prozone (Endpoint)	35	35
54	Expect. Value (L)	.....	.....
55	Expect. Value (H)	.....	.....
56	Instr. Fact. (a)	1.0	1.0
57	Instr. Fact. (b)	0.0	0.0
58	Key setting	.....	.....

..... Data entered by the operator

For further information please refer to the appropriate operator manual for the analyzer concerned, the respective application sheets, product information and the package inserts of all necessary components.

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## Literatur

- Goldstein DE, Little RR, Lorenz RA, Malone JL, Nathan D, Peterson CM. Tests of glycemia in diabetes. *Diabetes Care*. 1995;18:896-909.
- Goldstein DE, Little RR. More than you ever wanted to know (but need to know) about glycohemoglobin testing. *Diabetes Care*. 1994;17:938-939.
- The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl J Med*. 1993; 329:977-986.
- Santiago JV. Lessons from the diabetes control and complications trial. *Diabetes*. 1993;42:1549-1554.
- UK Prospective Diabetes Study (UKPDS) group. Intensive blood glucose control with sulfonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). *Lancet* 1998, 352, 837-853.
- Flückiger R, Mortensen HB. Review: glycated haemoglobins. *J Chromatogr*. 1988;429:279-292.
- Goldstein DE, Little RR, Wiedmeyer HM, England JD, McKenzie EM. Glycated hemoglobin: methodologies and clinical applications. *Clin Chem*. 1986;32:B64-B70.
- Rohlfing CL, Wiedmeyer HM, Little RR, England JD, Tennill A, Goldstein DE. Defining the relationship between plasma glucose and HbA1c. *Diabetes Care* 2002, 25, 275-278.

- Bunn HF, Gabbay KH, Gallop PM. The glycosylation of hemoglobin: relevance to diabetes mellitus. *Science*. 1978; 200:21-27.
- Karl J et al. Development and Standardization of a New Immunospectrometric HbA1c Assay. *Klin lab* 1993;39:991-996.
- Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR 1910.1030). *Federal Register*. July 1, 2001;17:260-273.
- Richtlinien des Rates (2000/54/EG). *Abt. d. Europ. Gem. Nr. L262 vom 17.10.2000*
- Guder WG, Narayanan S, Wisser H, Zawta B. The Quality of Diagnostic Samples. Broschüre in: *Samples: From the Patient to the Laboratory*, 2<sup>nd</sup> edition. Darmstadt: GIT Verlag 2001.
- Kobold U, Jeppsson JO, Duelffer T, Finke A, Hoelzel W, Miedema K. Candidate reference methods for hemoglobin A1c based on peptide mapping. *Clin Chem* 1997, 43, 1944-1951.
- Jeppsson JO, Kobold U, Finke A, Hoelzel W, Hoshino T, Miedema K, Mosca A, Mauri P, Paroni R, Thienpont L, Umemoto M, Weykamp C. Approved IFCC reference method for the measurement of HbA1c in human blood. *Clin Chem Lab Med* 2002, 40, 78-89.
- Dokumentation Roche Diagnostics.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986, 32, 470-474.
- Greiling H, Gressner AM, eds. *Lehrbuch der Klinischen Chemie und Pathobiochemie*, 3<sup>rd</sup> ed. Stuttgart/New York: Schattauer, 1995:279.
- American Diabetes Association. Standards of Medical Care for patients with diabetes mellitus. *Diabetes Care [Suppl]* 18/1, 1995:8-15.
- Sacks BW, Brunts DE, Goldstein DE, Maclaren NK, McDonald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. *Clin Chem* 2002,48, 436-472.
- Martina WV, Martijn EG, van der Molen M, Schermer JG, Muskiet FAJ. β-N-terminal glycohemoglobins in subjects with common hemoglobinopathies: relation with fructosamine and mean erythrocyte age. *Clin Chem*. 1993;39:2259-2265.
- Weykamp CW, Penders TJ, Muskiet FAJ, van der Slik W. Influence of hemoglobin variants and derivatives on glycohemoglobin determinations, as investigated by 102 laboratories using 16 methods. *Clin Chem*. 1993;39:1717-1723.
- Junge W, Wilke B, Halabi A et al. Determination of reference intervals in adults for hemoglobin A1c (HbA1c). Poster presentation 18<sup>th</sup> International Diabetes Federation Congress, Paris 2003.
- Bablok W et al. A General Regression Procedure for Method Transformation. *J Clin Chem Clin Biochem* 1988;26:783-790.

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## Bibliographie

- Goldstein DE, Little RR, Lorenz RA, Malone JL, Nathan D, Peterson CM. Tests of glycemia in diabetes. *Diabetes Care*. 1995;18:896-909.
- Goldstein DE, Little RR. More than you ever wanted to know (but need to know) about glycohemoglobin testing. *Diabetes Care*. 1994;17:938-939.
- The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl J Med*. 1993;329:977-986.
- Santiago JV. Lessons from the diabetes control and complications trial. *Diabetes*. 1993;42:1549-1554.
- UK Prospective Diabetes Study (UKPDS) group. Intensive blood glucose control with sulfonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). *Lancet* 1998, 352, 837-853.

- Flückiger R, Mortensen HB. Review: glycosylated haemoglobins. *J Chromatogr*. 1988;429:279-292.
- Goldstein DE, Little RR, Wiedmeyer HM, England JD, McKenzie EM. Glycosylated hemoglobin: methodologies and clinical applications. *Clin Chem*. 1986;32:B64-B70.
- Rohlfing CL, Wiedmeyer HM, Little RR, England JD, Tennill A, Goldstein DE. Defining the relationship between plasma glucose and HbA1c. *Diabetes Care* 2002;27:275-278.
- Bunn HF, Gabbay KH, Gallop PM. The glycosylation of hemoglobin: relevance to diabetes mellitus. *Science*. 1978;200:21-27.
- Karl J et al. Development and Standardization of a New Immunosorbent HbA1c Assay. *Klin Lab* 1993;39:991-996.
- Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR 1910.1030). Federal Register. July 1, 2001;17:260-273.
- Directive du Conseil de l'Europe (2000/54/CEE). Journal Officiel de la Communauté Européenne n° L374 du 17.09.2000.
- Guder WG, Narayanan S, Wisser H, Zawta B. The Quality of Diagnostic Samples. Brochure dans: *Samples: From the Patient to the Laboratory*, 2<sup>e</sup> édition. Darmstadt: GIT Verlag, 2001.
- Kobold U, Jeppsson JO, Duelffer T, Finke A, Hoelzel W, Miedema K. Candidate reference methods for hemoglobin A1c based on peptide mapping. *Clin Chem* 1997;43:1944-1951.
- Jeppsson JO, Kobold U, Finke A, Hoelzel W, Hoshino T, Miedema K, Mosca A, Mauri P, Paroni R, Thienpont L, Umemoto M, Weykamp C. Approved IFCC reference method for the measurement of HbA1c in human blood. *Clin Chem Lab Med*. 2002;40:78-89.
- Documentation de Roche Diagnostics
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-474.
- Greiling H, Gressner AM, eds. *Lehrbuch der Klinischen Chemie und Pathobiochemie*, 3<sup>e</sup> édition. Stuttgart/New York: Schattauer 1995:279.
- American Diabetes Association. Standards of Medical Care for patients with diabetes mellitus. *Diabetes Care [Suppl]* 18/1, 1995:8-15.
- Sacks DB, Bruns DE, Goldstein DE, Maclaren NK, McDonald JM, Parrott M. Guidelines and Recommendations for Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus. *Clin Chem* 2002;48, 436-472.
- Martina WV, Martijn EG, van der Molen M, Schermer JG, Muskiet FAJ.  $\beta$ -N-terminal glycohemoglobins in subjects with common hemoglobinopathies: relation with fructosamine and mean erythrocyte age. *Clin Chem*. 1993;39:2259-2265.
- Weykamp CW, Penders TJ, Muskiet FAJ, van der Sliik W. Influence of hemoglobin variants and derivatives on glycohemoglobin determinations, as investigated by 102 laboratories using 16 methods. *Clin Chem*. 1993;39:1717-1723.
- Junge W, Wilke B, Halabi A et al. Determination of reference intervals in adults for hemoglobin A1c (HbA1c). Poster presentation 18th International Diabetes Federation Congress, Paris, 2003.
- Bablok W et al. A General Regression Procedure for Method Transformation. *J Clin Chem Clin Biochem* 1988;26:783-790.

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1. Goldstein DE, Little RR, Lorenz RA, Malone JL, Nathan D, Peterson CM. Tests of glycemia in diabetes. *Diabetes Care*. 1995;18:896-909.
2. Goldstein DE, Little RR. More than you ever wanted to know (but need to know) about glycohemoglobin testing. *Diabetes Care*. 1994;17:938-939.
3. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl J Med*. 1993;329:977-986.
4. Santiago JV. Lessons from the diabetes control and complications trial. *Diabetes*. 1993;42:1549-1554.
5. UK Prospective Diabetes Study (UKPDS) group. Intensive blood glucose control with sulfonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). *Lancet* 1998;352:837-853.
6. Flückiger R, Mortensen HB. Review: glycated haemoglobins. *J Chromatogr*. 1988;429:279-292.
7. Goldstein DE, Little RR, Wiedmeyer HM, England JD, McKenzie EM. Glycated hemoglobin: methodologies and clinical applications. *Clin Chem*. 1986;32:B64-B70.

8. Rohlfing CL, Wiedmeyer HM, Little RR, England JD, Tennill A, Goldstein DE. Defining the relationship between plasma glucose and HbA1c. *Diabetes Care* 2002; 25, 275-278.
9. Bunn HF, Gabbay KH, Gallop PM. The glycosylation of hemoglobin: relevance to diabetes mellitus. *Science*. 1978; 200:21-27.
10. Karl J et al. Development and Standardization of a New Immunoturbidimetric HbA1c Assay. *Klin Lab* 1993;39:991-996.
11. Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR 1910.1030). *Federal Register*. July 1, 2001; 17:260-273.
12. Directiva del Consejo (2000/54/CE). *Diario oficial de las Comunidades Europeas*, N° L262 del 17/10/2000.
13. Guder WG, Narayanan S, Wissner H, Zawta B. The Quality of Diagnostic Samples. Folleto en: *Samples: From the patient to the Laboratory*, 2<sup>nd</sup> edition. Darmstadt: GIT Verlag 2001.
14. Kobold U, Jeppsson JO, Duelffer T, Finke A, Hoelzel W, Miedema K. Candidate reference methods for hemoglobin A1c based on peptide mapping. *Clin Chem* 1997;43:1944-1951.
15. Jeppsson JO, Kobold U, Finke A, Hoelzel W, Hoshino T, Miedema K, Mosca A, Mauri P, Paroni R, Thienpont L, Umemoto M, Weykamp C. Approved IFCC reference method for the measurement of HbA1c in human blood. *Clin Chem Lab Med* 2002; 40, 78-89.
16. Documentación de Roche Diagnostics.
17. Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-474.
18. Greiling H, Gressner AM, eds. *Lehrbuch der Klinischen Chemie und Pathobiochemie*, 3<sup>rd</sup> ed. Stuttgart/New York: Schattauer, 1995:279.
19. American Diabetes Association. Standards of Medical Care for patients with diabetes mellitus. *Diabetes Care* [Suppl] 18/1, 1995;8-15.
20. Sacks BW, Bruns DE, Goldstein DE, Maclaren NK, McDonald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. *Clin Chem* 2002;48, 436-472.
21. Martina WV, Martijn EG, van der Molen M, Schermer JG, Muskiet FAJ. β-N-terminal glycohemoglobins in subjects with common hemoglobinopathies: relation with fructosamine and mean erythrocyte age. *Clin Chem*. 1993;39:2259-2265.
22. Weykamp CW, Penders TJ, Muskiet FAJ, van der Slik W. Influence of hemoglobin variants and derivatives on glycohemoglobin determinations, as investigated by 102 laboratories using 16 methods. *Clin Chem*. 1993;39:1717-1723.
23. Junge W, Wilke B, Halabi A et al. Determination of reference intervals in adults for hemoglobin A1c (HbA1c). Poster presentation 18<sup>th</sup> International Diabetes Federation Congress, Paris 2003.
24. Bablok W et al. A General Regression Procedure for Method Transformation. *J Clin Chem Clin Biochem* 1988;26:783-790.

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## Bibliografia

- Goldstein DE, Little RR, Lorenz RA, Malone JI, Nathan D, Peterson CM. Tests of glycemia in diabetes. *Diabetes Care*. 1995;18:896-909.
- Goldstein DE, Little RR. More than you ever wanted to know (but need to know) about glycohemoglobin testing. *Diabetes Care*. 1994;17:938-939.
- The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl J Med*. 1993; 329:977-986.
- Santiago JV. Lessons from the diabetes control and complications trial. *Diabetes*. 1993;42:1549-1554.
- UK Prospective Diabetes Study (UKPDS) group. Intensive blood glucose control with sulfonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). *Lancet* 1998, 352, 837-853.
- Flückiger R, Mortensen HB. Review: glycated haemoglobins. *J Chromatogr*. 1988;429:279-292.
- Goldstein DE, Little RR, Wiedmeyer HM, England JD, McKenzie EM. Glycated hemoglobin: methodologies and clinical applications. *Clin Chem*. 1986;32:B64-B70.
- Rohlfing CL, Wiedmeyer HM, Little RR, England JD, Tennill A, Goldstein DE. Defining the relationship between plasma glucose and HbA1c. *Diabetes Care* 2002, 25, 275-278
- Bunn HF, Gabbay KH, Gallop PM. The glycosylation of hemoglobin: relevance to diabetes mellitus. *Science*. 1978; 200:21-27.
- Karl J et al. Development and Standardization of a New Immunoturbidimetric HbA1c Assay. *Klin Lab* 1993;39:991-996.
- Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR 1910.1030). Federal Register. July 1, 2001;17:260-273.
- Directiva do Conselho (2000/54/CEE). *Jornal Oficial das Comunidades Europeias* N° L262 de 17 de Out., 2000.
- Guder WG, Narayanan S, Wisser H, Zawta B. The Quality of Diagnostic Samples. Brochure in: *Samples: From the Patient to the Laboratory*, 2<sup>nd</sup> edition. Darmstadt: GIT Verlag, 2001.
- Kobold U, Jeppsson JO, Duellfer T, Finke A, Hoelzel W, Miedema K. Candidate reference methods for hemoglobin A1c based on peptide mapping. *Clin Chem* 1997, 43, 1944-1951.
- Jeppsson JO, Kobold U, Finke A, Hoelzel W, Hoshino T, Miedema K, Mosca A, Mauri P, Paroni R, Thienpont L, Umemoto M, Weykamp C. Approved IFCC reference method for the measurement of HbA1c in human blood. *Clin Chem Lab Med* 2002, 40, 78-89.
- Documentação da Roche Diagnostics.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986, 32, 470-474.
- Greiling H, Gressner AM, eds. *Lehrbuch der Klinischen Chemie und Pathobiochemie*, 3<sup>rd</sup> ed. Stuttgart/New York: Schattauer, 1995:279.
- American Diabetes Association. Standards of Medical Care for patients with diabetes mellitus. *Diabetes Care* [Suppl] 18/1, 1995;8-15.
- Sacks BW, Bruns DE, Goldstein DE, Maclaren NK, McDonald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. *Clin Chem* 2002;48, 436-472.
- Martina WV, Martijn EG, van der Molen M, Schermer JG, Muskiet FAJ. β-N-terminal glycohemoglobins in subjects with common hemoglobinopathies: relation with fructosamine and mean erythrocyte age. *Clin Chem*. 1993;39:2259-2265.
- Weykamp CW, Penders TJ, Muskiet FAJ, van der Slik W. Influence of hemoglobin variants and derivatives on glycohemoglobin determinations, as investigated by 102 laboratories using 16 methods. *Clin Chem*. 1993;39:1717-1723.
- Junge W, Wilke B, Halabi A et al. Determination of reference intervals in adults for hemoglobin A1c (HbA1c). Poster presentation 18<sup>th</sup> International Diabetes Federation Congress, Paris 2003.
- Bablok W et al. A General Regression Procedure for Method Transformation. *J Clin Chem Clin Biochem* 1988;26:783-790.

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