INSTRUCTION MANUAL

Cordless Pipette Controller

Cat. No: E4866-0021 (UK Plug)

Cat. No: E4866-1126 (Euro Plug)

www.starlabgroup.com

English	2 – 8
Français	9 – 15
Deutsch	16 – 22
Italiano	23 – 29

English

INTRODUCTION

The StarPet Pro has been carefully designed and manufactured as a lightweight, cordless pipetting aid for use with all plastic/glass pipettes in volumes from 1 to 100ml.

Constructed in recyclable materials, the StarPet Pro operates on two rechargeable, long-life, environmentally friendly, 1.2V nickel metal hydride batteries, and has uniquely formed finger triggers to minimize any high pressure points. These triggers also control the intake and dispense speeds.


A battery warning light (4) is triggered whenever the charge in the batteries is low to let you know the batteries should be recharged. Pipetting with a StarPet Pro can also take place while the unit is connected to the charger (5).

The StarPet Pro has a conveniently located slide switch (3) which allows you to select between three modes of operation:

Slide switch position 'High' (top) 	Suction and dispense at full speed (Blow-out)
Slide switch position 'Low' (middle) 	Suction and dispense at reduced speed (Blow-out)
Slide switch position 'G' (bottom)	Suction at reduced speed, dispense by gravity (To Deliver)

TECHNICAL DATA

Batteries	2 x 1.2V nickel metal hydride
Battery service Life	500–1000 cycles from empty
Charger	75mA/h for 14 hours plug-in; UK, European, USA, Indian, Japanese or Australian versions Installation (over voltage) category II Double insulation indicated on charger by symbol 
Overall weight	180g without stand
Filter (8)	Hydrophobic 25mm dia 0.45µm (fitted), plus two additional spare filters 1 x 0.45µm 1 x 0.20µm
Case material	ASA
Nosecone material (10)	ASA
Pipette holder (9)	Silicone
Pipette types	Plastic or glass 1–100ml

CE 2004/108/EEC Electromagnetic Compatibility
EN 61326-1:2006 Electrical Equipment for Measurement,
Control and Laboratory Use.

CONTENTS

The StarPet Pro kit comprises

- (a) StarPet Pro cordless pipette controller
- (b) Charger; either
 - 220V-240V 50/60Hz 3-pin (UK), 220V-240V 50/60Hz 2-pin (Europe),
 - 100V-120V 50/60Hz 2-pin (USA), 100V-120V 50/60Hz 2-pin (Japan),
 - 220V-240V 50/60Hz 3-pin (India), or 220V-240V 50/60Hz 2-pin (Australia)
- (c) Two spare hydrophobic filters, 1 x 0.45µm and 1 x 0.20µm
- (d) Bench stand / wall bracket
- (e) Instruction manual

SAFETY INSTRUCTIONS

The StarPet Pro has been designed with safety in mind. It is lightweight and balances safely in the hand. The triggers have been carefully designed to eliminate pressure points and control of intake and delivery is, therefore, made easier.

It is recommended when using deleterious liquids that the user wears protective clothing including goggles and gloves. If in doubt ask your safety officer. Please consult and establish appropriate health and safety practices applicable to the liquid being dispensed within your working environment.

WARNING!

- Only use this instrument within the recommended limits of its chemical resistance.
- Never use with fuming liquids which are incompatible with ASA, silicone or polypropylene.
- Never use with highly flammable liquids.
- Observe the reagent manufacturer's safety information.
- Never use the instrument in an atmosphere with danger of explosion.
- Always store on bench stand or wall bracket when not in use.
- Never use force when in use.
- Should the instrument fail to work to specification, stop pipetting and check instructions headed Trouble Shooting (see page 7).
- Only recharge with manufacturer's original charger.
- All repairs must be carried out by an authorised dealer.
- Only use the manufacturer's accessories and spare parts.
- Charger must only be used indoors and protected against moisture.

RECHARGING

Before starting to operate, carefully check mains voltage and recharge batteries.


To charge the StarPet Pro:

- 1 Recharging temperature +10°C to +35°C.
- 2 Plug the charger jack (6) into the base of the StarPet Pro.
- 3 Plug the charger (5) into a suitable wall socket. The red indicator on the charger will light up.
- 4 Place StarPet Pro on bench stand.
- 5 Allow unit to charge for 14 hours from empty to fully charged.
- 6 Once the StarPet Pro is fully charged, the charger automatically switches to a trickle charge and the red indicator will turn off. If the charge to the StarPet Pro is interrupted and then re-started, the charger starts its full 14-hour charging cycle again. This will not damage the battery.
- 7 Only use the charger indoors and protect against moisture.


FITTING THE WALL BRACKET

The StarPet Pro bench stand is supplied ready to use as a wall bracket. When used as a bench stand it is recommended that the wall fitting section is removed.

Bench Stand

Wall Bracket


CLEANING

Before cleaning disconnect from the mains charger. Do not clean the mains charger.

The StarPet Pro can be cleaned with most common laboratory disinfectants. The silicone pipette holder (9) can be autoclaved.

To remove the nosecone and pipette holder turn the nosecone anti-clockwise. When the nosecone is removed the silicone pipette holder can also be disengaged from the filter by gently turning.

The hydrophobic filter (8) will require occasional replacement when wetted by the liquid and whenever suction capacity is slowing.

PIPETTING

To aspirate

- 1 Carefully hold the pipette and insert it securely into the polypropylene nose cone and silicone holder.
- 2 Select mode of operation—High, Low, for blowout pipettes. 'G' for 'to deliver'.
- 3 Depress upper trigger (1) (The harder the pipetting triggers are pressed the faster the aspiration.)

To dispense

- 1 In BLOW-OUT mode, (mode switch on high or low) depress the lower trigger (2) (The harder the pipetting triggers are pressed, the faster the dispense).
- 2 If 'G' gravity mode is chosen, the upper aspirate trigger automatically operates at low speed while the lower trigger opens the pipette to atmosphere, allowing liquid to run out.

WARNING!

Please ensure that no liquid enters the instrument. The StarPet Pro should never be used without the filter (8). The filter also prevents liquid from entering the instrument. Where liquid has entered please clean and replace filter.

TROUBLE SHOOTING

The StarPet Pro is virtually service-free. If the recommendations above do not resolve the problem, please contact your authorised dealer.

PROBLEM	POSSIBLE CAUSE	REMEDY
Suction capacity slow	Filter Blocked, dirty or damaged	Clean or replace membrane filter
	Batteries low	Recharge
Pipette drips	Pipette not securely positioned in silicone pipette holder	Re-attach
	Pipette holder or pipette damaged	Replace

Note

When the meniscus level has been reached liquid should not drip out of the pipette. If the pipette does drip, check that the adaptor assembly and filter are correctly assembled and the pipette is properly located into the silicone pipette holder.

RETURN FOR REPAIR

In the unlikely event of repair, or where the damage to the unit necessitates return to your distributor or the manufacturer, please decontaminate the unit and complete the required decontamination certificate where appropriate. Please also give details of the fault.

VOLTAGES

UK	220V-240V 50/60Hz 3-pin
Europe	220V-240V 50/60Hz 2-pin
USA	100V-120V 50/60Hz 2-pin
Japan	100V-120V 50/60Hz 2-pin
India	220V-240V 50/60Hz 3-pin
Australia	220V-240V 50/60Hz 2-pin

ORDERING DATA

StarPet Pro comes complete with charger, bench stand/ wall bracket with two spare membrane filters.

UK Plug	Cat. No. E4866-0021
European Plug	Cat. No. E4866-1126

Replacement charger	Cat. No. E4266-0050
ASA nosecone (white)	Cat. No. P2131-0002
5 spare 0.20µm filters	Cat. No. E4266-0890
5 spare 0.45µm filters	Cat. No. E4266-0030
Silicone Pipette holder	Cat. No. E4266-0040

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STARLAB (UK), Ltd
STARLAB S.r.l.