EN TPLA REAGENT

INTENDED USE

For the quantitative determination of anti-Treponema *pallidum* antibodies in human serum

SUMMARY

Syphilis a chronic infection is caused by *Treponema pallidum* which can be transmitted congenitally or by sexual contact. It is characterized by episodes of active disease interrupted by periods of latency¹.

There are two types of serological test for syphilis, non-treponemal and treponemal. The most widely used non-treponemal antibody tests for syphilis are the rapid plasma reagin (RPR) and venereal disease research laboratory (VDRL) tests, both measuring antibodies against a cardiolipin-lecithin-cholesterol antigen complex. Treponemal tests measure antibodies to native (Nichols strain) or recombinant *T. pallidum* antigens.¹ This automated assay is based on a latex immunological agglutination test principle.

PRINCIPLE

The polystyrene latex, coated with antigen components derived from *Treponema pallidum* (Nichols strain), is exposed to the test sample under certain conditions, to induce formation of anti-Treponema antibody-latex aggregate. The elevation of turbidity due to the formation of this aggregate (the change in turbidity) relative to the pre-exposure level is measured, to determine the anti-Treponema antibody titer in the test sample.

REAGENTS

Composition

Component	Ingredients	Concentration
Reagent 1	Phosphate buffer	100 mmol/L pH 7.2-7.4
	Bovine serum albumin	
	Stabilizers Sodium Azide	<0.1%
Reagent 2	<i>Treponema pallidum</i> antigen-coated latex	1.5 - 3.5mg/mL
	Phosphate buffer	100 mmol/L pH 7.2-7.5
	Bovine serum albumin Sodium Azide	<0.1%

Precautions and Warnings

1. For In Vitro Diagnostic Use.
2. Do not use the reagents beyond the expiration date printed on the label.
3. **Warning** : All specimens used in the test should be considered potentially infectious. Universal precautions as they apply to your facility should be used for handling and disposal of materials during and after testing.²
4. TPLA Reagents must be used with the TPLA Calibrator Set.
5. **Caution** : Avoid freezing reagents.
6. **Caution** : Reagents 1 and 2 contain <0.1% sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form potentially explosive metal azide buildup. Flush with copious amounts of water when discarding material.
7. **Caution** : Do not pool reagents within a kit or between reagent kits.
8. Disposal of all waste material should be in accordance with local guidelines.

Preparation

Reagent 1 : Liquid, ready to use.

Reagent 2 : Liquid, ready to use.

Invert to mix before use. Avoid the formation of foam.

Storage and Stability

Unopened reagent is stable until the expiration date shown on the label when stored at 2 - 8°C.

Once opened, the reagent is stable up to 4 weeks at 2 - 8°C.

DO NOT FREEZE.

Onboard stability

Reagents are stable open on the Hitachi 7170 analyzer for 4 weeks at 2 - 8 °C.

Indications of Deterioration

Presence of turbidity in reagent 1 or microbial growth in either reagents may indicate deterioration.

Inability to recover control values.

SPECIMEN COLLECTION AND PREPARATION

Serum is the recommended collection media. Use standard sample collection and preparation methods.³


If not analyzed promptly, serum specimens may be stored at 2 - 8°C for 1 week or at 15 - 25°C for 1 day.⁴ If specimens need to be stored for more than 1 week, they may be preserved at -20°C or below for up to 4 weeks. ⁴

If the samples have been frozen, centrifuge at 15000xg for 10 minutes before measurement. Samples may be frozen and thawed once. ⁴

PROCEDURE

Assay

Below is a general example of the TPLA assay procedure for an automated analyzer. All analyzer applications should be validated.

Materials Provided

TPLA Reagents 1 and 2 are required for the measurement of anti-Treponema *pallidum* antibodies. The TPLA reagents are packaged and sold as a kit.

Description	Configuration	Catalog Number
TPLA Reagent 1	1 x 45 mL	486647
TPLA Reagent 2	1 x 10 mL	

Materials Required but not Provided

Description	Configuration	Catalog Number
TPLA Calibrator Set	5 levels x 1 mL	486654
TPLA Control Set	2 x 1mL level A 2 x 1mL level B	486661

- Analyzer capable of running two-reagent chemistries.

Calibration

Only the TPLA Calibrator set should be used to calibrate the TPLA assay. The assigned values of the TPLA Calibrators are traceable to an in-house standard.

Refer to the instrument operator's manual for analyzer specific calibration procedures and for guidance in determining calibration frequency.

Quality Control values should be within the expected ranges.

Quality Control

Reliability of test results should be monitored routinely with quality control materials or serum pools that reasonably represent performance with patient specimens. Controls or serum pools should be used to monitor that the reagents are functioning properly and that correct procedures are being followed. An acceptable range for each lot of control material should be established by the laboratory. If control values are not within the expected range, follow normal troubleshooting procedures. If assistance is required, please contact your local distributor.

Quality control requirements should be established in accordance with local, state and/or federal regulations, or accreditation requirements.

RESULTS

Results are expressed in titer units (T.U.). T.U. are based on the *Treponema pallidum* haemagglutination (TPHA) assay. 1280 T.U. is equivalent to a titer of 1:1280 TPHA.

To convert from T.U. (Titer Units) to mIU multiply T.U. units by 2.

Limitations/Interfering Substances

Criterion: Recovery within $\pm 15\%$ of initial value

Lipemia: No significant lipemic interference up to 2850 (Formazin degree). Centrifuge at 15000xg for 10 minutes before measurement.

Hemoglobin concentration of up to 480 mg/dL did not interfere (bias <15%) in samples with syphilitic anti-lipid antibodies levels of 106 T.U.

Conjugated Bilirubin concentration of up to 20 mg/dL did not interfere (bias <15%) in samples with syphilitic anti-lipid antibodies levels of 110 T.U.

Unconjugated Bilirubin concentration of up to 18 mg/dL did not interfere (bias <15%) in samples with syphilitic anti-lipid antibodies levels of 109 T.U.

Serum samples from patients in the early stage of antibody production due to compromised immune function contain a small amount of antibody and may test negative.

A non-specific immune response may occur in serum samples from patients with autoimmune diseases. The test result should be evaluated based on other test results and clinical symptoms.

Serum samples from patients receiving blood products containing immunoglobulin may test positive. Evaluate the test result carefully.

Due to the large range of specimen concentrations possible, sample probe washing must be adequate to prevent sample carryover.

Expected Values

A measurement of 10 T.U. or higher indicates that the sample is antibody positive.

A positive antibody test should be followed by subsequent tests and should be evaluated along with other test results and clinical symptoms. A final diagnosis of syphilis should be made by a physician.

Each laboratory should confirm the reference interval for the patient population it serves.


SPECIFIC PERFORMANCE CHARACTERISTICS

Method Comparison

Comparative performance studies were conducted using the TPLA Reagent on the Hitachi 7170 analyzer and a commercially available *Treponema pallidum* hemagglutination (TPHA) method. 171 serum samples, with positive TPLA sample concentrations between 27 and 585 T.U.

The regression analysis is provided below:

Sekisui Diagnostics TPLA vs. an existing TPHA test method (n = 171)	
Slope	1.06
Intercept (T.U.)	10.5
Correlation Coefficient (r)	0.828

Precision

Precision of the TPLA reagent was determined using 2 levels of controls and positive and negative human serum sample.

Within run studies used 10 samples and total precision used 21 samples which were run three times per day for 21 days on the Hitachi 7170 analyzer.

Within-Run Precision

Serum pool / Control	Mean Recovery (T.U.)	Standard Deviation (T.U.)	% CV
Human Serum (Syphilis Positive)	137.5	1.52	1.1
Human Serum (Syphilis Negative)	0.0	0.00	—
TPLA Control A	24.4	0.51	2.1
TPLA Control B	77.4	0.59	0.8

Total Precision

Serum pool/ Control	Mean Recovery (T.U.)	Standard Deviation (T.U.)	% CV
Human Serum	134.7	3.2	2.4
Human Serum (Syphilis Negative)	0.0	0.0	—
TPLA Control A	24.9	0.7	2.7
TPLA Control B	76.7	1.1	1.4

Limit of Detection (LoD)

The limit of detection is the actual concentration at which an observed test result is 2.6 SD above that of the lowest calibrator (calibrator 1). Limit of detection was established using the level 1 calibrator and 10 measurements with TPLA Reagent on the Hitachi 7170 clinical analyzer.

The limit of detection has been defined as: 7 T.U.

Specificity

A high dose hook effect was not observed at analyte concentrations up to 327 T.U.

Rheumatoid factor was tested to 550 IU/mL and did not affect performance.

A negative result that does not match the clinical signs may occur with patients with hyperglobulinemia. In such cases the sample may be diluted with 0.9% NaCl solution and rerun.

Linearity

The TPLA method is linear from 7 T.U. to 250 T.U.

Specimens above 250 T.U. may be diluted with antibody-negative serum or diluted with 0.9% NaCl solution. Multiply the result by the dilution factor to obtain the TPLA concentration for the sample.

Other Performance Studies

Cross Reactivity (Analytical specificity)

Samples containing potentially interfering substances were analysed for cross-reactivity. Samples tested with the Sekisui Diagnostics TPLA assay were:

- from patients with collagenosis and patients undergoing dialysis
- from pregnant women

Specimen	Number	TPLA Reactive	Rate
Collagenosis patients	28	0	0%
Pregnant women	26	0	0%
Dialysis patients	50	0	0%

No false-positive results were found.

Clinical sensitivity^{5,6,7}

A total of 268 confirmed syphilis-positive samples in various stages of the disease were tested with the Sekisui Diagnostics TPLA assay. The sensitivity was found to be 100%

Syphilis- positive samples	TPLA		
	Number	Positive	Negative
268	268	0	

Clinical specificity^{4,5,6}

A total of 3427 syphilis-negative samples were tested with the Sekisui Diagnostics TPLA assay. The specificity with these samples was 99.6%.

		TPLA	
Syphilis-negative samples	Number	Positive	Negative
	3427	13	3414

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1. Longo, DL, Fauci, AS, Kasper DL, Hauser SL, Jameson JL and Loscalzo J. Eds Harrison's Online-Harrison's Principles of Internal Medicine. 18e. Part 8. Infectious Diseases. McGraw-Hill
2. Richardson JH and Barkley WE, eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, HHS Publication No. (CDC) 84-8395, Washington, DC: 1984.
3. Clinical and Laboratory Standards Institute, Procedures for the Handling and Processing of Blood Specimens: Approved Guideline. CLSI Document H18-A, Wayne, PA: 1990. 1990.
4. Shibazaki M et al. An Automated Measurement of Anti-Treponema Antibody Titer by MEDIACE TPLA, a Latex Agglutination Test using Hitachi 7170 Automatic Analyzer. The Journal of Clinical Laboratory Instruments and Reagents 1996; 19 (4):635-639.
5. Osato K et al. Clinical Evaluation of Latex Agglutination Test Kits for Detecting Anti-syphilitic Lipoidal Antibodies and Anti-treponemal Antibodies. Japanese Journal of Sexually Transmitted Diseases 2002; 13 (1):124-130.
6. Osato K et al. Clinical Evaluation of Automated Latex Agglutination Test Kits (TPLA) for Syphilis Diagnosis. The Journal of Clinical Laboratory Instruments and Reagents 1991; 14 (4):739-743.
7. Kataniwa Y et al. Clinical Evaluation of Latex Reagent Sedia TPLA for Diagnosis of Syphilis. The Journal of Clinical Laboratory Instruments and Reagents 1991; 14 (4):735-738.
8. Data on file at SEKISUI

Definitions for Symbols/ Definition der Symbole/ Définition des symboles/ Definizione dei simboli/ Definiciones de los símbolos/ Symboldefinitioner/ Sembollerin Tanımları

REF

Catalog number
Bestellnummer
Référence du catalogue
Riferimento di Catalogo
Número de catálogo
Listnummer
Katalog numarası

Temperature limitation
Temperaturbegrenzung
Limites de temperature
Limiti di temperature
Limitación de temperature
Temperaturbegränsning
Sıcaklık sınırlaması

Use by
Verwendbar bis
Utiliser jusqu'au
Utilizzare entro
Fecha de caducidad
Används före
Bu tarihe kadar kullanın

LOT

Batch code
Loscode
Code du lot
Codice di lotto
Código de lote
Satsnummer
Seri kodu

Caution, consult accompanying documents
Achtung, Beilage beachten
Attention, voir la notice d'instructions
Attenzione, consultare i documenti allegati
Precaución, consultar documentación adjunta
Försiktighet, se medföljande document
Dikkat, birlikte gelen belgelere bakın

EC REP

Authorized Representative in the European Community
Bevollmächtigter in der Europäischen Gemeinschaft
Mandataire dans la Communauté européenne
Mandatario nella Comunità Europea
Representante autorizado en la Comunidad Europea
Auktoriserad representant inom Europeiska gemenskapen
Avrupa'da yetkili temsilci

Sekisui Medical Co. Ltd.
13-5, Nihombashi 3-chome, Chuo-ku
Tokyo, 103-0027, Japan

EC REP

Sekisui Diagnostics (UK) Limited
50 Gibson Drive
Kings Hill, West Malling
KENT ME19 4AF, UK
www.sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS

April 2012
K1486647.02

IVD

For in vitro diagnostic use
Für die In-vitro-Diagnostik bestimmt
Pour diagnostic in vitro
Per uso diagnostico in vitro
Para uso diagnóstico in vitro
Används vid in vitro-diagnostik
In vitro tanı amaçlı kullanım için

Manufactured by
Hergestellt von
Fabriqué par
Prodotto da
Fabricado por
Tillverkare
Üretici

Consult instructions for use
Gebrauchsanweisung beachten
Consulter la notice d'utilisation
Consultare le istruzioni per l'uso
Consúltense las instrucciones de uso
Läs bruksanvisningen
Kullanım talimatlarına bakın

CONT

Contents
Inhalt
Contenu
Contenuto
Contenido
Innehåll
İçindekiler