# Parker TrachView™

TV-2000

## *Intubating Videoscope* Operator's Manual

**Notes**

**IMPORTANT**

The Parker TrachView™ Intubating Videoscope has been designed to facilitate orotracheal and nasotracheal intubation and airway inspection. It should be used only for those procedures and only by medical personnel who are trained to perform them.

This manual describes the recommended procedures for preparing and operating the TrachView™ Videoscope. It does not describe how any medical procedure is to be performed on a patient with this instrument.

Read all instructions in this manual carefully before using the TrachView™ on patients.

Carefully follow all safety instructions to prevent fire, electrical shock, damage to this device, and injury to user and patient.

To maintain this device in optimal condition, follow all recommendations in this manual for its handling, cleaning, and storage.

If you have questions regarding the use of this device, please contact:

PARKER MEDICAL  
109 Inverness Dr. East, Suite J  
Englewood, CO 80112-5105  
USA  
Tel: 303-799-1990  
FAX: 303-799-1996  
www.parkermedical.com

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## 1.1

## LIST OF FIGURES

Figure		Page
1	TrachView™ video box, front view	6
2	TrachView™ video box, rear view	7
3	TrachView™ fiberscope	7
4a,b	Adjusting brightness and color	8
5	Aligning and connecting fiberscope	9
6	Focusing the fiberscope	9
7a,b,c	Connecting to external devices	10-11
8a,b	Switching to backup lamp	12
9	Refocusing the fiberscope	13
10	Lamp access	16
11	Lamp release	16
12	Lamp disconnect	17
13	Removing fuseholder	18
14	Removing fuses	18
15	Replacing fuseholder	18

## 6.1 Limited Warranty

The Parker TrachView™ Intubating Videoscope, when delivered to you in new condition in the original container, is warranted to be free from defects in material or workmanship for one year from the date of shipment to you by the distributor, upon proof of date. Within this time period, parts that are returned, freight prepaid, to Parker Medical (“Parker”) and are determined by Parker to be defective will be repaired or replaced by Parker without charge for parts, labor, or return ground shipping costs. Parker will make every effort to accomplish this repair or replacement within two weeks. After the warranty period, you must pay all charges for repair and replacement. This warranty does not cover light bulbs or other consumable items.

The above actions by Parker shall constitute your exclusive remedy and Parker’s sole obligation under this warranty. Parker shall not be responsible for warranty claims made after the warranty period.

To obtain warranty repair service, you must contact Parker to obtain a Return Material Authorization (“RMA”) number, then return the product, freight prepaid, to Parker or to a service facility authorized by Parker. The RMA number and a complete explanation of the problem must be included with the product being returned to Parker for warranty service. The product to be repaired must be returned in its original box and packaging, or a similar box and packaging affording an equivalent degree of protection. Upon completion of repairs, Parker will return the product to you, freight prepaid. The warranty period for replacement parts shall begin upon shipment of same, but shall in no event exceed the warranty period of the defective part.

Parker shall have no liability or obligation for a product that has been subjected to any of the following: improper use, abuse, negligent care or handling, accident, faulty installation, improper cleaning, improper maintenance, lightning or other indications of excess voltage. This warranty is also void if the product has been repaired or modified without prior written authorization from Parker, if the end-user has failed to follow the instructions or heed the warnings or specifications in the Operator’s Manual, or if the product’s serial number has been altered or removed.

**EXCEPT FOR THE FOREGOING WARRANTIES, PARKER HEREBY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY AND/OR ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. PARKER HEREBY DISCLAIMS ANY REPRESENTATIONS OR WARRANTY THAT THIS PRODUCT OR ANY OF ITS PARTS IS COMPATIBLE WITH NON-PARKER PRODUCTS OTHER THAN VIDEO EQUIPMENT ATTACHED TO ITS VIDEO OUTPUTS, AS DESCRIBED IN THE OPERATOR’S MANUAL.**

**THE LIABILITY OF PARKER, IF ANY, AND PURCHASER’S SOLE AND EXCLUSIVE REMEDY FOR DAMAGES FOR ANY CLAIM OF ANY KIND WHATSOEVER, REGARDLESS OF THE LEGAL THEORY, SHALL NOT BE GREATER IN AMOUNT THAN THE PURCHASE PRICE OF THE PRODUCT SOLD BY PARKER THAT CAUSED ANY ALLEGED DAMAGE. IN NO EVENT SHALL PARKER BE LIABLE TO PURCHASER FOR ANY SPECIAL, INDIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OF ANY KIND.**

## 5.1 System Specifications.

<b>Parker TrachView™</b>	
Model number	CE45
<b>Mechanical</b>	
Size	33 cm (13 in) long, 20.3 cm (8 in) wide, 12 cm (4.75 in) high
Weight	4.5 kg (10 lbs)
<b>Electrical</b>	
Input voltage	100-240~
Frequency range	50/60 Hz
Input current	1.2 A
Power connector	IEC 320
Classification	Class I equipment Type B equipment
Mode of operation	Continuous
Fuse	3.15A 250V Type F
<b>Video</b>	
Video camera	Internal
Video standard	NTSC (color)
Camera resolution	480 lines
Camera sensor	0.635 cm (0.25 in) CCD
External video output 1	NTSC video, 75Ω
External video output 2	S-VHS video, 75Ω
Integrated display	LCD 16.25 cm (6.4 in) TFT, active matrix (960Hx 234V)
<b>Light Source</b>	
Type	Halogen
Lamp	Osram HLX 64615
Power	75W, 12V
<b>Fiberscope</b>	
Length	122 cm (4 ft)
Diameter	5.5 mm (0.217 in) max diameter
Optical	10,000 pixel fused silica image fiber
Field of view	100 degrees
Depth of field	3-75 mm (0.12-2.95 in)
Illumination	16, 0.5 mm plastic light fibers
<b>Environmental</b>	
Operational altitude	-200 to 10,000 ft.
Non-operational altitude	-200 to 40,000 ft.
Temperature/humidity	
Operational environment	+10 deg C, 0-10% RH; +40 deg C, 95% RH
Storage environment	-5 deg C, 0% RH; +65 deg C, 95% RH

## CERTIFICATIONS:

## 1.2

## SAFETY INSTRUCTIONS

- All safety instructions and warnings in this manual and on the TrachView™ itself should be read, understood, and carefully followed.
- Notice and follow all warnings and instructions marked on the instrument itself.
- To avoid fire hazard and electrical shock:
  - Use only the power cord specified for this product.
  - Connect the grounding conductor to earth ground.
  - Use only the fuse type and rating specified for this device.
  - Unplug the power cord from the outlet before cleaning or changing lamps.
  - Do not operate this product with covers or panels removed.
  - Do not apply an input voltage that is outside the specified range.
  - Do not operate in wet/damp conditions.
  - Do not operate the video box in an explosive atmosphere.
  - Do not allow foreign objects inside the video box.
  - Do not disassemble the video box.
  - Never spill liquids of any kind on the video box or immerse it in liquids.
- Connect to hospital-grade outlet only.
- To avoid personal injury and damage to this device:
  - Do not attempt to replace lamps while they are hot.
  - Do not place on an unstable cart, stand, or table.
  - Do not operate with suspected failures. If you suspect there is damage to this product, have it inspected by qualified personnel.
  - Do not block ventilation slots or openings—provide proper ventilation.
- Unplug the device from the wall outlet if any of the following occur:
  - The power cord is damaged.
  - Liquid has been spilled on the video box, or it has been exposed to rain or other moisture.
  - The video box has been dropped.
  - The product displays a distinct change in performance.
  - The product does not operate properly when the operating instructions are followed.
- WARNING: POSSIBLE EXPLOSION IF USED IN THE PRESENCE OF FLAMMABLE ANESTHETICS.**
- If the device does not operate correctly, refer to Section 4.4, Troubleshooting. If the device still does not function correctly, or appears to require service other than lamp and fuse replacement, call Parker Medical for factory authorized inspection and repair procedure.
- For connection only to hospital-grade outlet.
- CAUTION!** Applied part is connected to ground. The patient may become grounded. Do not use with BF or CF type equipment.
- Inspect the fiberscope cable for possible damage before each use.
- If a stylet is used beside the fiberscope in the tracheal tube, the stylet must be inserted into the tracheal tube and withdrawn from the tracheal tube together with the fiberscope. Separate insertion and withdrawal of the stylet can scratch and damage the fiberscope.

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## 4.4 Troubleshooting

**Unit does not turn on.**

- a. Check the connection of the power cord to the video box and the wall outlet.
- b. Make sure the lamp access cover is completely closed (Fig. 1).
- c. Check input fuses (Fig. 2).

**Image is too bright or too dark.**

- a. Adjust the LCD monitor brightness control (Section 2.3).

**Color is unsatisfactory.**

- a. Adjust the LCD monitor color control (Section 2.3).

**Image is hazy or foggy.**

- a. Refocus by turning the focus adjustment knob (Fig. 9).
- b. Clean the distal end of the Fiberscope (Fig. 3) with isopropyl alcohol and a cotton swab.
- c. Clean the proximal window of the optical plug (Fig. 3) with isopropyl alcohol and a cotton swab.

## 4.3 Cleaning the video box

WARNING: To avoid the risk of electrical shock, never clean the video box with the power cord plugged into the wall outlet.

WARNING: Never immerse the video box in liquids or allow liquids to seep into the video box. Electrical shock and damage to the device may result.

WARNING: Make sure that all the parts and components of the video box are completely dry before plugging into the wall outlet.

- ▶ Do not use abrasive or caustic materials to clean the TrachView™.
- ▶ Clean the device only with a mild, warm soap solution, a 5% bleach solution, or multipurpose disinfectant cleaner.
- ▶ Use a soft cloth or sponge to clean the device surfaces, including the LCD monitor surface.

## 1.3

## REGULATORY COMPLIANCE STATEMENTS

This device has been tested and found to comply with the limits for a class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interferences in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio or television reception. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause interference to radio and television reception, which can be determined by turning the equipment on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that which the receiver is connected.

This digital apparatus does not exceed the Class B limits for radio noise emissions from digital apparatus as set out in the radio interference regulations of Industry Canada.

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

**4.2 Fuse replacement**

- a. Unplug the TrachView™ power cord from the wall. outlet. This will ensure that no electrical power is present when replacing the fuses.
- b. Use a standard screwdriver to pry the fuse holder free from the back of the unit as shown in Figure13.


Figure 13. Removing fuse holder

- c. Disengage the fuses as shown in Figure 14. Replace with new fuses.


Figure 14. Removing fuses

**WARNING: USE 3.15A, 250V, TYPE F FUSES ONLY.**

- d. Place the fuse holder back in the unit (Figure 15). Push firmly until it snaps into place.


Figure 15. Replacing fuse holder

f. Disconnect the lamp by pulling off the connector plug located on the back of the lamp as shown in Figure 12.


Figure 12. Lamp disconnect

g. Discard the burned out lamp.

h. Attach the connector to the rear of the new lamp. Use only Osram HLX 6461575 (75 watt) lamps. These can be obtained in the U.S.A., from OSRAM SYLVANIA Inc., 100 Endicott Street, Danvers, MA 01923, Tel.: 978-777-1900, Fax: 978-777-2152. For International suppliers, please visit their website: [www.osram.com](http://www.osram.com).

i. Reposition the lamp release lever to the upright position.

j. Slide the new lamp downward in the lamp holder until it touches the floor of the lamp holder. (Fig. 12)

k. Close the lamp access cover and from the bottom of the unit, turn the latch screw clockwise using a coin or flathead screw driver until it stops.

l. Reconnect the power cord.

m. Test the lamp to make certain it works.

## 2.1 Component description


Figure 1. TrachView™ video box, front view.


Figure 2. TrachView™ video box, rear view.


Figure 3. TrachView™ fiberscope.

4.1 Lamp replacement

- a. Unplug the TrachView™ power cord from the wall outlet.
- b. From the bottom of the unit, turn the latch screw counterclockwise using a coin or flathead screw driver until it releases the access cover (Figure 10a).
- c. Lift the lamp access cover as shown in Figure 10b to expose the lamps. Make certain that neither lamp is hot.


Figure 10. Lamp access

- d. Press the lamp release lever beside the lamp backward to release the lamp from the holder as shown in Figure 11.
- e. Slide the lamp upward out of the lamp holder.


Figure 11. Lamp release

3. After immersion in disinfectant solution for the recommended period of time, remove the fiberscope from the disinfectant solution. Rinse by successive one-minute immersions in large volumes of fresh rinse water, according to the specific recommendations of the disinfectant manufacturer. All toxic residue must be removed from the fiberscope. Discard the water after each rinse.

4. Thoroughly rinse the fiberscope under running water. Use sterile water if possible. If this is not feasible, use water which has been filtered to improve its microbiological quality.

5. If non-sterile water has been used, wipe the fiberscope using a cloth saturated with 70% ethyl or isopropyl alcohol. Alcohol is flammable. Handle with care.

6. Wipe fiberscope dry and inspect it for possible damage.

7. If no damage has been detected, the fiberscope is ready for use.

### 3.2 Storage of fiberscope

1. The fiberscope must be dried thoroughly prior to storing.

2. The storage area should be clean, dry, well-ventilated, and maintained at a normal temperature. Avoid direct sunlight, high temperature, high humidity, and X-ray exposure.

3. The fiberscope should be stored with the flexible fiberbundle as straight as possible. If it must be coiled for storage, do not coil the fiberbundle tighter than its coiled radius in the carrying case.

### 2.2 Power ON/OFF

Connect the power cord to the power input on the back of the video box. Plug the other end of the power cord into a properly grounded hospital-grade electrical outlet.

Connect the TrachView™ fiberscope to the video box.

To turn the instrument on, push the power ON/OFF button on the video box. The power indicator light should turn on.

Push the ON/OFF button again to turn the instrument off.

### 2.3 LCD monitor

The LCD monitor turns on when the power is turned on.

Use the recessed grip area on the top edge of the LCD monitor to pull the monitor up to the desired position.

Rotate the brightness control to adjust the brightness level of the picture on the monitor, as shown in Fig. 4a.


Figure 4a. Adjusting brightness

Rotate the color control to adjust the color intensity of the picture on the monitor, as shown in Fig. 4b.


Figure 4b. Adjusting color

### 2.4 Connecting the fiberscope

Before connecting the fiberscope to the video box, the fiberscope must be properly cleaned and disinfected as described in Section 3. For best image quality and to ensure patient safety, verify that both ends of the fiberscope are free of foreign particles.

Connecting the fiberscope to the video box:

- a. Align the colored dot on the handle of the fiberscope with the matching dot on the fiberscope input port of the video box, as shown in Figure 5.


Figure 5. Aligning and connecting the fiberscope

- b. Once the dots are aligned, insert the optical plug into the input port. Turn the focus adjustment knob clockwise two complete turns to secure the fiberscope to the input port.
- c. Turn the power on and hold the distal tip of the fiberscope approximately 2.0 cm from a surface with fine details.
- d. Continue turning the focus adjustment knob clockwise, as shown in Figure 6, until the fine-detailed surface at the distal end of the scope comes into optimal focus on the LCD video monitor.


Figure 6. Focusing the fiberscope

The TrachView™ is now focused and ready for use.

### 3.1 Cleaning and disinfecting the fiberscope

The Parker TrachView™ fiberscope is shipped to the user in a carrying case that cannot protect the fiberscope from exposure to potentially infectious material. The fiberscope must therefore be properly cleaned and disinfected before use. High-level disinfection is required.

High-level disinfection destroys all fungi, vegetative bacteria (including *Mycobacterium tuberculosis*), viruses (including the human immunodeficiency virus, or HIV), and some spores. High level disinfection is most commonly achieved with chemicals that are also sterilants if exposure time is extended. Examples of chemical sterilants that can achieve high-level disinfection at appropriate concentrations are glutaraldehyde, ortho-phthalaldehyde, hydrogen peroxide, and peracetic acid.

**IMPORTANT:** The Parker TrachView™ fiberscope is **NOT** designed to withstand high pressure, high temperature, ultrasonic cleaning techniques, autoclaving or ethylene oxide sterilization. *Use of these techniques will invalidate the warranty.*

#### PROCEDURE FOR CLEANING FIBERSCOPE

**DO NOT AUTOCLAVE!**

1. Clean the fiberscope immediately after the procedure.
2. Wipe the flexible portion of the fiberscope with a clean, lint-free cloth.
3. Immerse the entire fiberscope in an enzymatic pre-soak/holding solution (such as Powerzyme). This keeps debris soft and begins breaking down infectious material. Do not immerse for longer than 1 hour.
4. Remove and rinse the fiberscope.
5. Immerse the fiberscope in a large basin filled with a fresh, low-sudsing, neutral pH, medical instrument detergent (such as Enzol) at the temperature and concentration recommended by the manufacturer of the detergent. Do not use household detergents. Do not immerse for longer than 1 hour.
6. Scrub all external surfaces with a soft brush.
7. Remove the fiberscope. Rinse thoroughly under running water and dry the fiberscope.

#### PROCEDURE FOR HIGH-LEVEL DISINFECTION OF FIBERSCOPE

**DO NOT AUTOCLAVE!**

The fiberscope must be meticulously cleaned, as outlined above, prior to undergoing high-level disinfection.

Following such cleaning, it may be reprocessed in the Steris SYSTEM 1 Sterile Processing System, made by the STERIS Corporation, 5960 Heisley Road, Mentor, Ohio 44060-1834. [www.steris.com/products](http://www.steris.com/products)

Alternatively, the fiberscope may be cleaned in a disinfectant, as follows;

1. Immerse the fiberscope into a compatible, high-level disinfectant (such as CIDEX OPA Solution).
2. Allow the fiberscope to remain in the disinfectant solution for the period of time recommended by the disinfectant manufacturer to achieve high-level disinfection of the fiberscope.

c. Once the fiberscope distal end is illuminated, refocus the image using the focus adjustment knob, as shown in Figure 9.


Figure 9. Refocusing the fiberscope

d. Replace the burned-out lamp at earliest convenience, according to the instructions in Section 4.1.

2.5 Connecting external video devices

a. The TrachView™ can be connected to a standard video monitor or an S-VHS video monitor, or both. Consult your video monitor manual for the corresponding connections. Figure 7a

To make these connections, plug the connecting cable(s) into the video output(s) on the rear of the TrachView™ and the corresponding video input(s) of the external devices, as shown in Figure 7a, Figure 7b, and Figure 7c.

A BNC-to-RCA adapter is supplied with the TrachView™ for use with an RCA cable, so that such connections can be made with an RCA cable as well as a BNC cable. The adapter allows an RCA cable to be plugged into the BNC video output on the rear of the TrachView™ video box.

b. The TrachView™ can be connected to a standard VCR, consult your VCR manual for the corresponding connections. Figure 7b

c. The TrachView™ can be connected to a standard camcorder, consult your camcorder manual for the corresponding connections. Figure 7c

No cable is supplied with the TrachView™.


Figure 7a. Connecting to external monitors


Figure 7b. Connecting to VCR


Figure 7c. Connecting to camcorder

2.6 Switching to the backup lamp

When the fiberscope is connected to the video box and the power is turned on, a bright light should appear at the distal end of the fiberscope. If this light is not present or it ceases during use, you can quickly switch to a backup lamp.

To switch to the backup lamp:

- a. Push the fiberscope handle in toward the video box, as shown in Figure 8a.


Figure 8a. Switching to backup lamp

- b. While pushing, rotate the fiberscope handle 180 degrees in a counterclockwise direction, as shown in Figure 8b.


Figure 8b. Switching to backup lamp

This rotation will turn on the backup lamp, and a bright light will appear at the distal end of the fiberscope.