SafeDrive deluxe

Personal ALCOHOL DETECTOR

Operation manual

WARNING/AGREEMENT:

Due to fact that this alcohol sensor device can be affected by shock, temperature and humidity, the user of this product agrees not to hold the manufacturer, dealer, wholesaler, or distributor responsible for the consequences of any decision based on the use of this product, to use and operate any vehicle, boat, or aircraft or other equipment and to hold them harmless from the claims of others arising out of any such decision.

IF YOU DO NOT AGREE, YOU SHOULD PROMPTLY RETURN THIS PRODUCT TO THE PLACE OF PURCHASE, UNUSED, FOR A FULL REFUND.

- This unit is designed to detect the percentage of the alcohol in your breath, thus allowing you to personally analyse the level of alcohol units you have consumed!
- The unit uses the innovative MEMS (Micro Electro Mechanical System) Gas Sensor.

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1

IMPORTANT NOTICE

Because the sensor reference can change depending on temperature and humidity, it is recommended that you calibrate every 2 weeks at least, before using this breath analyser device.

After consuming an alcoholic drink you must then wait 30 minutes before using the breath analyser. This ensures the breath sample that is taken is from your lungs, and not from any remaining alcohol that could still reside inside your mouth, which may result in a higher reading.

Because the analyser is highly sensitive, the unit will not function correctly if it is wet or damaged. Keep the unit in a secure, clean and dry place at room temperature.

The alcohol sensor can be affected by temperature or humidity. Therefore, there could be a possibility of a margin of error.

Each individual's body responds uniquely to the consumption of alcohol. A Person's ability to drive may become legally impaired at a point below the legal blood alcohol limit.

Blood alcohol content changes as the time passes and measurement might be affected by environment conditions such as temperature and humidity.

The analyser is designed to be used as a reference measure only to aid in promoting reasonable driving: it should not be used in anticipation of the use or operation of vehicles, machinery, or equipment. The analyser makes no references whatsoever to blood alcohol contents or other such measures of intoxication levels that are set by government, medical or other agencies or bodies.

2

CALIBRATION

You need to calibrate the analyser before use.


1. Press the calibration (reset) button under LCD for 1-2 seconds and release it. The "C-A-L-1" appears on the display in sequence. Then numeral on the display counts down 5 to 0 with "WAIT" flashing on the display. While the 0 stays on the display, blow into the inhaler at a distance of 2cm for 3-4 seconds. Exhale a consistent amount of air each time and at the same angle. "BLOW" flashes. (1st calibration).

2. After 1st calibration, "C-A-L-2" appears in sequence on the display. Then, the analyser starts count down from 9 to 1 this time. Finally, 0 appears. Depending on the condition of the analyser, "A-D-J" might appear on the display in sequence and restart countdown from 5 or 9.


3. While the 0 stays on the display, blow into the inhaler at a distance of 2cm for 3-4 seconds. Exhale a consistent amount of air each time and at the same angle. "BLOW" flashes. (2nd calibration) The display goes off.

3


USING THE UNIT

<FOR % BAC>

<FOR % BAC and g/l>

4

USING THE UNIT

1. Press the power button under LCD for 1-2 seconds and release it. The display counts down from 5 to 0, with "WAIT" flashing. Depending on the condition of the analyser, "A-D-J" might appear on the display in sequence then the analyser starts count down from 5 or 9.

2. When 0 appears, a beep sounds. Blow into the inhaler for 3-4 seconds. Exhale a consistent amount of air each time and at the same angle. "BLOW" flashes.

3. The test result appears with beep sound. Then the power turns off automatically after 10 seconds. Press the power button one more, you can turn off the analyser when you want.

5

INSTALLING THE BATTERIES

Your analyser requires two AAA batteries. We recommend alkaline batteries.

1. Slide open the battery compartment cover on the back of the unit.

2. Install two AAA batteries, observing the polarity (+ and -) marked inside.

3. Close the battery compartment cover.

4. If the battery is discharged, the character "L" appears on the LCD display and blinks 10 times with beep sound, then should change the battery.

6

CARE

Keep the Breath Analyzer dry; if it gets wet, wipe it dry immediately.

Use and store the Breath Analyzer only in normal temperature environments. Handle Breath Analyzer carefully; do not drop it.

Keep the Breath Analyzer away from dust and dirt, and wipe it with a damp cloth occasionally to keep it looking new.

Modifying or tampering with the analyser's internal components can cause a malfunction and might invalidate its warranty.

SPECIFICATIONS

- Power source: Two AAA batteries
- Battery life: More than 150 measurement
- Sensor Type: MEMS Semiconductor Gas Sensor
- Dimensions: 90 x 30 x 17mm
- Weight: 43g including batteries.

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