
**Intended Use**

The IL Test™ Alkaline Phosphatase is intended for the quantitative *in vitro* diagnostic determination of alkaline phosphatase in human serum and plasma using the ILab Chemistry Systems. Increases in alkaline phosphatase can be a result of liver disease, bone disease and some carcinomas. Refer to the publication "Effects of Diseases on Clinical Laboratory Tests" for a summary of causes of increased and decreased alkaline phosphatase.<sup>1</sup>

**Principle/Test Summary**

Rate analysis. Based on the methodology first proposed by Bowers and McComb<sup>2,3</sup> using p-nitrophenylphosphate (pNPP) as the substrate and optimized by Tietz *et al.*<sup>4</sup> to include a metal ion buffer and a 2-amino-2-methyl-1-propanol (AMP) buffer. Conforms to IFCC recommendations.<sup>4</sup>

The rate of increase in absorbance due to production of p-nitrophenol (pNP) is directly proportional to sample alkaline phosphatase (ALP) activity.<sup>4</sup> Absorbance measurements are taken at 410 nm (ILab 900/1800) or 405 nm (ILab 600).

**Reactive Ingredients of Final Assay Mixture**

- p-nitrophenyl phosphate	16 mmol/L
- 2-amino-2-methyl-1-propanol	350 mmol/L
- Zn <sup>2+</sup>	1.0 mmol/L
- HEDTA	2.0 mmol/L
- Mg <sup>2+</sup>	2.0 mmol/L

**Precaution**

CAUTION: ALP SUBSTRATE AND ALP R1 REAGENTS CAN CAUSE IRRITATION. Avoid contact with eyes, skin and clothing. In case of contact, immediately flush eyes with plenty of water for at least 15 minutes. Call a physician. Flush skin with water.

CAUTION: ILab's waste stream contains approximately 0.02 mg of p-nitrophenol per test. Refer to specific Federal, State, & Local environmental regulations applicable to your establishment.

Handle all specimens as potentially infectious for HIV or Hepatitis B. Avoid contact with skin and eyes (S24/25). Do not empty into drains (S29). Wear suitable protective clothing (S36). This product is For *In vitro* Diagnostic Use.

**Reagent preparation**

Add the entire contents of a vial of ALP Substrate into a container of ALP R1. Reseal the ALP R1 container, invert to mix, and place into reagent tray (ILab 900/1800: place R1 into tray #1).

For smaller workload, use an empty vial and prepare the working reagent by mixing the two components proportionately i.e. 1 part ALP Substrate to 6.5 parts ALP R1.

**Storage of Unopened Reagents**

Refrigerate (2-8°C) until expired. See expiration date on carton.

**Working Reagent Storage and Stability**

7 days at 2-8°C and on board.

**Indicators of Reagent Deterioration**

**Physical appearance:** The working reagent should be pale yellow and clear. **Control Assay:** The proper range should be obtained during a controlled assay. **Absorbance Flange:** Changes in typical reagent blank absorbances may indicate reagent deterioration.

**Specimens**

Serum, lithium heparinized plasma. Analyze samples as soon as possible within 4 hours. Alkaline phosphatase activity increases slowly with storage. Increases of 5% to 10% after storage at 4°C for 24 hours are typical.<sup>4</sup>

**Calibration**

Enzyme activity is measured by the rate of production of a chromophore with well defined optical properties. The molar absorptivity constant of p-nitrophenol is 17.6 x 10<sup>4</sup> at 410 nm and 18.45 x 10<sup>4</sup> at 405 nm. Maintenance of optical properties, fluidics and thermal control is critical to the accuracy of this assay. The ILab continuously monitors fluidics, thermal control and some optical parameters. Automatic daily maintenance monitors additional optical properties. Maintenance procedures are detailed in the ILab Operator's manual. Any adjustments to the ILab should be immediately followed by performance of daily maintenance and quality control procedures. Any failures in the continuous monitoring, daily maintenance or quality control procedures indicate improper system performance. Reagent blanking should be performed daily or before patient assays if the reagent is used less frequently.

**Data Analysis**

The ILab calculates ALP activity as follows:

$$ACT_{\text{sample}} \text{ (U/L)} = EF \times \Delta A \quad EF = \frac{V_R}{V_S} \times \frac{1}{E} \times \frac{1}{PL} \times 10^6$$

Where:  
ACT = Activity  
EF = Enzyme factor  
V<sub>S</sub> = Total sample volume (mL)  
PL = Pathlength of cuvette (cm)  
ΔA = Change in absorbance per minute  
V<sub>R</sub> = Total reaction volume (mL)  
E = Molar absorptivity constant

\*The measurement in the ILab cuvette (0.5 cm for ILab 600 and 0.6 cm for ILab 900/1800) is automatically converted to 1 cm pathlength by the ILab.

$$\text{For ALP in ILab 900/1800: } EF = \frac{0.306}{0.006} \times \frac{1}{17.6 \times 10^4} \times \frac{1}{1} \times 10^6 = 2898$$

$$\text{For ALP in ILab 600: } EF = \frac{0.204}{0.004} \times \frac{1}{18.45 \times 10^4} \times \frac{1}{1} \times 10^6 = 2764$$

To convert U/L to nkat/L, multiply by 16.67.<sup>5</sup>

**Quality Control**

Use normal and abnormal controls in accordance with good laboratory practices. Analyze at least once per day. Use at least 2 levels of IL SeraChem<sup>®</sup> controls, such as SeraChem<sup>®</sup> Level 1 (Cat. No. 0018162412) and SeraChem<sup>®</sup> Level 2 (Cat. No. 0018162512). The target values and expected ranges can be found on the insert sheets. For identification and resolution of out-of-control situations, references such as Westgard *et al.*<sup>6</sup> are recommended.

**Interferences**

ILab 600: No interference from lipemia up to sample absorbance of 4.9/cm at 660 nm (1100 mg/dL triglycerides). No interference up to 125 mg/dL (0.075 mmol/L) hemoglobin. No interference up to 16 mg/dL (272 mmol/L) bilirubin.

ILab 900/1800: No interference from lipemia up to sample absorbance of 3.0/cm at 660 nm. No interference up to 100 mg/dL (0.06 mmol/L) hemoglobin. No interference up to 20 mg/dL (340 mmol/L) bilirubin. For a comprehensive review of interfering substances, refer to the publication by Young.<sup>7</sup>

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