• Indicates **cobas c** systems on which reagents can be used

Order information

Cholesterol Gen.2

400 tests	Cat. No. 03039773 190
Calibrator f.a.s. (12 x 3 mL)	Cat. No. 10759350 190
Calibrator f.a.s. (12 x 3 mL, for USA)	Cat. No. 10759350 360
Precinorm U plus (10 x 3 mL)	Cat. No. 12149435 122
Precinorm U plus (10 x 3 mL, for USA)	Cat. No. 12149435 160
Precipath U plus (10 x 3 mL)	Cat. No. 12149443 122
Precipath U plus (10 x 3 mL, for USA)	Cat. No. 12149443 160
Precinorm U (20 x 5 mL)	Cat. No. 10171743 122
Precipath U (20 x 5 mL)	Cat. No. 10171778 122
Precinorm L (4 x 3 mL)	Cat. No. 10781827 122
Precipath L (4 x 3 mL)	Cat. No. 11285874 122
Diluent NaCl 9 % (50 mL)	Cat. No. 04489357 190

	Roche/Hitachi cobas c systems	
	cobas c 311	cobas c 501/502
System-ID 07 6726 3	•	•
Code 401		
Code 401		
Code 300		
Code 300		
Code 301		
Code 301		
Code 300		
Code 301		
Code 304		
Code 305		
System-ID 07 6869 3		

English

System information

For **cobas c** 311/501 analyzers:

CHO2I: ACN 798: ID/MS Standardization

CHO2A: ACN 433: Abell/Kendall Standardization

For **cobas c** 502 analyzer:

CHO2I: ACN 8798: ID/MS Standardization

CHO2A: ACN 8433: Abell/Kendall Standardization

Intended use

In vitro test for the quantitative determination of cholesterol in human serum and plasma on Roche/Hitachi **cobas c** systems.

Summary

Cholesterol is a steroid with a secondary hydroxyl group in the C3 position. It is synthesized in many types of tissue, but particularly in the liver and intestinal wall. Approximately three quarters of cholesterol is newly synthesized and a quarter originates from dietary intake. Cholesterol assays are used for screening for atherosclerotic risk and in the diagnosis and treatment of disorders involving elevated cholesterol levels as well as lipid and lipoprotein metabolic disorders.

Cholesterol analysis was first reported by Liebermann in 1885 followed by Burchard in 1889. In the Liebermann-Burchard reaction, cholesterol forms a blue-green dye from polymeric unsaturated carbohydrates in an acetic acid/acetic anhydride/concentrated sulfuric acid medium. The Abell and Kendall method is specific for cholesterol, but is technically complex and requires the use of corrosive reagents. In 1974, Roeschlaub and Allain described the first fully enzymatic method. This method is based on the determination of Δ^4 -cholestenone after enzymatic cleavage of the cholesterol ester by cholesterol esterase, conversion of cholesterol by cholesterol oxidase, and subsequent measurement by the Trinder reaction of the hydrogen peroxide formed. Optimization of ester cleavage (> 99.5 %) allows standardization using primary and secondary standards and a direct comparison with the CDC and NIST reference methods.^{1,2,3,4,5,6,7,8,9} Nonfasting sample results may be slightly lower than fasting results.^{10,11,12}


The Roche cholesterol assay meets the 1992 National Institutes of Health (NIH) goal of less than or equal to 3 % for both precision and bias.¹²

The assay is optionally standardized against Abell/Kendall and isotope dilution/mass spectrometry. The performance claims and data presented here are independent of the standardization.

Test principle

Enzymatic, colorimetric method.

Cholesterol esters are cleaved by the action of cholesterol esterase to yield free cholesterol and fatty acids. Cholesterol oxidase then catalyzes the oxidation of cholesterol to cholest-4-en-3-one and hydrogen peroxide. In the presence of peroxidase, the hydrogen peroxide formed effects the oxidative coupling of phenol and 4-aminophenazone to form a red quinone-imine dye.

The color intensity of the dye formed is directly proportional to the cholesterol concentration. It is determined by measuring the increase in absorbance.

Reagents – working solutions

R1 PIPES buffer: 225 mmol/L, pH 6.8; Mg²⁺: 10 mmol/L; sodium cholate: 0.6 mmol/L; 4-aminophenazone: ≥ 0.45 mmol/L; phenol: ≥ 12.6 mmol/L; fatty alcohol polyglycol ether: 3 %; cholesterol esterase (*Pseudomonas* spec.): ≥ 25 $\mu\text{kat/L}$ (≥ 1.5 U/mL); cholesterol oxidase (*E. coli*): ≥ 7.5 $\mu\text{kat/L}$ (≥ 0.45 U/mL); peroxidase (horseradish): ≥ 12.5 $\mu\text{kat/L}$ (≥ 0.75 U/mL); stabilizers; preservative

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Safety data sheet available for professional user on request.

Disposal of all waste material should be in accordance with local guidelines.

Reagent handling

Ready for use.

Storage and stability

CHOL2

Shelf life at 2-8 °C:

See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer:

4 weeks

Diluent NaCl 9 %

Shelf life at 2-8 °C:

See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer:

12 weeks

Specimen collection and preparation

For specimen collection and preparation, only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K₂-EDTA plasma

Do not use citrate, oxalate or fluoride.¹³

Fasting and nonfasting samples can be used.¹¹

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer. Centrifuge samples containing precipitates before performing the assay.

Stability:^{14,15} 7 days at 15-25 °C

7 days at 2-8 °C

3 months at (-15)-(-25) °C

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

See "Order information" section.
General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Application for serum and plasma

cobas c 311 test definition

Assay type	1 Point		
Reaction time / Assay points	10 / 57		
Wavelength (sub/main)	700/505 nm		
Reaction direction	Increase		
Units	mmol/L (mg/dL, g/L)		
Reagent pipetting		Diluent (H ₂ O)	
R1	47 µL	93 µL	
<i>Sample volumes</i>	<i>Sample</i>	<i>Sample dilution</i>	
		<i>Sample</i>	<i>Diluent (NaCl)</i>
Normal	2 µL	–	–
Decreased	2 µL	15 µL	135 µL
Increased	4 µL	–	–

cobas c 501/502 test definition

Assay type	1 Point		
Reaction time / Assay points	10 / 70		
Wavelength (sub/main)	700/505 nm		
Reaction direction	Increase		
Units	mmol/L (mg/dL, g/L)		
Reagent pipetting		Diluent (H ₂ O)	
R1	47 µL	93 µL	
<i>Sample volumes</i>	<i>Sample</i>	<i>Sample dilution</i>	
		<i>Sample</i>	<i>Diluent (NaCl)</i>
Normal	2 µL	–	–
Decreased	2 µL	15 µL	135 µL
Increased	4 µL	–	–

Calibration

Calibrators	S1: H ₂ O S2: C.f.a.s.
Calibration mode	Linear
Calibration frequency	2-point calibration
	- after reagent lot change
	- and as required following quality control procedures

Traceability: This method has been standardized according to Abell/Kendall¹² and also by isotope dilution/mass spectrometry.¹⁶

Quality Control

For quality control, use control materials as listed in the "Order information" section.

Other suitable control material can be used in addition.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

Roche/Hitachi **cobas c** systems automatically calculate the analyte concentration of each sample.

Conversion factors: mmol/L x 38.66 = mg/dL
mmol/L x 0.3866 = g/L
mg/dL x 0.0259 = mmol/L

Limitations – interference¹⁷

Criterion: Recovery within ± 10 % of initial values at a cholesterol concentration of 5.2 mmol/L (200 mg/dL).

Icterus: No significant interference up to an I index of 16 for conjugated bilirubin and 14 for unconjugated bilirubin (approximate conjugated bilirubin concentration 274 µmol/L (16 mg/dL) and approximate unconjugated bilirubin concentration 239 µmol/L (14 mg/dL)).

Hemolysis: No significant interference up to an H index of 700 (approximate hemoglobin concentration: 435 µmol/L (700 mg/dL)).

Lipemia (Intralipid): No significant interference up to an L index of 2000. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{18,19}

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on Roche/Hitachi **cobas c** systems. The latest version of the Carry over evasion list can be found with the NaOHD/SMS/Multiclean/SCCS or the NaOHD/SMS/SmpCln1 + 2/SCCS Method Sheets. For further instructions refer to the operator manual.

cobas c 502 analyzer: All special wash programming necessary for avoiding carry over is available via the **cobas** link, manual input is not required.

Where required, special wash/carry over evasion programming must be implemented prior to reporting results with this test.

Limits and ranges

Measuring range

0.1-20.7 mmol/L (3.86-800 mg/dL)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:10 dilution. Results from samples diluted by the rerun function are automatically multiplied by a factor of 10.

Lower limits of measurement

Lower detection limit of the test

0.1 mmol/L (3.86 mg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying three standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

Expected values

Clinical interpretation according to the recommendations of the European Atherosclerosis Society:²⁰

	mmol/L	mg/dL	Lipid metabolic disorder
Cholesterol	< 5.2	(< 200)	No
Triglycerides	< 2.3	(< 200)	
Cholesterol	5.2-7.8	(200-300)	Yes, if HDL-cholesterol < 0.9 mmol/L (< 35 mg/dL)
Cholesterol	> 7.8	(> 300)	
Triglycerides	> 2.3	(> 200)	Yes

Recommendations of the NCEP Adult Treatment Panel for the following risk-cutoff thresholds for the US American population:²¹

Desirable cholesterol level	< 5.2 mmol/L	(< 200 mg/dL)
Borderline high cholesterol	5.2-6.2 mmol/L	(200-240 mg/dL)
High cholesterol	≥ 6.2 mmol/L	(≥ 240 mg/dL)

CHOL2

Cholesterol Gen.2

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using human samples and controls in an internal protocol. Repeatability* (n = 21), intermediate precision** (3 aliquots per run, 1 run per day, 21 days). The following results were obtained:

Repeatability*	Mean	SD	CV
	mmol/L (mg/dL)	mmol/L (mg/dL)	%
Precinorm U	2.29 (88.5)	0.02 (0.8)	1.1
Precipath U	4.74 (183)	0.04 (2)	0.9
Human serum 1	2.85 (110)	0.03 (1)	1.1
Human serum 2	7.39 (286)	0.05 (2)	0.7
Intermediate precision**	Mean	SD	CV
	mmol/L (mg/dL)	mmol/L (mg/dL)	%
Precinorm U	2.31 (89.3)	0.04 (1.6)	1.6
Precipath U	4.85 (188)	0.08 (3)	1.6
Human serum 3	1.97 (76.2)	0.03 (1.2)	1.6
Human serum 4	7.13 (276)	0.10 (4)	1.4

* repeatability = within-run precision

** intermediate precision = total precision / between run precision / between day precision

Method comparison

Cholesterol values for human serum and plasma samples obtained on a Roche/Hitachi **cobas c** 501 analyzer (y) were compared with those determined using the same reagent on a Roche/Hitachi 917 analyzer (x). Sample size (n) = 266

Passing/Bablok ²²	Linear regression
y = 1.002x + 0.045 mmol/L	y = 1.012x - 0.015 mmol/L
r = 0.953	r = 0.997

The sample concentrations were between 1.53 and 18.5 mmol/L (59.1 and 715 mg/dL).

References

- Greiling H, Gressner AM, eds. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer, 1995.
- Liebermann C. Ber Dtsch chem Ges 1885;18:1803.
- Burchard H. Beiträge zur Kenntnis der Cholesterine. Dissertation, Rostock 1889.
- Abell L et al. Standard Methods in Clinical Chemistry 1958;26:2.
- Allain CC et al. Clin Chem 1974;20:470.
- Roeschlau P et al. Z Klin Chem Klin Biochem 1974;12:226.
- Trinder P. Ann Clin Biochem 1969;6:24.
- Siedel J, Hägele EO, Ziegenhorn J et al. Clin Chem 1983;29:1075.
- Wiebe DA, Bernert JT. Clin Chem 1984;30:352.
- Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459.
- Pisani T, Gebiski CP, Leary ET et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995;119:1127.
- Recommendations for Improving Cholesterol Measurement: A Report from the Laboratory Standardization Panel of the National Cholesterol Education Program. NIH Publication No. 90-2964, February 1990.
- Nader R, Dufour DR, Cooper GR. Preanalytical Variation in Lipid, Lipoprotein, and Apolipoprotein Testing. In: Rifai N, Warnick GR, and Dominiczak MH, editors. Handbook of Lipoprotein Testing. 2nd ed. Washington: AACC press; p.176.
- Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders Company, 1995:130-131.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO publication WHO/DIL/LAB/99.1 Rev. 2. Jan. 2002.

- Siekman L, Hüskes KP, Breuer H. 1976: Determination of cholesterol in serum using mass fragmentography - a reference method in clinical chemistry. Z Anal Chem 279, 145-146.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Study Group, European Atherosclerosis Society. Strategies for the prevention of coronary heart disease: A policy statement of the European Atherosclerosis Society. European Heart Journal 1987;8:77.
- Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). NIH Publication No 01-3670; May 2001.
- Passing H, Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C, PRECINORM and PRECIPATH are trademarks of Roche. Other brand or product names are trademarks of their respective holders. Significant additions or changes are indicated by a change bar in the margin. © 2010, Roche Diagnostics

Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336

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Literatur

- Greiling H, Gressner AM, eds. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer, 1995.
- Liebermann C. Ber Dtsch chem Ges 1885;18:1803.
- Burchard H. Beiträge zur Kenntnis der Cholesterine. Dissertation, Rostock 1889.
- Abell L et al. Standard Methods in Clinical Chemistry 1958;26:2.
- Allain CC et al. Clin Chem 1974;20:470.
- Roeschlau P et al. Z Klin Chem Klin Biochem 1974;12:226.
- Trinder P. Ann Clin Biochem 1969;6:24.
- Siedel J, Hägele EO, Ziegenhorn J et al. Clin Chem 1983;29:1075.
- Wiebe DA, Bernert JT. Clin Chem 1984;30:352.
- Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459.
- Pisani T, Gebiski CP, Leary ET et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995;119:1127.
- Recommendations for Improving Cholesterol Measurement: A Report from the Laboratory Standardization Panel of the National Cholesterol Education Program. NIH Publication No. 90-2964, February 1990.
- Nader R, Dufour DR, Cooper GR. Preatalytical Variation in Lipid, Lipoprotein, and Apolipoprotein Testing. In: Rifai N, Warnick GR, and Dominiczak MH, editors. Handbook of Lipoprotein Testing. 2nd ed. Washington: AACC press; p.176.
- Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders Company; 1995:130-131.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO publication WHO/DIL/LAB/99.1 Rev. 2. Jan. 2002.
- Siekman L, Hüskes KP, Breuer H. 1976: Determination of cholesterol in serum using mass fragmentography - a reference method in clinical chemistry. Z Anal Chem 279, 145-146.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Study Group, European Atherosclerosis Society. Strategies for the prevention of coronary heart disease: A policy statement of the European Atherosclerosis Society. European Heart Journal 1987;8:77.
- Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). NIH Publication No 01-3670; May 2001.
- Passing H, Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.

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Références

- Greiling H, Gressner AM, éd. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3^e édition. Stuttgart/New York: Schattauer, 1995.
- Liebermann C. Ber Dtsch chem Ges 1885;18:1803.
- Burchard H. Beiträge zur Kenntnis der Cholesterine. Dissertation, Rostock 1889.
- Abell L et al. Standard Methods in Clinical Chemistry 1958;26:2.
- Allain CC et al. Clin Chem 1974;20:470.
- Roeschlau P et al. Z Klin Chem Klin Biochem 1974;12:226.
- Trinder P. Ann Clin Biochem 1969;6:24.
- Siedel J, Hägele EO, Ziegenhorn J et al. Clin Chem 1983;29:1075.
- Wiebe DA, Bernert JT. Clin Chem 1984;30:352.
- Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459.
- Pisani T, GebSKI CP, Leary ET, et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995;119:1127.
- Recommendations for Improving Cholesterol Measurement: A Report from the Laboratory Standardization Panel of the National Cholesterol Education Program. NIH Publication No. 90-2964, février 1990.
- Nader R, Dufour DR, Cooper GR. Preanalytical Variation in Lipid, Lipoprotein, and Apolipoprotein Testing. Dans: Rifai N, Warnick GR, and Dominiczak MH, éditeurs. Handbook of Lipoprotein Testing. 2^e édition. Washington: AACC press; p.176.
- Tietz NW, éditeur. Clinical Guide to Laboratory Tests, 3^e édition. Philadelphie, PA: WB Saunders Company, 1995:130-131.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO publication WHO/DIL/LAB/99.1 Rev. 2. Janvier 2002.
- Siekman L, Hüskes KP, Breuer H. 1976: Determination of cholesterol in serum using mass fragmentography - a reference method in clinical chemistry. Z Anal Chem 279, 145-146.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Study Group, European Atherosclerosis Society. Strategies for the prevention of coronary heart disease: A policy statement of the European Atherosclerosis Society. European Heart Journal 1987;8:77.
- Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). NIH Publication No. 01-3670; mai 2001.
- Passing H, Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.

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- Greiling H, Gressner AM, eds. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer, 1995.
- Liebermann C. Ber Dtsch chem Ges 1885;18:1803.
- Burchard H. Beiträge zur Kenntnis der Cholesterine. Dissertation, Rostock 1889.
- Abell L et al. Standard Methods in Clinical Chemistry 1958;26:2.
- Allain CC et al. Clin Chem 1974;20:470.
- Roeschlau P et al. Z Klin Chem Klin Biochem 1974;12:226.
- Trinder P. Ann Clin Biochem 1969;6:24.
- Siedel J, Hägele EO, Ziegenhorn J et al. Clin Chem 1983;29:1075.
- Wiebe DA, Bernert JT. Clin Chem 1984;30:352.
- Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459.
- Pisani T, Gebiski CP, Leary ET et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995;119:1127.
- Recommendations for Improving Cholesterol Measurement: A Report from the Laboratory Standardization Panel of the National Cholesterol Education Program. NIH Publication No. 90-2964, February 1990.
- Nader R, Dufour DR, Cooper GR. Preanalytical Variation in Lipid, Lipoprotein, and Apolipoprotein Testing. En: Rifai N, Warnick GR, and Dominiczak MH, editors. Handbook of Lipoprotein Testing. 2nd ed. Washington: AACC press; p.176.
- Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders Company, 1995:130-131.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO publication WHO/DIL/LAB/99.1 Rev. 2. Jan. 2002.
- Siekmann L, Hüskes KP, Breuer H. 1976: Determination of cholesterol in serum using mass fragmentography - a reference method in clinical chemistry. Z Anal Chem 279, 145-146.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Study Group, European Atherosclerosis Society. Strategies for the prevention of coronary heart disease: A policy statement of the European Atherosclerosis Society. European Heart Journal 1987;8:77.
- Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). NIH Publication No 01-3670; May 2001.
- Passing H, Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.

La barra del margen indica cambios o suplementos significativos.
© 2010, Roche Diagnostics

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Letteratura

- Greiling H, Gressner AM, eds. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3^a ed. Stuttgart/New York: Schattauer, 1995.
- Liebermann C. Ber Dtsch chem Ges 1885;18:1803.
- Burchard H. Beiträge zur Kenntnis der Cholesterine. Dissertation, Rostock 1889.
- Abell L et al. Standard Methods in Clinical Chemistry 1958;26:2.
- Allain CC et al. Clin Chem 1974;20:470.
- Roeschlau P et al. Z Klin Chem Klin Biochem 1974;12:226.
- Trinder P. Ann Clin Biochem 1969;6:24.
- Siedel J, Hägele EO, Ziegenhorn J et al. Clin Chem 1983;29:1075.
- Wiebe DA, Bernert JT. Clin Chem 1984;30:352.
- Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456–2459.
- Pisani T, Gebiski CP, Leary ET et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995;119:1127.
- Recommendations for Improving Cholesterol Measurement: A Report from the Laboratory Standardization Panel of the National Cholesterol Education Program. NIH Publication No. 90-2964, febbraio 1990.
- Nader R, Dufour DR, Cooper GR. Preanalytical Variation in Lipid, Lipoprotein, and Apolipoprotein Testing. In: Rifai N, Warnick GR e Dominiczak MH, editori. Handbook of Lipoprotein Testing. 2^a ed. Washington: AACC press; p. 176.
- Tietz NW, ed. Clinical Guide to Laboratory Tests, 3^a ed. Philadelphia, PA: WB Saunders Company, 1995:130–131.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO publication WHO/DIL/LAB/99.1 Rev. 2. Gen. 2002.
- Siekman L, Hüskes KP, Breuer H. 1976: Determination of cholesterol in serum using mass fragmentography – a reference method in clinical chemistry. Z Anal Chem 279, 145–146.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470–475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385–386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376–385.
- Study Group, European Atherosclerosis Society. Strategies for the prevention of coronary heart disease: A policy statement of the European Atherosclerosis Society. European Heart Journal 1987;8:77.
- Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). NIH Publication No. 01-3670; maggio 2001.
- Passing H, Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783–790.

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Bibliografia

- Greiling H, Gressner AM, eds. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer, 1995:
- Liebermann C. Ber Dtsch chem Ges 1885;18:1803.
- Burchard H. Beiträge zur Kenntnis der Cholesterine. Dissertation, Rostock 1889.
- Abell L et al. Standard Methods in Clinical Chemistry 1958;26:2.
- Allain CC et al. Clin Chem 1974; 20:470.
- Roeschlau P et al. Z Klin Chem Klin Biochem 1974;12:226.
- Trinder P. Ann Clin Biochem 1969;6:24.
- Siedel J, Hägele EO, Ziegenhorn J et al. Clin Chem 1983; 29:1075.
- Wiebe DA, Bernert JT. Clin Chem 1984; 30:352.
- Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459.
- Pisani T, Gebiski CP, Leary ET et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995;119:1127.
- Recommendations for Improving Cholesterol Measurement: A Report from the Laboratory Standardization Panel of the National Cholesterol Education Program. NIH Publication No. 90-2964, February 1990.
- Nader R, Dufour DR, Cooper GR. Preanalytical Variation in Lipid, Lipoprotein, and Apolipoprotein Testing. In: Rifai N, Warnick GR, and Dominiczak MH, editors. Handbook of Lipoprotein Testing. 2nd ed. Washington: AACC press; p.176.
- Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders Company, 1995:130-131.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO publication WHO/DIL/LAB/99.1 Rev. 2. Jan. 2002.
- Siekman L, Hüskes KP, Breuer H. 1976: Determination of cholesterol in serum using mass fragmentography - a reference method in clinical chemistry. Z Anal Chem 279, 145-146.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug Effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Study Group, European Atherosclerosis Society. Strategies for the prevention of coronary heart disease: A policy statement of the European Atherosclerosis Society. European Heart Journal 1987;8:77.
- Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). NIH Publication No 01-3670; May 2001.
- Passing H, Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.

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