CYTORICH® RED PRESERVATIVE FLUID

REF 05NG000018

3600 ml

15 °C 30 °C

CYTORICH® Red Preservative Fluid is a hemolytic fixative. It is for cytology and small biopsy processing. It is packaged in leak-resistant plastic containers. Its active ingredients include alcohols, formaldehyde, and non-toxic demulcents, emollients, and buffers.

INTENDED USE

1. CYTORICH® Red Preservative Fluid is intended to preserve cells and small tissue fragments in suspension.
2. It is used to prepare cells and small tissue fragments for cytological and histological examination.
3. CYTORICH® Red Preservative Fluid lyses red blood cells, and solubilizes proteins.
4. The preserved samples are compatible with immuno-histochemical staining and results are similar to those achieved with neutral buffered formalin.

SUMMARY AND EXPLANATION

1. CYTORICH® Red Preservative Fluid is not sold as a fixative for screening cervicovaginal cytology.
2. CYTORICH® Red is used for diagnostic non-cervicovaginal cytology. In this context, it is a fixative for cells and small tissue samples that have been collected from human body sources including, but not limited to, sputum, washings, brushings, scrapings, body fluids, exudates, and fine needle aspirates.
3. Cells and small tissue fragments that have been fixed in CYTORICH® Red Preservative Fluid may be processed using the PREPSTAIN® slide processor, cytocentrifugation, or cell block methods.

MECHANISMS OF HEMOLYSIS AND FIXATION

CYTORICH® Red Preservative Fluid lyses red blood cells, cross links soluble globulin proteins, and alcohol fixes cells and small tissue fragments. One milliliter of CYTORICH® Red Preservative Fluid can lyse red blood cells and stabilize soluble proteins from 25 – 50 microliters of whole blood.

1. Red blood cell membranes are selectively emulsified and further disrupted by a combination of alcohol and osmotic lysis. The red blood cells of rare individuals with hemoglobinopathy may lyse slowly or resist complete lysis. Senescent red blood cells resist lysis.
2. Globulin proteins that are released from red blood cells and collected with blood plasma and interstitial tissue fluids are cross-linked by formaldehyde in solution before they are

denatured by alcohols and are thus prevented from forming protein super-aggregates (flocculent precipitates).

3. Other cell membranes act as selective barriers that allow alcohols to pass more rapidly than formaldehyde into cells' cytoplasm and nucleoplasm. This causes cells to appear alcohol fixed.
4. Subsequently, intercellular junctions and intracellular proteins are cross-link fixed by formaldehyde, which accounts for the long-term stability of the cells and tissue fragments in suspension.

STORAGE, STABILITY, AND SPECIMEN TRANSPORT REQUIREMENTS

1. CYTORICH® Red Preservative Fluid should be stored at room temperature (15° - 30° C).
2. The shelf life, extended from its manufacture date, is two years.
3. Specimens fixed in CYTORICH® Red Preservative Fluid are stable for at least 30 days.
4. Morphological stability has been demonstrated in specimens stored up to six months at room temperature.
5. Specimens fixed in CYTORICH® Red Preservative Fluid can withstand shipment at temperatures ranging from about -5° to about +45° C.

GENERAL PRECAUTIONS

1. Do not pipette by mouth.
2. Do not allow any reagents to come in contact with open wounds.
3. Wear powder-free gloves, lab coat, and eye protection during laboratory procedures.
4. Follow appropriate biohazard precautions when handling cell and tissue specimens.
5. Do not ingest (contains alcohols other than ethanol and formaldehyde).

WARNING: Cytologic specimens may contain infectious agents. Wear gloves and follow appropriate biohazard precautions when handling samples.

GENERIC SPECIMEN PREPARATION, BACKGROUND

One milliliter of CYTORICH® Red Preservative Fluid can lyse and stabilize 25 – 50 microliters of whole blood.

Using this as a guideline, one drop (50 microliters) of aqueous cell slurry may be cleared of its red cells and protein and fixed by one milliliter of CYTORICH® Red Preservative Fluid.

GENERIC SPECIMEN PREPARATION, METHOD

1. Collect cells in their native fluid, in a balanced salt solution, or in physiological saline. Fluids may be anti-coagulated (e.g., heparinized) in order to inhibit the formation of fibrin clots.

2. Concentrate the cell suspension by conventional centrifugation (600 x gravity for 5 minutes).
3. Decant and properly discard the supernatant solution.
4. Vortex the cell button to form a homogeneous cell slurry with about 50% cell-crit (50% cell-crit = 50% v/v cells + supernatant).
5. Transfer between 5 – 10 drops of the cell slurry to 10 milliliters of CYTORICH® Red Preservative Fluid and mix to an even cell suspension.
6. The sample is now stable for up to several weeks. It should be allowed to set for at least 30 minutes before it is processed.
7. Process the sample with:
 - PREPSTAIN® slide processor (see PREPSTAIN® Operator's Manual).
 - CYTORICH® Yellow and cytocentrifugation (see CYTORICH® Yellow package insert).
 - Cell block

RECOMMENDED SPECIMEN PREPARATION PROTOCOLS

BODY CAVITY FLUIDS

1. Concentrate the anti-coagulated sample by conventional centrifugation (600 x gravity for 5 minutes).
2. Decant and properly discard the supernatant solution.
3. Vortex the cell button to form a homogeneous cell slurry with about 50% cell-crit.
4. Transfer between 5 – 10 drops of the cell slurry to 10 milliliters of CYTORICH® Red Preservative Fluid and mix to an even cell suspension. Very hypocellular samples such as cerebrospinal fluids can be fixed in as little as 1 milliliter of CYTORICH® Red Preservative Fluid.
5. The sample is now stable for up to several weeks. It should be allowed to set for at least 30 minutes before it is processed.

WASHING, BRUSHINGS, AND SCRAPINGS

1. Collect cell washings into anti-coagulated aqueous salt solution. If the brush or other abrasive device is to be used several times, it may be intermittently washed in an anti-coagulated aqueous salt solution.
2. Concentrate the sample by conventional centrifugation (600 x gravity for 5 minutes).
3. Decant and properly discard the supernatant solution.
4. Vortex the cell button to form a homogeneous cell slurry with about 50% cell-crit.
5. Transfer between 5 – 10 drops of the cell slurry to 10 milliliters of CYTORICH® Red Preservative Fluid and mix to an even cell suspension.

- The sample is now stable for up to several weeks. It should be allowed to set for at least 30 minutes before it is processed.
- If the brush or other abrasive device is to be used only once (i.e., will not re-enter the patient after the collection has been performed), it may be placed directly into an adequate amount of CYTORICH® Red Preservative Fluid.
- An alternative to Step #7, is to rinse the brush off in CYTORICH® Red Preservative Fluid and then rinse it in sterile saline solution before it re-enters the patient.
- Endometrial brushes, endocervical brushes, and wooden spatulas generally require about 10 – 15 milliliters of preservative.
- Cells are readily dislodged from brushes and wooden spatulas by vigorous agitation or vortex mixing.
- The sample should be allowed to fix for at least 30 minutes before it is processed. The sample is stable for up to several weeks.
- For processing a sample that is admixed with collection devices, cover the orifice of the original collection vessel with tulle (bridal veil lace) and pour the liquid sample into a conical test tube.

FINE NEEDLE ASPIRATION BIOPSIES, METHOD 1 (SLIDES ARE NOT MADE ON SITE)

- Pre-fill a syringe with 2 milliliters of anti-coagulated salt solution.
- Pre-fill a test tube or other capped container with 10 milliliters of CYTORICH® Red Preservative Fluid.
- Perform an aspiration biopsy allowing the cell sample to just rise above the hub of the needle.
- Eject the contents of the needle and syringe into CYTORICH® Red Preservative Fluid.
- Visible small fragments and clouding of the preservative fluid are generally (> 95% of the time) associated with an adequate sampling.
- The sample should be allowed to fix for at least 30 minutes before it is processed. The sample is stable for up to several weeks.
- For processing a sample that is admixed with tissue fragments, post-fix the tissue fragments in formalin and process with cell block.

FINE NEEDLE ASPIRATION BIOPSIES, METHOD 2 (SLIDES ARE MADE ON SITE)

- Pre-fill a syringe with 2 milliliters of air. Heparinization of the needle-syringe assembly is optional but recommended.
- Perform an aspiration biopsy allowing the cell sample to just rise above the hub of the needle.

- Expel the contents of the needle onto glass slide(s) and make either air-dried or fixed smear preparations.
- Backwash the needle and syringe with 5 – 10 milliliters of CYTORICH® Red Preservative Fluid in order to clean the needle-syringe assembly.
- The sample should be allowed to fix for at least 30 minutes before it is processed. The washing from the needle-syringe assembly is now stable for up to several weeks.
- For processing a sample that is admixed with tissue fragments, cover the orifice of the original collection vessel with tulle (bridal veil lace) and pour the liquid sample into a conical test tube. Post-fix the tissue fragments in formalin and process with cell block.

SPUTUM

Unlike Saccomanno's solution, CYTORICH® Red Preservative Fluid does not harden mucus. It is not strictly mucolytic, but demulcents and emollients moisten and soften mucus. Cells can be separated from softened liquid mucus by vigorous agitation. Hardened nasal secretions, mucus blobs and small tissue fragments can be separated from the liquid sample and examined by cell block histology.

- Mix one volume of sputum with about 5 volumes of CYTORICH® Red Preservative Fluid (usually 25 milliliters) in a capped specimen container (e.g., 150 milliliters specimen jar).
- Add a stirring bar to the mixture.
- Mix vigorously on a magnetic stirrer for 15 – 30 minutes.
- In a biohazard hood, remove the container's lid, cover the container's orifice with tulle, and filter the liquid portion of the specimen into a conical test tube.
- The sample should be allowed to fix for at least 30 minutes before it is processed. The sample is stable for up to several weeks.
- Hardened nasal secretions, mucus blobs and small tissue fragments can be separated from the liquid sample, post-fixed in formalin and examined by cell block histology.

LIMITATIONS OF THE PROCEDURE

- Cytologic samples should be fixed in CYTORICH® Red Preservative Fluid as soon as possible after collection.
- A sample that has become degraded prior to fixation will be unsatisfactory for examination.
- CYTORICH® Red Preservative Fluid should not be used to fix tissue fragments larger than 5 millimeters in average diameter.

Catalog number

For *in vitro* diagnostic use

Consult instructions for use

Contains / Sufficient for <n> tests

Caution, consult accompanying document

Storage Temperature Limitations

Batch Code (Lot)

Use by

Manufacturer

TECHNICAL SUPPORT

USA

Telephone: 1-877-822-7771

Fax: 1-336-290-8333

Email: tripathtechsupport@bd.com

Europe

Telephone: +32 (0)53 720 673

Fax: +32 (0)53 720 678

TriPath Imaging, Inc.

780 Plantation Drive

Burlington, NC 27215 USA

Benex Limited

Rineanna House

Shannon Free Zone

Shannon, County Clare

Ireland

SurePath, SurePath PreCoat slides, CytoRich and PrepStain are products and registered trademarks of TriPath Imaging, Inc.

© 2011 TriPath Imaging, Inc. All Rights Reserved.

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©2011 TriPath Imaging, Inc. Wszelkie prawa zastrzeżone.