DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your metallurgical microscope Zuzi 101. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the video microscope Zuzi 101, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the video microscope avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the video microscope to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the video microscope.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Zuzi Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the video microscope warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the video microscope warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

parts.

NOTE: Never disassemble the objective lenses

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INSTRUCTIONS ON ENVIRONMENT PROTECTION

over at a collection point for the recycling of electrical and electronic appliances. The symbol on the product, the instructions for use or the packing will inform about the methods for disposal.

The materials are recyclable as mentioned in its marking. By recycling, material recycling or other forms or re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The microscope is supplied with a standard wire and it should be connected to a current wire provided with an earth wire, it should be handy to be disconnected in case of emergency.
- ◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the microscope (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Zuzi Technical Assistance Department.

Changing the lamp

IMPORTANT: before changing the lamp bulb or the fuse, make sure that the microscope is disconnected from the net. Check you are using the proper lamps, as other types can provoke malfunction.

- ◆ Turn off and unplug the microscope.
- ◆ Remove the lid of the lamp compartment (4.8).
- ◆ Take the lamp and pull firmly to remove it. If the lamp is hot, take care and do not touch the lamp to avoid burns.
- ◆ Replace the lamp for another one, making sure it is a 20 W Zuzi original one. Please take care and do not touch the lamp with bare hands. Do not use lamps with a higher power, as they could provoke over heating or any malfunction.
- ◆ Put again the lamp compartment lid (4.8).
- ◆ If necessary, center again the lamp as previously indicated.

Cleaning

NOTE: Zuzi recommends using the microscope cleaner kit (code 89000001).

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the microscope and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.
- ◆ Lenses must not be disassembled by the user. Were there any dust or dirt to be cleaned, clean it with a natural horse hair brush or a smooth piece of cloth, fluff-free, dampened with a bit of xylene or toluene.
- ◆ To remove dust from lenses, blow them with a plastic bulb or clean them with a natural horsehair brush.
- ◆ Use non-corrosive lubricants in metallic parts, being careful of not touching optical

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish 2-11
 English 12-21
 French..... 22-32

INDEX OF CONTENTS


1. USES OF THE INSTRUMENT 13
 2. DESCRIPTION 13
 3. TECHNICAL SPECIFICATIONS 14
 4. INSTALLATION / SETTING UP..... 15
 5. MAINTENANCE AND CLEANING 19
 ANNEX I: CE CERTIFICATE 32

1. USES OF THE INSTRUMENT

Zuzi metallurgical microscope model 101 is aimed to the microscopic study of any kind of opaque samples covering a wide range of applications in industry and research. It has been especially designed to require a minimum optical and mechanical maintenance assuring years of high quality and reliable service.

2. DESCRIPCIÓN

1.1 Filter holder	1.13 Eyepieces
1.2 Field diaphragm centering screw	1.14 Binocular head holding screw
1.3 Polarizer adjustment lever	1.15 Field diaphragm adjustment lever
1.4 Trinocular tube	1.16 Aperture diaphragm adjustment lever
1.5 Sliding rod to direct light	1.17 Socket
1.6 Analyzer	1.18 Light intensity control
1.7 Objectives	1.19 On/off switch
1.8 Stage	1.20 Stage cross movement knob
1.9 Upper focus limit control	1.21 Stage longitudinal movement knob
1.10 Coarse adjustment knob	1.22 Friction control ring
1.11 Fine adjustment knob	
1.12 Stage quick adjustment knob	

3. TECHNICAL SPECIFICATIONS

Code	50170101
Model	101
Head	Trinocular, Siedentopf type
Eyepeace	WF10x
Nosepeace	Quintuple
Objectives	Plan achromatic: 5x (A.N.: 0.12) 10x (A.N.: 0.25) 20x (A.N.: 0.40) 40x (A.N.: 0.60) 80x (A.N.: 0.80)
Stage	Mechanical, triple layer, 250x230 mm.Movement range154x154 mm
Diaphragms	Field and aperture
Filters	Matt, blue, green and yellow
Lamp	Halogen 6 V, 20 W
Polarizing system	Polarizer and analyzer can be slid in/out the light path

Using the polarization system

- ◆ The analyzer is used together with the polarizer for studying the optical properties of substances:
- ◆ Push the corresponding lever (1.6) to put the analyzer into the light path.
- ◆ Push the polarizer lever (1.3) to put this into the light path.
- ◆ Turn the polarizer lever (1.3), while observing through the eyepieces until analyzer and polarizer are orthogonally positioned.

Photography

The trinocular tube allows coupling a photograph camera to take images of the observed samples. For this purpose, the microscope presents a rod to direct light (1.5) with two positions:

- ◆ Outwards: the 100% of the light reaches the eyepieces for sample observation.
- ◆ Inwards: the 100% of the light reaches the trinocular tube for taking images.

If you are not using the microscope for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Security

- ◆ The microscope must be used by previously qualified staff that knows how to operate it thanks to the user manual.
- ◆ Put the microscope in a horizontal plane stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the microscope near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During operation dangerous materials such as flammable or pathological substances must be out of the working area.
- ◆ When you are not using the microscope for a long period of time please make sure it is unplugged from the net in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the microscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.


5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

Illumination adjustment (Figure 4)

- | | |
|--|---------------------------|
| 4.1. Trinocular tube | 4.7. Lamp |
| 4.2. Field diaphragm adjustment lever | 4.8. Lamp compartment lid |
| 4.3. Field diaphragm centering screw | 4.9. Lamp centering screw |
| 4.4. Filter holder | |
| 4.5. Aperture diaphragm adjustment lever | |
| 4.6. Focusing lens adjustment lever | |

The aperture diaphragm can be adjusted to change the contrast.

For centering the lamp:

- Put a piece of paper on the stage. Remove one of the objectives from the nosepiece and turn the nosepiece until the hole of the removed objective stays over the paper.
- Turn the microscope on and adjust the light intensity to have enough light.
- Focus the lamp filament on the paper. If the filament image is not centered into the field of view, reduce the aperture of the aperture diaphragm (4.5), remove the lid of the lamp compartment (4.8) and adjust the lamp centering screw (4.9) until the filament image is perfectly centered.
- Adjust the focusing lens lever (4.6) until a sharp image of the filament is obtained.

Put again the objective in the nosepiece.

Turn the focus knobs until obtaining a sharp image. If light is not uniform, re-adjust slightly the lamp centering screw (4.9) and the focusing lens lever (4.6) to make the light intensity be the same in all the field of view. Then put the lid compartment again (4.8).

Reduce the aperture of the field diaphragm (4.2). Center it by using the centering screw (4.3); to check it is well centered, open and close the diaphragm.

The filter holder (4.4) is provided with green, yellow, blue and matt glass filters. Push the filter holder until the desired filter is in the light path.

**4. INSTALLATION / SETTING UP**

Unwrap the microscope, take off the involving plastic and take off the polyspan protection in which it comes fitted. Take off all the protective items and, without connecting the microscope to the net, make sure that it does not present any damage because of the shipment. In case the microscope presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

- Binocular head
- Objectives 5x, 10x, 20x, 40x and 80x
- Eyepieces WF10x
- Filter holder
- Standard net wire
- Halogen lamp, 20 W
- Fuse
- User's manual
- Warranty

We will only accept the equipment return within 15 days after delivery provided it comes in its original wrapping.

Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.

- ◆ Please put the microscope on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the microscope near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ Loosen the binocular head holding screw (1.14). Insert the binocular head into the microscope and hold it with the screw.

NOTE: Do not release the binocular head until it has been perfectly fastened.

- ◆ Install the filter holder (1.1) and both levers for centering the field diaphragm (1.2) in the structure.
- ◆ Install the objectives (1.7) in the nosepiece.
- ◆ Remove the caps that protect the eyepiece tubes and insert both eyepieces into the binocular head.
- ◆ The microscope is supplied with a standard net wire. Please insert the wire that feeds

the AC electric current in the base of current 220 V, 50 Hz $\pm 10\%$ provided with earth wire and the other end to the microscope connector. The tension should be 220 V, 50 Hz $\pm 10\%$.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz $\pm 10\%$.

Setting up

Light controls

The on/off switch (1.19) and the light intensity control (1.18) are located at the lower right side of the microscope. The electrical system is protected by fuses; the fuse box is located near the socket.

Turn the light on (1.19) and adjust the light intensity by using the control knob (1.18). If it seems that light does not work please check the adjustment of the light intensity control (1.18).

Focusing


◆ Focusing adjustment is accomplished by using the large coarse adjustment knobs (1.10) located comfortably on each side of the frame. Fine adjustment is accomplished using the smaller knobs (1.11) located on the same focus shaft. This coaxial arrangement allows for easy, precise adjustment without drift or discomfort.

◆ Focus control: using the objective with lesser magnification, raise the stage by using the coarse adjustment knob until the image of the sample appears. Then, move the fine adjustment knob until the image is focused.

◆ Focus tension adjustment: depending on the observant, the tension of the coarse adjustment knob can be modified, since it is provided with a friction control ring (2.1) that is placed next to the right coarse adjustment knob of the microscope (figure 2). Turning the friction control ring toward the rear of the microscope increases the tension, and toward the front of the microscope loosens it. Avoid loosening it excessively since this could produce the fall of the stage together with problems on focusing accuracy.

Figure 2

- 2.1 Friction control ring
- 2.2 Upper focus limit control


◆ Pre-focusing or focus limit control: There is also a device for the control of the upper focus limit (2.2) located next to the left coarse adjustment knob (figure 2). This device avoids that the plate raises beyond a point, that can be determined once the sample has been focussed, and this fact can avoid possible accidental collisions during the observation. This mechanism also allows the possibility of getting a quick focusing, accelerating the observation process when working with different samples.

After focusing the sample with the low magnification objective, rotation of the lever toward the front of the microscope will set an upper limit on the coarse adjustment movement. After changing samples or objectives, focusing is easily accomplished by rotating the coarse adjustment knob to reach the pre-focused position, then making fine adjustments with the fine adjustment knob. The fine adjustment is not affected by using the upper focus limit control.

Diopter compensation and interpupillary distance adjustment (Figure 3)

◆ Diopter compensation: a more or less marked difference in the focusing is appreciated by closing alternatively one eye, then the other; the dioptre compensation corrects this effect. The microscope is provided with a dioptre compensation ring at the left eyepiece holder (3.1). First of all, adjust the dioptre ring to zero and looking through the right eyepiece, focus the sample specimen using the coarse and fine adjustment knobs. Once the sample is focused, look through the left eyepiece and adjust the dioptre compensation ring until obtaining a clear, sharp image.

◆ Interpupillary distance: Adjust the interpupillary distance by separating or joining the eyepieces until there is a total fusion of the two images.

Note: once the values of interpupillary distance and dioptre compensation are found it will be very useful to learn them by heart, mostly if the microscope is shared with other users. This is highly advisable to avoid repeating the adjustment each time the microscope is used.

Mechanical stage controls

The triple layer stage presents the following knobs:

- ◆ Cross movement knob (1.20)
- ◆ Longitudinal movement knob (1.21)
- ◆ Stage quick adjustment knob (1.12), to move quickly the stage cross and longitudinally.