# DADE BEHRING

## vWF Ag\*

Reagents for the immunoturbidimetric determination of von Willebrand factor antigen (vWF Ag\*)

### Intended Use

The Dade Behring vWF Ag\* test kit is intended for in vitro diagnostic use with Dade Behring coagulation analyzers for the quantitative determination of vWF Ag\* in human plasma by immuno-turbidimetry.

### Summary and Explanation

von Willebrand factor (vWF) is a glycoprotein secreted into plasma from endothelial cells and megacaryocytes<sup>1</sup>. It has a multimeric form and can have a molecular weight up to 15,000 kDa. After secretion it is vulnerable to digestion by plasmin unless it is not stored in the  $\alpha$ -granules of thrombocytes after endocytotic uptake.

The vWF plays an important role both in primary hemostasis by the formation of the hemostatic plug, due to its function in platelet adhesion and aggregation, and in the coagulation process by the stabilization of factor VIII<sup>2</sup>.

Qualitative and/or quantitative deficiency in vWF is the cause for the von Willebrand disease (vWD), which is the most common hereditary bleeding disorder, having a frequency estimated to be as high as one person in every hundred<sup>3</sup>. It has an autosomal dominant mode of inheritance. Three principal types of vWD were defined<sup>4</sup>.

Whereas the BC von Willebrand reagent (Dade Behring) can be used for determination of von Willebrand ristocetin cofactor activity (vWF:RCo) the vWF Ag\* assay is determining the level of specific antigen. Both are necessary for differential diagnosis of vWD subtypes.

vWD type 1 is the most frequent subtype with 70 to 80% of all vWD. It is characterized by a decrease of vWF Ag\* and vWF activity. Type 2 refers to a qualitative deficiency of vWF activity, whereby the antigen level still can be normal. vWD type 3 is characterized by a total absence of vWF both in plasma and in the cellular compartments.

Test	Type 1	Type 2 A	Type 2 B	Type 2 N	Type 3
vWF Ag*	↓	↓	±↓	Normal	Absent
vWF:RCo	↓	↓↓	↓ or ↓↓	Normal	Absent
Factor VIII	Normal or ±↓	Normal or ±↓	Normal or ±↓	↓↓	↓↓

Further, it should be mentioned that vWF deficiencies may be associated with various other diseases such as Systemic Lupus Erythematosus (SLE), myeloma, lymphoma, etc.. These cases may be referred to as acquired von Willebrand diseases. vWF Ag\* and vWF activity can be increased in all kinds of inflammatory processes as a result of damages of the vascular endothelium<sup>4</sup>. Besides this patients with myocardial infarction may also show high vWF levels<sup>5</sup>.

### Principle of the Method

Small polystyrene particles to which specific antibodies have been attached by covalent bonding are agglutinated when mixing with samples containing von Willebrand antigen. This agglutination is then detected turbidimetrically via the increase in turbidity, which is proportional to the antigen level present in the test sample.

### Reagents

#### Materials provided

Code No. OPAB: 4 x 6 mL (2 mL Latex Reagent + 4 mL Diluent) of vWF Ag\* reagent.

#### Composition

The vWF Ag\* kit contains the following reagents:

Buffer 4 x 5 mL (Glycine buffer)

Latex Reagent 4 x 2 mL (green screw cap, suspension of small polystyrene particles coated with rabbit anti-human vWF antibodies).

Diluent for Latex Reagent 4 x 4 mL (green screw cap, solution containing glycine intended for dilution of the Latex Reagent)

#### Warnings and Precautions

- For *in vitro* diagnostic use
- Reagents containing sodium azide must be handled with due caution.

Do not ingest or allow to contact skin or mucous membranes! If discarded into drain, flush with a large volume of water to prevent azide build-up. Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

#### Preparation of the Reagents

##### Buffer:

Before opening, allow the buffer vial to stand at room temperature (+15 to +25°C) for 15 minutes. Swirl the vial gently.

##### Latex Reagent:

Pour the whole volume of a vial Diluent for Latex Reagent (green screw cap) into a vial of Latex Reagent (also green screw cap) of the same kit. Ensure that the dilution liquid is completely transferred into the Latex Reagent vial. Then, swirl the mixture in the Reagent vial without creating any bubbles. Allow the final solution to stand at room temperature (+15 to +25°C) for 15 minutes and swirl gently afterwards.

##### Diluent for Latex Reagent:

This diluent is ready for use.

#### Storage and Stability

When stored at +2 to +8°C the reagents in intact vials are stable until the expiration date indicated on the box label.

Stability after reconstitution / once opened (closed vials)

+2 to +8 °C 15 days

+15 to +25 °C 7 days

Various analyzer do have individual stability data.

### Specimens

To separate the plasma, mix 1 part sodium citrate solution (0.109 mol/L, i.e., 3.2%) with 9 parts venous blood, taking care to avoid the formation of foam. Centrifuge at 2500 x g for at least 10 minutes, remove the supernatant plasma, and keep at +15 to +25°C until required in the assay, but not longer than 8 hours. Plasma samples can be kept at +2 to +8°C for up to 24 hours or alternatively stored deep frozen up to 1 month at -20 to -30°C. Frozen plasmas must be thawed directly at +37°C for at least 15 minutes and mixed before testing.

### Procedure

Always refer to the appropriate chapters of the Operator's Manual, particularly those for calibration, loading and quality control.

#### Materials required but not provided

Buffer: Dade® Owren's Veronal Buffer (Code No. B 4234-25)

For calibration: Standard Human Plasma (Code No. ORKL)

For quality control: Control Plasma N (Code No. ORKE) and P (Code No. OUPZ)

Dade Behring Coagulation Analyzer having an application for vWF Ag\*.

### Coagulation Analyzer

Dade Behring's vWF Ag\* can be used with a variety of coagulation analyzers. Please follow the Operator's Manual.

### Calibration

Assay calibration is performed with Dade Behring's Standard Human Plasma (calibrated against WHO Standard 97/586). The standard dilutions are automatically prepared by the Dade Behring coagulation analyzers by dilution with Dade® Owren's Veronal Buffer, according to the parameters entered in the analyzers for this assay. The buffer alone represents the 0%-point. The calibration curve can be examined on the analyzers respectively the PC screen in the Calibration menu (see appropriate Operator's Manual).

### Patient samples

Patient plasmas are used undiluted. Dilutions with Dade® Owren's Veronal Buffer are automatically prepared by the Dade Behring coagulation analyzers. Patients plasmas are identified by the Dade Behring coagulation analyzers either by scanning their bar codes across the bar code reader or by manually entering their identifications on the keyboard.

### Internal Quality Control

It is necessary to run controls to ensure accuracy and reproducibility of the results. Two different levels of controls should be included in each work-shift. Use Control Plasma N and P. Prepare the controls and transfer to the Dade Behring coagulation analyzers the information contained in the control plasma kits. These controls are used undiluted. Dilutions are being prepared automatically by the Dade Behring coagulation analyzers.

If the result of the controls is outside the confidence interval, the determination must be repeated. If the repeated determination confirms the deviation, a new reference curve should be established. Do not release patient results until the cause of deviation has been identified and corrected.

### Results

The vWF Ag\* assay is automatically carried out by the Dade Behring analyzers as soon as the samples have been loaded and the results are shown on the screen.

### Limitations of the Procedure and High-Dose-Hook Effect

Suitable coagulation analyzers are those for which assay protocols are provided by Dade Behring. The use of vWF Ag\* on other coagulation analyzers must be validated by the test laboratory. In such cases, the performance characteristics may differ from those stated.

Please refer to the Operator's Manual of your Dade Behring Coagulation Analyzer or to the application sheet of vWF Ag\* with respect to the possible influence of the high-dose-hook effect, hemoglobin, bilirubin, lipaemia or heparin.

The extremely rare presence of anti-bovine albumin and/or anti rabbit antibodies in certain subjects may lead to an over-estimation of vWF Ag\*. Furthermore, the presence of rheumatoid factor may lead to an overestimation of vWF Ag\*<sup>6</sup>.

### Expected Values

Reference values vary from laboratory to laboratory depending on technique used, therefore, each laboratory should establish its own reference interval. (For more information on establishing reference intervals see NCCLS document C28-A, How to define and determine Reference intervals in the Clinical Laboratory, Approved Guideline).

The vWF Ag\* plasma level increases with age, postmenopausal, physical exercises, and stress<sup>7</sup>.

The normal plasma vWF level in the adult population is usually in the range of 50 – 160%<sup>8,9</sup>.

### Specific Performance Characteristics

#### Measuring Range

Please refer to the Operator's Manual of your Dade Behring Coagulation Analyzer or to the application sheet of vWF Ag\*.

#### Specificity

The vWF Ag\* determines specifically the vWF antigen level in the test plasma sample.

#### Precision

Precision studies were performed with Dade Behring coagulation analyzers by evaluation of Control Plasma P and Control Plasma N.

The coefficient of variation (CV) within run (intra-assay precision) ranged from 1.4 % to 4.2 %, while the coefficient of variation from day to day (inter-assay precision) ranged from 0.9 % to 4.2 %.

### Bibliography

- Verweij CL: Biosynthesis of human von Willebrand factor. *Haemostasis* 1988; 18: 224-245.
- Ruggeri ZM, Ware J: The structure and function of von Willebrand factor. *Thromb. Haemostasis* 1992; 67: 594-599.
- Werner EJ, Broxson EH, Tucker EL et al.: Prevalence of von Willebrand disease in children: a multiethnic study. *J. Pediatr.* 1993; 123: 893-898.
- Blann A: Von Willebrand factor and the endothelium in vascular disease. *Br. J. Biomed. Sci.* 1993; 50: 125-134.
- Lip GYH, Blann AD: Von Willebrand factor and its relevance to cardiovascular disorders. *Br. Heart J.* 1995; 74: 580-583.
- Favaloro EJ, Aboud M, Arthur C: Possibility of potential vWD misdiagnosis or misclassification using LIA technology and due to presence of rheumatoid factor. *Am. J. Hematol.* 2001; 66: 53-56.
- Wamala S, Murray MA, Horsten M et al.: Socioeconomic status and determinants of hemostatic function in healthy women. *Arterioscler. Thromb. Vasc. Biol.* 1999; 19: 485-492.
- Holmberg L, Berntorp E, Donner M, Nilsson IM: von Willebrand's disease characterized by increased ristocetin sensitivity and the presence of all von Willebrand factor multimers in plasma. *Blood* 1986; 68 (3): 668-672.
- Nishino M, Girma JP, Rothschild C, Fressinaud E, Meyer D: New variant of von Willebrand disease with defective binding to factor VIII. *Blood* 1989; 74(5): 1591-1599.

Dade is a registered trademark of Dade Behring Behring Inc. in the USA, Germany and other countries.

\* vWF Ag is a registered trademark of Dade Behring Marburg GmbH in Germany and other countries. vWF Ag is not available in the USA.

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Dade é uma marca registada da Dade Behring Inc. nos EUA, Alemanha e noutros países.