Reflotron Hemoglobin

REF		SYSTEM
10744964	30	Reflotron

English

Intended use

Test for the quantitative determination of hemoglobin in blood with Reflotron systems.

Summary

Determination of the total hemoglobin in the blood, in combination with other parameters like the haematocrit and the red cell count, provides valuable information about functions of the circulatory system which are extremely important for the organism: oxygen transport, carbon dioxide transport and regulation of the acid-base balance.

Determination of hemoglobin is thus one of the fundamental components of any clinical-chemical evaluation. Determination of hemoglobin alone can already indicate the presence of anaemia, the type of which (reduced hemoglobin content or reduced numbers of erythrocytes) must then be established in further tests.

Test principle

The Reflotron Hemoglobin test is based on the oxidation of hemoglobin to methemoglobin by potassium hexacyanoferrate (III) ($\text{Fe}^{2+} \rightarrow \text{Fe}^{3+}$). In order to speed up the reaction, the erythrocytes are lysed by saponin:¹

The hemoglobin concentration (proportional to the color intensity) is measured at a wavelength of 567 nm and 37 °C, and is displayed after about 120 seconds in g/dL, g/L or mmol/L.

Reagents

Components per test:

$\text{K}_3[\text{Fe}(\text{CN})_6]$: 324 µg; saponin: 324 µg; buffer.

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

Avoid any contact with the application zone of a test strip (e.g., during pipetting of sample).

Reagent handling

Test strips are ready-for-use.

Storage and stability

Store at 2-30 °C.

Do not use the test strip after the specified expiration date.

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Use fresh venous blood immediately after collection.

EDTA- or heparinized blood kept in closed containers must be used within 24 hours. Shake the sample well before performing the test in order to ensure the homogeneous distribution of the cellular components.

Do not use fluoride as an anticoagulant.

Sample volume: 30 µL

Materials provided

1 container with 30 test strips

Materials required (but not provided)

Reflotron instrument

Reflotron pipette

Reflotron pipette tips

- Reflotron Precinorm HB

- 0.9 % NaCl

General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document. Refer to the appropriate operator's manual for instrument-specific instructions.

- Remove a test strip from the container. Tightly recap the container immediately after removing a test strip.
- All Reflotron tests require a sample volume of 30 µL
- Apply the required volume of sample onto the centre of the red application zone using a pipette (e.g., Reflotron pipette) – being careful not to touch the application zone. Avoid air-bubbles.
- Open the flap or sliding cover. Within 15 seconds of applying the sample, place the test strip onto the guide, and slide it forward horizontally until it locks into place. Close the sliding cover or flap.
- The test parameter abbreviation is shown on the display, if the test strip has been correctly inserted and the magnetic code has been read. The result is displayed depending on the setting of the instrument.

Calibration

The function curve for the Reflotron Hemoglobin assay for converting reflectance values into concentrations is defined for each lot using the Sysmex SLS-hemoglobin method.

The parameters of the curve are automatically transferred to the instrument via the magnetic strip during testing.

Quality control

For quality control, use Precinorm HB.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

The hemoglobin concentration is calculated automatically from the measurements taken, as well as function and conversion factors read from the magnetic strip on the lower face of each test strip. The concentration is displayed in g/dL, g/L or mmol/L depending on whether the instrument has been set to show conventional or SI units.

Conversion factors: $g/\text{dL} \times 0.621 = \text{mmol/L}$; $g/\text{dL} \times 10 = \text{g/L}$

Limitations - interference

Incomplete filling of the microcuvette due to improper dosing can lead to incorrect values. The test must be repeated. No influence of endogenous or exogenous substances on the result of the test has been found to date. Haematocrit values up to 55 % have no influence on the result of the test.^{2,3}

Capillary blood samples cannot be used as sample material with the Reflotron Hemoglobin test. Only venous blood is a suited sample material.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Measuring range

5.0-20.0 g/dL, 3.10-12.4 mmol/L, or 50-200 g/L.

If the measured hemoglobin value is above the measuring range of the Reflotron Hemoglobin assay, the sample may be diluted 1 + 1 with physiological saline solution. Multiply the result by a factor of 2.

Expected values

Men:⁴ 14-17.5 g/dL, 8.7-10.9 mmol/L or 140-175 g/L.

Women:⁴ 12.3-15.3 g/dL, 7.6-9.5 mmol/L or 123-153 g/L.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

The data for the Reflotron Hemoglobin assay were determined in evaluation studies. The majority of the test results were within the given ranges.

Precision

Repeatability (within-run precision):

CV (coefficient of variation) in the normal and in the pathological range: 3.98 %; sample material: EDTA or heparinized blood.

Intermediate precision (between-day precision):

CV in the normal and in the pathological range: 3.64 %; sample material: control materials.

Method comparison

A comparison of the Reflotron Hemoglobin assay (y) with the Sysmex SLS-hemoglobin method (x) using EDTA, heparinized blood gave the following correlations (n = 75):

EDTA venous blood: $y = 1.043x - 0.87$, $r = 0.86$

Heparinized blood: $y = 1.163x - 1.79$, $r = 0.83$

For further information, please refer to the appropriate operator's manual for the analyzer concerned, and the Method Sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Significant additions or changes are indicated by a change bar in the margin.

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