
MODELS
865XXXXX
866XXXXX

OPERATING MANUAL

MINISTAND

STAND UP LIFT

MINISTAND

2001, Tanguay Street
Magog (Québec)
Canada
J1X 5Y5
Tel.: (819) 868-0441
Toll Free 1-800-868-0441
Fax : (819) 868-2249
www.bhm-medical.com

TABLE OF CONTENTS

SYMBOLS	4
SAFETY INSTRUCTIONS & WARNINGS	5
GENERAL	5
SHOCK PREVENTION	7
FIRE AND EXPLOSION	7
EQUIPMENT WARNING LABELS	7
OPERATING THE MINISTAND.....	8
INTRODUCTION.....	8
ADJUSTING THE WIDTH OF THE LEGS.....	8
POWER BASE MODELS	9
CORRECT USE OF THE BRAKES	9
EMERGENCY SHUT-OFF.....	10
BATTERY PACK.....	10
REMOVABLE BATTERY PACK.....	11
CHARGING THE BATTERY.....	11
RAISING AND LOWERING THE PATIENT	12
DAILY VISUAL CHECKLIST.....	13
BAND SLING INSTALLATION	13
Band Sling sizes and product numbers	14
BAND SLING INSTALLATION (PATIENT)	14
Fixing the band sling to the Ministand.....	15
BAND SLING POSITIONS	15
TRANSFERRING THE PATIENT	16
EMERGENCY LOWERING.....	17
MONTHLY INSPECTION DETAILS.....	18
DIAGRAM.....	19
TROUBLESHOOTING GUIDE	20
TRANSFER SLING SAFETY	21
HOW TO CONDUCT A VISUAL SAFETY INSPECTION.....	22
LAUNDERING INSTRUCTIONS.....	24
VISUAL SAFETY INSPECTION RECORD.....	24
WARRANTY	25

© 2004 **WARNING!**

All rights reserved. CONFIDENTIAL. The reproduction of this document or the transmittal in any form of any information, contained herein, without the authority in writing of an officer of the manufacturer is prohibited.

TECHNICAL SPECIFICATIONS AND DESIGN SUBJECT TO CHANGE WITHOUT NOTICE.

SYMBOLS

WARNING: this symbol is intended to alert the user to hazards or unsafe practices which could result in serious bodily harm.

CAUTION: this symbol is intended to alert the user of the presence of important operating and maintenance instructions which could prevent product damage or possible personal injury.

NOTE: this symbol offers helpful information concerning certain operating procedures.

CONSULT ACCOMPANYING DOCUMENTS

BOOM "UP"


BOOM "DOWN"

BASE "CLOSED"

BASE "OPENED"

INDICATES CHARGING MODE

INDICATES CHARGER POWER IS "ON"

DO NOT ATTEMPT TO USE THIS EQUIPMENT WITHOUT UNDERSTANDING THIS MANUAL.

To ensure safe operation, read carefully the entire manual, especially the section on "Safety Instructions and Warnings", before installing, operating or servicing this equipment.

If anything is not completely understood, please contact your supplier for more details. Failure to comply with warnings in this manual may result in injury.

Keep this manual with lift and refer to it as required. Contents of this manual are subjected to change without prior notice to users.

SAFETY INSTRUCTIONS & WARNINGS


GENERAL

- IMPORTANT – CAREFULLY READ THESE INSTRUCTIONS OR SERIOUS INJURY MAY RESULT.
- KEEP THESE INSTRUCTIONS WITH THE LIFT AT ALL TIMES.
- READ OPERATING MAINTENANCE INSTRUCTIONS IN THIS MANUAL BEFORE INSTALLING, OPERATING OR SERVICING THIS EQUIPMENT.
- YOUR LIFT is for transferring patients only. Do not use the lift for any other purpose.
- ALWAYS carry out the daily checklist before using the lift.
- BHM Medical patient lifts are designed specifically for use with slings and accessories manufactured by BHM Medical. We cannot ensure the safety of a transfer, and therefore cannot be held responsible, for any incident that may occur as a result of the use of slings or accessories produced by another manufacturer. In order to ensure the safety of the patient and that of the caregiver use only products manufactured by BHM Medical with this patient lift.
- This mobile floor lift is intended to be used for patients within the specified weight limit indicated for the lift. Do not attempt to lift more than the weight limit indicated.
- BEFORE attempting to transfer, the patient must be assessed by a qualified professional.
- This mobile patient lift must be used by a caregiver with proper training to work with the patient to be transferred. This patient lift **should never be used** by a patient on their own.
- ONLY trained and qualified caregivers should transfer a patient. DO NOT attempt to use the lift if you have not been properly trained to do so.
- ALWAYS be prepared before attempting to transfer a patient.
- DO NOT use a sling that is not recommended for the lift.
- NEVER use a damaged, torn or frayed sling.
- ALWAYS place the sling around patient according to the instructions enclosed.
- FOLLOW lifting procedures outlined in this manual.
- USE all controls and safety features only according to the rules specified in this manual. Never attempt to force a control or button on the lift.
- DO NOT store the lift in a shower, bath or other area with high humidity.
- IMPORTANT: keep all components of the lift clean and dry, and have electrical and mechanical safety checkpoints done as instructed in the “Maintenance” section of this manual.
- Replace any precautionary or instruction labels that cannot be easily read.
- Do not attempt to manoeuvre the lift by pushing on the mast, motor shaft, boom or patient.
- ALWAYS manoeuvre the lift with the handles provided.
- DO NOT push a loaded lift at speeds which exceed a slow walking pace (3 km/hour or 0.8 meter/second).

SAFETY INSTRUCTIONS & WARNINGS

- DO NOT push the lift over uneven or rough ground, particularly if loaded.
- DO NOT attempt to push or pull a loaded lift over a floor obstruction which the casters are unable to ride over easily.
- DO NOT BUMP the lift down steps, loaded or unloaded.
- DO NOT park a loaded lift on any sloping surface.
- DO NOT use an electric lift in a shower
- This mobile floor lift is not intended to be a transport device or to be used to transport a patient over any long distance.
- BASE should not be opened as much as possible for optimum safety.
- DO NOT USE, in any way, the actuator as a handle to push or pull the lift. If the actuator is used as a handle, the caregiver and the patient safety is greatly compromised.
- Do not put fingers, hands or feet where space is limited (see diagram below). This could pinch, cut or seriously injure any parts exposed.

* All models

SAFETY INSTRUCTIONS & WARNINGS

SHOCK PREVENTION

- AVOID violent shock during transportation.
- DO NOT touch or use a lift with bare conductors or a damaged power cord. Electrically live equipment can electrocute a patient. If the lift or charger has any exposed or damaged wires, contact your local dealer immediately.
- DO NOT splash or expose electric parts of the device to water or moisture.
- CHECK nameplate for voltage and cycle requirements. These requirements differ by country. Do not attempt to use the lift in an area that has a different voltage and cycle requirements.
- READ the battery and charger instructions thoroughly before using or storing them.

FIRE AND EXPLOSION

- Batteries may explode, leak and cause personal injury if not disposed of properly.
- Do not place or store the battery under direct sunlight or near a heat source.
- Do not dispose of in fire.
- Do not short the battery terminals.
- Do not incinerate.
- Flush with water if electrolyte (Acid) comes in contact with skin or eyes.
- Batteries must be recycled or disposed of according to local law regulations. When returning batteries, insulate their terminals with adhesive tape, etc. otherwise the residual electricity in use batteries may cause fire or explosion.

IMPORTANT :

BATTERIES NEED TO BE CHARGED FOR A MINIMUM OF 8 HOURS PRIOR TO USING THIS LIFT FOR THE FIRST TIME.

EQUIPMENT WARNING LABELS

- Inspect all precautionary labels on the equipment. Order and replace all labels that cannot be easily read.

OPERATING THE MINISTAND

INTRODUCTION

CAREFULLY READ THE "SAFETY INSTRUCTIONS AND WARNINGS" SECTION BEFORE ATTEMPTING TO USE THE LIFT.

Always carry out the daily checklist before each lift use.

Operating the Ministand is simple and requires little instruction. Using the lift, an operator is able to raise, transfer, and lower a patient safely.

IMPORTANT!

The Ministand should not be used in conjunction with any patient slings other than the Manufacturer slings which are expressly designed for the system. The manufacturer disclaims any responsibility for improper use employing any other sling.

NOTE:

Instructions provided here are described from the operator's point of view, ie: when positioned behind the lift and ready to transfer the patient.

ADJUSTING THE WIDTH OF THE LEGS

The base of the lift can be set with the legs opened or closed. Adjustments of the base varied according to the model:

- **Shifter base:** The adjustments are made by changing the position of the shifter: Disengage the shifter plate by pulling the shifter arm then push it to the left will close the legs. Pushing the shifter arm to the right will open the legs (Fig.1).
- **Pedal base :** The adjustments are made by changing the position of the pedal. Press down the left side of the pedal will close the legs. Press down the right side of the pedal will open the legs (Fig. 2).
- **Power base :** The adjustments are made by using the control buttons OPEN and CLOSE on the control box or on the hand control (Fig. 3).

The position of the legs of the mobile floor lift must be fully set in either the opened or closed position. Do not position the legs between these two points.


Fig.1


Fig. 2


Fig. 3

Since the residents centre of gravity is positioned over the centre of the chassis, in most situations the resident can be comfortably transported with the base closed. Keeping the resident at the lowest point within the base, lowers the centre of gravity, increasing resident security. The residents lateral movement is also reduced by the design of the arms.

OPERATING THE MINISTAND

POWER BASE MODELS

The base of the lift can be set with the legs opened or closed by using the two buttons on the hand controller or by utilizing the auxiliary buttons located on the panel. (Fig. 2)


Fig. 2

CORRECT USE OF THE BRAKES

Foot operated brakes are fitted on both rear casters.

Brakes should only be used in the following situations:

- When lifter and resident are momentarily at rest, for example, during preparation for transferring to a bed or a chair.
- Whenever movement of the lift has to be halted while transporting a resident.
- Travelling across or stopping on slightly inclined floors may call for temporary use of the brakes.

Brakes should not be used in the following situation:


- When raising or lowering a resident in the lifter over a bed or chair. Leaving the brakes off allows the lifter (with resident) to maintain its centre of gravity throughout the transfer and reduces the tendency for the patient to move laterally outside the base.

OPERATING THE MINISTAND

CORRECT USE OF THE BRAKES


Brakes application and release:

- To apply brakes, step on the back portion of pad. (Fig 3a)
- To release brakes, step on the front portion of pad or lift with toe of shoe. (Fig. 3b)

EMERGENCY SHUT-OFF

The operator can shut off the electrical power at any time by pressing the red mushroom shaped emergency button on the control box. The button is reset by a half-turn clockwise in the direction of the arrows. First time users should practice the shut-off manoeuvre before operating the lift with a resident. (Fig. 4)

BATTERY PACK

The battery pack is medically approved according to EN 60601-1, CAN/CSA-C22-2, No. 601-1 M90 and UL 2601-1.

The battery pack assembly is two 12 volt, 5 Ah (24 volt capacity) fuse protected batteries, delivering up to 100 lifts per charge (200 lbs or 90 kg).

Battery life is variable (2-3 years) and is influenced by proper charging practices and load exertion.

OPERATING THE MINISTAND

REMOVABLE BATTERY PACK

The removable battery pack reduces the time your lift is out of service because of a discharged battery. To remove discharged battery pack, use both hands, grasp the battery pack by the handles and at the same time, push the 2 latch buttons located on each side of the battery pack, and pull straight out towards you. (Replace the pack with a fully charged one (second battery pack not included) from the wall mounted charging unit (not included). (Fig. 5)


Fig. 5

CHARGING THE BATTERY

Plug the charger into a wall socket compatible with universal voltage input of 90 V up to 240 V. (Fig. 6a)


Fig. 6a

The charger has two LED indicators; yellow and green. Their significance is:

- Green light "ON" with no battery pack inserted: indicates the charger is on and ready for use.
- Green & Yellow light "ON": indicated the battery pack is being charged.
- Green light "ON" with battery pack inserted: indicates pack in the charger is fully charged and ready to use. (Fig. 6b)


Fig. 6b

DO NOT OPERATE THE CHARGER WITH A DAMAGED CORD OR IF THE UNITS HOUSING HAS BEEN DAMAGED.

DO NOT forcibly bend the power cord or place a heavy object on it. This will damage the cord and may cause a fire or electrical shock.

OPERATING THE MINISTAND

RAISING AND LOWERING THE PATIENT

The electrical actuator operated in both directions allows the operator to raise or lower the patient without any physical effort. The hand control can be used to operate the lift while at the same time the caregiver moves to a position where they can be with the patient continuously. (Fig. 7)


Fig. 7

To raise the patient:

Touch the “UP” button on the hand control. Keep your finger on the button until the patient is at the desired height, release finger pressure and the unit stops.

To lower the patient:

Touch the “DOWN” button on the hand control. Keep your finger on the button until the patient is at the desired height, release finger pressure and the unit stops.

Auxiliary “UP” and “DOWN” buttons are located on the control box attached to the mast. These controls serve two functions:

They are secondary controls in the event the hand controller stops functioning or if two attendants are involved in the transfer, one can use the auxiliary controls while the other moves into a position where they are with the patient continuously. (Fig. 7b)

NOTE: FIRST FUNCTION BUTTON PUSHED WILL OVERRIDE ALL OTHER FUNCTIONS.

OPERATING THE MINISTAND

DAILY VISUAL CHECKLIST

Proceed to the following verification daily prior using the Ministand :

- Make sure that split rings are correctly fixed on the actuator pivot shaft.
- Make sure the mast is fixed to the base.
- Verify the good functioning of the rear casters brakes.
- Make sure that knee pad is correctly fixed.

Before lifting the patient:

- Make sure that all straps are attached to the arms.
- Make sure the patient is comfortable.
- Make sure the sling is not caught on an obstruction (wheelchair brake or arm of chair).

**If any of the above occurs – DO NOT LIFT THE PATIENT
and correct the problem IMMEDIATELY.**

BAND SLING INSTALLATION

This sling is specifically for mobile stand assist lift only. The use of this sling in conjunction with the Ministand assists in changing incontinence pads or to transfer onto a toilet as it provides an open area from the middle of the back to the feet. The soft polyester nylon material is padded and quilted to add comfort. Strips sewn into the fabric in the back of the sling covers a wide area and distributes the weight evenly on the patient's back and reduce pressure under the arms. PVC mesh material on the inside of the sling provides a grip that is soft to the touch, yet holds to the clothing of the patient. Additional padding under the arms reduces pressure and makes the sling very comfortable. Waist belt with velcro closure and belts makes the sling easy to adjust and secure around the patient In order to use this sling, the patient must have adequate muscle-tone in their shoulders and lower body. Available in sizes S, M, L, XL. Weight capacity 180 kg (400 lbs.)

Do not use the band sling on a patient who does not have good muscle-tone in their shoulders and neck. Check with the patient's physician, nurse or medical specialist before using this sling with a patient.


OPERATING THE MINISTAND

Band Sling sizes and product numbers

	SMALL	MEDIUM	LARGE	X-LARGE
<i>CHEST SIZE</i>	65-105 cm (26-42 po)	75-110 cm (30-44 po)	80-115 cm (32-46 po)	85-120 cm (34-48 po)
<i>PRODUCT #</i>	TST-S	TST-M	TST-S	TST-XL
<i>COLOR CODE</i>	Red	Yellow	Green	Blue

Band sling description

Sling has two sides, the white PVC material indicates the side that is intended to be next to the patient, the label and colored handles are on the outside. The top of the sling is rounded and the padding on the arm bands should be upward. The sling has a waist section with velcro and belts.

BAND SLING INSTALLATION (PATIENT)

Whenever possible have someone assist you to ensure the patient's safety at all time.

1. Lean patient forward in the chair and make sure the patient is supported.
2. With the white PVC material of the sling facing the patient, place the sling approximately halfway down the back, or so the sling is below the shoulder blades. Lean the patient against the back of the chair.
3. Fold the fabric waist band with the velcro and belts attached across the mid-section of the patient. Velcro the belt into place, then clip the belts closed. Make sure the buckle is secured.

4. Tighten the belt until snug by pulling gently on the end of it.

OPERATING THE MINISTAND

Fixing the band sling to the Ministand

- Make sure the sling is correctly installed around the patient waist (as previously shown).
- Choose the appropriate loop colour according to the patient needs. Install it into the hooks of the Ministand arm.
- Make sure the colour of both loops is the same on both hooks.
- Make sure the loops of the sling are correctly installed in the hook – see Picture A. If not, do not lift the patient and install loop properly.
- Ensure the patient's arms are over the sling and that the sling is not caught on any obstructions (eg. Wheelchair brake or handle). See picture B.

A

B

BAND SLING POSITIONS

Sling position	
	FRONT
	BLUE
	GREEN
	YELLOW
	RED

OPERATING THE MINISTAND

TRANSFERRING THE PATIENT

HAVE SOMEONE ASSIST YOU WHEN ATTEMPTING TO TRANSFER A PATIENT

THIS MOBILE FLOOR LIFT IS NOT INTENDED TO BE A TRANSPORT DEVICE. PATIENT IN THE LIFT SHOULD NOT BE MOVED MORE THAN A FEW METERS (YARDS)

DO NOT ATTEMPT TO MANŒUVRE THE LIFT BY PULLING IN THE MAST, BOOM, ACTUATOR OR PATIENT.

BEFORE LIFTING THE PATIENT :

1. Lower the boom using DOWN button so it is as low as possible.
2. Bring the lift towards the patient so that the Ministand foot rest is right in front of the patient's feet.
3. Gently place the patient's feet on the foot rest and push the lift until the patient's shins are in the pads.
4. Lock the rear wheels on the Ministand.
5. Place the sling loops on the Ministand arms hooks according to the instructions shown in the sling installation section of this manual. The patient's arms must be over top of the sling, and if the patient is able, holding on to the Ministand arms.
6. Press the UP button to raise the patient slightly. Check to make sure that there are no obstacles in the way. If you find that the sling has loosen a bit around the waist of the patient, tighten the belts around the patient's waist. Once the patient is secure, continue to raise the patient until he/she is above the chair or toilet before moving to the next transfer point.
10. Move the patient directly over the point of transfer, lock the wheels on the wheelchair or bed. When the patient is positioned against the chair or bed, lock the rear wheels on the Ministand.
11. Slowly lower the patient by pressing the DOWN button. As you lower the patient, use the handles on the back of the sling to position the patient in the chair.
12. Lower the arms of the Ministand so they are as far down as possible.
13. Remove the straps of the sling from the arms of the Ministand. Lift the patient's feet off the foot rest and move the lift a few inches away from the patient.
14. Sling may be removed from around the patient.

OPERATING THE MINISTAND

EMERGENCY LOWERING

The Stand-up™ is equipped with an emergency lowering function, which is located just above the emergency stop button.

To operate the emergency lowering function:

Press and hold the button located on the front of the control box until patient is safely lowered to desired height. (Fig 7).


FIG 7

MONTHLY INSPECTION DETAILS

Refer to diagram on the following page for the location of each inspection area.


These procedures are supplied as a guideline for personnel who are responsible for the maintenance and safety inspections of patient handling devices. The manufacturer accepts no responsibility for ongoing inspections and checks. These are guidelines only and are not intended as a substitute for a regular, thorough equipment maintenance program.

Note : if during your inspection any bolt, nut or structural fastener is removed, it is imperative that you reinstall all associated safety devices such as cotter pins, split rings and lock nuts. If one of these devices is defective, you must replace it with an original new component.

For all Stand-up lifts models (with shifter, with pedal or power base)

ITEM #	PROCEDURE	√
1	ENSURE THE PIVOT CLEVIS PIN IS WELL SECURED BY THE SPLIT RING AT BOTH END OF THE ACTUATOR.	
2	ENSURE THAT THIS PIVOT SHOULDER BOLT IS SECURELY FASTENED WITH THE LOCKNUT.	
3	CHECK FRONT AND REAR CASTERS REGULARLY FOR HAIR AND DEBRIS, CLEAN AS NECESSARY, CHECK BRAKE FUNCTION ON REAR CASTERS.	
4	INSPECT ALL WELD SITES FOR CRACKING OR SEPARATION.	
5	ENSURE ALL NUTS ARE SECURE ON THE STUDS, ENSURE JAM NUTS ARE SECURE ON SHAFT	
6	MAKE SURE THE COTTER PIN IS IN PLACE	
7	RUN ACTUATOR FROM LIMIT TO LIMIT TO ENSURE LIMIT STOPS ARE FUNCTIONING.	
8	PRESS THE EMERGENCY BUTTON TO ENSURE THAT ALL ELECTRICAL POWER IS CUT OFF.	
9	CHECK ALL THE FUNCTIONS ON THE HAND CONTROL – CHECK TO ENSURE WIRES ARE NOT PULLED OUT OF CONNECTOR ON HAND CONTROL CASE.	
10	CHECK FOR THE SMOOTH FUNCTION OF THE UP/DOWN SWITCHES.	
11	CHECK BATTERY PACKS FOR DAMAGE – CHECK THAT THEY WILL ACCEPT A CHARGE BY PUTTING THEM IN THE CHARGER.	
12	CHECK THE FUNCTION OF THE CHARGER LIGHTS (ON THE LIFT OR ON WALL CHARGER) BY INSERTING A BATTERY PACK – ENSURE THAT POWER LIGHT IS FUNCTIONAL.	
13	CHECK THAT THE SAFETY LATCHES ARE IN PLACE AND FUNCTIONAL.	
14	MAKE SURE BOTH LEGS ARE PARALLEL	

DIAGRAM

TROUBLESHOOTING GUIDE

LIFT TROUBLE	CHECK
HAND CONTROLLER DOES NOT RESPOND.	CHECK LARGE RED EMERGENCY STOP BUTTON ON CONTROL BOX CHECK CONNECTOR ON HAND CONTROLLER CORD. CHECK BATTERIES CONDITION (REPLACE WITH A FULLY CHARGED BATTERY PACK).
UP & DOWN BUTTONS ON CONTROL BOX DO NOT RESPOND.	CHECK LARGE RED EMERGENCY STOP BUTTON ON CONTROL BOX. CHECK BATTERIES CONDITION (REPLACE WITH A FULLY CHARGED BATTERY PACK).
ACTUATOR DOES NOT RESPOND.	CHECK LARGE RED EMERGENCY STOP BUTTON ON CONTROL BOX. CHECK BATTERIES ARE INSTALLED CORRECTLY AND FULLY CHARGED. TEST WITH A NEW FULLY CHARGED BATTERY PACK. CHECK HAND CONTROLLER IS CONNECTED. CHECK ACTUATOR IS CONNECTED TO CONTROL BOX.
AUDIBLE "BEEP" IS HEARD FROM CONTROL BOX	BATTERIES ARE LOW, REPLACE WITH A FRESHLY CHARGED BATTERIES
ACTUATOR "STALLS" DURING LIFT	BATTERIES ARE LOW, REPLACE WITH A FRESHLY CHARGED BATTERIES. MAKE SURE TO NOT EXCEED THE LIFTING CAPACITY.
CHARGER TROUBLE	CHECK
"POWER ON" LIGHT ON CHARGER IS NOT LIT	CHECK CHARGER IS PLUGGED INTO WALL RECEPTABLE
CHARGER IS PLUGGED IN BUT "POWER ON" LIGHT IS NOT LIT	CHECK THERE IS POWER TO WALL OUTLET*
YELLOW INDICATOR DOES NOT LIGHT WHEN BATTERY IS INSERTED IN CHARGER AND THE GREEN LIGHT IS "ON"	CHECK BATTERIES ARE PROPERLY SEATED IN CHARGE.
BATTERY TROUBLE	CHECK
BATTERIES ARE PROPERLY SEATED BUT NO LIGHTS ARE VISIBLE.	CALL FOR SERVICE (CHARGER MAY BE FAULTY)
YELLOW INDICATOR LIGHT DOES NOT GO OFF AFTER SEVERAL HOURS OF CHARGING TIME.	INTERNAL BATTERIES NEED REPLACING CALL FOR SERVICE OR CONTACT YOUR MAINTENANCE DEPARTMENT.
BATTERY PACK INDICATES IT'S FULLY CHARGED WHEN IN THE CHARGER BUT WHEN PLACED IN THE LIFT WILL ONLY DO A FEW LIFTS.	AMPERAGE IN THE BATTERIES IS TOO LOW AND WILL NOT TAKE A FULL CHARGE, REPLACE INTERNAL BATTERY IN THE BATTERY PACK SYSTEM**.

(*) SOME WALL OUTLETS ARE CONTROLLED BY WALL/LIGHT SWITCHES. ENSURE POWER TO WALL OUTLET CONTINUES AFTER WALL/LIGHT SWITCH IS TURNED OFF.

(**) GENERALLY WITH LOW AMPERAGE THE ACTUATOR WILL EMIT A HUMMING NOISE INDICATING INSUFFICIENT BATTERY POWER.

TRANSFER SLING SAFETY

Patient Transfer Sling Safety Inspection & Care Guidelines

Sling and patient lift products are not interchangeable.
To ensure safe patient transfers, use only original equipment slings with this patient lift.

Due to the nature of their use, it is imperative that a patient transfer sling be inspected prior to each use. A documented monthly inspection program should be established to formally inspect all slings to ensure the safest possible transfer of a patient.

It is the manufacturer position that the number of factors influencing the life span of a patient transfer sling are so varied that a sling should be taken out of service after 2 years. This 2 years span is a guideline for the useful life of a sling and, in fact, it may be shorter or longer depending on how the slings condition is affected by the number of washing, washing temperature, detergents, disinfectants, frequency of use, patient weight and/or numerous other factors. Currently there is no method to measure the strength of a sling once it is put into service and has been laundered multiple times, without damaging the sling itself. All of our sling models are manufactured to the highest standards and under ideal circumstances it will provide many years of service.

With this in mind, the manufacturer has developed a set of visual guidelines to assess the safety of a sling currently in use. Any visual inspection, as these guidelines are, is a subjective evaluation and therefore can never be considered a guarantee of a slings safety. It will however dramatically reduce the risk of failure.

If there is ever any question concerning the safe condition of a sling, take it out of service IMMEDIATELY!


TRANSFER SLING SAFETY

HOW TO CONDUCT A VISUAL SAFETY INSPECTION


How to conduct a visual safety inspection:

Lay sling out on a flat surface so that all areas of the sling are visible. (Fig. 10)

1. Check all loops at their connection/stress points. Twists these with your fingers and look for any signs of fraying. Depending on the sling model there may be up to 12 loop points in total. See the diagram of a common sling below to assist you in locating loop points and other key areas.
2. Check the stitching of the entire slings, look for any fraying or loose stitching.
3. Check the sling for heat damage. This may be detected as an over all shrinking of the sling or may be noticed on the padded section and be identified by a shrinking or scrunching of the leg portion. Additionally, heat damage may be found on other areas by noticing a brittle or ridged/stiff feel to the fabric.
4. Check the body of the sling for any rips or holes.
5. Check the sling for signs of exposure to bleach. This may be suspected if there is fading of the sling ID and/or caution labels. Reject any sling laundered with bleach.
6. Check the sling for excessive staining. While some staining may occur through use by an incontinent patient other staining may indicate exposure to chemicals.

Failure to pass your inspection in any one of the above areas, demands for safety reasons, you remove the sling from service. Remember, if there is ever any question concerning the safe condition of a sling, TAKE IT OUT OF SERVICE!

TRANSFER SLING SAFETY

Chemical Attack

Oil, grease, creosote or paint stains are harmless but other forms of chemical attack of a sufficient degree may lead to deterioration or extreme softening of the webbing, which would lead to the fibres being rubbed off (almost as powder in extreme cases). It is best to avoid fumes, spray or mists of acid and alkalis or organic solvents. If contamination is suspected, wash out well in warm water. Avoid contact with excessive heat which is likely to affect the product.

Cleaning

Certain chemical substances mentioned above may be removed by a diluted heavy duty cleaner and warm water. Ensure that any cleaning agents are removed by rinsing thoroughly in warm water, maximum 60°C (140°F), and allowed to dry naturally away from any source of direct heat.

For normal washing and cleaning, follow the care instructions on the label of the product.

Storage

Slings should be stored away from direct sunlight, stress, pressure, excessive heat or humidity. The sling should be kept away from contact with sharp objects, corrosives or other possible damaging material.

Take note of inspection results in the log book at the end of this manual.

LAUNDERING INSTRUCTIONS

Your laundry staff or service must be made aware of these care instructions as their handling of this sling will have a direct impact on its condition.

- Machine wash in mild soap solution.
- Never use bleach.
- Rinse thoroughly.
- Tumble dry on a “cool” setting.
- Do not place in contact with a heat source.
- Do not dry clean.
- Do not iron.

This care label be found on all slings manufactured by the Manufacturer

MAGOG, CANADA
1-800-868-0441
www.bhm-medical.com

CE CONFORMS TO
ISO 10535

MODEL
SIZE
PART #

HAMMOCK-6
SMALL
THA6-S

100% POLYESTER SHELL / 100% NYLON STRAP / UHMW REINFORCEMENT

SWL 272 Kg (600 lbs)

READ INSTRUCTION FOR CORRECT COMBINATION OF
LIFT, CARRY BAR, SLING AND PATIENT.

ALWAYS INSPECT THIS SLING BEFORE EACH USE
DO NOT USE IF TORN OR DAMAGED

VISUAL SAFETY INSPECTION RECORD

- Use this record each time you conduct safety inspection.
- Retain this record on file so it can continue to be completed.

	Inspection Date		Inspected by		Conditions Notes	
	Year 1	Year 2	Year 1	Year 2	Year 1	Year 2
January						
February						
March						
April						
May						
June						
July						
August						
September						
October						
November						
December						

Serial Number:		In-Service Date:	
Model Number:			

WARRANTY

This warranty is extended only to the original purchaser/user of BHM products.

BHM Medical Inc. warrants its products to be free from defects in material under normal use and service, within the periods stated below from the date of purchase. If within such warranty period any such product shall be proven to be defective, such product shall be repaired or replaced at BHM Medical option. This warranty does not include any labor or shipping charges incurred in replacement part installation or repair of any such product. BHM Medical sole obligation and your exclusive remedy under this warranty shall be limited to such repair and/or replacement.

Patient Lifter	1 year
Tracks and installation	Life time warranty*
Weighing Devices	1 year
Accessories on Lifter	1 year
Slings	1 year
Batteries - Voyager Portable	3 months
Batteries - All other lifts	1 year
Easytrack System	1 year

For warranty service, please contact the dealer from whom you purchased the BHM Medical product. [In the event that you do not receive satisfactory warranty service, please contact BHM Medical (see contact information in Table of Contents).]

Do not return products to our factory without prior authorization. BHM Medical will issue a Return Merchandise Authorization (RMA) Number. C.O.D. shipments will be refused; all shipments to BHM Medical must be prepaid. For this warranty to be valid, the purchaser must present its original proof of purchase at the moment of the claim. The defective unit, assembly or part must be returned to BHM Medical for inspection. The part or components repaired or replaced are guaranteed for the remaining period of the initial warranty.

Limitations and Exclusions:

The warranty above does not apply to serial numbered products if the serial number has been removed or defaced.

No warranty claim shall apply where the product or any other part thereof has been altered, varied, modified, or damaged; either accidentally or through improper or negligent use and storage. Warranty does not apply to products modified without BHM Medical express written consent (including but not limited to products modified with unauthorized parts or attachments), products damaged by reason of repairs made to any component without the specific consent of BHM Medical, or to products damaged by circumstances beyond BHM Medical control. BHM Medical will solely determine evaluation of warranty claim. The warranty does not apply to problems arising from normal wear or failure to adhere to the instructions in this manual. BHM Medical Inc. slings are void of warranty if not laundered as per instructions on the Sling Label.

BHM Medical Inc. shall not be liable for damages losses or inconveniences caused by a carrier.

This warranty replaces any other warranty expressed or implicit and constitutes BHM Medical Inc. only obligation towards the purchaser. BHM Medical shall not be liable for any consequential or incidental damages whatsoever.


* Valid only if BHM Medical did the original installation. Warranty voids if tracks/installation have been modified.

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