PRESSA STABIL PLANT


STABIL-PLANT PRESS

PRESSE STABIL-PLANT

HANDBUCH FÜR
STABIL-PLANT PRESSE

PRESSA STABIL-PLANT

- MANUALE ISTRUZIONI
- OPERATOR'S HANDBOOK
- MANUEL D'UTILISATION
- GEBRAUCHSANWEISUNG
- MANUAL DE INSTRUCCIONES

CATTANI S.p.A.

PRESSA STABIL-PLANT

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Modello: *Pressa Stabil-Plant*

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Type of appliance: *B*

Use: *intermitt. service*

Protection against liquids: *common*

Level of protection against direct or indirect contact type: *B* 

| | | |
|---|----------------------------|----------------------|
|  | <i>Alternating current</i> | <i>IEC 417-5032</i> |
|  | <i>Earthing</i> | <i>IEC 417-5019</i> |
|  | <i>Type B appliance</i> | <i>IEC 878-02-02</i> |
|  | <i>Off</i> | <i>IEC 417-5008</i> |
|  | <i>On</i> | <i>IEC 417-5007</i> |

This appliance cannot work in contact with a flammable anaesthetic mixture with air, oxygen or nitrous oxide

Motor thermal protector 120° C + 1-5° C

Output power: 220 VA

Maximum flow: 14 l min

Maximum operating head for continuous service: 600 mm H₂O

Sound pressure level Stabil Press Plant 60.5 dB(A)*

**Sound pressure level tested (with canalized air) according to ISO 3746-1979 (E) regulation.*

Parameters: r = 1 - background noise ≤ 34 dB(A) - instruments: Brüel & Kjær Type 2232.

Manufactured by CATTANI S.P.A. - PARMA - ITALY

Introduction

Our Stabil-Plant Vacuum press (front page) enables the podiatrist to adapt a thermo-plastic plantar to a customized cast in a vacuum, in a practical and safe way (fig. 1).

Installation

The press must be installed in compliance with CEI 62-5/2° 1991 regulation for electro-medical appliances (the Stabil-Plant press has been designed accordingly). Prior to plugging in the press, check the specifications on the label and make sure that the mains are compatible with the appliance and protected against overcurrent according to CEI 64-8 regulations.

The press must be protected against direct and indirect contacts for class I appliances according to CEI 64-8 and 64-4 reg. The plug is equipped with earthing protective cable: do not remove the cable in any case and make sure that the socket is compatible.

Functioning (illus. 2)

Once you have plugged the press into a compatible socket, the appliance is ready to start operating. Open lid C and place the customized cast in the centre of the upper level of the press; make sure you fill all the gaps around the cast with plasticine, so that the cast can be pulled out once the plantar has been adapted.

Soften the thermoplastic sheet by using an electrical resistance or a suitable flame; place the softened sheet over the cast, close the lid and lock fastener B; press switch A to create a vacuum (in a few seconds manometer D will indicate - 8 bar); keep the press in a vacuum for a few seconds or until the thermoplastic sheet has hardened. Turn on tap E and re-establish atmospheric pressure inside the press, pull out the cast, cut out and finish the plantar.

Maintenance

Keep the press clean, with special care for the rims of the lid; when the rubber is worn out replace rubber seals, as dry rubber crumbles and allows air penetration, slowing down the creation of a vacuum. In case of failure address to an authorised technician who can supply original spare parts.

Important notice

The appliance is guaranteed for one year from date of sale, provided that guarantee card addressed to manufacturer is returned to the manufacturer reporting date of sale, retailer stamp and customer's name.

Guarantee and manufacturer liability cease in case appliances and/or plants are treated with unsuitable products, not suggested by the manufacturer and improperly used and in case they are found tampered by any kind of action performed by unable and thus unauthorised people.

The manufacturer, together with distributors, agents, authorised retailers and service engineers, is willing to supply spare parts, literature, instructions and any other useful information.

Technical requirements, improvement of the product, legislation and logistic problems, availability of products or semi-manufactured products can lead to modifications without notice.

Transport and storage

Packed equipment can be transported and stored at a temperature range of -10° C +60° C. Packages must be kept away from water and splashing and cannot tolerate humidity 70%.

Packages with the same weight can be stored in piles of three only.

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ESAM S.p.A.