# CoaguChek® XS PT Controls

REF 04696522 190 **4 tests**

## English

### Intended use

CoaguChek XS PT Controls are used for system checks and quality control testing of prothrombin time\* with the CoaguChek XS Plus or CoaguChek XS Pro meter and CoaguChek XS PT Test Strips.  
The CoaguChek XS PT Controls are intended for use for professional users.  
a) Also known as Quick's test or Quick's value or under the abbreviation PT/INR.

**Reagents - working solutions**  
Each bottle of control contains 307 mg of anticoagulated rabbit plasma (before lyophilization) and 3 mg of additives.  
Each diluent-filled dropper contains 0.84 mg of calcium chloride, 0.1 µg of preservatives in water.

### Target values and ranges

The attached code chip contains all relevant lot-specific information to perform a quality control test. The CoaguChek XS Plus or XS Pro meter displays the control range, the reading obtained and an evaluation of the result (inside or outside of the control range). The result is automatically saved as a control reading in the memory of the meter.  
Each laboratory or doctors office should establish corrective measures to be taken if values fall outside the limits. Follow the applicable local guidelines and regulations for quality control.

### Precautions and warnings

For in vitro diagnostic use.  
All components of the pack can be discarded along with household waste or, when used in a laboratory or doctors office, according to your local guidelines. Safety data sheet available for professional user on request.

### Handling

Dissolve carefully the contents of one bottle by adding the contents of 1 dropper and allow to stand closed for at least 1 minute to reconstitute. Mix carefully by gently swirling the bottle.

### Storage and stability

Store at 2-8 °C. The lyophilized control plasma is stable up to the stated expiration date.  
The reconstituted control material is stable for 30 minutes.

### Materials provided

4 bottles of lyophilized control plasma Level 1  
4 diluent-filled droppers  
1 code chip

### Materials required (but not provided)

CoaguChek XS Plus or CoaguChek XS Pro meter  
CoaguChek XS PT Test strips

### Assay

1. Prepare the meter and a test strip according to the instructions given in the operator's manual and in the instructions for use for the test strips.
2. Insert the quality control code chip in the meter.
3. Use the dropper to apply one hanging drop of the control material to the sample application area of the test strip. Save the remaining control solution until the test result is obtained, and for repeating the quality control check if necessary. Use the control material within 30 minutes of reconstitution.

### System checks

The CoaguChek XS Plus or CoaguChek XS Pro meter has a number of inbuilt system checks. Please refer to the operator's manual for details. The CoaguChek XS PT Test strip has an inbuilt quality control function.

### Reproducibility

Precision 2.8 % CV for an INR of 1.8.

For further information please refer to the appropriate operator's manual for the meter concerned, product information and the package inserts of all necessary components.

### Symbols

	Catalogue number
	Batch code/lot number
	In vitro diagnostic medical device
	Manufacturer
	Use by
	Temperature limitation (Store at)
	Consult instructions for use
	Contains sufficient for <math>n</math>- tests
	This product fulfills the requirements of the European Directive 98/79/EC for in vitro diagnostic medical devices

Significant additions or changes are indicated by a change bar in the margin.

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