YOU-Q

VR2 (PG DT)

GEbruIKSHANDLEIDING BEDIENINGSKASTEN
MANUEL D'UTILISATION - BOÎTIER DE COMMANDE
BEDIENUNGSANLEITUNG STEUERKASTEN
USER MANUAL CONTROLS
BRUGSANVISNING STYRESYSTEMER
KÄYTTÖOHJE OHJAINLAITTEET
BRUKSANVISNING STYRENHETER
BRUKSANVISNING KONTROLLE

Nederlands	4
English	14
Français	24
Deutsch.....	34
Dansk	44
Suomi	54
Svenska	64
Norsk.....	74

Nederlands

© 2009 YOU-Q

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6. Technische specificaties

Supply Voltage:	24Vdc
Operating Voltage:	16Vdc to 35Vdc
Peak Voltage:	35Vdc
Reverse Battery Voltage:	40Vdc
PWM Frequency:	20kHz \pm 0.5%
Brake Voltage:	12Vdc or 24Vdc
Brake Current:	100mA min. 1A max.
Charger Connector:	Use only Neutrik NC3MX
Batt. Charging Current:	12Arms max.
Actuator Current:	12A max.
Maximum Drive Current:	VR2 50 50A VR2 60 60A VR2 70 70A VR2 90 90A
Moisture Resistance:	IPx4
Operating Temperature:	-25°C to +50°C
Storage Temperature:	-40°C to +65°C


EMC tested on sample wheelchair:

Susceptibility:	Tested at 30V/m to EN12184 (1999) and ANSI/RESNA requirements
Emissions:	To EN55022 Class B
ESD:	IEC801 part 2

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
© 2009 YOU-Q

All rights reserved.

The information provided herein may not be reproduced and/or published in any form, by print, photoprint, microfilm or any other means whatsoever (electronically or mechanically) without the prior written authorisation of YOU-Q.

The information provided is based on general data concerning the constructions known at the time of the publication of this manual. YOU-Q executes a policy of continuous improvement and reserves the right to changes and modifications.

The information provided is valid for the product in its standard version. YOU-Q cannot be held liable for possible damage resulting from specifications of the product deviating from the standard configuration.

The available information has been prepared with all possible diligence, but YOU-Q cannot be held liable for possible errors in the information or the consequences thereof. YOU-Q accepts no liability for loss resulting from work executed by third parties.

Names, trade names, etc. used by YOU-Q may not, as per the legislation concerning the protection of trade names, be considered as being available.

2009-01

Preface.....	16
1. The controller	17
1.1 VR2 controller	17
2. Driving the wheelchair with the controller.....	19
2.1 Switching the controller on or off	19
2.2 Driving the wheelchair.....	19
2.3 Speed.....	19
3. Operating the electrical adjustment options	19
3.1 Selecting the desired adjustment option	20
4. Lights	20
5. Troubleshooting.....	21
6. Technical specifications	21
7. Technical diagrams	22
7.1 Technical diagram.....	22
7.2 Technical diagram battery charger.....	23

Preface

This manual

This user manual describes the use of the electric wheelchair via the controller. The total user manual for this electric wheelchair consists of three booklets. Read the entire user manual carefully before taking the product into use. The information provided by this manual is essential for the safe use and proper care (cleaning) of the wheelchair.

- The general user manual
- The user manual for the seating system (QLASS)
- The user manual for the controller (this booklet)

If one of these user manuals was not included with your wheelchair, please contact your dealer immediately. In addition to this user manual, there is also a service manual for qualified specialists.

This general user manual will refer, where necessary, to one of the other user manuals. This is indicated as follows:

 **GENERAL:** Refers to the general user manual.

 **CONTROLLER:** Refers to the user manual for the controllers.

 **QLASS:** Refers to the user manual for the seating system.

1. The controller

1.1 VR2 controller

A controller will usually have three basic functions:

- Driving and steering a wheelchair
- Operating electrical seat adjustments
- Charging the wheelchairs batteries

There are many different control systems for wheelchairs on the market. If the controller on your wheelchair does not resemble the one in figure 1.1, contact your dealer.

Figure 1.1 From the left to the right ; VR2 drive only, VR2 Lights, VR2 Seats&Lights

VR2 by PG Driving Technologies is a collective name for the entire control system of the wheelchair. The wheelchair is operated by means of a controller, which includes the following components:

The controller


Figure 1.2 From the left to the right ; VR2 drive only, VR2 Lights, VR2 Seats&Lights

Part	Function
A. Joystick	In the 'drive' mode: driving and steering In the 'adjustment options' mode: <ul style="list-style-type: none"> • Left/right to select the adjustment options • Front/back to select the adjustment options
B. On/off switch	Switching the controller on or off
C. Horn	Warning signal with sound
D1: Battery indicator	Displays the power level of the battery
D2: Maximum speed	Displays the maximum speed limit as set by the user
D3: Adjustment option	Displays the selected adjustment options
E. "Mode" button	Changing between the 'driving' and the 'adjustment options' mode
F1: Speed regulator	Reduce driving speed (slower)
F2: Speed regulator	Increase driving speed (faster)
G. Charge connector	Input for the battery charger
H. Charge plug for the battery charger	Connector for the battery charger
I. Lights button	Switching the lights on or off
J. Hazard lights	Warning signal with lights
K1: Direction indicator left	Switches the left direction indicator on or off
K2: Direction indicator right	Switches the right direction indicator on or off

2. Driving the wheelchair with the controller

2.1 Switching the controller on or off

To be able to drive or operate the electronic adjustment options of the wheelchair, the controller must be switched on. Press the on/off button (B in figure 1.2).

2.2 Driving the wheelchair

Driving an electric wheelchair is done by operating a joystick. Move the joystick forwards and the wheelchair will also move forwards. Steer left and right and the wheelchair will turn.

2.3 Speed

The maximum speed can be controlled by the speed regulator on the controller (F1 en F2 in figure 1.2). De maximum speed is displayed by the LED lights above the speed regulators (D2 figure 1.2). Speed can be controlled with the joystick while driving. If the joystick is moved a little, the wheelchair will move more slowly.

3. Operating the electrical adjustment options

Not every wheelchair has been equipped with electronic adjustment options. We make a distinction between four different adjustment options to the seating system:

- Tilt adjustments
- High/low adjustments
- Backrest adjustments
- Legrest adjustments

To use electrical adjustment options, the VR2 Seats&Lights controller is required (maximum of two options).

Figure 3.1 Electronic adjustment options

3.1 Selecting the desired adjustment option

1. Switch on the controller
2. Press on the 'mode' button (E in figure 1.2) to select the 'adjustment options' mode. The controller will now be in the 'adjustment options' mode. It is also possible to follow the status.

The joystick is used to select and operate the electronic adjustment options.

3. Move the joystick to the left or the right to select the desired adjustment option. The selected adjustment option will be visible on the display screen.
4. Moving the joystick forwards and/or backwards will activate the selected electronic adjustment option (see table 3.1). Move the joystick forwards or backwards until the desired adjustment option has been attained.

	Move joystick backwards	Move joystick forwards
Tilt adjustment	The entire chair will tilt backwards	The entire chair will tilt forwards
High/low adjustments	The entire chair will be raised	The entire chair will be lowered
Backrest adjustments	The backrest will tilt backwards	The backrest will tilt forwards
Legrest adjustments	The legrest angle will increase, the footplate will be raised	The legrest angle will decrease, the footplate will be lowered

Table 3.1 Electrical adjustments

5. To return to the 'drive' mode: Press on the 'mode' button to select the 'drive' mode.

Note: If you move the chair while adjusting the high/low option, the speed will be reduced.

4. Lights

Not every wheelchair has been equipped with lights. You can activate the following functions only by using VR2 Lights and/or VR2 Seat&Lights:

- Lights (I in figure 1.2)
- Hazard Lights (J in figure 1.2)
- Direction indicators (K1 & K2 in figure 1.2)

5. Troubleshooting

If the wheelchair will not function while the batteries are fully charged, check the following points before consulting your dealer:

- Switch the controller off and then switch it on again. Check to see if the malfunction has been solved.
- Check if the free wheel switch was switched to Drive.
- Check if the joystick was in the 0 position when the controller was switched on. In other words, the joystick must not be moved when the controller is being switched on or off.

Malfunctions list

An extended malfunction list can be found in the appendix of the service manual (for qualified specialists only). The service manual can be found on www.YOU-Q.com

6. Technical specifications

Supply Voltage:	24Vdc
Operating Voltage:	16Vdc to 35Vdc
Peak Voltage:	35Vdc
Reverse Battery Voltage:	40Vdc
PWM Frequency:	20kHz \pm 0.5%
Brake Voltage:	12Vdc or 24Vdc
Brake Current:	100mA min. 1A max.
Charger Connector:	Use only Neutrik NC3MX
Batt. Charging Current:	12Arms max.
Actuator Current:	12A max.
Maximum Drive Current:	VR2 50 50A VR2 60 60A VR2 70 70A VR2 90 90A
Moisture Resistance:	IPx4
Operating Temperature:	-25°C to +50°C
Storage Temperature:	-40°C to +65°C
EMC tested on sample wheelchair:	
Susceptibility:	Tested at 30V/m to EN12184 (1999) and ANSI/ RESNA requirements
Emissions:	To EN55022 Class B
ESD:	IEC801 part 2

7. Technical diagrams

7.1 Technical diagram

The technical diagram can also be found in the battery box.


Figure 7.1:

7.2 Technical diagram battery charger

The controller's standard configuration includes a '3-pin connection'. Ensure that the battery charger is properly connected so that the 'negative pole' and the 'inhibit' are connected, enabling the system to prevent the wheelchair from moving when the battery is being charged.


Figure 7.2

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YOU-Q B.V.