 HENRY SCHEIN®

900-3250

HS-CU 1000™

HALOGEN CURING LIGHT

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Users manual for *HS-CU 1000*

Table of Contents

SETTING UP	1
OPERATION	1
HIGH POWER MODE	1
SOFT-UP MODE	2
BUILT-IN ILLUMINOMETER/LIGHT METER.....	2
MAINTENANCE	3
TECHNICAL DATA	3
TROUBLE SHOOTING.....	3

SETTING UP :

- (1) Unpack the curing light. Place on firm flat surface or wall mount if preferred.
- (2) Insert the fiber optic probe into the handpiece until firmly seated. Attach the protective shield onto the probe if required.
- (3) Securely place the handpiece on the rest.
- (4) Plug the unit into an AC outlet within a range from 100V to 240V. (**The machine operates with a universal power supply)
- (5) The HS-CU 1000 light curing system is ready for use.

OPERATION :

- (1) Turn on the power switch located on the left side of base unit.
- (2) The LEDs lights will run through a power up test. When the test is completed, a “beep” sound will confirm. The unit will retain previous curing mode set up in memory.
- (3) Apply composite in mm’s as indicated by composite manufacturer. Select curing mode to begin curing process.

HIGH POWER MODE :

By pressing the “HIGH POWER” button, you will have 3 choices of curing time- **5 sec.**, **10 sec.**, and **20 sec.** – all working at full strength.

**The curing cycles are activated by pressing the trigger switch on the handpiece, which simultaneously activates the magnetic fan for cooling.

- 5 sec.** – The curing cycle ends with a “beeping” signal at 5 seconds, but the fan running will continue its automatic cooling phase (varies from mode to mode).
- 10 sec.** – A “beeping” signal is sounded at 5 second intervals before the auto shutdown.
- 20 sec.** – A “beeping” signal is sounded at 10 second intervals before the auto shutdown.

SOFT-UP MODE :

The button “SOFT-UP” works the curing cycle with increasing intensity from low to high power light output.

- 10 sec.** – This curing cycle begins with 5 seconds of climbing intensity until reaching the maximum and stays at maximum curing power for 5 seconds.
 (“Beeping” signal sounds at 5 second intervals.)
- 20 sec.** – This curing cycle begins with 7 seconds of climbing intensity until reaching the maximum and stays at maximum curing power for 13 seconds.
 (“Beeping” signal sounds at 10 second intervals.)
- 30 sec.** – This curing cycle begins with 10 seconds of climbing intensity until reaching the maximum and stays at maximum curing power for 20 seconds.
 (“Beeping” signal sounds once at 10 seconds, twice at 20 seconds and again at 30 seconds.)

NOTE: The mode is designed to prevent composite shrinkage during polymerization. It is also very beneficial to use when curing dual-cure resin cement, all-ceramic and poly-ceramic systems.

****NOTE:** pressing the trigger switch on the hand-piece will **start or stop** the curing light at anytime.

BUILT-IN ILLUMINOMETER/LIGHT METER :

This feature is designed to measure intensity of the lamp output and use these readings to auto adjust the curing time to compensate for any light deterioration. **All of the** aforementioned operations are **FULLY COMPUTERIZED** and **AUTOMATED** wherever necessary. Any requirement to adjust the curing time will be shown by LED indicators as follows:

- GREEN** – This colour indicates that the light intensity is within “**Accepted**” curing range at **over the 600mw/cm2**.
- YELLOW** – This colour indicates that the light intensity has fallen into a weaker-than-normal range **between 200~600mw/ cm2**. At this point, the curing-time adjuster will **EXTEND** the curing time relative to HOW MUCH the halogen bulb has deteriorated, using the following guidelines.
 - (a) If the light intensity reading shows between **400 MW/ cm2 to 600mw/ cm2**,
The selected normal curing time will be extended by **3 seconds**
 - (b) If the light intensity reading shows between **200 MW/ cm2 to 400mw/ cm2**,
The selected normal curing time will be extended by **6 seconds**

** Any changes will be stored in memory data throughout the same day.

****This memory data can be erased by completely switching off the main power located on the side of the curing unit, should the user wish to remeasure the intensity for subsequent operation.**

- RED** – Indicates that the light intensity has fallen **below 200mw/ cm2** . At this point, a continuous short “beeping” signal will sound as a reminder to replace your halogen bulb.

MAINTENANCE :

1 - The fibre-optic tip

After each patient use, the fibre-optic tip should be processed by disinfection and/or sterilization. It is autoclaveable up to 135°C.

2 - Any dirt or stains on the control unit or surface of the handpiece use a cloth wipe or neutral cleansing/disinfection agent to clean.

NOTE: Never use a solvent or volatile oil to for cleaning.
This will cause surface label printing to fade.

TECHNICAL DATA :

Main supply voltage:

100-240Vac, 50/60Hz

Effect:

A

Classification:

Class I

Protection against electric shock:

Type B applied part ⚡

Protection against harmful ingress of water:

Ordinary equipment, IPX0

Mode of operation:

Continuous operation or 30seconds

TROUBLE SHOOTING: (Guideline for fault diagnosis) :

PROBLEM	CAUSE	METHOD
No power even when switch is on	Cord or plug connections	Inspect the connection of power Cable plug and co-axial cable
No power even when switch is on	Fuse	Inspect fuse and change if needs, If it fails again, please send for repair
No light from Handpiece	Halogen bulb not working	Replace bulb
Unexpected "beeping"	Some buttons are pressed	Check all buttons
No light after a few uses (if bulb okay)	Cooling fan not working	Send for check and repair

**** In case of more serious malfunctions, please return to manufacture for repair.**

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