FIDIS™ anti-CCP

REF **MX 104**

DEFINITION

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Pratesi F., Tommasi C., Anzilotti C., Puxeddu I., Sardano E., Di Colo G., Migliorini P. Antibodies to a new viral citrullinated peptide, VCP2: fine specificity and correlation with anti-cyclic citrullinated peptide (CCP) and anti-VCP1 antibodies. *Clin. Exp. Immunol.* 2011, 164, 337.

Anzilotti C, Merlini G, Pratesi F, Tommasi C, Chimenti D, Migliorini P. Antibodies to viral citrullinated peptide in rheumatoid arthritis. *J Rheumatol.* 2006 Apr;33(4):647-51.

Arnett FC, Edworthy SM, Bloch DA, et al. The American Rheumatism Association 1987 revised criteria for the classification of rheumatoid arthritis. *Arthritis Rheum* 1988;31(3):315-24.

Van Venrooij WJ, Hazes JM, Visser H. Anticitrullinated protein/peptide antibody and its role in the diagnosis and prognosis of early rheumatoid arthritis. *Neth J Med* 2002;60(10):383-8.

Nienhuis RL, Mandema E, Smids C. A new serum factor in patients with rheumatoid arthritis. The antiperinuclear factor. *Ann Rheum Dis* 1964;23:302-05.

Young BJ, Mallya RK, Leslie RD, et al. Anti-keratin antibodies in rheumatoid arthritis. *Br Med J* 1979;2:97-9.

Hoet RM, Boerbooms AM, Arends M, et al. Antiperinuclear factor, a marker autoantibody for rheumatoid arthritis: colocalisation of the perinuclear factor and profilaggrin. *Ann Rheum Dis* 1991;50:611-8.

Sebbag M, Simon M, Vincent C, et al. The antiperinuclear factor and the so-called antikeratin antibodies are the same rheumatoid arthritis-specific autoantibodies. *J Clin Invest* 1995;95:2672-9

Schellekens GA, de Jong BA, van den Hoogen FH, et al. Citrulline is an essential constituent of antigenic determinants recognized by rheumatoid arthritis-specific autoantibodies. *J Clin Invest* 1998;101(1):273-81

Vossenaar ER, van Venrooij WJ. Anti-CCP antibodies, a highly specific marker for (early) rheumatoid arthritis. *Clin Applied Immunol Rev* 2004;4:239-62.

Pruijn GJ, Vossenaar ER, Drijfhout JW, et al. Anti-CCP antibody detection facilitates early diagnosis and prognosis of rheumatoid arthritis. *Current Rheumatology Reviews* 2005;1(1):1-7.

Rantapaa-Dahlqvist S, de Jong BA, Berglin E, et al. Antibodies against cyclic citullinated peptide and IgA rheumatoid factor predict the development of rheumatoid arthritis. *Arthritis Rheum* 2003; 48(10):2741-49

Meyer O, Labarre C, Dougados M, et al. Anticitrullinated protein / peptide antibody assays in early Rheumatoid Arthritis for predicting five year radiographic damage. *Ann Rheum Dis* 2003; 62:120-26

Aletaha D, Neogi T, Silman AJ, et al. Rheumatoid Arthritis Classification Criteria. An American College of Rheumatology/European League against Rheumatism Collaborative Initiative. *Arthritis Rheum* 2010; 62 (9):2569-81

Susan S. Copple, Thomas B. Martins, C. Masterson, E. Joly, Harry R. Hill. Comparison of Three Multiplex Immunoassays for Detection of Antibodies to Extractable Nuclear Antigen Using Clinically Defined Sera. *Ann. N.Y Acad. Sci.* 1109: 464-472 (2007).






Buliard A, Fortenfant F, Ghillani-Dalbin P, Musset L, Oksman F, Olsson N.O. Related Articles, Links [Analysis of nine autoantibodies associated with systemic autoimmune diseases using the Luminex technology. Results of a multicenter study]. *Ann Biol Clin (Paris)*, 2005; Jan-Feb, 63 (1): 51-8.

Rouquette A-M, Desgruelles C, Laroche P. Evaluation of the new multiplexed immunoassay, FIDIS, for simultaneous quantitative determination of antinuclear antibodies and comparison with conventional methods. *Hematology Laboratory, Tenon Hospital, 75020 Paris, France*
Am J Clin Pathol., 2003; 120 (5): 676-81.

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REF **MX 104**

INTENDED USE

FIDIS™ anti-CCP () is a semi-quantitative homogeneous fluorescent-based microparticles immunoassay using flow cytometry readings. It is designed for the simultaneous detection of 4 autoantibodies (IgG isotype) directed against 4 synthetic citrullinated peptides : HCP1, HCP2, VCP1 and VCP2. HCP1 and HCP2 are from human protein origin whereas VCP1 and VCP2 are from viral protein origin. These autoantibodies are able to detect cyclic citrullinated peptides (CCP) and deaminated fibrinogen and are associated to anti-citrullinated peptide/protein antibodies (ACPA).

SUMMARY AND DIAGNOSTIC VALUES

Rheumatoid arthritis is one of the most widespread systemic auto-immune diseases. The prevalence of the disease reaches 0.5 to 1% of adult population. The diagnostic of rheumatoid arthritis is essentially based on clinical, radiological and immunological characteristics.

The measurement of Rheumatoid Factor (RF) is widely used. Although RF test is sensitive, the specificity is low. Indeed, RF are found in patients with rheumatoid diseases, inflammatory or auto-immune diseases, as well as patients subject to chronic infections and in healthy population.

In 1998, Schellekens et al. showed that auto-antibodies directed against synthetic citrullinated peptides were very specific for rheumatoid arthritis. Further studies shown that auto-antibodies directed against cyclic citrullinated peptides were specific for rheumatoid arthritis and were more sensitive than linear peptide.

It has been shown that anti-CCP antibodies can be detected at the onset of the disease, before any kind of clinical symptoms. Several publications show that a high level of anti-CCP antibodies is correlated to a very active disease. These results demonstrate the importance of the anti-CCP antibodies for the diagnosis of rheumatoid arthritis in the early state of the disease.

In 2008, the search for anti-CCP antibodies was introduced to the Nomenclature of the Medical Biology's acts in France.

In 2010, the American College of Rheumatology (ACR) and the European League Against Rheumatism (EULAR), published criteria for the classification of rheumatoid arthritis. For the first time, serological criteria included the measurement of anti-citrullinated peptide/protein antibodies (ACPA), such as anti-CCP antibodies for example.

ASSAY PRINCIPLE

FIDIS™ anti-CCP () is based on the use of distinct uniform size color-coded microspheres and a benchtop flow cytometer interfaced to digital signal processing hardware and software. A red diode laser beam in the flow cytometer classifies each set of microspheres on the basis of its unique fluorescence intensity (red to orange) thus identifying which analyte is being tested. At the same time, a green laser beam illuminates the external second molecule fluorescence to quantify the reaction related to the specific analyte.

Each antigen (synthetic peptides : HCP1, HCP2, VCP1 and VCP2) required for the assay is coupled to an individual set of microspheres through its surface functional groups.

The different peptides coupled microspheres are mixed together, constituting the final microsphere reagent.

The **FIDIS™ anti-CCP** () allows the detection of 4 autoantibody specificities: anti-HCP1, anti-HCP2, anti-VCP1 and anti-VCP2.

The test is performed in a 96 well microplate with a filtering membrane.

- In the first step, the sample is distributed in each well containing the microspheres mixture. If this sample contains one or more of the suspected antibodies, this(ese) antibody(ies) bind to the corresponding antigen(s) on the various sets of microspheres.
- After incubation, a wash step using a filtration process removes the unbound antibodies.
- A phycoerythrin labelled **anti-human IgG** conjugate is then added that binds to the previously bound antibodies.
- A final wash step allows to stop the reaction.
- The reaction is then directly measured by the flow cytometer, which differentiates each set of microspheres according to its fluorescence color while simultaneously measuring the average fluorescence emitted by the conjugate.
- A calibration system allows the determination of the titer (AU/mL) of each sample by interpolation for each following antigenic specificity : HCP1, HCP2, VCP1 and VCP2.

The **FIDIS™ anti-CCP** kit could be used with **CARIS™** system (diluting and dispensing device).

CAUTION

- Reagents in solution contain less of 0.1% (w/v) sodium azide and 0.02% (w/v) of Proclin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.
- Avoid to use reagents if signs of contamination or other visible changes occur.
- The **FIDIS™ anti-CCP** (Theradag) has been developed according to CE Directives 67/548/EEC and 1999/45/EC regarding classification, packaging and labelling of dangerous preparations.
- The **FIDIS™ anti-CCP** (Theradag) has been optimized for use as described in this procedure. Do not substitute other manufacturer reagents. Dilution or alteration of these reagents may also alter the performance of the test. Follow thoroughly the test procedure to insure optimal performance.
- *Calibrators and controls are from human origin. The human sera used in the preparation of these products were tested and found non-reactive for antibodies to HIV-1, HIV-2, anti-HCV and Hepatitis B antigen. Because no test method can offer complete assurance that infectious agents are absent, handle as if capable of transmitting infectious diseases.*

SPECIMEN COLLECTION AND HANDLING

- The test is performed on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- If the test is not run immediately, the samples should be stored at +2°C/+8°C for a maximum of 5 days, for longer storage, frozen undiluted samples at -20°C.
- To avoid non-specific binding, it is recommended to centrifuge and filter the cloudy samples or samples frozen for more than 6 months.

REAGENTS SUPPLIED

96 wells microplate with filtering membrane and lid. MP	1 plate
Vial (A) of color-coded microsphere set of 4 sensitized by HCP1, HCP2, VCP1 and VCP2 synthetic peptides. <u>Lyophilized</u> (to be reconstituted with the buffer named D) MICROSPHERES	qs 6mL
Vial (B) of sample dilution buffer (white vial) <u>Ready to use</u> DIL SPE	2 x 115mL
Vial of calibrator* <u>Ready to use</u> Each titer is printed on the vial label CAL	1 x 1.5mL
Vial of positive control* concentrated. This control has a standard reactivity, which provides evidence of the proper reagents activity and proper assay performance. <u>To be diluted</u> Expected values are printed on the vial label. CONTROL +	1 x 250 µL
Vial of negative control concentrate <u>To be diluted</u> CONTROL -	1 x 250µL
Vial of anti-human IgG coupled to phycoerythrin <u>Ready to use</u> CONJ IgG	1 x 12mL
Vial (C) of washing buffer (black vial) <u>Ready to use</u> BUF WASH	1 x 100mL
Vial (D) of reconstitution buffer for the microsphere set <u>Ready to use</u> BUF MICROSPHERES	1 x 6mL

Calibrator titers* are expressed in arbitrary units per ml (AU/ml).

ADDITIONAL MATERIAL – NOT SUPPLIED

- Precision pipettes capable of accurately delivering 5µL to 1000µL
- Multichannel Pipettes or dispensers capable of accurately delivering 40µL to 300µL or 5µL to 2 mL
- Vortex mixer
- Laboratory timer to monitor incubation steps
- Ultrasonic bath
- Absorbent towel
- Serological pipettes
- Microplate sealing films

STABILITY AND STORAGE

- Store reagents in their original packaging at +2°C to +8°C.
- Do not freeze reagents.
- Do not use kits beyond the expiration date.
- After use, store all components immediately back at +2°C to +8°C.

TEST SET-UP

FIDIS™ Washer

☞ Check the tubing connecting the pump to the stand and the manometer setting (wheel completely closed).

Using FIDIS™ Analyzer and MLX-BOOSTER™ Software

☞ See the User's Manual provided with the **FIDIS™ Instrument** for detailed instructions on running the equipment and for calculation. For additional information and/or trouble shooting problems, please contact  subsidiaries or distributors.

1. Start the run as described in the User's Manual.
NOTE: The FIDIS™ takes 30 minutes to warm up after being turned on. A new warm up is necessary after 4 hours of system inactivity.
2. Calibration and controls are described in the User's Manual. Calibration and controls should be routinely performed once per month and each time a new lot of Sheath fluid reagent is used in order to insure optimal instrument performance.
Calibration should also be performed when the temperature is out of the determined range shown on the run batch screen.
3. Programming a batch or a multibatch of test protocol on the FIDIS™ as described in the User's Manual.
4. Load the microplate into the plate holder of the FIDIS™ as described in the User's Manual.
5. Analyze the results according the User's Manual.
6. When finished for the day, perform the sanitizing and soak operations, prior to turning the Analyzer off according to the shutting down procedure described in the User's Manual.

Reagents Preparation

☞ Remove the individual components from storage, allow them to warm up to **room temperature** (+18°C to +25°C), and mix them well.

1. Microsphere preparation

- Reconstitute the microsphere solution by adding the vial **D** of reconstitution buffer to the vial **A** of lyophilized microspheres. Wait 5 minutes and vortex.
- Stable 2 months at +2°C to +8°C, after reconstitution.

2. Preparation of samples and controls.

- Dilute the samples and controls at **1:201** in dilution buffer (**B**).
Ex.: 10µL sample in 2000µL dilution buffer (B).
- Vortex vigorously.

3. Assay configuration

Use the work-sheet included in the kit to identify the location of the samples.

a. *When setting up the test, systematically take into account the following well requirements:* ⇒ See examples below

- ✗ 1 "reagent-blank" well
- ✗ 1 well for negative control
- ✗ 1 well for positive control
- ✗ 2 "calibrator" wells

b. *Calculation of the correct number of wells necessary and their location.*

In the following examples, 2 different configurations are described, according to the number of wells needed, and the availability of unused wells. The sample dispensing must be systematically carried out in a column. Do not leave empty wells.

The same microplate can be used for more than one serie of tests if unused wells are protected by microplate sealing film. Any single well can only be used once.

To avoid any error in using a well more than once, identify the used wells with a marker (see example 2 below).

c. *Programming of the test protocol on the FIDIS™*

☞ Refer to the FIDIS™ user's manual: "Creating a batch or a multi-batch".

Example 1:

- *The microplate is totally blank (no well has yet been used): Use A1 as the 1st well.*

	1	2	3	4	5	6	7	8	9	10	11	12
A	Blank	S4	S12	S20								
B	Ctrl-	S5	S13	S21								
C	Ctrl+	S6	S14	S22								
D	Cal	S7	S15	S23								
E	Cal	S8	S16	S24								
F	S1	S9	S17	S25								
G	S2	S10	S18									
H	S3	S11	S19									

Example 2:

- 16 wells have already been used: Use A3 as the 1st well.

	1	2	3	4	5	6	7	8	9	10	11	12
A			Blank	S4	S12	S20	S28	S36	S44	S52		
B			Ctrl-	S5	S13	S21	S29	S37	S45	S53		
C			Ctrl+	S6	S14	S22	S30	S38	S46			
D			Cal	S7	S15	S23	S31	S39	S47			
E			Cal	S8	S16	S24	S32	S40	S48			
F				S1	S9	S17	S25	S33	S41	S49		
G				S2	S10	S18	S26	S34	S42	S50		
H				S3	S11	S19	S27	S35	S43	S51		

ASSAY PROCEDURE

1. Microspheres dispensing

- Before to start the assay and if it is necessary, protect the unused wells with a microplate sealing film.
- **Vortex the microspheres reagent vigorously for 20s** and dispense 50µL in each well to be used.

PLEASE NOTE:

- If **CARIS™** is used, the microspheres should be sonicated for 5 minutes.
- If a dispensing/diluting device is used (like **CARIS™**), the procedure steps 1 and 2 should be reversed:
⇒ dispense diluted sera in first and then, dispense microspheres

2. Addition and Incubation of the samples

- Dispense 100µL of sample dilution buffer (**B**) for reagent blank.
- Dispense 100µL of prediluted controls, of ready to use calibrator and of prediluted samples.
- Cover the plate and incubate 30 minutes at room temperature, away from direct sunlight and without shaking.

3. Wash step 1

Wash the plate 2 times with washing buffer (**C**) using filtration unit:

- Remove the microplate lid (verify if the vacuum break valve is closed) and place it on the washer.
- Start the pump. Stop it when the whole incubation liquid has been removed in totality and open the vacuum break valve. Close the vacuum break valve.
- Dispense 300µL of washing buffer (**C**) in all the used wells.
- Start the pump. Stop it when the whole buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- Dispense 300µL of washing buffer (**C**) in all the used wells.
- Start the pump. Stop it after 5 additional seconds when all the buffer has been removed and open the vacuum break valve. Close the vacuum break valve.

- Remove the microplate from the washer and remove the residual buffer by tapping vigorously the plate 10 times on an absorbent towel.
- Place again the microplate on the washer and start the pump for 5 seconds. Stop it and quickly open the vacuum break valve.
- Remove the microplate from the washer and remove the residual buffer by blotting vigorously the plate on an absorbent towel.
- Place the plate on a totally dry surface before starting the conjugate incubation.

4. Incubation of the conjugate

- Dispense 100µL of ready to use conjugate in each well used.
- Cover the plate and incubate 30 minutes at room temperature away from direct sunlight and without shaking. The incubation time starts after the conjugate has been added to all wells. If this timing is not followed, the results might be erroneous.

Note : the conjugate is photosensitive.
⇒ After using, close the vial carefully.

5. Wash step 2

Wash the plate using the filtration unit on 1 cycle with dilution buffer (**B**).

Use two different reagent reservoirs for the 2 washing steps.

*Either clean and dry correctly the reservoir between both wash steps (avoid to mix both **B** and **C** buffers).*

- Remove the microplate lid (verify that the vacuum break valve is closed) and place it on the washer.
- Start the pump. Stop it when the whole incubation liquid has been removed and open the vacuum break valve. Close the vacuum break valve.
- Dispense 100µL of dilution buffer (**B**) in all the used wells.
- Start the pump. Stop it when the whole buffer has been removed and open the vacuum break valve. Close the vacuum break valve.
- Remove the microplate from the washer and remove the residual buffer by blotting the plate on an absorbent towel.
- Dispense 100µL of dilution buffer (**B**) in all the used wells. Place the plate on a totally dry surface at room temperature.

The reading must be doing in the hour following adding the solution (B). During this time, keep the plate at room temperature away from direct sunlight.

6. Reading of the test

• Follow the **MLX-BOOSTER™** and **FIDIS™** user's manual: Processing batches.

VALIDATION CRITERIA OF THE RESULTS

Calibrator, negative and positive controls have to be run with every batch of samples to ensure that all reagents and procedures performed properly.

In order to validate the results, all the criteria listed below must be met. Otherwise, the test is invalid and must be repeated.

- a. The positive controls should show a value within the limits printed on the corresponding vial labels.
- b. The negative controls should be less than 25 AU/mL.

The microsphere set contains Internal standard beads allowing to verify the presence of serum in wells, as well as the good respect of the test procedure. Wells presenting a not corresponding fluorescent signal will be no valid with MLX BOOSTER™ software.

INTERPRETATION OF RESULTS

Determine for each specificity, the presence or the absence of antibodies, according to the following table:

Arbitrary Units (AU/mL)	< 25 AU/mL	25 - 30 AU/mL	> 30 AU/mL
Interpretation	Negative	Equivocal(*)	Positive

(*) Borderline results should be controlled on a second sample and the interpretation of the results should be done in the frame of additional testing and taking into account the clinical status of the patient.

Each laboratory should establish and maintain its own references (normal) range values, based on the patient population and other local factors.

Determine the final status of the sample:

- a. A sample is considered « **positive** » when at least one result from HCP1, HCP2, VCP1 and VCP2 antibody specificity is found « positive ».

(In other terms, a result is considered “positive” if at least one result from each antibody specificity is higher than 30 AU/mL)

- b. A sample is considered « **negative** » when the result from each specificity is less than 25 AU/mL.

- c. A sample is considered « **equivocal** » when at least one of the result is comprised between 25 AU/mL and 30 AU/mL.

CALCUL OF RESULTS

The results are automatically calculated by MLX-BOOSTER™ software and can be printed for each assay.

LIMITATION

Hemolytic, lipemic, icteric samples or samples with abnormal concentration of IgG and/or complement levels or samples with rheumatoid factor may confound the results of this assay.

The serums from Immunodeficient patients will give a no valid result.

QUALITY CONTROL

It is recommended to use internally and externally sourced control material for the different specificities.

REPRODUCIBILITY

Antigenic Specificity	Within-run (30 tests in the same run)		Between-run (2 tests in 20 different runs)	
	Mean value AU/mL	CV (%)	Mean value AU/mL	CV (%)
HCP1	48	6.5	55	9
HCP2	37	8.5	40	8
VCP1	36	6.5	40	9
VCP2	42	7	41	7

BIBLIOGRAPHY

WO2004087747 (2004-10-14): Citrullinated synthetic peptides and uses thereof, Inv: Migliorini Paola.

WO/2011/061720 (2011.05. 26): Viral citrullinated peptides and uses thereof, Inv: PRATESI, Federico; MIGLIORINI, Paola.

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Pratesi F., Tommasi C., Anzilotti C., Puxeddu I., Sardano E., Di Colo G., Migliorini P. Antibodies to a new viral citrullinated peptide, VCP2: fine specificity and correlation with anti-cyclic citrullinated peptide (CCP) and anti-VCP1 antibodies. *Clin. Exp. Immunol.* 2011, 164, 337.

Anzilotti C, Merlini G, Pratesi F, Tommasi C, Chimenti D, Migliorini P. Antibodies to viral citrullinated peptide in rheumatoid arthritis. *J Rheumatol.* 2006 Apr;33(4):647-51.

Arnett FC, Edworthy SM, Bloch DA, et al. The American Rheumatism Association 1987 revised criteria for the classification of rheumatoid arthritis. *Arthritis Rheum* 1988;31(3):315-24.

Van Venrooij WJ, Hazes JM, Visser H. Anticitrullinated protein/peptide antibody and its role in the diagnosis and prognosis of early rheumatoid arthritis. *Neth J Med* 2002;60(10):383-8.

Nienhuis RL, Mandema E, Smids C. A new serum factor in patients with rheumatoid arthritis. The antiperinuclear factor. *Ann Rheum Dis* 1964;23:302-05.

Young BJ, Mallya RK, Leslie RD, et al. Anti-keratin antibodies in rheumatoid arthritis. *Br Med J* 1979;2:97-9.

Hoet RM, Boerbooms AM, Arends M, et al. Antiperinuclear factor, a marker autoantibody for rheumatoid arthritis: colocalisation of the perinuclear factor and profilaggrin. *Ann Rheum Dis* 1991;50:611-8.

Sebbag M, Simon M, Vincent C, et al. The antiperinuclear factor and the so-called antikeratin antibodies are the same rheumatoid arthritis-specific autoantibodies. *J Clin Invest* 1995;95:2672-9

Schellekens GA, de Jong BA, van den Hoogen FH, et al. Citrulline is an essential constituent of antigenic determinants recognized by rheumatoid arthritis-specific autoantibodies. *J Clin Invest* 1998;101(1):273-81

Vossenaar ER, van Venrooij WJ. Anti-CCP antibodies, a highly specific marker for (early) rheumatoid arthritis. *Clin Applied Immunol Rev* 2004;4:239-62.

Pruijn GJ, Vossenaar ER, Drijfhout JW, et al. Anti-CCP antibody detection facilitates early diagnosis and prognosis of rheumatoid arthritis. *Current Rheumatology Reviews* 2005;1(1):1-7.

Rantapaa-Dahlqvist S, de Jong BA, Berglin E, et al. Antibodies against cyclic citullinated peptide and IgA rheumatoid factor predict the development of rheumatoid arthritis. *Arthritis Rheum* 2003; 48(10):2741-49

Meyer O, Labarre C, Dougados M, et al. Anticitrullinated protein / peptide antibody assays in early Rheumatoid Arthritis for predicting five year radiographic damage. *Ann Rheum Dis* 2003; 62:120-26

Aletaha D, Neogi T, Silman AJ, et al. Rheumatoid Arthritis Classification Criteria. An American College of Rheumatology/European League against Rheumatism Collaborative Initiative. *Arthritis Rheum* 2010; 62 (9):2569-81

Susan S. Copple, Thomas B. Martins, C. Masterson, E. Joly, Harry R. Hill. Comparison of Three Multiplex Immunoassays for Detection of Antibodies to Extractable Nuclear Antigen Using Clinically Defined Sera. *Ann. N.Y Acad. Sci.* 1109: 464-472 (2007).

Buliard A, Fortenfant F, Ghillani-Dalbin P, Musset L, Oksman F, Olsson N.O. Related Articles, Links [Analysis of nine autoantibodies associated with systemic autoimmune diseases using the Luminex technology. Results of a multicenter study].




Ann Biol Clin (Paris), 2005; Jan-Feb, 63 (1): 51-8.

Rouquette A-M, Desgruelles C, Laroche P. Evaluation of the new multiplexed immunoassay, FIDIS, for simultaneous quantitative determination of antinuclear antibodies and comparison with conventional methods. Hematology Laboratory, Tenon Hospital, 75020 Paris, France
Am J Clin Pathol., 2003; 120 (5): 676-81.

SUMMARY OF THE TEST PROCEDURE

	VOLUME TO BE DISTRIBUTED	REAGENTS	INCUBATION CONDITIONS
Dispense 50µL of microspheres reagent in each well to be used			
Incubation of the samples	100µL	Dilution buffer (B): reagent blank	30 min. Room temperature
	100µL	Diluted negative control in dilution buffer (B)	
	100µL	Diluted positive control in dilution buffer (B)	
	100µL	Ready-to-use Calibrator (in duplicate)	
	100µL	Diluted samples in dilution buffer (B)	
Washing 1	Wash twice in washing buffer (C) (300µL/well)		
Incubation of the conjugate	100µL	Ready to use Anti-IgG conjugate	30 min Room temperature
Washing 2	Wash 1 time in Dilution buffer (B) (100µL/well)		
Reading	Add 100µL/well of Dilution buffer (B) and proceed immediately with the reading in FIDIS™ Instrument The reading must be doing in the hour following adding the solution (B). During this time, keep the plate at room temperature away from direct sunlight		

SYMBOLS USED

	Biological risk		Temperature limitation		Catalog Number
	Read instructions for use		In Vitro Diagnostic Use		Lot Number
	Number of tests		Use by		EC Declaration of Conformity
	Microplate		Negative control		Microspheres
	Specimen diluent		Positive control		IgG Conjugate
	Calibrator		Wash Buffer		Microsphere reconstitution buffer
	FIDIS test		Manufacturer		

4-6 bld de Beaubourg,
CS 90138 CROISSY BEAUBOURG
77435 MARNE LA VALLEE CX2
France

Tel : +33 (0)1 64 62 10 12
Fax : +33 (0)1 64 62 09 66

E-mail : support@theradiag.com
Internet : www.theradiag.com