FLEX
Monoclonal Rabbit
Anti-Human
Postmeiotic Segregation Increased 2
 Clone EP51
Ready-to-Use
 (Link)

English
Code IR087

Intended use	For in vitro diagnostic use. FLEX Monoclonal Rabbit Anti-Human Postmeiotic Segregation Increased 2, Clone EP51, Ready-to-Use (Link) is intended for use in immunohistochemistry together with Autostainer Link instruments. Antibodies to postmeiotic segregation increased 2 (PMS2) may be useful for the differential identification of colorectal carcinomas. When deficient, PMS2 is associated with the onset of hereditary non-polyposis colorectal cancer (HNPCC) (1, 2). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Synonyms for antigen	hPMS2, mismatch repair endonuclease PMS2, PMS1 protein homolog 2
Summary and explanation	<i>PMS2</i> is a ubiquitously expressed gene encoding the protein PMS2, an ~100 kDa nuclear protein involved in the DNA mismatch repair (MMR) pathway. The MMR pathway is utilized by normal proliferating cells to repair mutations that may occur during DNA replication. PMS2 forms a heterodimer with the mismatch repair protein MutL protein homolog 1 (MLH1), and the MLH1-PMS2 heterodimeric complex is recruited to the mismatch DNA sequence by a MutS MMR heterodimer consisting of MSH2 and MSH6, which binds directly to the base mismatch. The MLH1-PMS2 complex then initiates downstream functions including the excision of the mismatched DNA strands and repair through the recruitment of nucleases, polymerases, and other proteins (1). PMS2 deficiency is often the result of germline mutations in MMR deficient individuals (1). PMS2 loss of function results in an increased mutation rate and is associated with the development of HNPCC (also called Lynch syndrome). HNPCC, which accounts for 1-5% of colorectal cancers, is an autosomal dominant condition that results from germline mutations in any of at least four MMR genes including <i>MSH2</i> (~40%), <i>MLH1</i> (~50%), <i>MSH6</i> (~5%) and <i>PMS2</i> (<5%) and is associated with increased susceptibility to colon and endometrial cancer and to a lesser degree other types of neoplasms (1-3). A hallmark feature of cancers with defects in the MMR pathway is microsatellite instability (MSI), alterations in the number of microsatellite repeats in the DNA as a result of the faulty repair process (4, 5). A germline mutation in a single MMR gene copy predisposes patients to HNPCC, however, cancer development requires loss of function of the remaining wildtype copy. Mutations in <i>PMS2</i> account for approximately 5% of all cases of HNPCC (6), however <i>PMS2</i> is dependent upon <i>MLH1</i> for stable expression, thus at the protein level, both sporadic and hereditary colorectal cancers deficient in <i>MLH1</i> generally exhibit concomitant loss of <i>PMS2</i> (7). Antibodies to PMS2 are useful for identifying mismatch repair deficiencies in tumors of the gastrointestinal tract including HNPCC and associated extracolonic cancers by immunohistochemistry (IHC). PMS2 deficiency as determined by IHC has been reported in endometrial, ovarian, breast, esophageal, and ampullary carcinomas, and melanomas (8-13). Refer to Dako's <i>General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures.
Reagent provided	Ready-to-use monoclonal rabbit antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L sodium azide. <u>Clone:</u> EP51. <u>Isotype:</u> Rabbit IgG.
Immunogen	Synthetic peptide corresponding to residues in human PMS2 protein.
Specificity	In Western blotting of SHSY5Y cells, anti-PMS2, clone EP51, recognizes a major band at 110 kDa corresponding to the expected molecular weight of PMS2.
Precautions	- For professional users. - This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. - As with any product derived from biological sources, proper handling procedures should be used. - Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.

5. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Quick guide

Step		Comments
Fixation	Formalin	
Pre-treatment	EnVision™ FLEX, High pH (Code K8004)	20 min HIER, 3-in-1 using PT Link and PT Link Rinse Station
Dilution	Ready-to-use	30 min incubation
Dilution buffer	Pre-diluted	
Negative Control	FLEX Negative Control, Rabbit (Code IR600)	30 min incubation
Visualization	EnVision™ FLEX, High pH (Code K8000)	30 min incubation, 2x5 min DAB+
Counterstain	EnVision™ FLEX Hematoxylin (Code K8008)	5 min incubation
Control Tissue	Appendix, colon	Nuclear staining
Slides	FLEX IHC Microscope Slides (Code K8020)	Recommended for greater adherence of tissue sections to glass slides
Mounting	Non-aqueous, permanent mounting required	After staining, the sections must be dehydrated, cleared and mounted using permanent mounting medium.
Instrumentation	Autostainer Link 48	Use instrument-specific vials (Code SK201-SK203)

**The user must always read the package insert for detailed instructions of the staining procedure and handling of the product.*

Specimen preparation

Paraffin sections: The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm.

Pre-treatment: Pre-treatment of formalin-fixed, paraffin-embedded tissue sections with heat-induced epitope retrieval (HIER) is required. Pretreating tissues with HIER using diluted EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8004) is recommended. Deparaffinization, rehydration and epitope retrieval can be performed in Dako PT Link (Code PT100/PT101). For details, please refer to PT Link User Guide. The following parameters should be used for PT Link: Pre-heat temperature: 65 °C; epitope retrieval temperature and time: 97 °C for 20 (±1) minutes; cool down to 65 °C. Remove slide rack from PT tank and immediately dip slides in jar/tank (e.g., PT Link Rinse Station (Code PT109)) containing diluted room temperature EnVision™ FLEX Wash Buffer (20x) (Code K8007). Leave slides in Wash Buffer for 1-5 minutes.

The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended. After staining the sections must be dehydrated, cleared and mounted using permanent mounting medium.

Staining procedure

Visualization: The recommended visualization system is EnVision™ FLEX, High pH (Link) (Code K8000).

Program: The staining steps and incubation times are pre-programmed into the Autostainer Link software. The recommended reagent application volume is 1 x 200 µL or 2 x 150 µL per slide. Please refer to the proper Autostainer Link User Guide for detailed instructions on loading slides and reagents. If the protocols are not available on the Autostainer Link platform, please contact Dako Technical Services. An Autostainer Link software installer can be found at www.dako.com/Installer. This installer will update DakoLink software with protocol and reagent information for Anti-PMS2, Clone EP51. The short name for the antibody protocol is PMS2 EP51. All incubation steps should be performed at room temperature.

Counterstaining: The recommended counterstain is EnVision™ FLEX Hematoxylin (Link) (Code K8008). Non-aqueous, permanent mounting medium is recommended.

Controls: Positive and negative controls should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include appendix or colon and the cells/structures should display reaction patterns as described for this tissue in the "Performance characteristics" section. The recommended negative control reagent is FLEX Negative Control, Rabbit (Link) (Code IR600).

Staining interpretation

The cellular staining pattern is predominantly nuclear.

Performance characteristics

Normal tissues: Monoclonal Rabbit Anti-Human Postmeiotic Segregation Increased 2, Clone EP51 positively labels the nuclei of a variety of cell types in normal tissue. Although not noted in the table below, positive nuclear staining was observed in endothelial cells in many tissue types. (14)

Tissue Type (# tested)	Positive Tissue Elements	Tissue Type (# tested)	Positive Tissue Elements
Adrenal (3)	(3/3) Cortex (10-50%), nuclear	Parathyroid (3)	(3/3) Epithelium (20-80%), nuclear
	(2/3) Medulla (10-50%), nuclear	Pituitary (3)	(3/3) Endocrine cells (50-90%), nuclear

Appendix (1)	(1/1) Epithelium (80%), nuclear	Prostate (3)	(3/3) Glandular and ductal epithelium (40-90%), nuclear
	(1/1) Extra-follicular lymphocytes (100%), nuclear		(3/3) Stromal cells (50-70%), nuclear
	(1/1) Smooth muscle (50%), nuclear	Salivary gland (3)	(3/3) Ductal epithelium (60-80%), nuclear
Bone marrow (2)	(2/2) Blast cells (10-20%), nuclear		(3/3) Glandular epithelium (20-80%), nuclear
Breast (3)	(1/1) Glandular epithelium (80%), nuclear		(1/1) Smooth muscle (30%), nuclear
	(2/2) Ductal epithelium (80%), nuclear		(1/1) Lymphoid cells (50%), nuclear
	(1/1) Connective Tissue (30%), nuclear		(2/2) Connective tissue cells (20%), nuclear
Cerebellum (3)	(3/3) Purkinje cells (80-90%), nuclear	Skeletal muscle (3)	(0/3) Muscle cells
	(3/3) Glial cells (50%), nuclear		(3/3) Connective tissue (10-20%), nuclear
Cerebrum (3)	(3/3) Glial cells (20-60%), nuclear	Skin (3)	(2/2) Basal epithelium (50-80%), nuclear
Cervix (3)	(3/3) Basal epithelium (80-100%), nuclear		(3/3) Sweat gland (30-90%), nuclear
	(3/3) Stromal cells (50%), nuclear	Small intestine (3)	(3/3) Epithelium (60-80%), nuclear
Colon (3)	(3/3) Epithelium (50-80%), nuclear		(3/3) Smooth muscle (30-40%), nuclear
	(3/3) Lymphocytes (20-40%), nuclear	Spleen (3)	(3/3) Resting lymphocytes (20-70%), nuclear
	(3/3) Connective tissue (10-40%), nuclear		(1/1) Reactive lymphocytes (80%), nuclear
Esophagus (3)	(3/3) Squamous epithelium (10-50%), nuclear		(3/3) Sinus lining (red pulp) (30-50%), nuclear
	(1/1) Connective tissue (30%), nuclear	Stomach (3)	(1/1) Connective tissue (30%), nuclear
	(0/1) Lymphocytes		(1/2) Epithelium (10%), nuclear
Heart (3)	(3/3) Cardiac cells (50-60%), nuclear		(1/1) Gastric gland (10%), nuclear
	(3/3) Connective tissue (30-40%), nuclear		(2/2) Smooth muscle (20-60%), nuclear
Kidney (3)	(3/3) Tubule (50-60%), nuclear		(1/1) Lymphocytes (60%), nuclear
	(3/3) Bowman's capsule (10-20%), nuclear	Testis (3)	(3/3) Leydig cells (60-80%), nuclear
	(3/3) Glomerulus (10%), nuclear		(3/3) Tubular cells (10-20%), nuclear
Liver (3)	(0/3) Hepatocytes	Thymus (3)	(3/3) Cortex (10-30%), nuclear
	(1/3) Bile duct (10%), nuclear		(3/3) Medulla (10%), nuclear
Lung (3)	(2/2) Lymphocytes (80%), nuclear	Thyroid (3)	(2/3) Follicular epithelium (50%), nuclear
	(3/3) Alveolar epithelium (50-60%), nuclear		(3/3) Parafoallicular C cells (10-30%), nuclear
Mesothelial cells (3)	(2/2) Mesothelium (20-50%), nuclear	Tonsil (3)	(3/3) Follicular lymphocytes (40-60%), nuclear
	(3/3) Connective tissue (50%), nuclear		(3/3) Extra-follicular lymphocytes (40-50%), nuclear
Nerve, peripheral (3)	(3/3) Schwann cells (10-30%), nuclear		(3/3) Epithelium (60-80%), nuclear
	(3/3) Connective tissue (30-50%), nuclear	Uterus (3)	(3/3) Glandular epithelium (30-50%), nuclear
Ovary (3)	(3/3) Stromal cells (20-40%), nuclear		(3/3) Stromal cells (50-60%), nuclear
	(3/3) Follicular epithelium (10-50%), nuclear		(3/3) Myometrium (50%), nuclear
Pancreas (3)	(3/3) Acinar cells (50-60%), nuclear		
	(3/3) Islet cells (10-30%), nuclear		
	(3/3) Pancreatic ductal cells (10-40%), nuclear		

Abnormal tissues: In PMS2 deficient tumors, no nuclear staining should be observed in the neoplastic cells. A population of benign cells should demonstrate nuclear staining and serve as an internal tissue control.

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Literaturangaben

1. Peltomäki P. Role of DNA mismatch repair defects in the pathogenesis of human cancer. *J Clin Oncol* 2003;21:1174-79.
2. Lynch H, Smyrk T. Hereditary nonpolyposis colorectal cancer (Lynch syndrome). *An updated review.* *Cancer* 1996;78:1149-67.
3. Marcus VA, Madlensky L, Gryfe R, Kim H, So K, Millar A, et al. Immunohistochemistry for hMLH1 and hMSH2: A practical test for DNA mismatch repair-deficient tumors. *Am J Surg Pathol* 1999;23:1248-55.
4. Lanza G, Gafà R, Maestri I, Santini A, Matteuzzi M, Cavazzini L. Immunohistochemical pattern of MLH1/MSH2 expression is related to clinical and pathological features in colorectal adenocarcinomas with microsatellite instability. *Mod Pathol* 2002;15:741-9.
5. Okuda T, Sekizawa A, Purwosunu Y, Nagatsuka M, Morioka M, Hayashi M, Okai T. Genetics of endometrial cancers. *Obstet Gynecol Int* 2010;984013.
6. Burgart L.J. Testing for defective DNA mismatch repair in colorectal carcinoma. A practical guide. *Arch Pathol Lab Med* 2005;129:1385-9.
7. Lynch HT, Lynch PM, Lanspa SJ, Snyder CL, Lynch JF, Boland CR. Review of the Lynch syndrome: history, molecular genetics, screening, differential diagnosis, and medicolegal ramifications. *Clin Genet* 2009; 76:1-18.
8. Garg K, Shih K, Barakat R, Zhou Q, Iasonos A, Soslow RA. Endometrial carcinomas in women aged 40 years and younger: tumors associated with loss of DNA mismatch repair proteins comprise a distinct clinicopathologic subset. *Am J Surg Pathol* 2009;33:1869-77.
9. Pal T, Permuth-Wey J, Sellers TA. A review of the clinical relevance of mismatch-repair deficiency in ovarian cancer. *Cancer* 2008;113:733-42.
10. Buerki N, Gautier L, Kovac M, Marra G, Buser M, Heinimann K. Evidence for breast cancer as an integral part of Lynch syndrome. *Gene Chrom Cancer* 2012;51:83-91.
11. Agaram NP, Shia J, Tang LH, Klimstra DS. DNA mismatch repair deficiency in ampullary carcinoma: a morphologic and immunohistochemical study of 54 cases. *Am J Clin Pathol* 2010;133:772-80.
12. Farris AB 3rd, Demicco EG, Le LP, Finberg KE, Miller J, Mandal R, et al. Clinicopathologic and molecular profiles of microsatellite unstable Barrett Esophagus-associated adenocarcinoma. *Am J Surg Pathol* 2011;35:647-55.
13. Alvino E, Marra G, Pagani E, Falcinelli S, Pepponi R, Perrera C, et al. High-frequency microsatellite instability is associated with defective DNA mismatch repair in human melanoma. *J Invest Dermatol* 2002;118:79-86.
14. PMS2 clone EP51 M3647.IR087 30 Normal Report D08987

Monoclonal Rabbit Anti-Human Postmeiotic Segregation Increased 2, Clone EP51 has been created by Epitomics Inc., using Epitomics' proprietary rabbit monoclonal antibody technology covered under Patent Nos 5,675,063 and 7,402,409.

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Dako North America, Inc.
Dako Denmark A/S