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ENGLISH

INDEX

	Page
— GENERAL RUNNING DATA OF TURBO-SMART 50/60 HZ	16
— INTRODUCTION	17
— SIGNALS AND WARNINGS	17
— GENERAL FEATURES	17
— TECHNICAL ROOM	18
— THE TECHNICAL ROOM MUST BE FITTED WITH	18
— PIPING LAYOUT	19
— INSTALLATION AND STARTING	20
— INSTRUCTIONS FOR THE USE OF THE "AC80" CIRCUIT PROGRAMMING KEYBOARD	21
— TESTING	22
— WORKING	22
— GENERAL WARNINGS AND BIOLOGICAL DANGER	23
— ROUTINE MAINTENANCE OF THE ASPIRATION SYSTEM FITTED WITH TURBO-SMART MEANT TO SPECIALLY TRAINED AND INSTRUCTED TO SURGERY STAFF	24
— EXTRAORDINARY MAINTENANCE OF THE ASPIRATION SYSTEM FITTED WITH TURBO-SMART	25
— IMPORTANT NOTICES	26
— TRANSPORT AND STORAGE	26
— TRANSPORT OF SECOND-HAND APPLIANCES	26
— ROUTINE MAINTENANCE CARD	27
— EXTRAORDINARY MAINTENANCE CARD	28
— PICTURES	71

**General running data of Turbo-Smart 50/60 Hz
Dental aspirator**

<i>Model</i>	<i>Turbo-Smart "A"</i>	<i>Turbo-Smart "B"</i>
<i>Rated Voltage</i>	<i>230V</i>	<i>230V</i>
<i>Rated frequency</i>	<i>50/60Hz</i>	<i>50/60Hz</i>
<i>Rated current</i>	<i>7A</i>	<i>9A</i>
<i>Isolation class</i>	<i>I</i>	<i>I</i>
<i>Type of appliance</i>	<i>B</i>	<i>B</i>
<i>Use</i>	<i>S1</i>	<i>S1</i>
<i>Protection against humidity</i>	<i>Common</i>	<i>Common</i>
<i>Level of protection against director indirect contacts type</i>		
<i>Max absorbed power</i>	<i>1,450 kW</i>	<i>1,870 kW</i>
<i>Max capacity flow</i>	<i>86 m³/h</i>	<i>105 m³/h</i>
<i>Max head for continuous service</i>	<i>2000 mmH2O</i>	<i>2000 mmH2O</i>
<i>Sound pressureversion without box **</i>	<i>/</i>	<i>73 db (A) (•)</i>
<i>Sound pressureversion with box</i>	<i>/</i>	<i>/</i>

<i>.....,</i>	<i>Alternating current</i>	<i>IEC 417-5032</i>
	<i>Earthing</i>	<i>IEC 417-5019</i>
<i>O</i>	<i>Open (disconnected from main electrical supply)</i>	<i>IEC 417-5008</i>
<i>I</i>	<i>Closed (connected to the main electrical supply)</i>	<i>IE C 417-5007</i>

- *Sound pressure level tested according to the ISO 3746-1979 (E) standard.*
 - *Parameters: r or d=1,5 – background noise ; 51 db (A) – instrument Bruel & Kjaer type 2232*
- ** Sound pressure level registered at maximum speed.*

Introduction

The following presentation is meant to illustrate the equipment and the appliance in order to clarify their working and maintenance and to inform about the dangers and precautions which are necessary for accident prevention.

Signals and warnings

- **Electrical shock risk: also 230 V ~ can be lethal.** 
- **Biological danger, risk of infections from epidemic diseases.** 
- **General danger sign.** 
- **High temperature.** 
- **Compulsory direction of flow or of rotation.** 

Signs cannot always fully express danger warnings, therefore it is necessary that the user reads the warnings and keeps them in due consideration. Failure to respect a danger sign or warning may harm operators or patients. Do not remove protections; do not tamper with machines or their operation.

Despite all our efforts, it is still possible that danger warnings are not exhaustive: we apologise with the users and kindly request them to care for all danger sources that might have pass unnoticed and to inform us accordingly.

General features (Fig. 1-2 and 3 pages 71-72-73)

Turbo-Smart is a dental aspirator basically made up of:

- an aspirating unit, (det. A1)
- a centrifugal separator (det. A2) separating air from liquid,
- an electric circuit (detail A3) controlling all the aspirator functions.

Turbo-Smart can be supplied complete with the following optional extras:

- Hydrocyclone (det. A2 and B2) and,
- plastic box, (det. A3) with handle (det. B3).

The centrifugal separator (Fig. 4 and 5 pages 74-75) is fitted with:

- thermostatic valve, (det. A4) reducing head values in case of high room temperature thus allowing for a gradual return to a proper operating temperature;
- overflow valve (det. A5) draining the liquids that might go into the centrifugal separator while the aspirator is off;

- recycling valve (det. C4) allowing for the gradual drainage of an extraordinary quantity of liquids entering the centrifugal separator or the Hydrocyclone;
- one-way valve (det. D1) is a device for installing in parallel two or more Turbo-Smart units.

Installing the valve inside the centrifugal separator unit has enabled us to build a complete appliance with small dimensions.

Turbo-Smart:

- allows for head regulation from zero to two meters in water column;
- standard version operates two dental chairs;
- by means of password (B) can be adjusted to operate three/four dental units;
- normally powered by 230 V single-phase (AC) - other voltages are available on demand.

The features, adjustments and regulation possible on Turbo-Smart make it for all dental practices: surgery, prosthetics and endodontics. Just regulate the head and you will have the correct aspiration for any procedure.

Introduction

Following notes are meant to offer some advice, which we deem useful in designing and organizing the technological room and for the aspiration and compressed air distribution systems maintenance, as well.

Apart from dimensions of surgery, there are some regulations to comply with, first of all for safety of operators and patients, but also useful for a long life of machines and plants.

Also a lobby before a bathroom, a room without windows or even a balcony can be transformed into a small technical room fitted with all safety devices required by regulations.

Technical room

- Aspirator and compressor must be installed in a room where access to customers, to extraneous people, and even to the surgery staff who has not been purposely trained is forbidden.
- If such a room is not available, machines must be protected by a special box, which must not be easy to take off, so that to avoid accidental contacts, dangerous because of risk of electrical shocks, of machine running temperature, of the possibility (unlikely, but not impossible) of fire, of bursting and of contaminating air and/or liquid leakage.
- Use boxes designed and produced by the manufacturer of the machines only.

The technical room must be fitted with:

A smoke-detector or a temperature sensor, for fire prevention, connected to a non-stop surveillance service, for example to secretary's office; install a suitable fire extinguisher outside the technical room and check regularly its efficiency; see that routine and extraordinary maintenance staff are properly protected against infected sprinkles and accidental contacts; ensure that machines are not modified, that safety protections are not tampered with, and that no repair is made on machines in operation or on ones connected to the mains.

- In the technical room temperature can range from a minimum of +5 °C to +35 °C max.
- Electrical supply must be sufficient for the loads specified on machine rate plates; three

phases + neutral + earth must be in place - grounding is an indispensable condition also for single-phase machines.

- Supply lines, control panels and machines must be protected against direct and indirect electrical contacts, against overloads and overcurrents, in compliance with C.E.I. 64-8 regulations concerning I class appliances (I.E.C. correspondents available on request).
- It is advisable to protect the feeding line against extra-current.
- Protect electrical control panels and machines from accidental sprinkles.
- Provide a floor liquid-drainage connected to the sewage system.
- Keep the technical room clear of anything that is not related to the machines herein contained, taking special attention of flammable materials; make sure that there is no risk of corrosive, flammable or explosive mixes.

Piping layout (Fig. 6 and 7, page 76)

The diameter of the main duct or collecting pipe (A) must be proportionate to the flow of the aspiration unit.

Branches (B) from the surgeries will be connected to the main duct by means of 45° connections (C) bent towards the surgeries. Where necessary, inspection holes (D) should be included in the plan to allow for inspection inside aspiration piping.

As a rule, **the location of the piping must comply with the laws of physics on the dynamics of fluids into a pipe.**

A slope of the piping towards the pump will facilitate its functioning, pipe on the same level will give no problems, but counter-slopes compromise functioning proportionately to the gradient of the slopes.

In liquid and semi-wet installations, siphons and counterslopes in the tubes cause stoppages in suction.

To avoid such problems it is necessary to insert a "supplementary aspiration draining valve" into the farthest point. When the piping is clogged by liquid the supplementary valve opens and lets air in, which drives liquids into centrifugal separator thus clearing the piping.

Once the piping is installed, it is advisable to carry out a tightness test: all outlets should be closed up, the piping should be put under slight pressure and checked with a suitable instrument; pressure must remain steady. An electrical cable for low-tension controls will run along the piping.

An A version of Turbo-Smart installed in a system fitted with a 40 mm main pipe and two 30 mm branches (Fig. 6) services two dental chairs. Aspiration values would remain optimal even in case of a 30 mm main duct, provided that the distance between the farthest dental chair and the plant room is less than 12 meters.

The version B of Turbo-Smart (with password) on a main 50 mm main pipe fits three dental chairs at a distance of 15 metres whereas on a main 50 mm cross-section pipe it can operate four dental chairs at the same distance (Fig. 7, page 76).

Our technical department is ready to be called on for the design of personalized, bigger and more complex systems; we are obviously willing to put our experience at our customers' disposal.

Installation and starting (wiring diagram Fig. 8, page 77)

- **Unpack the appliance following the instructions shown on the package.**
- **Dispose of the package in compliance with regulations.**
- **Verify that the appliance has not been damaged during transport.**
- **Do not connect damaged appliances to the mains.**
- Do not use extension leads, multiple plugs or sockets.
- Ascertain that the feeding line is adequate to feed the machine.
- **Ascertain that the appliance (or Blok-Jet) is installed in a clean spot, far from heat sources and from depots of contaminating substances, and that it aspirates clean air, free from dust, gases, and humidity.**
- Any machine installation must be carried out by a specialist, duly trained and equipped with necessary tools.
The installer must refer to the appliance manual, carry out the machine final test and instruct users on use and routine maintenance of new not used and therefore not yet contaminated machines.
- Before starting the compressor and/or the aspirator ascertain that piping is clean; heavy debris could damage the appliances.
- Danger signals and warnings (light and beepers) such as failing-machine or overheating warning must be removed from technical room and installed in room under frequent control to be signalled by a sound or visual warning.
- Air exhausted from aspirators should be filtered, through a suitable certified absolute filter and then conveyed outside.
- After finishing the installation, carry out functional and regulation tests; check the rotation sense of motors, the mains tension, and electrical absorption.
- Check the functioning of drainages, aspiration flow rate and head, compressor charging time, functioning of silica-gel regenerating device, and running of ventilators.
- Make sure that periodical checks are carried out on equipments: such an inspection is not only a means to protect the surgery from forced stops and to check its operation conditions, but it is also a proper way to prevent accidents.
- Our **updated** manual are in our web site **www.cattani.it**. We recommend to consult them especially for **safety** up-dating.

Install and connect the aspirator as per wiring diagram, (Fig. 8) select the ON position on the switch, which will light up. At this point aspiration will start.

The control panel is fitted with a keypad which allows access to the program, to regulate or to modify the head value and with a display showing the current working parameters of the system.

It is allowed to enter the program to interrogate the operation log, especially those concerning alarms. In case of alarms it is also possible to place a remote flashing signal, in the reception office or in another room.

Instructions for the use of the “AC80” circuit programming keyboard

1.Introduction

- The keyboard allows to regulate and control the “AC 80” inverter system.
- This document is an introduction to the basics for the regulation and the interpretation of the display menu.

2.Use of keyboard and display

- The display shows all the information about the aspirator working.
- When it is switched-on the system shows the opening menu. To scroll the menu use the keys with the up and down indicators.


Fig. 10
Draw. 10

ENGLISH

3.Menu consultation

Use the 6 available consultation keys to enter the settings.

Function	Consultation mode	Values change mode
	Menu scrolling upwards	Head value increase
	Menu scrolling downwards	Head value decrease
	Access to the head control under – menu	
	Return to the main menu	
	Head value change start mode	Confirmation of head change and return to consultation
	Alarms cancellation	Change cancellation and return to consultation

- For basic instructions **scroll** the menu with the keys  and 
- Once selected the **menu**, enter **I** with the key 
- To regulate the parameters enter the **“editing”** mode pushing the key 
- The cursor appearance shows that the control has been accepted.
- To change the values select the number with  and 
- The value can be **modified** through  and 
- Push again  and leave the **“editing”** function **to save** the change or push  in order **to cancel** the change.

Testing (Fig. 10, page 21)

Checking the “AC80” program and using the keyboard (Fig. 10) it is possible to check: voltage, working frequency, power consumption and head. Such values must correspond to the plate data of the machine and to the general working data quoted at page 2 of this manual; in case of discrepancies contact the manufacturer or the local engineer. The control program supplies other information which is mentioned in the chapter about maintenance.

Make sure that the piping, the coupling to dentist's units and the dentist's units themselves have no aspiration losses.

We recommend carrying out a working test using aspirated water to check the piping tightness and a visual check of the cabling, the absorption and loudness in presence of liquids. The weephole must be completely clear and possibly raked; to avoid bad smells we suggest placing a small trap on the waste pipe of the room provided that the trap does not obstruct the draining.

Working (Fig. 1-2-3-4 and 5 from page 71 to 75)

The aspirated fluid must be filtered through a filter with a particle size not larger than 0,8 mm placed on the dental chair, since together with the liquids solid particles are aspirated as well. Passing through the piping, the fluid enters the centrifugal separator (detail B1) through a second filter (detail H4); the liquid arrives at the bottom of the centrifugal separator where the impeller (detail C5) forces the liquid inside the Hydrocyclon while the air is aspirated by the aspirating unit. The air coming from the aspirator must be at first silenced (the sound absorbing filter is standard), then filtered by means of a bacteriological filter (available on demand) and pumped out.

The small particles of amalgam not removed by the filters are collected by the container (detail B2) at the base of the Hydrocyclon (see the ISO manual of the Hydrocyclon and its use) The resulting amalgam-free liquid from the Hydrocyclon is discharged to waste, the liquid exceeding the Hydrocyclon capacity is conveyed back inside the centrifugal separator through the inlet E4. The outlet F4 can be used to drain the liquids which are inside the Hydrocyclon and the centrifugal separator before a possible repair or extended maintenance.

The aspirator reacts quickly to the arrival of a high quantity of liquid as follows: the aspirating unit (detail A1) slows down the rotation speed thus limiting the incoming liquid,

whereas the centrifugal separator drains at maximum capacity, the C4 valve opens, thus allowing the recirculating of the excess liquid inside the centrifugal separator while the total electric absorption remains constant.

The auto-protection function of Turbo-Smart starts working in presence of high temperature as well. The printed circuit board is fitted with three thermocouples: one is connected to the motor group of the aspirating unit, the second one is recommended to check the room temperature of the room or (the plant room) where the machine is placed, the third one can detect the temperature of the exhausted air.

In case of temperature increase, which is dangerous for the integrity of the machine, the aspirating unit slows down its rotation speed to reduce the head and to allow a gradual temperature decrease.

When the temperature becomes dangerously high, the thermostatic valve opens (detail A4) which has the same safety function as the above described electric control.

The A4 valve opens also when the vacuum exceeds 250 mbar in water column fixed by the ISO standards. In both cases similar functions come into operation to improve safety.

General warnings and biological danger (Fig. 12, page 80)

- **Electrical shock risk: also 230 V can be lethal.**

- **Biological danger, risk of infections from epidemic diseases.**

- **General danger sign.**

- **High temperature.**

Before starting any servicing operation on appliances that have been used, clean with Puli-Jet plus or with other suitable substance. **Disconnect the appliance from the mains** and padlock main switch, if provided.

For any servicing operation wear disposable gloves, goggles, mask and overall. Aspirated liquids are contaminated and infected, and they can spread infectious diseases.

The centrifugal separator (Fig. 12 detail D12), puts the aspirated liquids under pressure, so, **in case of leakages from the pipe or from appliances under pressure, the biological hazard is higher.**

Routine maintenance of the aspiration system fitted with Turbo-Smart meant to specially trained and instructed to surgery staff

- **Before any maintenance operation, wear disposable gloves, mask, overalls, and goggles.**

The aspirated liquid, which is infected and contaminating, can be under pressure inside the machines and piping; in case of break-down of an under-pressure component, biological hazard can be more serious.

- **After surgical operations, or with high risk of infection,** aspirate a solution of Puli-Jet plus by means of Pulse Cleaner.

In the middle of the day (before lunch time) aspirate a solution of Puli-Jet plus by means of Pulse Cleaner on each dental chair.

Directions for use are printed on the bottle-label of Puli-Jet plus; in case of installations servicing two or more chairs, prepare Pulse Cleaner with the solution for each dental unit, place aspiration hose on the bucket terminals. Aspiration do not require any assistance by the operator. Pour the solution left in the bucket into the spittoon bowl.

- **We recommend the use of disinfecting Puli-Jet plus** after any operation on infection risk patients and after the last treatment in the evening; do not rinse, Puli-Jet plus will take effect during the night.
- **Also filters must be replaced daily.** To avoid spilling contaminating liquids, before taking off the filters from the dental unit, start the aspiration without aspirating any liquid, take off the filters used during the day, dip them in a **disinfecting**-detergent and put in place clean filters with **disinfecting** antifoam tablets.
- **Every day:** check the display showing the work setting of the aspirator. Contact the engineers in case of alarms such as : temperature exceeding the allowed limit, High Current Warning, failure of the condensers capacity, saturation of the amalgam container. Also all other filters of spittoon draining units and on aspiration units must be replaced in the same way. Isolate from the mains before maintenance on the aspirators (even if only to replace the filters).
Any amalgam trapped in the filters must be recovered in compliance with ISO 11143 regulation and conveyed to institutes authorised for disposal, being a toxic, harmful and contaminated waste.
- Check that machines installed in gardens, in spaces under staircases or in other facilities outside the surgery, are not choked by leaves, dust, rags or similar objects likely to hindrance the machine ventilation.
- **Periodically check the loudness, inform the engineer in case of alterations.**
- **Each operator will increase or decrease** the frequency of described operations according to his experience and work intensity.
- **Fill in the routine maintenance form/log book.**

Non routine-maintained machines block like the high-temperature-running ones as well; in the first case wet dirt particles are carried by aspirated air into the most sensitive point and create frictions; in the second case overheating and the following night cooling make the bearing grease to dry, create dimensional unbalancing to mechanical components, and destroy insulating materials and electrical components.

- **Every five years replace the bearings placed on the aspirating unit and on the centrifugal separator.**

For any anomaly that the user cannot understand address the manufacturer, do not intervene with spare parts which are not original.

Maintenance which is carried out carefully according to the recommended methods and with the use of suitable products allows the regular working of the machines and can prevent premature failure.

Fill in the guarantee slip and send it to the manufacturer.

- * *The recommended frequency for the various interventions is only indicative and it can be reduced for surgeries with little activity or increased for professional surgeries with intense activity.*

Extraordinary maintenance of the aspiration system fitted with Turbo-Smart (Fig. 1-2-7-11 pages 71-72-76-79)

- **Electrical shock risk: also 230 V can be lethal.**

- **Biological danger, risk of infections from epidemic diseases.**

- **High temperature.**

- The extraordinary maintenance must be entrusted to an engineer in possession of original spare parts, the work must be carried out periodically with a frequency based on the work load of the surgery.
- The recommended engineer must train the surgery staff (according to the instructions of this manual) and make sure that the ordinary maintenance is carried out carefully and methodically using Magnolia products (the manufacturer will not consider under warranty machines which have been treated with other products).
- The first intervention of extraordinary maintenance must be carried out after three months from the installation of the machine; afterwards an extraordinary maintenance every six months can be enough provided that ordinary maintenance is carried out carefully.
- On the occasion of every maintenance intervention it is necessary to wear: gloves, glasses, half mask respirator and disposable apron. The aspirated liquid is infected and possibly carrying disease pathogens. Please note that the liquid inside hoses A12 and B12 is under pressure: this particularly increases the biological risk.
- Any intervention on the machine must be preceded by a series of washes with Puli-Jet and hot water using Pulse Cleaner.
- Isolate from the mains supply before maintenance.
- First control: check the tightening of the terminals inside the electric control panel.
- First control: according to the instructions of the "AC80" program (page 21), visualize on the display all information about working, causes of working anomalies such as: maximum registered temperatures, maximum registered amperage, irregular supply voltage of the condensers, high amalgam level, will have to be eliminated.
- On occasion of any intervention: check that the filters on the dentist's units and the H4 filter on Turbo-Smart are clean.
- Every six months: drain the liquids which are inside the Hydrocyclon as follows: loosen screw E12, turn on valve D11, take out the upper part A11, check that the inner chamber B11 is clean and check that the 5 holes are clear.
- Every six months: drain the liquids which are inside the centrifugal separator C5 as follows: disconnect the hose from the hole holder F4, collect with a small bowl the residual liquids (these liquids, which have not been treated by the amalgam separator, probably still contain amalgam particles and must be aspirated when the machine is reassembled), remove the base B4, check the conditions of the fan C5 of the centrifugal separator, replace the parts which are frayed or damaged. If these parts are very dirty or crusted it is recommended to check the cleaning regime, the times and the products used for ordinary maintenance.
- Once a year check the shaft of the recircle valve C4 and replace it in case of wear. We recommend to replace also the rubber part of the same valve and to replace once a year the draining valves D4 and A5 together with all the sealing "O" rings.
- We recommend replacing once a year the whole piping, especially hoses under pressure A12 and B12.
- Once you have checked the machine and reassembled the parts, test the working of the machine.
- Every five years replace the bearings placed on the aspirating unit and on the centrifugal separator.

ENGLISH

*For any anomaly that the user cannot understand address the manufacturer, do not intervene with spare parts which are not original.
Maintenance which is carried out carefully according to the recommended methods and with the use of suitable products allows the regular working of the machines and can prevent premature failure.
Fill in the guarantee slip and send it to the manufacturer.*

Important notices

- *Appliances are guaranteed for a year from the date of sale, provided that the guarantee card addressed to the factory is returned with date of sale and retailer's and customer's name.*
- *Guarantee and manufacturer's liability cease when the appliances: are treated with products unsuitable or different from those recommended by the manufacturer; are improperly used; are tampered with operations of any kind done by persons unauthorised by the manufacturer.*
- *The manufacturer, concessionaires, agents, and authorised engineers are at customers' disposal for advice and assistance and to supply literature, spare parts and anything useful.*
- *The manufacturers reserves the right to modify the products for improvements, for technical, normative and functional reasons or for problems due to the availability of semi finished products, without prior notice.*
- *On the web-site: www.cattani.it you can find our **updated** manuals.
We suggest to consult them especially as for **safety** updating.*

Transport and storage

- *Packed equipment can be transported and stored at a temperature ranging from -10 °C to + 60 °C.*
- *Packages must be kept away from water and splashing and cannot tolerate humidity >70%.*
- *Packages with the same weight can be stored in piles of three only.*

Transport of second-hand appliances

- *Before packing, cleanse and sanitize with Puli-Jet plus (see paragraphs "General warnings and biological danger" and "Maintenance to be carried out by the surgery staff").*
- *Close with polyethylene plugs all machine inlets and outlets.*
- *Place the unit into a polyethylene bag, seal and pack it in 3-layer corrugated board.*

Signals and warnings

- **Electrical shock risk: also 230 V ~ can be lethal.** 
- **Biological danger, risk of infections from epidemic diseases.** 
- **General danger sign.** 
- **High temperature.** 
- **Compulsory direction of flow or of rotation.** 

**ASPIRATION SYSTEM EQUIPPED WITH TURBO-SMART
ROUTINE MAINTENANCE CARD**

Ordinary maintenance must be entrusted to specially instructed and trained staff.

Photocopy this file and once filled in, store it.

Maintenance date User signature		
Every day: check the possible alarm signal on the display*	YES	NO
Every day: evening aspiration of the solution made up of Puli-Jet plus Disinfectant, hot water and Pulse Cleaner	YES	NO
Occasionally: check the aspirator noisiness*	YES	NO
Disconnect the machine before any maintenance	YES	NO
Every day: replace the dirty filters with clean and disinfected filters, place correctly the antifoaming disinfectant tablets on the filters of the dentist's unit	YES	NO
At every maintenance: make sure that the aspirator ventilation is not obstructed (by dust, paper or other things)	YES	NO
Keep the plant room free from anything not concerning the machines, especially from flammable material.	YES	NO
Make sure that there is no possibility for the formation of corrosive, flammable and explosive mixtures.	YES	NO

* Contact the engineer in case of alarms or increased noisiness.

*This maintenance sheet is a useful reference, but must be used in conjunction with the full user manual
On the website www.cattani.it it is possible to consult our **updated** manuals. We recommend the consultation especially for updates about **safety**.*

ENGLISH

Signals and warnings

- **Electrical shock risk: also 230 V ~ can be lethal.** 
- **Biological danger, risk of infections from epidemic diseases.** 
- **General danger sign.** 
- **High temperature.** 
- **Compulsory direction of flow or of rotation.** 

**ASPIRATION SYSTEM EQUIPPED WITH TURBO-SMART
EXTRAORDINARY MAINTENANCE CARD**

Extraordinary maintenance must be entrusted to a trained engineer in possession of original spare parts.

The first intervention must be carried out after three months from the installation, the other interventions should take place every six months.

Photocopy this file and once filled in, store it.

Maintenance date	User signature	
At every intervention: check on the display the maximum temperatures registered and the alarms*	YES	NO
At every intervention: check the aspirator noisiness*	YES	NO
At every intervention: check the state of ordinary maintenance**	YES	NO
At every intervention: make sure that Magnolia products and methods are used**	YES	NO
Before any intervention disconnect the machine and lock the isolator if necessary	YES	NO
Every six months: replace the absolute filter (if present) on the exhausted air	YES	NO
With the first intervention: make sure that the terminals inside the electric control panel are fastened	YES	NO
Every six months: check the Hydrocyclone conditions from the inside***	YES	NO
Ogni sei mesi: controllare internamente le condizioni del separatore centrifugo***	YES	NO
Once a year check the shaft of the recirculating valve and replace the rubber parts	YES	NO
Once a year replace the rubber parts of the draining valves	YES	NO
Every 5 years or 10.000 hours replace the bearings	YES	NO

* In case of alarms : eliminate the causes which have provoked it.

** In case the maintenance is not carried out or in case it is not carried out properly retrain the staff if necessary and inform the person in charge.

***At every intervention replace the sealing O rings.

*This maintenance file is useful to keep in mind some practical operations but it does not let the users off reading the manual.
On the website www.cattani.it it is possible to consult our **updated** manuals. We recommend the consultation especially for the updating about **safety**.*

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