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## Important information for audiologists: How to adjust the priority circuit

The priority circuit serves to enhance the system's signal-to-noise ratio and to improve voice intelligibility. Although the EK 2013 PLL receiver comes with the priority circuit switched on and preset (preadjustment: attenuation = 15 dB for an input impedance of the audio socket of 2 kΩ), the "System 2013 PLL" should be specially matched to the customer's hearing aid.

Make the "System 2013 PLL" ready for operation. Switch it on (position "M"). Connect the hearing aid to the audio socket of the EK 2013 PLL. Use one of the cables listed in Sennheiser's "Instructions for use", paragraph "Accessories". Make sure that the audio signal is not attenuated, e.g. through decoupling resistors in the audio shoe.

Set controls P1 and P2 on the EK 2013 PLL receiver to "0", which serves to mute signals from the transmitter and the in-built microphone. This is the only way to obtain exact measurements!

Place hearing aid in test box and apply sound pressure level of about 70 dB(N). Neither AGC (I) nor AGC (O) should be activated! Adjust arbitrary amplification and write down output level of hearing aid.

Switch on both transmitter and receiver. Speak into the SK 2013 PLL transmitter to cause the green LED on the EK 2013 PLL to light permanently.

Set automatic volume phasing facility on EK 2013 PLL to value that reduces the hearing aid's sound pressure level by about 15 dB when speaking into the SK 2013 PLL. The said adjustment is effected with the help of control E.

Activate the priority circuit. Use controls P1 and P2 to set the output level of the hearing aid to the original value.

The priority circuit is now perfectly matched to the hearing aid.

Subject to alterations and corrections.

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