9. EFECTOS COLATERALES

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ENGLISH

1. BLEED STP COMPOSITION

- Natural vegetable polysaccharide.

2. FUNCTION AND INDICATION

Bleed STP has as function the blood natural hemostasy through hydrophilic particles of natural vegetable polysaccharide. The product is indicated for general surgical procedures.

There are not specific studies regarding the use of the product in the ophthalmological area.

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OBSERVATION

The product is absorbed by the organism in a period between 24 - 48 hours.

3. ACTION MECHANISM

Bleed STP, in contact with the blood, absorbs the most fluid part of this (plasm), forming a gel. In the extremities of this gel, platelets, red globules and the other blood constituent are concentrated, speeding up the natural process of hemostasy

4. PRESENTATION FORMS

BLEED STP is commercialized in the following Kits:

- **Kit 1:**
 - 1 applicator (wind type) with 1 or 3 or 5 grams of Bleed STP;
 - 1 User Manual.
- **Kit 2:**
 - 3 applicators (wind type) with 1 or 3 or 5 grams of Bleed STP each;
 - 1 User Manual.

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6. STERILE PRODUCT

The product is supplied sterile, therefore the packing must be sealed up before the use. Bleed STP is sterilized with gamma ray.

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
7. USE INSTRUCTIONS

- 1) Open the external packing (quaternary) and verify if the tertiary and secondary packings (surgical degree paper) are unbroken;
- 2) Open the tertiary packing (surgical degree paper) at the moment of the surgical table preparation;
- 3) Open the secondary packing (surgical degree paper) only at the moment of the product application;
- 4) Apply the product directly on the place that occurs the bleed and pressure it lightly with a gauze. If the bleed persists, remove the product excess and repeat the application;

- 5) After 30 seconds, remove the gauze and verify if the bleed was interrupted. Use water or serum to humidify the gauze at the moment of removing it.

WARNING

Bleed STP increases its volume when in contact with liquids. In this way, the material excess must be removed after the hemostasy process.

8. SECURITY GENERAL INSTRUCTIONS

- Restricted use of the professional Dentist Surgeon /Surgeon and professionals of the rescue area;
- The application of this hemostatic must be restricted to the treatment area;
- This product must be manipulated by the operator with protection gloves;
- Once the hemostasy is gotten, the excess of Bleed STP must be removed from the application place irrigating and inhaling, mainly when used in osseous foramen, regions of osseous confines, the spinal marrow, and/or the nerve and optic chiasmus. The product increases its volume when in contact with liquids. The dry and white Bleed STP that remains after the application must be removed, therefore the removal of the dry excess material reduces the possibilities of the product to intervene in the normal physiological functions and/or to cause necrosis for compression of adjacent tissues;
- The security and the effectiveness of Bleed STP had not been tested clinically in children or in pregnant women. It was observed that the amylase activity is lower in babies born and babies until the first ten months of life; therefore, Bleed STP can take more than 48 hours to be absorbed in this population;
- Bleed STP must be used with care if there is infection or if used in contaminated areas of the body. If there is any signal of abscess or infection after the application of Bleed STP, a new intervention for draining can be necessary;
- Safety and effectiveness in ophthalmological procedures are not known.
- Bleed STP must not be used to control hemorrhage after-childbirth or menorrhagy.

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10. COUNTERINDICATION

- Do not mix this substance with vase-constrictors: epinephrine, racemic epinefrine, iron chlorides, ferric sulphate, aluminum chloride, etc;
- Do not inject, nor use Bleed STP in sanguineous vases, therefore it can occur embolization or death;
- Bleed STP is not recommended in the primary treatment of coagulation riots;
- Tests with the use of Bleed STP in osseous surfaces before the rank of prosthetic materials with adhesives had not been carried through; therefore, such use is not recommended;
- Do not apply more than 50 g of Bleed STP in diabetic patients, therefore, according to calculations, superior amounts than this can influence the glucose levels.

11. CARES

- Use Bleed STP until the validity date;
- After the rupture of the surgical degree paper packing (sterile packing), use Bleed STP immediately;
- Bleed STP was projected to be used in the solid state, and its hemostatic effectiveness diminishes when in contact with saline solutions or antibiotics before the hemostasy;
- In urologic procedures, must not leave Bleed STP in renal pelvis or ureter, in order to eliminate possible focus of calculus formation.

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ENGLISH

1. BLEED STP+ COMPOSITION

- Natural vegetable polysaccharide and carboxymethylcellulose (CMC).

2. FUNCTION AND INDICATION

The Bleed STP + has the function of natural hemostasis of the blood through hydrophilic particles of natural vegetable polysaccharide and adhesion control through the carboxymethylcellulose (CMC). The product is indicated for general surgical procedures and dressings on the outer part of the body.

There is no scientific study regarding the use of the product in the ophthalmologic area.

**WARNING**

This product is of restricted use to the surgeons dentists/doctors and professionals of the rescue area.

**OBSERVATION**

The product is absorbed by the organism in a period between 28 - 40 days.

3. ACTION MECHANISM


Bleed STP+, in contact with the blood, absorbs the most fluid part of this (plasm), forming a gel. In the extremities of this gel, platelets, red globules and the other blood constituent are concentrated, speeding up the natural process of hemostasy. The carboxymethylcellulose helps protect the affected region, functioning as a protective film to prevent adhesion formation.

4. PRESENTATION FORMS

Bleed STP+ is commercialized in the following Kits:

- **Kit 1:**
 - 1 applicator (wind type) with 1 or 3 or 5 g;
 - 1 cannula;
 - 1 User Manual.

- **Kit 2:**
 - 3 applicators (wind type - as picture aside) with 1 or 3 or 5 g each;
 - 3 cannulas;
 - 1 User Manual.

5. STORAGE AND TRANSPORT

Keep in place protected of direct light, heat and far from the children reach. Store the product in ambient temperature.

6. STERILE PRODUCT

The product is supplied sterile, therefore the packing must be sealed up before the use. Bleed STP+ is sterilized with gamma ray.

**WARNING**

The material cannot be re-sterilized. In case that it has material surplus after the use, this will have to be discarded.

7. USE INSTRUCTIONS


7.1 Open the external packing (quaternary) and verify if the tertiary (blister)

and secondary packings (surgical degree paper) are unbroken;


7.2 Open the tertiary packing (blister) at the moment of the surgical table preparation;

7.3 Open the secondary packing (surgical degree paper) only at the moment of the product application;

7.4 For nasal surgery and laparoscopy, attach the cannula to the bellows-type applicator, as pictures below:


Remove the applicator tip by unscrewing it in a **counterclockwise direction**.

Screw the cannula in the applicator tip, screwing it in **clockwise direction**.

Hemostatic Function

- Apply the product directly on the place that occurs the bleed and pressure it lightly with a gauze.
- After 30 seconds, remove the gauze and verify if the bleed was interrupted. If the bleed persists, remove the product excess and repeat the application. Use water or serum to humidify the gauze at the moment of removing it.

Control of Adhesion Formation

- At the end of surgery, apply a thin film of product directly on the site where will occur possible adhesion formation. PRODUCT SHALL NOT BE REMOVED.
- Finish the surgical procedures as needed.

8. SECURITY GENERAL INSTRUCTIONS

- Restricted use of the professional Dentist Surgeon /Surgeon and professionals of the rescue area;
- The application of this hemostatic must be restricted to the treatment area;
- This product must be manipulated by the operator with protection gloves;
- Once the hemostasy is gotten, the excess of Bleed STP+ must be removed from the application place irrigating and inhaling, mainly when used in osseous foramen, regions of osseous confines, the spinal marrow, and/or the nerve and optic chiasmus. The product increases its volume when in contact with liquids. The dry and white Bleed STP+ that remains after the application must be removed, therefore the removal of the dry excess material reduces the possibilities of the product to intervene in the normal physiological functions and/or to cause necrosis for compression of adjacent tissues;
- The security and the effectiveness of Bleed STP+ had not been tested clinically in children or in pregnant women. It was observed that the amylase activity is lower in babies born and babies until the first ten months of life; therefore, Bleed STP+ can take more than 28-40 days to be absorbed in this population;
- Bleed STP+ must be used with care if there is infection or if used in contaminated areas of the body. If there is any signal of abscess or infection after the application of Bleed STP+, a new intervention for draining can be necessary;
- The security and effectiveness in ophtalmologic procedures had not been established;
- Bleed STP+ must not be used to control hemorrhage after-childbirth or menorraghy.

9. COLLATERAL EFFECT

There aren't stories about collateral effect of this product.

10. COUNTERINDICATION

- Do not mix this substance with vase-constrictors: epinephrine, racemic epinefrine, iron chlorides, ferric sulphate, aluminum chloride, etc;
- Do not inject, nor use Bleed STP+ in sanguineous vases, therefore it can occur embolization or death;
- Bleed STP+ is not recommended in the primary treatment of coagulation riots;
- Tests with the use of Bleed STP+ in osseous surfaces before the rank of prosthetic materials with adhesives had not been carried through; therefore, such use is not recommended;
- Do not apply more than 50 g of Bleed STP+ in diabetic patients, therefore, according to calculations, superior amounts than this can influence the glucose levels.

11. CARES

- Use Bleed STP+ until the validity date;
- After the rupture of the surgical degree paper packing (barren packing), use Bleed STP+ immediately;
- Bleed STP+ was projected to be used in the solid state, and its hemostatic effectiveness diminishes when in contact with saline solutions or antibiotics before the hemostasy;
- In urologic procedures, must not leave Bleed STP+ in renal pelvis or ureter, in order to eliminate possible focus of calculus formation.

12. VALIDITY

The product possesses validity of 2 years from the manufacture date. After the opening of the surgical degree paper packing (tertiary and secondary), the same must be used immediately.

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