



RESmart<sup>®</sup> BPAP System

# **Service Manual**

## Table of Contents

1. Introduction	1
1.1 Intended Use	1
1.2 Implement	1
2. Structure	1
2.1 Outer Illumination	1
2.2 Inner Construction	2
2.3 Spare Part	2
3. Software	2
3.1 Instruction for Upgrading Device Software	2
3.2 Hidden Function	3
4. Malfunction and Countermeasure	4
4.1 Error Code Table and Information	4
4.2 Malfunction Determination	5
Appendices	6
Appendix A:	6
Device Configuration	6
1. Main device	6
2. Humidifier	8
Components Illumination	10
1. Device Inside	10
2. Device Outside	13
3. Humidifier	14
Appendix B:	15
Spare Part Configuration	15
Spare Part List	17

# IMPORTANT

All data in this manual have been verified correctly. It is sufficient for servicing and repairing this device. If applying this manual on other purpose, the validation should be confirmed by BMC Medical Co., Ltd. Other wise, BMC Medical Co., Ltd. has no responsibility for the result. All information in this manual is protected by law and regulation. All copyrights belong to BMC Medical Co., Ltd.

IMPORTANT, CAUTION AND WARNING in this manual are to emphasize dangers to service people.

**WARNING:** If do not operate properly, may cause damage to people and environment.

**CAUTION:** If do not operate properly, may cause damage to instrument.

**IMPORTANT:** Important information for servicing and repairing.

## Safety Notice

### Electric Shock

There is above 100VAC voltage inside this device, please be sure to repair device after power off.

### Chemical Safety

There may be risk of virus after touching by patient. Please clean device or wear protective glove before servicing and repairing. Please deal with the waste according to regulation. Wash hands by disinfectant after operating this device.

# 1. Introduction

## 1.1 Intended Use

The RESmart® BPAP system is a BIPAP (Bi-level Positive Airway Pressure) device. It is intended to provide non-invasive ventilation for patients with respiratory insufficiency or obstructive sleep apnoea (OSA), in the hospital or home.

The RESmart® BPAP is to be used only on the instruction of a licensed health care professional. Your home care provider will make the correct pressure settings according to your health care professional's prescription.

This device is not intended for life support.

This service manual is used to help service people maintenance RESmart® BPAP system of BMC Medical Co., Ltd. more efficiently. Instructions in this manual may help device work best. Engineer can find needed instruction from this manual quickly.

### **IMPORTANT!**

Read and understand instruction in this manual before operating.

### **IMPORTANT!**

RESmart® BPAP system can not be used for life support.

## 1.2 Implement

### 1. Implement for repairing

Type '+' screw driver


Pincers


### 2. Implement for measuring

Pressure meter

# 2. Structure

## 2.1 Outer Illumination





## 2.2 Inner Construction

See appendix A

## 2.3 Spare Part

See appendix B


## 3. Software

### 3.1 Instruction for Upgrading Device Software

Note: Once the RESmart® BPAP device is powered on, the LED display screen will show the software version of the device as the following picture. Please remember the software version in the red frame.



1. Power off the RESmart® BPAP device, and then connect PC and the device via data cable.
2. Run the software “BMC RESmart® Tools for Upgrade” from PC and display as below:




- a) Select the upgrade file which has same suffix letter as the machine version. (For example Version 1.88-TN should choose file name 1.98-TN)
  - b) Select the right serial port, default is COM1.
3. Click “OK”, display as below:



Power on the BPAP device. If the device is powered on from beginning, please power off the device and power on again after 5 seconds or more.

4. After the device is turned on, the following will display:



Click “Upgrade” to start. If success, the right version number will be displayed on BPAP device screen when every power on.

**Note:** If failed during upgrading, please repeat from above 1-4 steps.

If you select an incorrect upgrade file (for example: apply a CPAP upgrade file to a BPAP device), the device will display “Error 08” and cannot work. In this case, apply correct software and upgrade again.

## 3.2 Hidden Function

### 1. Humidifier matching

Assemble humidifier on main device, without water chamber. Enter the maintenance menu, and set Ramp=40, then press Heated Humidifier Button when the Init P appears. Thus the humidifier is matched. When ‘OK’ appears on screen, humidifier matching succeeded.

### 2. Reset (Available on version 1.58 and later)

Press and hold the Ramp Button until the user menu appears (about 3 seconds) when the device is on standby, press +/- User Buttons to access the Date Setting, set ‘YYYY/MM/DD’=‘2097/07/01’, then press ‘Pressure Start/Stop Button’ when ‘Hour’ setting appears, thus all device settings and therapy records (doesn’t include Use Days and Use Time) will be reset.

If set ‘Minute’=‘16’ during the above procedure, then not only device settings and therapy records, but also patient information will be reset.

If set ‘Minute’=‘26’ during the above procedure, then not only device settings, therapy records and

patient information, but also Use Days and Use Time will be reset.

### 3. Pressure calibration (Available on version 1.58 and later)


Connect pressure meter to the air outlet of the device, in the Date and Time Setting menu, set 'YYYY/MM/DD'='2098/08/18', press 'Pressure Start/Stop Button' when Hour Setting appears to start the Pressure Calibration. The device output pressure is based on 20 hPa (cmH<sub>2</sub>O), and a number will appear on the screen, when the output pressure is stable, if the value on the pressure meter is not 20 hPa (cmH<sub>2</sub>O), press '+/-' User Buttons to adjust the output pressure, finally press 'Pressure Start/Stop Button' to finish the Pressure Calibration.

## 4. Malfunction and Countermeasure

### 4.1 Error Code Table and Information

Error Index	Problem	Causing
Error1	Over temperature inside or Temperature sensor fault	Malfunction on main board
Error2	Motor stop	Malfunction on motor or main board
Error3	Motor low speed	Malfunction on motor or main board
Error4	Low temperature inside	Malfunction on main board
Error7	Pressure sensor fault	Malfunction on main board
Error8	Wrong software version	Software
Error9	Host parameter error	Malfunction on main board
Error12	Over temperature inside or flowrate sensor fault	Malfunction on main board
Error13	Temperature below 0°C or flowrate sensor fault	Malfunction on main board
Error14	The motor is over temperature or motor temperature sensor fault	Malfunction on motor or Motor temperature sensor
Error15	Motor temperature sensor fault	Malfunction on motor temperature sensor

## 4.2 Malfunction Determination





# Appendices

## Appendix A:

### Device Configuration

#### 1. Main device

Fig. A - Main board (backside)





Fig. B – Outlet (leftside)



Hole #1, 2 and 3 in Fig. B, connect to sensor YL and LL via silicon rubber canal. Specified in Table 1:

Fig. B	Fig. A	Specification of canal
1	Hole #1 on sensor YL	φ 2mm, Φ 4mm, L12cm
2	Hole #2 on sensor LL	φ 2mm, Φ 4mm, L12cm
3	Hole #1 on sensor LL	φ 2mm, Φ 4mm, L8cm

**Table 1**

Fig. C – Power PCBA (front side) Sockets figuration




Fig. D – Wire connecting on power PCBA




Fig. E – Power supply connecting





Fig. F – Power supply to humidifier



Connecting relationship:

Fig. \ Part	1	2	3	4
D	Back of device	Fig. C: CN1	/	
E	Fig. C: CN2	Fig. A: P24	Fig. F: Left	Fig. F: Right
Blower 8-pin plug	Fig. A: PM			

Table 2

## 2. Humidifier

Fig. G – Humidifier PCBA (back side)




Fig. H – Power supply to humidifier





Fig. I – Heater PCBA









**Connecting relationship:**





Fig. \ Part	1	2	3	4
I	Fig. G: P2	Fig. G: P1	Fig. H: Left	Fig. H: Right




**Table 3**

## Components Illumination






### 1. Device Inside

S/N	Assembly	Figure
1	Insert humidifier locker on the shield bottom.	
2	Insert humidifier power supply reed on the shield bottom as per Fig. E and F.	
3	Assembly the humidifier locker and button spring.	
4	Put fixing on the spring and fix by screw (2×3mm*8mm).	
5	Put power wire.	
6	Assembly the shield back on the bottom, fix by screw (2×3mm*8mm).	




<p>7</p>	<p>Assembly the foam pad (non-woven side up), the front and back foam blocks of the blower unit in the shield bottom.</p>	 <p>Front foam block</p> <p>Foam pad</p> <p>Back foam block</p>
<p>8</p>	<p>Align the location column of blower unit to the location port in the front block, Place the blower unit right between the front and back foam blocks in the shield bottom and lay it down till the bottom.</p>	 <p>Location column</p>
<p>9</p>	<p>Clean the air outlet of the blower and insert it into the fixing hole, and set the protruding tube of the blower unit into the air outlet. It is required that the top of blower tube stand up against the bottom of air outlet.</p>	 <p>Blower tube standing up against outlet bottom</p>
<p>10</p>	<p>Fit silicon rubber cap onto the top surface of blower unit.</p>	

<p>11</p>	<p>1. Connect the air hole of the pressure sensor YL on the main board with the bottom hole(hole#1) in the front of the blower outlet with a canal of 12cm in length and 2mm in inner-diameter. 2. Crosswise connect the two air holes of the pressure sensor LL on the main board with the two top holes (hole#2 &amp; hole#3) around the blower outlet via two canals of 8cm in length and 2mm in inner-diameter.</p>	 <p>Canal to sensor LL 8cm in length</p> <p>Canal to sensor LL 8cm in length</p> <p>Canal to sensor YL 12cm in length</p>
<p>12</p>	<p>Clean up the blower wires and wind them with a tube. Insert the plugs into sockets PM1 and PM2. Put the lead of thermistor into socket PRT on main board. Connect the main board to 24V power supply.</p>	 <p>24V power supply</p> <p>Thermistor lead to PRT</p> <p>Blower cable to PM2</p> <p>Blower driving lead to PM1</p> <p>Blower wire and winding tube</p>
<p>13</p>	<p>Draw the shield back backwards and put the 9-pin plug into the corresponding socket on its back. Lay the main board onto the 4 supporting poles on the shield bottom and set it into the power board via the fixing pole on PCB and get the main board fixed. Put blower wires on the side into the gap between blower unit and shield bottom.</p>	


2. Device Outside







S/N	Assembly	Figure
14	Put handle in the socket on shield cover.	
15	Assembly handle fixing to handle and fix by screw (4×3mm*8mm).	
16	Assembly and fix the panel.	
17	Put on the key button.	
18	Assembly the infrared window, and put the shield cover on the bottom, adjust the position of outlet and make sure it is clamped in the right socket.	



19	Fix the shield by screws (4×M3*10mm).	 <p>Screw hole</p>
20	Stick four device pads.	 <p>Device pad</p>
21	Put in filter and assembly the cover.	
22	Put on all plugs if necessary (without humidifier).	 <p>Plug (up)</p> <p>Power plug</p> <p>Plug (down)</p>



### 3. Humidifier




S/N	Assembly	Figure
1	As per Fig. H, assembly the heater plane and wire. Fix them by screw (3mm*8mm).	 <p>Black</p> <p>Orange</p>

2	Put spring on each of three poles on heater plane.	
3	Assembly the scaleboard on the heater plane and fix by screw (3mm*8mm).	
4	As per Table 3, assembly heater PCBA and wires.	
5	Put on light window.	
6	Assembly the inner and outer shield.	
7	Fix on screw (M3*10mm).	

## Appendix B:

### Spare Part Configuration

S/N	Name	Qty	Figure
1	Shield Cover	1	
2	Handle	1	
3	Handle Fixing	2	

4	Shield Bottom	1	
5	Humidifier Clip	1	
6	Humidifier Locker	1	
7	Humidifier Fixing	1	
8	Filter Cover	1	
9	Shield Back	1	
10	Outlet	1	
11	PCBA Fixing	2	
12	Infrared Window	2	
13	Panel	1	
14	Humidifier Outer Shield	1	

15	Humidifier Inner Shield	1	
16	Light Window	1	
17	Pole Platelet	1	
18	Humidifier Power Plug	1	
19	Humidifier Plug (left-up)	1	
20	Humidifier Plug (left-down)	1	
21	Humidifier Plug (right-up)	1	
22	Humidifier Plug (right-down)	1	
23	Device Pad	6	
24	Humidifier Inner Connector	1	
25	Key Button	1	
26	Humidifier Scaleboard	1	
27	Heater Plane Spring	3	
28	Button Spring	1	

## Spare Part List

Part Name	S/N	Qty	Unit	Classify	Use
Main Board PCBA	140022	1	PC	PCBA	Main Device
Heater PCBA	140017	1	PC	PCBA	Humidifier

Power Supply PCBA	270031	1	PC	PCBA	Main Device
Blower Unit	120069	1	PC	Assembly	Main Device
Heater Plane	120006	1	Set	Assembly	Humidifier
Power Supply PCBA Wire	120002	1	Set	Wire	Main Device
Main Board PCBA Wire	120003	1	Set	Wire	Main Device
Silicon Rubber Cap	220029	1	PC	Assembly	Main Device
Silicon Rubber Canal	290001	3	PC	Assembly	Main Device
Foam (front/back)	240048	1	PC	Foam	Main Device
Shield Cover	210067	1	PC	Plastic	Main Device
Handle	210068	1	PC	Plastic	Main Device
Handle Fixing	210011	2	PC	Plastic	Main Device
Shield Bottom	210069	1	PC	Plastic	Main Device
Humidifier Clip	210070	1	PC	Plastic	Main Device
Humidifier Locker	210012	1	PC	Plastic	Main Device
Humidifier Fixing	210013	1	PC	Plastic	Main Device
Filter Cover	210071	1	PC	Plastic	Main Device
Shield Back	210072	1	PC	Plastic	Main Device
Outlet	210080	1	Set	Plastic	Main Device
PCBA Fixing	210014	2	PC	Plastic	Main Device
Infrared Window	210010	2	PC	Plastic	Main Device
Panel	210073	1	PC	Plastic	Main Device
Humidifier Outer Shield	210074	1	PC	Plastic	Humidifier
Humidifier Inner Shield	210075	1	PC	Plastic	Humidifier
Pole Platelet	210025	1	PC	Plastic	Humidifier
Light Window	210023	1	PC	Plastic	Humidifier
Button Spring	230001	1	PC	Metal	Main Device
Humidifier Scaleboard	230004	1	PC	Metal	Humidifier
Heater Plane Spring	230005	3	PC	Metal	Humidifier
Humidifier Power Plug	220023	1	PC	Rubber	Main Device
Humidifier Plug (left-up)	220024	1	PC	Rubber	Main Device
Humidifier Plug (left-down)	220025	1	PC	Rubber	Main Device
Humidifier Plug (right-up)	220026	1	PC	Rubber	Main Device
Humidifier Plug (right-down)	220027	1	PC	Rubber	Main Device
Device Pad	220001	6	PC	Rubber	All
Key Button	220002	1	PC	Rubber	Main Device
Humidifier Inner Connector	220009	1	PC	Rubber	Humidifier
Screw-M3*6	240008	4	PC	Screw	All
Screw-M3*10	240009	4	PC	Screw	All
Screw-M3*8	240010	8	PC	Screw	Main Device
Filter	240014	1	PC	Accessory	All
Power Cord	110099	1	PC	Accessory	All